CERTEST H. pylori+ TRANSFERRIN

ONE STEP H. pylori +Transferrin
COMBO CARD TEST

CERTEST BIOTEC S.L.

CERTEST *H. pylori*+Transferrin

One Step *H. pylori* and Transferrin combo card test /Prueba combo de un solo paso para *H. pylori* y transferrina en cassette

ENGLISH

INTENDED USE

The CerTest *H. pylori*+Transferrin one step combo card is a colored chromatographic immunoassay for the qualitative detection of *H. pylori* and transferrin in stool samples.

CerTest *H. pylori*+Transferrin combo card offers a highly sensitive, simple and non-invasive assay to determine, on the one hand, the presence of *H. pylori* that causes gastric and duodenal ulcers and on the other, the presence of gastrointestinal bleeding and moreover, to monitor patients treatment (*H. pylori* eradication).

INTRODUCTION

Helicobacter pylori (H. pylori) is a spiral-shaped bacterium that is found in the gastric mucous layer or adherent to the epithelial lining of the stomach. *H. pylori* causes more than 90% of duodenal ulcers and up to 80% of gastric ulcers. The importance of *Helicobacter pylori* testing has increased due to the strong correlation between the presence of bacteria and confirmed gastrointestinal diseases (stomach and duodenum) like gastritis, peptic ulcer disease and gastric carcinoma. Patients suffering from gastro-duodenal ulcer may develop complications such as upper gastrointestinal bleeding or perforation. A precise diagnosis of *H. pylori* after eradication reduces significantly the risk of ulcer reactivation, and thus prevents new episodes of upper gastrointestinal bleeding.

Transferrin is a blood-derived component. In gastrointestinal bleeding diseases, it may be leaked into gastrointestinal tract and then discharged with the feces. Transferrin is stable in feces and a good marker to detect gastrointestinal bleeding.

This immunochromatographic assay detects simultaneously *H. pylori* and human transferrin in stool samples, obtaining more accurate testing results regarding the *H. pylori* infection.

TEST PRINCIPLE

The **CerTest *H. pylori*+Transferrin** combo card is based on the principle of a qualitative immunochromatographic assay for the determination of *H. pylori* and transferrin in stool samples. Each strip consists of a membrane pre-coated with antibodies on the test band (result region), against antigens, as follows:

Strip **A** consists of a membrane pre-coated with mouse monoclonal antibodies on the result region against *H. pylori* and colored conjugate (anti-*H. pylori* antibodies-red microspheres) pad pre-dried on the end of the membrane.

Strip **B** consists of a membrane pre-coated with mouse monoclonal antibodies on the test band (test line) against transferrin and colored conjugate (anti-transferrin antibodies-red microspheres) pad pre-dried on the end of this membrane.

During testing in each strip, the sample is allowed to react with the conjugate and the mixture then moves upward on the membrane, by capillary action. As the sample flows through the test membrane, the colored particles migrate. In the case of a positive result, the specific antibodies present on the membrane will capture the colored conjugate. A colored band will be visible, A and/or B depends on the antigen sample content. These bands are used to interpret the result.

The mixture continues to move across each membrane to the immobilized antibody placed in the control band region; this green colored band always appears. The presence of this green band serves as 1) verification that sufficient volume is added, 2) proper flow obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The assay should be performed immediately after collection. The samples can be stored in the refrigerator (2-8 °C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing.


Please take note that there are two different stool collection tubes:

1. “*H. pylori*” stool collection tubes (without reducer piece)

Specimen preparation (see illustration):


1. Take out the cap of the stool collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick 4 times, collect large-quantity fecal sample (approx. 250mg) and add to the stool collection tube.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion.

2. “Transferrin” stool collection tube (with a reducer piece):

Specimen preparation (see illustration):

- (1) Unscrew the cap and use the stick by introducing the fecal specimen four times to pick up the sample (approx. 15mg). Add the sample into the stool collection tube.
- (2) Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion.


Precautions: Patients suffering from menstrual period, bleeding hemorrhoids, blood in urine or strain during bowel movement should not collect samples.

MATERIALS

MATERIALS PROVIDED

- CerTest *H. pylori*+Transferrin combo card
- Instructions for use
- “*H. pylori*” stool collection tubes with diluent
- “Transferrin” stool collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE


Allow the tests, stool samples and controls to reach to room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay. Please take note that there are two procedures by using different stool collection tubes:

1. Proceed to shake each stool collection tube in order to assure good sample dispersion.
2. Remove the CerTest *H. pylori*+Transferrin combo device from its sealed bag just before using.
3. Take the *H. pylori* stool collection tube, cut the end of the cap (3) and dispense 4 drops or 100µL in the circular window marked with the letter A (4). Avoid adding solid particles with the liquid.
4. Take the Transferrin stool collection tube, cut the end of the cap (5) and dispense 4 drops or 100µL in the circular window marked with the letter B (6).
5. Read the results at **10 minutes** (the colored bands appear).


CERTEST H. pylori+ Transferrin

One Step H. pylori and Transferrin combo card test /Prueba combo de un solo paso para H. pylori y transferrina en cassette


In case the tests did not run due to solid particles fallen into the round window, stir the sample added or dispense a drop of extraction buffer until seeing the liquid running through the reaction zone.

a. H. pylori strip-A procedure


b. Transferrin strip-B procedure


INTERPRETATION OF RESULTS (please refer to the illustration below)

1. NEGATIVE


A: Green → Negative H. pylori
 B: Green → Negative Transferrin

2. POSITIVE/NEGATIVE


A: Green-Red → Positive H. pylori
 B: Green → Negative Transferrin

3. POSITIVE/POSITIVE


A: Green-Red → Positive H. pylori
 B: Green-Red → Positive Transferrin

4. NEGATIVE/POSITIVE

A: Green → Negative H. pylori
 B: Green-Red → Positive Transferrin

INVALID

	A (H. pylori)	B (Transferrin)	Interpretation of the results
1.	- GREEN	- GREEN	No infection caused by <i>H. pylori</i> , no gastrointestinal bleeding problem.
2.	+ GREEN-RED	- GREEN	<i>H. pylori</i> infection without gastrointestinal bleeding.
3.	+ GREEN-RED	+ GREEN-RED	<i>H. pylori</i> infection with gastrointestinal bleeding presumably caused by <i>H. pylori</i> infection.
4.	- GREEN	+ GREEN-RED	No <i>H. pylori</i> infection but gastrointestinal bleeding problem.
5.	Any other result	Any other result	Invalid result: Invalid result either A or B, we recommend repeating the assay using the same sample with another combo device.

CERTEST *H. pylori*+Transferrin

One Step *H. pylori* and Transferrin combo card test /Prueba combo de un solo paso para *H. pylori* y transferrina en cassette

INVALID: A total absence of the control colored band regardless the appearance or not of the result line. Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test. If the problem still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red colored band in the result line region (T) will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value, nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. The GREEN line appearing in the control region (C) is an internal control. It confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours of opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
3. Some stool samples can decrease the intensity of the control line.
4. This test provides a presumptive diagnosis of *Helicobacter pylori* infections and a possible gastrointestinal bleeding problem. A confirmed infection or bleeding diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated must be based in the correlation of the results with further clinical observations.
5. Patients suffering from menstrual period, bleeding hemorrhoids, blood in urine or strain during bowel movement should not collect samples to detect transferrin (false positive result).
6. Positive results confirm the presence of human transferrin in fecal samples; nevertheless, it can be also due to several causes besides colorectal bleeding, such as hemorrhoids, blood in urine or stomach irritations. A positive result should be followed up with additional diagnostic procedures to determine the exact cause and source of the blood in the stool.
7. Transferrin negative results do not exclude bleeding, as some polyps and colorectal cancers may bleed intermittently or not during certain stages of the disease. Moreover, blood may not be uniformly distributed in stool samples.

EXPECTED VALUES

Peptic ulcer is mostly the common cause of acute intestinal hemorrhage. *H. pylori* is a pathogenetic agent in chronic active gastritis and is strongly associated with ulcer disease.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit)

A sample containing human transferrin at concentration equal to or higher than 0.4µg/g of feces produces positive results when using CerTest *H. pylori*+Transferrin test.

Different transferrin dilutions were tested directly in the extraction buffer or spiked in a negative stool sample in accordance with the kit instructions to determinate the detection limit of the test. Detection limit values: 4ng/mL of human transferrin.

Clinical Sensitivity and Specificity

An evaluation was performed comparing the immunochromatographic test to EIA assay to detect *H. pylori* infection. The samples were directly taken from patients presumably suffering from *H. pylori* infection. The results were as follows:

	EIA Test	
	Sensitivity	Specificity
IC (immuno-chromatographic assay)	>94%	>99%

Cross reaction

The CerTest *H. pylori*+Transferrin (transferrin strip) test is specific for transferrin, showing no cross-reaction with transferrin from bovine and pig.

The use of a mouse monoclonal antibody in CerTest *H. pylori*+Transferrin assures high degree of specificity for the detection of *H. pylori* and transferrin as marker. The antibodies used to elaborate this combo test recognize epitopes present in the antigen found in the patients stool, as well as the preparations from the bacteria (*H. pylori*) cultures *in vitro*.

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