CERTEST

Calprotectin

ONE STEP human calprotectin
CARD TEST

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ENGLISH

INTENDED USE

CerTest Calprotectin one step card test is a coloured chromatographic immunoassay (cut-off qualitative test) for the qualitative determination of human calprotectin (hCp) in stool samples that may reflect gastrointestinal inflammation caused by several pathologies (inflammatory bowel disease, colorectal cancer and some enteropathies).

CerTest Calprotectin card test offers a simple, a highly sensitivity and non-invasive screening assay to determine intestinal inflammatory activity, to monitor treatment response and to predict risk of relapse.

INTRODUCTION

Calprotectin is a neutrophil cytosolic protein with antimicrobial properties, which is present at increased concentration in stool samples during bowel inflammation. The stability of the protein to degradation keeps it stable in faeces for up to 7 days at room temperature, making it an ideal analyte. Calprotectin is released by activation of leukocytes, giving increased levels in plasma, cerebral spinal fluid, synovial fluid, urine or stools as a consequence of disease in the relevant organ(s). Calprotectin inhibits zinc-dependent enzyme systems, as a result kills microbes and induces apoptosis in normal and cancer cells. In the presence of calcium, calprotectin is remarkably resistant to proteolytic degradation and so is stable in stools kept at room temperature for 7 days.

This is a non-invasive marker of intestinal inflammation (for example in Ulcerative Colitis (UC) and Crohn's Disease (CD)).

PRINCIPLE OF THE TEST

CerTest Calprotectin is based on the principle of a cut-off qualitative immunochromatographic assay for the determination of human calprotectin in stool samples.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against human calprotectin and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-human calprotectin) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-human calprotectin monoclonal antibodies-red polystyrene microspheres), which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-human calprotectin antibodies present on the membrane (test line) will capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no human calprotectin presence and yet, the antigens may be present in a concentration lower than the cut-off value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-human calprotectin antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red test line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appears. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until through the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contaminations errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contains preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the CerTest Calprotectin test. Do not use any other commercial kit component.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.


SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. For longer storage, maximum 6 months, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

1. Take out the cap 1 of the stool collection tube and use the stick to pick up sufficient sample quantity. Introduce the stick once into 4 different parts of the stool sample (1); making sure that, in each insertion, the stick's screw is covered with the sample to take the appropriate amount of stool. Then, introduce the sample in the stool collection tube (2). For liquid samples, add approx. 20µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample (2). Shake the tube in order to assure good sample dispersion (3). The stool collection tube with dilute sample can be stored for 7 days in the refrigerator (2-8°C) prior to testing.

MATERIALS

MATERIALS PROVIDED

- CerTest Calprotectin card tests
- Instructions for use
- Stool collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE


Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove the CerTest Calprotectin card test from its sealed bag just before using it.
3. Take the stool collection tube, open the cap 2 (4) and dispense 3 drops in the circular window marked with the letter S (5). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF THE RESULTS (please refer to the illustration below)

	Calprotectin (hCp)	Interpretation of the results
1.	- GREEN	There is no human calprotectin presence which might mean neither active gastrointestinal inflammation, nor risk of relapse (CD or UC relapse).
2.	+ GREEN-RED	There is human calprotectin present which might mean active gastrointestinal inflammation or risk of relapse in clinical remission.
3.	Any other result	Invalid result, we recommend repeating the assay using the same sample with another test.

NEGATIVE: Only the GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the result window will vary depending on the concentration of human calprotectin present in the specimen.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

- The test must be carried out within 2 hours after opening the sealed bag.
- An excess of stool sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
- The intensity of test line may vary from very strong at high human calprotectin concentration to faint when the human calprotectin concentration is close to the cut-off value of the test.
- CerTest Calprotectin should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
- Positive results determine the presence of human calprotectin in faecal samples; nevertheless, it can be due to several causes (Inflammatory bowel disease, colorectal cancer and some other enteropathies). A positive result should be followed up with additional diagnostic invasive procedures, a colonoscopy and a biopsy in order to confirm the diagnosis and to establish the inflammation extent.
- A negative result is not meaningful because of it is possible the human calprotectin concentration in the stool sample is lower than the cut-off value. If symptoms or situation still persist, calprotectin determination should be carried out invasive techniques. Negative results do not exclude inflammation, some diseases such as celiac sprue and microscopic colitis polyps that mainly involve mononuclear inflammation.
- Stool samples from patients with non-steroidal anti-inflammatory drug (NSAID) treatment could show positive result.
- Neonatal faecal calprotectin levels have been reported higher than those in normal children with a mean of 167µg/g (range 22-860µg/g).
- Stool from patients with active inflammatory bowel disease that usually involve significant neutrophilic inflammation of the intestine, such as Crohn's disease and ulcerative colitis, would be positive for faecal calprotectin (the level of calprotectin in Crohn's disease is higher than in ulcerative colitis) and would be positive for faecal lactoferrin. CerTest Calprotectin could be sensitive for this diagnosis in patients with chronic diarrhea.

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10. Bloody stool samples can contain components that may cause non-specific reactions in the test. Every bloody stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
11. Mucous stool samples may cause non-specific reactions in the test. Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.

EXPECTED VALUES

Some studies established equal or higher 50µg hCp/g faeces as cut-off value to allow detect adult patients with gastrointestinal inflammatory pathology that will require to diagnosis additional diagnostic invasive procedures.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

Cut-off value of CerTest Calprotectin is 500 ng/mL (50 µg hCp/g faeces) for human calprotectin.

Clinical sensitivity and specificity

An evaluation was performed comparing an immunochromatographic assay (CerTest Calprotectin, CerTest) and another commercial immunoassay (Calprest®, Eurospital). The results were as follows:

Data for cut-off value 50 µg hCp/g

		ELISA test: Calprest®		
		+	-	Total
IC test: CerTest Calprotectin	+	31	7	38
	-	1	41	42
	Total	32	48	80

CerTest Calprotectin vs Calprest®		
	Mean Value	95% CI (confidence interval)
Sensitivity	96.9%	83.8 – 99.9%
Specificity	85.4%	72.2 – 93.9%
PPV	81.6%	65.7 - 92.3%
NPV	97.6%	87.4 –99.9%

Sensitivity for a cut-off value of 100 µg hCp/g: 100.0% (mean value); 85.2-100.0% (95% confidence interval).

Specificity for a cut-off value of 25 µg hCp/g: 94.9% (mean value); 82.7-99.4% (95% confidence interval).

The results showed a high sensitivity and specificity to detect human calprotectin using CerTest Calprotectin.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest Calprotectin; no cross reactivity against other gastrointestinal pathogens, other organism, substances and/or faecal markers occasionally present in faeces:

<i>Adenovirus</i>	<i>Clostridium perfringens</i>	<i>Human Lactoferrin</i>	<i>Salmonella typhi</i>
<i>Astrovirus</i>	<i>Cryptosporidium</i>	<i>Human Transferrin</i>	<i>Salmonella typhimurium</i>
<i>Bovine Haemoglobin</i>	<i>Entamoeba dispar</i>	<i>Legionella pneumophila</i>	<i>Shigella boydii</i>
<i>Bovine Lactoferrin</i>	<i>Entamoeba histolytica</i>	<i>Listeria monocytogenes</i>	<i>Shigella dysenteriae</i>
<i>Bovine Transferrin</i>	<i>Escherichia coli O111</i>	<i>Norovirus GI</i>	<i>Shigella flexneri</i>
<i>Campylobacter coli</i>	<i>Escherichia coli O26</i>	<i>Norovirus GII</i>	<i>Shigella sonnei</i>
<i>Campylobacter jejuni</i>	<i>Escherichia coli O157</i>	<i>Pig Haemoglobin</i>	<i>Streptococcus pneumococcal</i>
<i>Clostridium difficile antigen GDH</i>	<i>Giardia</i>	<i>Rotavirus</i>	<i>Streptococcus pyogenes</i>
<i>Clostridium difficile Toxin A</i>	<i>Helicobacter pylori</i>	<i>Salmonella enteritidis</i>	<i>Yersinia Enterocolitica O:3</i>
<i>Clostridium difficile Toxin B</i>	<i>Human Haemoglobin</i>	<i>Salmonella paratyphi A</i>	<i>Yersinia Enterocolitica O:9</i>

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1. Angriman I. et al. Enzymes in feces: Useful markers of chronic inflammatory bowel disease. *Clinica Chimica Acta* 381 Feb 2007, p. 63-68.
2. Quail, M.A. et al. Fecal Calprotectin Complements Routine Laboratory Investigations in Diagnosing Childhood Inflammatory Bowel Disease. *Inflamm Bowel Dis*, Vol 15 No 5; May 2009, p. 756-759.
3. Gaya D.R., et al. Faecal calprotectin in the assessment of Crohn's disease activity. *Q J Med* 2005, Vol 98, May 2005, p. 435-441.
4. Langhorst, M.D. et al. Noninvasive Markers in the Assessment of Intestinal Inflammation in Inflammatory Bowel Diseases: Performance of Fecal Lactoferrin, Calprotectin and PMN-Elastase, CRP, and Clinical Indices. *Am. J. Gastroenterol.* 2008; Vol 103, p. 162-169.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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