---

# CERTEST FOB

---

ONE STEP human haemoglobin  
CARD TEST

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### ENGLISH

#### INTENDED USE

CerTest FOB one step card test is a coloured chromatographic immunoassay (cut-off qualitative test) for the qualitative determination of human haemoglobin (hHb) in stool samples.

CerTest FOB card test offers a simple, a highly sensitive and non-invasive screening assay to make a presumptive diagnosis of gastrointestinal bleeding.

#### INTRODUCTION

Colorectal cancer is the second leading cause of illness and death in Western world. The screening with faecal occult blood tests is based on the concept that important target colonic neoplasm, such as early-stage cancer and large adenomatous polyps, will bleed, for which may be detected by an occult blood test. Colorectal cancer is also associated with local acute inflammatory reaction being visualized, in some cases, by white cell neutrophil scanning.

Haemoglobin is the iron-containing oxygen-transport protein in the red blood cells of all vertebrates that may be leaked into gastrointestinal tract and then discharged with the feces in gastrointestinal bleeding diseases.

When gastrointestinal blood is lost, the stool will contain a combination of intact or nearly intact haemoglobin, intact heme and heme-derived porphyrins in amounts that depend on the site and amount of bleeding and the transit time through the gut. Immunochemical tests detect intact or nearly intact human haemoglobin, being a very specific technique for detecting loss of blood from the lower intestine, because blood from lower sites is less degraded during transit.

#### TEST PRINCIPLE

CerTest FOB is based on the principle of a cut-off qualitative immunochromatographic assay for the determination of human haemoglobin in stool samples.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against human haemoglobin and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-human haemoglobin) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-human haemoglobin monoclonal antibodies-red polystyrene microspheres), which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-human haemoglobin antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no human haemoglobin presence and yet, the antigens may be present in a concentration lower than the cut-off value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-human haemoglobin antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appears. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

#### STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

#### PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest FOB test. Do not use any other commercial kit component.


- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

## SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. For longer storage, maximum 6 months, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.

### Specimen preparation (see illustration):

1. Take out the cap 1 of the stool collection tube and use the stick to pick up sufficient sample quantity. Introduce the stick once into 4 different parts of the stool sample (1), making sure that, in each insertion, the stick's screw is covered with the sample to take the appropriate amount of stool. Then, introduce the sample in the stool collection tube (2). For liquid samples, add approx. 20µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample (2). Shake the tube in order to assure good sample dispersion (3). The stool collection tube with dilute sample can be stored until 2 days at room temperature or until 5 days in the refrigerator (2-8°C) prior to testing.

**Precautions:** Patients should not collect samples during their menstrual period, if they have bleeding haemorrhoids, blood in urine or if they have strained during bowel movement.

## MATERIALS

### MATERIALS PROVIDED

- CerTest FOB card tests
- Instructions for use
- Stool collection tubes with diluent

### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

## TEST PROCEDURE


**Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.**

1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove the CerTest FOB card test from its sealed bag just before using it.
3. Take the stool collection tube, open the cap 2 (4) and dispense 3 drops in the circular window marked with the letter S (5). Avoid adding solid particles with the liquid.
4. **Read the result at 10 minutes.** Do not read the test result later than 10 minutes.


If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

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### INTERPRETATION OF RESULTS (please refer to the illustration below)

|    | FOB (hHb)        | Interpretation of the results                                                                                                                                                 |
|----|------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1. | -<br>GREEN       | There is no human haemoglobin presence which might mean no faecal occult blood and no gastrointestinal bleeding.                                                              |
| 2. | +<br>GREEN-RED   | There is human haemoglobin presence which might mean gastrointestinal bleeding caused by ulcers, a colonic neoplasm, such as early-stage cancer and large adenomatous polyps. |
| 3. | Any other result | Invalid result, we recommend repeating the assay using the same sample with another test.                                                                                     |

**NEGATIVE:** Only the GREEN line appears across the results window in the site marked with the letter C (control line).

**POSITIVE:** In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

**INVALID:** Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

### NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of human haemoglobin present in the specimen.

### QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

### LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. An excess of stool sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
3. The intensity of test line may vary from very strong at high human haemoglobin concentration to faint when the human haemoglobin concentration is close to the cut-off value of the test.

4. CerTest FOB should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
5. Positive results determine the presence of human haemoglobin in faecal samples; nevertheless, it can be due to several causes, besides colorectal bleeding, such as haemorrhoids, blood in urine or stomach irritations. A positive result should be followed up with additional diagnostic procedures to determine the exact cause and source of the blood in the stool.
6. A negative result is not meaningful because of it is possible the human haemoglobin concentration in the stool sample is lower than the cut-off value. Negative results do not exclude bleeding, as some polyps and colorectal cancers may bleed intermittently or not during certain stages of the disease. Moreover, blood may not be uniformly distributed in stool samples.
7. This test may be less sensitive for detecting upper gastrointestinal bleeding because blood degrades as it passes through the gastrointestinal tract.
8. Patients should not collect samples during their menstrual period, if they have bleeding haemorrhoids, blood in urine or if they have strained during bowel movement.
9. Mucous stool samples may cause non-specific reactions in the test. Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
10. Faecal samples whose antigen concentration is near to the cut-off value of the CerTest FOB test could be lost if the sample has been diluted and preserved until 2 days at room temperature or until 5 days at 2-8°C until to be tested. For faecal samples near to the detection limit value is better to dilute the sample and make immediately the test.

## EXPECTED VALUES

Colorectal cancer affects both men and women of all racial and ethnic groups, being most often found in people aged 50 years or older. Regarding men, colorectal cancer is the third most common cancer after prostate and lung cancers. Regarding women, colorectal cancer is the third most common cancer after breast and lung cancers.

Some studies established equal or higher 5.1µg hHb /g faeces as cut-off values to allow detect adult patients with gastrointestinal bleeding that will require to diagnosis additional diagnostic invasive procedures.

## PERFORMANCE CHARACTERISTICS

### Analytical sensitivity

Cut-off value of CerTest FOB is 50ng/mL (5.1µg hHb/g faeces) for human haemoglobin.

### Clinical sensitivity and specificity

An evaluation with faecal samples was performed comparing the results obtained by an immunochromatographic test (CerTest FOB, CerTest) and an ELISA test (Human Hemoglobin ELISA kit, ICL).

The results were as follows:

#### Data for cut-off value 5 µg hHb/g

| IC test: CerTest FOB |   | ELISA test: Human Hemoglobin ELISA kit |    |       |
|----------------------|---|----------------------------------------|----|-------|
|                      |   | +                                      | -  | Total |
|                      | + | 38                                     | 11 | 49    |
|                      | - | 0                                      | 49 | 49    |
| Total                |   | 38                                     | 60 | 98    |

| CerTest FOB vs Human Hemoglobin ELISA kit |            |                         |
|-------------------------------------------|------------|-------------------------|
|                                           | Mean Value | 95% confidence interval |
| <b>Sensitivity</b>                        | 100.0%     | 90.7– 100.0%            |
| <b>Specificity</b>                        | 81.7%      | 69.6 – 90.5%            |
| <b>PPV</b>                                | 77.6%      | 63.4 - 88.2%            |
| <b>NPV</b>                                | 100.0%     | 92.7 – 100.0%           |

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The results showed a high sensitivity and specificity to detect human haemoglobin using CerTest FOB.

### Cross-Reactivity

An evaluation was performed to determine the cross reactivity of CerTest FOB; no cross reactivity against other gastrointestinal pathogens, other organism, substances and/or faecal markers occasionally present in faeces:

|                                          |                                |                               |                                    |
|------------------------------------------|--------------------------------|-------------------------------|------------------------------------|
| <i>Adenovirus</i>                        | <i>Clostridium perfringens</i> | Human lactoferrin             | <i>Salmonella typhi</i>            |
| <i>Astrovirus</i>                        | <i>Cryptosporidium</i>         | Human transferrin             | <i>Salmonella typhimurium</i>      |
| Bovine haemoglobin                       | <i>Entamoeba dispar</i>        | <i>Legionella pneumophila</i> | <i>Shigella boydii</i>             |
| Bovine transferrin                       | <i>Entamoeba histolytica</i>   | <i>Listeria monocytogenes</i> | <i>Shigella dysenteriae</i>        |
| Bovine lactoferrin                       | <i>Escherichia coli</i> O111   | <i>Norovirus GI</i>           | <i>Shigella flexneri</i>           |
| <i>Campylobacter coli</i>                | <i>Escherichia coli</i> O26    | <i>Norovirus GII</i>          | <i>Shigella sonnei</i>             |
| <i>Campylobacter jejuni</i>              | <i>Escherichia coli</i> O157   | Pig haemoglobin               | <i>Streptococcus pneumococcal</i>  |
| <i>Clostridium difficile</i> antigen GDH | <i>Giardia</i>                 | <i>Rotavirus</i>              | <i>Streptococcus pyogenes</i>      |
| <i>Clostridium difficile</i> Toxin A     | <i>Helicobacter pylori</i>     | <i>Salmonella enteritidis</i> | <i>Yersinia enterocolitica</i> O:3 |
| <i>Clostridium difficile</i> Toxin B     | Human calprotectin             | <i>Salmonella paratyphi</i> A | <i>Yersinia enterocolitica</i> O:9 |

No special diet is recommended prior to testing. There are not interferences with any foods (Vitamin C, broccoli, carrots...).

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1. Towler BP, Irwig L, Glasziou P, Weller D, Kewenter J. Screening for colorectal cancer using the faecal occult blood test, Hemoccult. Cochrane Database Syst Rev. 2000;(2): CD001216.
2. Ransohoff DF and Lang CA. Screening for colorectal cancer with the Fecal Occult Blood Test: a background paper. Ann Intern Med. 1997; 126: 811-822.
3. Ransohoff DF and Lang CA. Suggested technique for Faecal Occult Blood testing and interpretation in colorectal cancer screening. Ann Intern Med. 1997; 126: 808-810.

## SYMBOLS FOR IVD COMPONENTS AND REAGENTS/ SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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