CERTEST Lactoferrin


ONE STEP human lactoferrin
CARD TEST

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ENGLISH

INTENDED USE

CerTest Lactoferrin one step card test is a coloured chromatographic immunoassay (cut-off qualitative test) for the quantitative determination of human lactoferrin (hLf) in stool samples that may reflect intestinal inflammation caused by several pathologies (inflammatory bowel disease, colorectal cancer and some enteropathies).

CerTest Lactoferrin card test offers a simple, a highly sensitivity and non-invasive screening assay to determine intestinal inflammatory activity, to monitor treatment response in these patients and to predict risk of relapse.

INTRODUCTION

Lactoferrin is a glycoprotein component of neutrophil secondary granules, a primary component of the acute inflammatory response and is released from fecal leukocytes. This protein is resistant to proteolysis in the faeces and may serve as a marker of inflammation in the intestine. The major cause of fecal neutrophils in patients with chronic diarrhoea is chronic inflammatory bowel disease of the colon (i.e., Crohn's Disease and ulcerative colitis). Lactoferrin has been also studied as a predictor of infection with invasive enteropathogens in children with diarrhoea. Bacterial inflammatory diarrhea may be caused by *Shigella*, *Salmonella*, *Campylobacter* and *Clostridium difficile*.

TEST PRINCIPLE

CerTest Lactoferrin is based on the principle of a cut-off qualitative immunochromatographic assay for the determination of human lactoferrin in stool samples.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against human lactoferrin and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-human lactoferrin) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-human lactoferrin monoclonal antibodies-red polystyrene microspheres), which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-human lactoferrin antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no human lactoferrin presence and yet, the antigens may be present in a concentration lower than the cut-off value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-human lactoferrin antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red test line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appear. The presence of this green line serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest Lactoferrin test. Do not use any other commercial kit component.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.


SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. For longer storage, maximum 6 months, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

1. Take out the cap 1 of the stool collection tube and use the stick to pick up sufficient sample quantity. Introduce the stick once into 4 different parts of the stool sample (1), making sure that, in each insertion, the stick's screw is covered with the sample to take the appropriate amount of stool. Then, introduce the sample in the stool collection tube (2). For liquid samples, add approx. 20µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample (2). Shake the tube in order to assure good sample dispersion (3). The stool collection tube with dilute sample can be stored for 7 days in the refrigerator (2-8°C) prior to testing.

MATERIALS

MATERIALS PROVIDED

- CerTest Lactoferrin card tests
- Instructions for use
- Stool collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE


Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove CerTest Lactoferrin card test from its sealed bag just before using it.
3. Take the stool collection tube, open the cap 2 (4) and dispense 3 drops in the circular window marked with the letter S (5). Avoid adding solid particles with the liquid.
4. **Read the result at 10 minutes.** Do not read the test result later than 10 minutes.


If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

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INTERPRETATION OF THE RESULTS (please refer to the illustration below)

	Lactoferrin (hLf)	Interpretation of the results
1.	- GREEN	There is no human lactoferrin presence which might mean neither active gastrointestinal inflammation.
2.	+ GREEN-RED	There is human lactoferrin presence which might mean an inflammatory intestinal disease.
3.	Any other result	Invalid result, we recommend repeating the assay using the same sample with another test.

NEGATIVE: Only a GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of human lactoferrin present in the specimen.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. An excess of stool sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
3. The intensity of test line may vary from very strong at high human lactoferrin concentration to faint when the human lactoferrin concentration is close to the cut-off value of the test.
4. CerTest Lactoferrin should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.

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5. Positive results determine the presence of human lactoferrin in faecal samples; nevertheless, it can be also due to several causes besides inflammatory bowel disease (IBD). A positive result should be followed up with additional diagnostic procedures. Endoscopy and histology on biopsy specimens are the methods for detecting and quantifying bowel inflammation.
6. A negative result is not meaningful because of it is possible the human lactoferrin concentration in the stool sample is lower than the cut-off value. Negative results do not exclude inflammation, some diseases such as celiac sprue and microscopic colitis polyps that involve mainly monocuclear inflammation.
7. Lactoferrin is a component of breast milk; the test will be positive in breast fed children and should not be used to evaluate neonates receiving breast milk.
8. Stool from patients with active inflammatory bowel diseases that usually involve significant neutrophilic inflammation of the intestine, such as Crohn's disease and ulcerative colitis, would be positive for faecal lactoferrin. CerTest Lactoferrin could be sensitive for this diagnosis in patients with chronic diarrhea.
9. Bloody stool samples can contain components that may cause non-specific reactions in the test. Every bloody stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
10. Mucous stool samples may cause non-specific reactions in the test. Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.

EXPECTED VALUES

Some studies established equal or higher 10µg hLf/g faeces as cut-off value to allow detect adult patients with gastrointestinal inflammatory pathology that will require to diagnosis additional diagnostic invasive procedures.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

Cut-off value of CerTest Lactoferrin is 100ng/mL (10µg hLf/g faeces) for human lactoferrin.

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest Lactoferrin, CerTest) and other commercial ELISA tests (Human Lactoferrin ELISA kit, ICL). The results were as follows:

Data for cut-off value 10 µg hLf/g

		ELISA test: Human Lactoferrin ELISA kit		
		+	-	Total
IC test: CerTest Lactoferrin	+	27	13	40
	-	3	47	50
Total		30	60	90

CerTest Lactoferrin vs Human Lactoferrin ELISA kit		
	Mean Value	95% CI (confidence interval)
Sensitivity	90.0%	73.5 – 97.9%
Specificity	78.3%	65.8 – 87.9%
PPV	67.5%	50.9 - 81.4%
NPV	94.0%	83.5 –98.7%

Specificity for a cut-off value of 2.5 µg hLf/g: 100.0% (mean value); 83.9-100.0% (95% confidence interval).

The results showed a high sensitivity and specificity to detect human lactoferrin using CerTest Lactoferrin.

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Cross reactivity


An evaluation was performed to determine the cross reactivity of CerTest Lactoferrin; no cross reactivity against other faecal markers occasionally present in faeces:

<i>Adenovirus</i>	<i>Clostridium perfringens</i>	<i>Human Haemoglobin</i>	<i>Salmonella typhi</i>
<i>Astrovirus</i>	<i>Cryptosporidium</i>	<i>Human Transferrin</i>	<i>Salmonella typhimurium</i>
<i>Bovine Haemoglobin</i>	<i>Entamoeba dispar</i>	<i>Legionella pneumophila</i>	<i>Shigella boydii</i>
<i>Bovine Lactoferrin</i>	<i>Entamoeba histolytica</i>	<i>Listeria monocytogenes</i>	<i>Shigella dysenteriae</i>
<i>Bovine Transferrin</i>	<i>Escherichia coli O111</i>	<i>Norovirus GI</i>	<i>Shigella flexneri</i>
<i>Campylobacter coli</i>	<i>Escherichia coli O26</i>	<i>Norovirus GII</i>	<i>Shigella sonnei</i>
<i>Campylobacter jejuni</i>	<i>Escherichia coli O157</i>	<i>Pig Haemoglobin</i>	<i>Streptococcus pneumococcal</i>
<i>Clostridium difficile antigen GDH</i>	<i>Giardia</i>	<i>Rotavirus</i>	<i>Streptococcus pyogenes</i>
<i>Clostridium difficile Toxin A</i>	<i>Helicobacter pylori</i>	<i>Salmonella enteritidis</i>	<i>Yersinia enterocolitica O:3</i>
<i>Clostridium difficile Toxin B</i>	<i>Human Calprotectin</i>	<i>Salmonella paratyphi A</i>	<i>Yersinia enterocolitica O:9</i>

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		ELISA test: Human Lactoferrin ELISA kit		
		+	-	Total
IC test: CerTest Lactoferrin	+	27	13	40
	-	3	47	50
	Total	30	60	90

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1. Angriman I. et al. Enzymes in feces: Useful markers of chronic inflammatory bowel disease. Clinica Chimica Acta 381 Feb 2007, p. 63-68
2. Guerrant R. et al. Measurement of Fecal Lactoferrin as a Marker of Fecal Leukocytes. Journal of Clinical Microbiology, Vol 30 No 5; May 1992, p. 1238-1242
3. Langhorst, M.D. et al. Noninvasive Markers in the Assessment of Intestinal Inflammation in Inflammatory Bowel Diseases: Performance of Fecal Lactoferrin, Calprotectin and PMN-Elastase, CRP, and Clinical Indices. Am. J. Gastroenterol. 2008; Vol 103, p. 162-169.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/

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