CERTEST

Clostridium difficile

GDH+Toxin A+B+

Lactoferrin

ONE STEP

Clostridium difficile antigen GDH,
Toxin A, Toxin B and *Lactoferrin*


COMBO CARD TEST

CERTEST BIOTEC S.L.

For information purposes only

For information purposes only

CERTEST *Clostridium difficile* GDH+Toxin A+B+Lactoferrin

One Step test to detect *Clostridium difficile* antigen GDH, Toxin A, Toxin B and Lactoferrin combo card format/
Prueba combo de un solo paso para detección del antígeno GDH de *Clostridium difficile*, Toxina A, Toxina B y Lactoferrin en formato cassette

ENGLISH

INTENDED USE

CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection of *Clostridium difficile* Glutamate Dehydrogenase (GDH), Toxin A and Toxin B in stool samples and for the semi-quantitative detection of human lactoferrin (hLf) in stool samples.

CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin combo card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of *Clostridium difficile* infection and intestinal inflammation presence.

INTRODUCTION

Clostridium difficile (*C. difficile*), a Gram-positive spore bearing anaerobic bacterium is the major aetiological agent of diarrhoea and colitis associated with antibiotics. *C. difficile* is the most common cause of health care-associated diarrhoea in developed countries and is a major source of nosocomial morbidity and mortality worldwide.

Disease due to *C. difficile* develops when the organism is allowed to proliferate in the colon, most commonly after antibiotic use has eliminated competing flora. *C. difficile* can release two high-molecular-weight toxins, toxin A and toxin B, which are responsible for the clinical manifestations, which range from mild, self-limited watery diarrhoea to fulminant pseudomembranous colitis, toxic megacolon, and death.

Clostridium difficile Glutamate Dehydrogenase (GDH) is an enzyme produced in large quantities by all toxigenic and non-toxigenic strains, making it an excellent marker for the organism.

The toxigenic culture (TC) is used as the gold standard technique to determine *Clostridium difficile* infection. This method consists in culture and isolation of *C. difficile* from faeces, followed by toxin testing of the isolate, a labour-intensive assay to obtain a result.

Lactoferrin is a glycoprotein component of neutrophil secondary granules, a primary component of the acute inflammatory response and is released from faecal leukocytes. This protein is resistant to proteolysis in the faeces and may serve as a marker of inflammation in the intestine. The major cause of faecal neutrophils in patients with chronic diarrhoea is chronic inflammatory bowel disease of the colon (i.e., Crohn's Disease and ulcerative colitis). Lactoferrin has been also studied as a predictor of infection with invasive enteropathogens in children with diarrhoea. Bacterial inflammatory diarrhoea may be caused by *Shigella*, *Salmonella*, *Campylobacter* and *Clostridium difficile*.

TEST PRINCIPLE

CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin is based on the principle of a qualitative immunochromatographic assay for the determination of *Clostridium difficile* Glutamate Dehydrogenase (GDH), Toxin A, Toxin B and human lactoferrin in stool samples.

Strip A consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against GDH and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-GDH) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

Strip B consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against Toxin A and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-Toxin A) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

Strip C consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against Toxin B and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-Toxin B) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

Strip D consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against human lactoferrin and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-human lactoferrin) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

If the sample is GDH positive, the antigen of the diluted sample reacts with the red-coloured conjugates complex (anti-GDH monoclonal antibodies-red polystyrene microspheres) in the strip A, if the sample is Toxin A positive, the antigens of the diluted sample react with the red-coloured conjugates complex (anti-Toxin A monoclonal antibodies-red polystyrene microspheres) in the strip B, if the sample is Toxin B positive, the antigens of the diluted sample react with the red-coloured conjugated complex (anti-Toxin B monoclonal antibodies-red polystyrene microspheres) in strip C and if the sample is human lactoferrin positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-human lactoferrin monoclonal antibodies-red polystyrene microspheres) in strip D, which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the

binding conjugate complexes migrate. The anti-GDH antibodies present on the membrane of strip A (test line), the anti-Toxin A antibodies present on the membrane of strip B (test line), the anti-Toxin B antibodies present on the membrane of strip C (test line) and the anti-human lactoferrin antibodies present on the membrane of strip D (test line) capture the coloured conjugate and the red line will be visible in the strips. These bands are used to interpret the result.

If the sample is negative, there is no GDH, Toxin A, Toxin B and human lactoferrin presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-GDH, the anti-Toxin A, the anti-Toxin B and the anti-human lactoferrin antibodies present on the membranes (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in the four strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on four membranes will capture control green-conjugate complex and four control lines will always appear. The presence of these green lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin combo card test. Do not use any other commercial kit component.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 24 hours prior to testing. For longer storage, the specimen must be kept frozen at -20°C. Freezing and thawing cycles are not recommended. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.


Please take note that there are two different stool collection tubes.

- a. "GDH+Toxin A+Toxin B" stool collection tubes (without reducer piece).

Specimen preparation (see illustration):

1. Take out the cap of the collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2) making sure that, at each insertion, only the stick's screw is covered with the sample to take the appropriate amount of faecal sample and add it to the collection tube. Any additional sample that exceeds the stick's screw could cause wrong results. For liquid samples, add approx. 125µL in the collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).


b. "Lactoferrin" stool collection tube (with a reducer piece).

Specimen preparation (see illustration):


1. Take out the cap 1 of the stool collection tube and use the stick to pick up sufficient sample quantity. Introduce the stick once into 4 different parts of the stool sample (1), making sure that, in each insertion, the stick's screw is covered with the sample to take the appropriate amount of stool. Then, introduce the sample in the stool collection tube (2). For liquid samples, add approx. 20µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample (2). Shake the tube in order to assure good sample dispersion (3). The stool collection tube with dilute sample can be stored for 7 days in the refrigerator (2-8°C) prior to testing.

MATERIALS

MATERIALS PROVIDED

- CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin combo card tests
- Instructions for use
- "GDH+Toxin A+B" stool collection tubes with diluent
- "Lactoferrin" stool collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE


Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay. Please take note that there are two procedures by using different stool collection tubes.

1. Proceed to shake both stool collection tubes in order to assure good sample dispersion.
2. Remove the CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin combo card test from its sealed bag just before using it.
3. Take the "GDH+Toxin A+B" stool collection tube, cut the end of the cap (4) and, using each time the same tube, dispense 3 drops in the circular window marked with the letter A (5), 3 drops in the circular window marked with the letter B (6) and 3 drops in the circular window marked with the letter C (7). Avoid adding solid particles with the liquid.
4. Take the "Lactoferrin" stool collection tube, open the cap 2 (8) and dispense 3 drops in the circular window marked with the letter D (9). Avoid adding solid particles with the liquid.
5. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.


CERTEST *Clostridium difficile* GDH+Toxin A+B+Lactoferrin

One Step test to detect *Clostridium difficile* antigen GDH, Toxin A, Toxin B and Lactoferrin combo card format/
 Prueba combo de un solo paso para detección del antígeno GDH de *Clostridium difficile*, Toxina A, Toxina B y Lactoferrin en formato cassette


GDH Strip A-procedure

Add 3 drops in A-circular window


Toxin A Strip B-procedure

Add 3 drops in B-circular window


Toxin B Strip C-procedure

Add 3 drops in C-circular window


Lactoferrin Strip D-procedure

Add 3 drops in D-circular window


INTERPRETATION OF THE RESULTS (please refer to the illustration below)


VALID RESULTS


A: Green→Negative GDH
 B: Green→Negative Toxin A
 C: Green→Negative Toxin B
 D: Green→Negative Lactoferrin


A: Green/Red→Positive GDH
 B: Green/Red→Positive Toxin A
 C: Green/Red→Positive Toxin B
 D: Green/Red→Positive Lactoferrin


A: Green/Red→Positive GDH
 B: Green/Red→Positive Toxin A
 C: Green/Red→Positive Toxin B
 D: Green→Negative Lactoferrin

A: Green/Red→Positive GDH
 B: Green/Red→Positive Toxin A
 C: Green→Negative Toxin B
 D: Green/Red→Positive Lactoferrin


A: Green/Red→Positive GDH
 B: Green→Negative Toxin A
 C: Green/Red→Positive Toxin B
 D: Green/Red→Positive Lactoferrin


A: Green/Red→Positive GDH
 B: Green/Red→Positive Toxin A
 C: Green→Negative Toxin B
 D: Green→Negative Lactoferrin

CERTEST *Clostridium difficile* GDH+Toxin A+B+Lactoferrin


One Step test to detect *Clostridium difficile* antigen GDH, Toxin A, Toxin B and Lactoferrin combo card format/
Prueba combo de un solo paso para detección del antígeno GDH de *Clostridium difficile*, Toxina A, Toxina B y Lactoferrin en formato cassette


A: Green/Red→Positive GDH
B: Green→ Negative Toxin A
C: Green/Red→ Positive Toxin B
D: Green→Negative Lactoferrin

A: Green/Red→ Positive GDH
B: Green→Negative Toxin A
C: Green→ Negative Toxin B
D: Green/Red→Positive Lactoferrin


A: Green/Red→ Positive GDH
B: Green→Negative Toxin A
C: Green→ Negative Toxin B
D: Green→ Negative Lactoferrin


A: Green→Negative GDH
B: Green→ Negative Toxin A
C: Green→ Negative Toxin B
D: Green/Red→Positive Lactoferrin


OTHER RESULTS


A: Green→Negative GDH
B: Green/Red→Positive Toxin A
C: Green/Red→Positive Toxin B
D: Green/Red→Positive Lactoferrin


A: Green→Negative GDH
B: Green/Red→Positive Toxin A
C: Green→Negative Toxin B
D: Green/Red→Positive Lactoferrin


A: Green→Negative GDH
B: Green→Negative Toxin A
C: Green/Red→Positive Toxin B
D: Green/Red→ Positive Lactoferrin

A: Green→Negative GDH
B: Green/Red→Positive Toxin A
C: Green/Red→Positive Toxin B
D: Green→ Negative Lactoferrin

A: Green→Negative GDH
B: Green/Red→Positive Toxin A
C: Green→Negative Toxin B
D: Green→Negative Lactoferrin

A: Green→Negative GDH
B: Green→Negative Toxin A
C: Green/Red→Positive Toxin B
D: Green→Negative Lactoferrin

	A (GDH)	B (Toxin A)	C (Toxin B)	D (Lactoferrin)	Interpretation of the results
1.	-	-	-	-	There are not GDH, Toxin A and Toxin B of <i>Clostridium difficile</i> presence and neither human lactoferrin. No infection caused by <i>Clostridium difficile</i> .
2.	+	+	+	+	There are GDH, Toxin A and Toxin B of <i>Clostridium difficile</i> and human lactoferrin presence. Infection caused by <i>Clostridium difficile</i> and there is an intestinal inflammation associated with this or other infection.
3.	+	+	+	-	There are GDH, Toxin A and Toxin B of <i>Clostridium difficile</i> presence. Infection caused by <i>Clostridium difficile</i> .
4.	+	+	-	+	There is GDH, Toxin A and human Lactoferrin presence. Infection caused by <i>Clostridium difficile</i> and there is an intestinal inflammation associated with this or other infection.
5.	+	-	+	+	There is GDH, Toxin B and human Lactoferrin presence. Infection caused by <i>Clostridium difficile</i> and there is an intestinal inflammation associated with this or other infection.
6.	+	+	-	-	There is GDH and Toxin A presence. Infection caused by <i>Clostridium difficile</i> .
7.	+	-	+	-	There is GDH and Toxin B presence. Infection caused by <i>Clostridium difficile</i> .

CERTEST *Clostridium difficile* GDH+Toxin A+B+Lactoferrin

One Step test to detect *Clostridium difficile* antigen GDH, Toxin A, Toxin B and Lactoferrin combo card format/
Prueba combo de un solo paso para detección del antígeno GDH de *Clostridium difficile*, Toxina A, Toxina B y Lactoferrin en formato cassette

8	+	-	-	+	There is GDH and human Lactoferrin presence. Infection caused by <i>Clostridium difficile</i> and there is an intestinal inflammation associated with this or other infection.
	GREEN-RED	GREEN	GREEN	GREEN-RED	
9.	+	-	-	-	There is GDH presence. Infection caused by <i>Clostridium difficile</i> .
	GREEN-RED	GREEN	GREEN	GREEN	
10.	-	-	-	+	There is human lactoferrin presence which might mean an intestinal inflammation associated with other infections or non-infectious causes of intestinal inflammation.
	GREEN	GREEN	GREEN	GREEN-RED	
11,14.	-	+	+	+	If this result appears it must be repeat the test using a fresh sample. If results are again positive for Toxin A and B and negative for GDH, the sample should be considered positive for Toxin A and B.
	GREEN	GREEN-RED	GREEN-RED	GREEN-RED	
12,15.	-	+	-	+	If this result appears it must be repeat the test using a fresh sample. If the result is again positive for Toxin A and negative for GDH, the sample should be considered positive for Toxin A.
	GREEN	GREEN-RED	GREEN	GREEN-RED	
13,16.	-	+	+	+	If this result appears it must be repeat the test using a fresh sample. If the result is again positive for Toxin B and negative for GDH, the sample should be considered positive for Toxin B.
	GREEN	GREEN-RED	GREEN-RED	GREEN-RED	
17.	Any other result				Invalid result: either A, B, C or D, we recommend repeating the assay using the same sample with another test.

11, 12, 13, 14, 15 and 16, a very low percentage of specimens might result negative for GDH but positive for toxins.

INVALID: Total absence of any control coloured line (GREEN) regardless the appearance or not of the test lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of GDH, Toxin A, Toxin B and human lactoferrin present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
3. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit of the test.
4. CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
5. Positive results determine the presence of GDH, Toxin A, Toxin B of *Clostridium difficile* and human lactoferrin in faecal samples. A positive result should be followed up with additional laboratory techniques (toxigenic culture) to determine the strain of *Clostridium difficile*. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
6. A negative result is not meaningful because of it is possible GDH, Toxin A, Toxin B of *Clostridium difficile* and/or human lactoferrin concentration in the stool samples is lower than the detection limit or cut-off value. If the symptoms or situation still persist, a *Clostridium difficile* determination should be carried out on a sample from an enrichment culture.
7. Lactoferrin is a component of breast milk; the test will be positive in breast fed children and should not be used to evaluate neonates receiving breast milk.
8. Stool from patients with active inflammatory bowel diseases that usually involve significant neutrophilic inflammation of the intestine, such as Crohn's disease and ulcerative colitis, would be positive for faecal lactoferrin.

CERTEST *Clostridium difficile* GDH+Toxin A+B+Lactoferrin

One Step test to detect *Clostridium difficile* antigen GDH, Toxin A, Toxin B and Lactoferrin combo card format/
Prueba combo de un solo paso para detección del antígeno GDH de *Clostridium difficile*, Toxina A, Toxina B y Lactoferrin en formato cassette

9. Bloody stool samples can contain components that may cause non-specific reactions in the test. Every bloody stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
10. Mucous stool samples may cause non-specific creations in the test. Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.

EXPECTED VALUES

Clostridium difficile is associated with 95-100% cases of pseudomembranous colitis, 60-75% cases of antibiotic-associated colitis and 35% of cases of antibiotic-associated diarrhoea cases.

Some studies established equal or higher 10µg hLf/g faeces as cut-off value to allow detect adult patients with gastrointestinal inflammatory pathology that will require to diagnosis additional diagnostic invasive procedures.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

Detection limit values of CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin (GDH strip) is: 0.78-0.10 ng/mL of *Clostridium difficile* GDH recombinant protein.

Detection limit values of CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin (Toxin A strip) is: 2-0.25 ng/mL of *Clostridium difficile* Toxin A

Detection limit values of CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin (Toxin B strip) is: 3.12-0.19 ng/mL of *Clostridium difficile* Toxin B.

Cut-off value of CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin (Lactoferrin strip) is: 100ng/mL (10µg hLf/g faeces) for human lactoferrin.

Clinical sensitivity and specificity (GDH+Toxin A+Toxin B)

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin, CerTest) and an Evaluation Criteria.

For strip A: GDH, the Evaluation Criteria consisted on the evaluation in parallel of two rapid test (CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin) strip A, vs other rapid test from the competitor; the discrepant results were confirmed by qPCR technique.

For strip B and C: Toxin A and B, the Evaluation Criteria consisted on the evaluation in parallel of two rapid test (CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin) strip B and C, vs other rapid test from the competitor; the discrepant results were confirmed by qPCR technique.

The results were as follows:

Clinical sensitivity and specificity (GDH)

IC test: CerTest <i>Clostridium difficile</i> GDH+Toxin A+B+Lactoferrin (Strip A: GDH)	Evaluation Criteria (GDH)		
		+	-
+	51	0	51
-	0	75	75
Total	51	75	126

CerTest <i>Clostridium difficile</i> GDH+Toxin A+B+Lactoferrin (GDH) vs Evaluation Criteria (GDH)		
	Mean Value	95% CI (confidence interval)
Sensitivity	100.0%	93.0 – 100.0%
Specificity	100.0%	95.2 – 100.0%
PPV	100.0%	93.0 - 100.0%
NPV	100.0%	95.2 –100.0%

CERTEST *Clostridium difficile* GDH+Toxin A+B+Lactoferrin

One Step test to detect *Clostridium difficile* antigen GDH, Toxin A, Toxin B and Lactoferrin combo card format/
Prueba combo de un solo paso para detección del antígeno GDH de *Clostridium difficile*, Toxina A, Toxina B y Lactoferrin en formato cassette

Clinical sensitivity and specificity (Toxin A)

IC test: CerTest <i>Clostridium difficile</i> GDH+Toxin A+B+Lactoferrin (Strip B: Toxin A)	Evaluation Criteria (Toxin A)			Total
		+	-	
+	28	0	28	
-	1	97	98	
Total	29	97	126	

CerTest <i>Clostridium difficile</i> GDH+Toxin A+B+Lactoferrin (Toxin A) vs Evaluation Criteria (Toxin A)		
	Mean Value	95% CI (confidence interval)
Sensitivity	96.6%	92.2 – 99.9%
Specificity	100.0%	96.2 – 100.0%
PPV	100.0%	87.7 - 100.0%
NPV	99.0%	94.4 –100.0%

Clinical sensitivity and specificity (Toxin B)

IC test: CerTest <i>Clostridium difficile</i> GDH+Toxin A+B+Lactoferrin (Strip C: Toxin B)	Evaluation Criteria (Toxin B)			Total
		+	-	
+	36	1	37	
-	0	89	89	
Total	36	90	126	

CerTest <i>Clostridium difficile</i> GDH+Toxin A+B+Lactoferrin (Toxin B) vs Evaluation Criteria (Toxin B)		
	Mean Value	95% CI (confidence interval)
Sensitivity	100.0%	90.3 – 100.0%
Specificity	98.9%	94.0 – 100.0%
PPV	97.3%	85.8 – 99.9%
NPV	100.0%	95.9 –100.0%

Clinical sensitivity and specificity (Lactoferrin)

CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin (strip D: Lactoferrin) was evaluated with a commercial ELISA (Human Lactoferrin ELISA kit, ICL). The results were as follows:

Data for cut-off value 10 µg hLf/g

IC test: CerTest <i>Clostridium difficile</i> GDH+Toxin A+B+Lactoferrin (Strip D: Lactoferrin)	ELISA test: Human Lactoferrin ELISA kit			Total
		-	-	
+	27	13	40	
-	3	47	50	
Total	30	60	90	

CERTEST *Clostridium difficile* GDH+Toxin A+B+Lactoferrin

One Step test to detect *Clostridium difficile* antigen GDH, Toxin A, Toxin B and Lactoferrin combo card format/
Prueba combo de un solo paso para detección del antígeno GDH de *Clostridium difficile*, Toxina A, Toxina B y Lactoferrin en formato cassette

CerTest <i>Clostridium difficile</i> GDH+Toxin A+B+Lactoferrin (Lactoferrin) vs Human Lactoferrin ELISA kit		
	Mean Value	95% CI (confidence interval)
Sensitivity	90.0%	73.5 – 97.9%
Specificity	78.3%	65.8 – 87.9%
PPV	67.5%	50.9 - 81.4%
NPV	94.0%	83.5 –98.7%

Specificity for a cut-off value of 2.5 µg hLf/g: 100.0% (mean value); 83.9-100.0% (95% confidence interval).

The results showed a high sensitivity and specificity to detect GDH, Toxin A, Toxin B of *Clostridium difficile* and human lactoferrin using CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *Clostridium difficile* GDH+Toxin A+B+Lactoferrin; no cross reactivity against gastrointestinal pathogens, other organism, substances and/or faecal markers occasionally present in faeces:

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1. WREN, M.W.D, et al. "Laboratory diagnosis of *Clostridium difficile* infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory". British Journal of Biomedical Science, 66 (1), 2009.
2. VAISHNAVI, Ch., "Clinical spectrum & pathogenesis of *Clostridium difficile* associated diseases". Indican J. Med. Res. 131, April 2010, pp 487-499.
3. Angriman I. et al. Enzymes in feces: Useful markers of chronic inflammatory bowel disease. Clinica Chimica Acta 381 Feb 2007, p. 63-68.
4. Guerrant R. et al. Measurement of Fecal Lactoferrin as a Marker of Fecal Leukocytes. Journal of Clinical Microbiology, Vol 30 No 5; May 1992, p. 1238-1242.
5. Langhorst, M.D. et al. Noninvasive Markers in the Assessment of Intestinal Inflammation in Inflammatory Bowel Diseases: Performance of Fecal Lactoferrin, Calprotectin and PMN-Elastase, CRP, and Clinical Indices. Am. J. Gastroenterol. 2008; Vol 103, p. 162-169.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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