# REAL TIME PCR DETECTION KIT

## *Zika, Dengue & Chikungunya*

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## ENGLISH

## 1. Intended of use

VIASURE *Zika, Dengue & Chikungunya* Real Time PCR Detection Kit is designed for specific detection and differentiation of Zika, Dengue and/or Chikungunya viruses in clinical samples from patients with signs and symptoms of Zika, Dengue and/or Chikungunya viruses infection. This test is intended for use as an aid in the diagnosis of the Zika, Dengue and/or Chikungunya viruses in combination with clinical and epidemiological risk factors. RNA is extracted from specimens, amplified using RT-PCR and detected using fluorescent reporter dye probes specific for Zika and Dengue and Chikungunya viruses.

## 2. Summary and Explanation

Zika (ZIKV), Dengue (DENV) and Chikungunya (CHIKV) are arthropod-borne viruses (arboviruses) that share mosquitos of the *Aedes* genus as common vector, specifically *A. aegypti* and *A. albopictus*. Environmental changes, urbanization and the globalization of travel have enhanced the spread of ZIKV, DENV and CHIKV resulting in epidemics and the co-circulation in overlapping geographic areas as well as the possibility of viral co-infection within a single host. Currently the major endemic regions include Africa, Southeast Asia, the Pacific Islands and the Americas.

Besides, these arboviruses cause similar clinical presentations, especially in the initial stages of infection (fever, headache, skin rash, muscle and joint pain, and gastrointestinal symptoms). In fact, due to neither virus possesses any specific distinguishing clinical features and the outcomes and management strategies for these three viruses are vastly different, an early and accurate diagnosis is imperative.

Therefore, ZIKV, DENV and CHIKV must be initially included in the differential diagnosis for a patient with suspicious clinical symptoms who is living in or returning from travel to an endemic area. Laboratory diagnostic tests are based on the detection of the virus, viral components (antigens or nucleic acid), or the host immunologic response to the virus. However, due to the cross-reactivity of the antibodies of these arboviruses limits the use of serology, Real Time RT-PCR is a detection method commonly used during the acute phase of the infection. In order to do that, clinical specimens that can be tested including blood, serum, plasma, urine and others.

## 3. Principle of the procedure

VIASURE *Zika, Dengue & Chikungunya* Real Time PCR Detection Kit is designed for the diagnosis of the Zika, Dengue and/or Chikungunya viruses in clinical samples. The detection is done in one step real time RT format where the reverse transcription and the subsequent amplification of specific target sequence occur in the same reaction well. The isolated RNA target is transcribed generating complementary DNA by reverse transcriptase which is followed by the amplification of a conserved region of the *envelope* gene (ZIKV), 3'Non coding region (DENV) and *NSP1* gen (CHIKV), using specific primers and a fluorescent-labelled probe.

VIASURE *Zika, Dengue & Chikungunya* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE *Zika, Dengue & Chikungunya* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase, retrotranscriptase) in an stabilized format, as well as an internal control to monitor PCR inhibition. Zika virus RNA targets are amplified and detected in Cy5 channel, Dengue Virus RNA targets are amplified and detected in FAM channel, Chikungunya virus RNA targets are amplified and detected in ROX channel and the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

#### 4. Reagents provided

VIASURE *Zika, Dengue & Chikungunya* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table 1 and Table 2:

Reference	Reagent/Material	Description	Colour	Amount
VS-ZDC1SL/ VS-ZDC1SH	<i>Zika, Dengue &amp; Chikungunya Virus</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPS, stabilizers and Internal control in stabilized format	White	6/12 X 8-well strip
VS-RB02	Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
VS-ZDC1C	<i>Zika, Dengue &amp; Chikungunya Virus</i> Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
VS-NC1	Negative control	Non template control	Violet	1 vial x 1 mL
VS-H2O	Water RNAse/DNAse free	Water RNAse/DNAse free	White	1 vial x 1 mL
VS-OCS	Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE *Zika, Dengue & Chikungunya* Real Time PCR Detection Kit with Ref. VS-ZDC106L, VS-ZDC106H, VS-ZDC112L and VS-ZDC112H.

Reference	Reagent/Material	Description	Colour	Amount
VS-ZDC1PL/ VS-ZDC1PH	<i>Zika, Dengue &amp; Chikungunya Virus</i> 96-well plate	A mix of enzymes, primers probes, buffer, dNTPS, stabilizers and Internal control in stabilized format	White	1 plate
VS-RB02	Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
VS-ZDC1C	<i>Zika, Dengue &amp; Chikungunya Virus</i> Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
VS-NC1	Negative control	Non template control	Violet	1 vial x 1 mL
VS-H2O	Water RNAse/DNAse free	Water RNAse/DNAse free	White	1 vial x 1 mL
VS-OCS	Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 X 8-cap strip

Table 2. Reagents and materials provided in VIASURE *Zika, Dengue & Chikungunya* Real Time PCR Detection Kit with Ref VS-ZDC113L and VS-ZDC113H.

## 5. Reagents and equipment to be supplied by the user

The following list includes materials that are required for using but not included in the VIASURE *Zika, Dengue & Chikungunya* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler) (to check compatibility see Annex I).
- RNA extraction kit.
- Centrifuge for 1.5 mL tubes.
- Vortex.
- Micropipettes (0.5-20  $\mu$ L, 20-200  $\mu$ L).
- Filter tips.
- Powder-free disposal gloves.

## 6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until expiration date stated in the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate in aliquots to minimize freeze and thaw cycles.
- Keep components away from sunlight.

## 7. Precautions for users

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious as well as all reagents and materials that have been exposed to the samples and must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

## 8. Test procedure

### 8.1. RNA EXTRACTION

Perform the sample preparation according to the recommendations appearing in the instructions for use of extraction kit used.

For RNA extraction from clinical samples (blood, serum, plasma, urine, tissues and others) you can use your manually or automatic routine optimized system. Also, you can use any commercially available RNA extraction kit and follow the manufacturer’s instructions for use. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- Maxwell® 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® 16 instrument (Promega).
- QIAamp Viral RNA Mini Kit, using QIAcube instrument (Qiagen).
- High Pure Viral Nucleic Acid Kit (Roche).

8.2. LYOPHILIZED POSITIVE CONTROL

*Zika, Dengue & Chikungunya Virus* Positive Control contains high copies template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Zika, Dengue & Chikungunya Virus* Positive Control (red vial) adding 100 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate in aliquots to minimize freeze and thaw cycles.

8.3. PCR PROTOCOL

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run. Peel off protective aluminium seal from plates or strips.

1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

2) Adding samples and controls.

Add 5 µL of RNA sample, reconstituted *Zika, Dengue & Chikungunya Virus* Positive Control (red vial) or Negative Control (violet vial) in different wells and close the wells with the caps provided. Centrifuge briefly.

Load the plate or the strips in the thermocycler.

3) Set up your thermocycler.

Program your thermocycler following the conditions below and start the run:

Cycles	Step	Time	Temperature
1	Reverse transcription	15 min	45°C
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 3. PCR protocol

Fluorogenic data should be collected during the extension step (\*) through the Cy5 (*Zika virus*), FAM (*Dengue*), ROX (*Chikungunya*) and HEX, JOE or VIC channels (Internal Control (IC)). Depending on the equipment used select the proper detection channel (see Annex 2). In Applied Biosystems 7500 Fast Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none.

## 9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in negative control well and the presence of signal for Zika, Dengue and Chikungunya viruses positive control well. Check Internal Control signal to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software itself of the used real time PCR equipment according to manufacturer’s instructions.

Using the following table read and analyze the results:

Zika Virus	Dengue Virus	Chikungunya Virus	Internal control	Negative control	Positive control	Interpretation
+	+	+	+/-	-	+	Zika , Dengue and Chikungunya viruses Positive
-	-	-	+/-	-	+	Zika, Dengue and Chikungunya viruses Negative
+	-	-	+/-	-	+	Zika Virus Positive, Dengue and Chikungunya Viruses Negative
+	+	-	+/-	-	+	Zika and Dengue Viruses Positive, and Chikungunya Virus Negative
+	-	+	+/-	-	+	Zika and Chikungunya Viruses Positive, and Dengue Virus Negative
-	+	-	+/-	-	+	Dengue Virus Positive, Zika and Chikungunya Viruses Negative
-	+	+	+/-	-	+	Dengue and Chikungunya Viruses Positive, Zika Virus Negative
-	-	+	+/-	-	+	Chikungunya Virus Positive, Zika and Dengue Viruses Negative
+	+	+	+	+	+	Experiment fail
-	-	-	-	-	-	Experiment fail

Table 4. Sample interpretation


+: Amplification curve

-: No amplification curve

A sample is considered positive if the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System.

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

## 10. Limitations of the test

- The result of the test should be evaluated by a health care professional and evaluated in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated with serum and urine samples.
- The quality of the test depends on the quality of the sample; proper RNA from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by Zika, Dengue and/or Chikungunya viruses, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

## 11. Quality control

VIASURE *Zika, Dengue & Chikungunya* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

## 12. Performance characteristics

### 12.1. ANALYTICAL SENSITIVITY

VIASURE *Zika, Dengue & Chikungunya* Real Time PCR Detection Kit has a detection limit of  $\geq 10$  RNA copies per reaction for ZIKV, DENV and CHIKV (Figure 2, 3 and 4).

Figure 2. Dilution series of Zika Virus ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Cy5 channel).


Figure 3. Dilution series of Dengue Virus ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (FAM channel).



Figure 4. Dilution series of Chikungunya Virus ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (ROX channel).

## 12.2. ANALYTICAL SPECIFICITY

The specificity of the Zika, Dengue and Chikungunya viruses assays were confirmed by testing a panel consisting of 10 different microorganisms representing the most common arbovirus.

No cross-reactivity of Zika Virus was detected against any of the following microorganisms tested.

Cross-reactivity testing			
Chikungunya virus strain S27 Petersfield	-	St Louis Encephalitis virus strain 17D	-
Dengue 1 virus strain Hawaii	-	West Nile virus strain H160/99	-
Dengue 2 virus strain New Guinea C	-	West Nile virus Heja	-
Dengue 3 virus strain H87	-	West Nile virus Ug37	-
Dengue 4 virus strain H241	-	Yellow Fever virus strain 17D	-

Table 6. Reference pathogenic microorganisms used in this study.

No cross-reactivity of Dengue Virus was detected against any of the following microorganisms tested.

Cross-reactivity testing			
Chikungunya virus strain S27 Petersfield	-	West Nile virus Heja	-
Zika virus strain MR 766	-	West Nile virus Ug37	-
St Louis Encephalitis virus strain 17D	-	Yellow Fever virus strain 17D	-
West Nile virus strain H160/99	-		

Table 7. Reference pathogenic microorganisms used in this study.

No cross-reactivity of Chikungunya Virus was detected against any of the following microorganisms tested.

Cross-reactivity testing			
Zika virus strain MR 766	-	St Louis Encephalitis virus strain 17D	-
Dengue 1 virus strain Hawaii	-	West Nile virus strain H160/99	-
Dengue 2 virus strain New Guinea C	-	West Nile virus Heja	-
Dengue 3 virus strain H87	-	West Nile virus Ug37	-
Dengue 4 virus strain H241	-	Yellow Fever virus strain 17D	-

Table 8. Reference pathogenic microorganisms used in this study.

## 12.3. ANALYTICAL REACTIVITY

The reactivity of VIASURE *Zika, Dengue & Chikungunya* Real Time PCR Detection Kit for Zika Virus was evaluated against Zika virus strain MR 766 showing positive results.

The reactivity of VIASURE *Zika, Dengue & Chikungunya* Real Time PCR Detection Kit for Dengue Virus was evaluated against Dengue 1 virus strain Hawaii, Dengue 2 virus strain New Guinea C, Dengue 3 virus strain H87 and Dengue 4 virus strain H241 showing positive results.

The reactivity of VIASURE *Zika, Dengue & Chikungunya* Real Time PCR Detection Kit for Chikungunya was evaluated against Chikungunya virus S27 Petersfield showing positive results.

ANNEX 1:

**COMPATIBILITY OF THE MOST COMMON REAL TIME PCR EQUIPMENT**

Low profile strips can be used in all PCR thermocyclers equipped with low profile block, like systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with high or regular profile block, like systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Applied Biosystems	7500 Fast Real-Time PCR System
Applied Biosystems	7500 Fast Dx Real-Time PCR System
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
Bio-Rad	CFX96 Touch™ Real-Time PCR Detection System
Roche	LightCycler® 480 Real-Time PCR System
Roche	LightCycler® 96 Real-Time PCR System
Agilent Technologies	AriaMx Real-Time PCR System
DNA-Technology	DTlite Real-Time PCR System
DNA-Technology	DTprime Real-time Detection Thermal Cycler
Qiagen	Rotor-Gene® Q*
Cepheid	SmartCycler®

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Bio-Rad	CFX96 Touch™ Deep Well Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ 5 Real-Time PCR Detection System
Eppendorf	Mastercycler™ ep <i>realplex</i>
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
Abbott	Abbott m2000 RealTime System
BIONEER	Exicycler™ 96
DNA-Technology	DTlite Real-Time PCR System
DNA-Technology	DTprime Real-time Detection Thermal Cycler
Qiagen	Rotor-Gene® Q*
Cepheid	SmartCycler®

Table A1/A2. Compatible low and high profile Real Time PCR systems.

\* The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into specific Rotor-Gene® Q or SmartCycler® tubes.

ANNEX 2:

**DETECTION CHANNELS OF MOST COMMON REAL TIME PCR EQUIPMENT**

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/60	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

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
### 13. Bibliography/Bibliografía

1. R.S. Lanciotti *et al.* Genetic and serologic properties of Zika virus associated with an epidemic, Yap State, Micronesia, 2007. *Emerging Infectious Diseases* 2008; 14(8): 1232-1239.
2. O. Faye *et al.* Quantitative real-time PCR detection of Zika virus and evaluation with field-caught mosquitoes. *Virology Journal* 2013; 10: 311.
3. A.C. Gourinat *et al.* Detection of Zika virus in urine. *Emerging Infectious Diseases* 2015; 21(1): 84-86.
4. D. Musso *et al.* Detection of Zika virus in saliva. *Journal of Clinical Virology* 2015; 68: 53-55.
5. D. Musi *et al.* Potential Sexual Transmission of Zika Virus. *Emerging Infectious Diseases* 2015 21(2): 359-361.
6. J. Mlakar. Zika Virus Associated with Microcephaly. *New England Journal of Medicine* 2016.
8. C. Drosten *et al.* Rapid detection and quantification of RNA of Ebola and Marburg viruses, Lassa virus, Crimean-Congo hemorrhagic fever virus, Rift Valley fever virus, dengue virus, and yellow fever virus by real-time reverse transcription-PCR. *Journal of Clinical Microbiology* 2002; 40(7): 2323-2330.
9. A. Dumoulin *et al.* Pan-dengue virus detection by PCR for travelers returning from the tropics. *Journal of Clinical Microbiology* 2008; 46(9): 3104-3106.
10. J. Liu *et al.* Development of a TaqMan Array Card for Acute-Febrile-Illness Outbreak Investigation and Surveillance of Emerging Pathogens, Including Ebola Virus. *Journal of Clinical Microbiology* 2016; 54(1): 49-58.
11. S.K. Mardekian and A.L. Roberts. Diagnostic Options and Challenges for Dengue and Chikungunya Viruses. *BioMed Research International* 2015: 834371.
12. K.A. Tsetsarkin *et al.* Interspecies transmission and chikungunya virus emergence. *Current Opinion in Virology* 2016; 16: 143-150.
13. R.S. Lanciotti *et al.* Chikungunya virus in US travelers returning from India, 2006. *Emerging Infectious Diseases journal* 2007; 13(5): 764-767.
14. Centers for Disease Control and Prevention. (<http://www.cdc.gov/>).
15. World Health Organization. (<http://www.who.int/>).

### 14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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VIASURE *Zika, Dengue & Chikungunya* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96 Touch™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, Rotor-Gene®Q and SmartCycler®. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™ is registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

CERTEST BIOTEC S.L.  
Pol. Industrial Río Gállego II, Calle J, Nº 1,  
50840, San Mateo de Gállego, Zaragoza (SPAIN)

[www.certest.es](http://www.certest.es) 