CERTEST

S. pneumoniae+
Legionella

ONE STEP

S. pneumoniae+*Legionella*

COMBO CARD TEST

CERTEST BIOTEC S.L.

For information purposes only

ENGLISH**INTENDED USE**

CerTest *S. pneumoniae*+*Legionella* one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection of *Streptococcus pneumoniae* and *Legionella pneumophila* serogroup 1 in human urine samples.

CerTest *S. pneumoniae*+*Legionella* combo card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of *pneumoniae* and/or legionellosis in infected humans.

INTRODUCTION

Gram positive bacteria, *Streptococcus pneumoniae* is one of the most important pathogen that affects mainly children and elderly and can cause life-threatening diseases. *Streptococcus pneumoniae* infection leads to many clinical manifestations including meningitis, septicaemia, bacteraemia, pneumonia, acute otitis media, and sinusitis. It has caused great morbidity and mortality worldwide, especially in children. Pneumococcal infection annually has caused approximately 14.5 million cases of invasive pneumococcal disease (IPD) and 0.7–1 million deaths in children under five years old, mostly in developing and underdeveloped countries. Based on differences in the polysaccharide capsules of the pneumococcal cell wall, *Streptococcus pneumoniae* is classified into over 90 serotypes, which present different antigenic properties and induce different inflammatory responses.

Legionnaires' Disease, named after the outbreak in 1976 at the American Legion convention in Philadelphia, is caused by *Legionella pneumophila* and is characterized as an acute febrile respiratory illness ranging in severity from mild illness to fatal pneumonia. The disease occurs in both epidemic and endemic forms and sporadic cases are not easily differentiated from other respiratory infections by clinical symptoms. An estimated 25,000 to 100,000 cases of *Legionella* infection occur in the United States annually. The resulting mortality rate, ranging from 25% to 40%, can be lowered if the disease is diagnosed rapidly and appropriate antimicrobial therapy is instituted early. Known risk factors include immunosuppression, cigarette smoking, alcohol consumption and concomitant pulmonary disease. The young and the elderly are particularly susceptible. *Legionella pneumophila* is responsible for 80-90% of reported cases of Legionella infection with serogroup 1 accounting for greater than 70% of all legionellosis. Current methods for the laboratory detection of pneumonia caused by *Legionella pneumophila* require a respiratory specimen (e.g. expectorated sputum, bronchial washing, transtracheal aspirate, lung biopsy) or paired sera (acute and convalescent) for an accurate diagnosis.

The CerTest *S. pneumoniae*+*Legionella* allow for early diagnosis of *S. pneumoniae* and/or *Legionella pneumophila* serogroup 1 infection through detection of some specific soluble antigens present in the urine of patients with *S. pneumoniae* and/or *L. pneumophila*. The test is rapid, giving a result at 15 minutes, and utilizes a urine specimen which is convenient for collection, transport, and subsequent detection of early, as well as later, stages of disease.

TEST PRINCIPLE

CerTest *S. pneumoniae*+*Legionella* is based on the principle of a qualitative immunochromatographic assay for the determination of *Streptococcus pneumoniae* (*S. pneumoniae*) and *Legionella pneumophila* (*L. pneumophila*) in urine samples.

Strip A consists of a nitrocellulose membrane pre-coated with monoclonal antibodies on the test line (T), in the results window, against *S. pneumoniae* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (monoclonal antibodies anti-*S. pneumoniae*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to blue polystyrene latex, forming two coloured conjugate complexes.

Strip B consists of a nitrocellulose membrane pre-coated with polyclonal antibodies on the test line (T), in the results window, against *L. pneumophila* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (polyclonal antibodies anti-*L. pneumophila* conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to blue polystyrene latex, forming two coloured conjugate complexes.

If the sample is *Streptococcus pneumoniae* positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*S. pneumoniae* monoclonal antibodies-red polystyrene microspheres) in the strip A, and if the sample is *L. pneumophila* positive, the antigens of the diluted sample react with the red-coloured conjugated complex (anti-*L. pneumophila* polyclonal antibodies-red polystyrene microspheres) in strip B, which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*S. pneumoniae* antibodies present on the membrane of strip A (test line) and the anti-*L. pneumophila* antibodies present on the membrane of strip B (test line) capture the coloured conjugate and the red line will be visible in both strips. These bands are used to interpret the result.

If the sample is negative, there is no *S. pneumoniae* and *L. pneumophila* presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-*S. pneumoniae*

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and anti-*L. pneumophila* antibodies present on the membranes (test lines) will not capture the antigen-red -coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in both strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on both membranes will capture control blue-conjugate complex and both control lines will always appear. The presence of these blue lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the CerTest *S. pneumoniae*+*Legionella* combo card test. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

SPECIMEN COLLECTION AND PREPARATION

Urine specimens should be collected in standard containers. The samples can be stored at room temperature (15-30°C) if assayed within 24 hours of collection. Alternatively, specimens may be stored at 2-8°C for up to 14 days or at -10°C to -20°C for longer periods before testing. When necessary, urine specimens should be shipped in leak-proof containers at 2-8°C or frozen. Allow all specimens to equilibrate to room temperature before testing.

MATERIALS

MATERIALS PROVIDED

- CerTest *S. pneumoniae*+*Legionella* combo card tests
- Plastic pipettes
- Testing tubes
- Reagent (sample and controls diluent)
- Positive Control Swab: Inactivated *S. pneumoniae* swab + testing tube+ pipette
- Positive Control Swab: Inactivated *L. pneumophila* swab + testing tube + pipette
- Negative Control Swab: negative swab + testing tube + pipette
- Instructions for use

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE

Allow tests, urine samples, diluent and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay

Procedure A: Patient samples


1. Using a separate testing tube or vial for each sample. Add 6 drops of urine sample (1).
2. Add 2 drop of Reagent into the testing tube or vial and mix (2). Homogenize the sample.
3. Remove CerTest *S. pneumoniae*+*Legionella* card test from its sealed bag just before using it.

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- Use a separate pipette and device for each sample or control. Dispense 3 drops from the testing tube in the circular window marked with the letter A (3) and 3 drops, using the same testing tube, in the circular window marked with the letter B (4).
- Read the results at 15 minutes.** Do not read the test result later than 15 minutes.


(1)
Add 6 drops of Urine Sample

(2)
Add 2 drops of Reagent and mix

S. pneumoniae Strip A-procedure


Add 3 drops in A-circular window

(3)

Legionella Strip B-procedure

Add 3 drops in B-circular window


(4)

Procedure B: Positive Control Swabs.


- Hold Reagent vertically. Add slowly 13 free falling drops of Reagent into the testing tube (1a).
- Immediately remove *S. pneumoniae* Positive Control swab from the pouch and put the swab into the testing tube with reagents, mix 1 minute and extract as much liquid possible from the swab, squeezing the sides of the tube as the swab is withdrawn. Discard the swab (2a).
- Remove the CerTest *S. pneumoniae*+*Legionella* device from its sealed bag just before using.
- Use a separate pipette and device for each sample or control. Dispense 3 drops from the testing tube, into the circular window marked with the letter A (3a).
- Read the result at 15 minutes.** Do not read the test result later than 15 minutes.

Repeat the procedure for Positive Control Swabs using CerTest *Legionella* Positive Control swab (1b, 2b), to dispense 3 drops from the testing tube, into the circular window marked with the letter B (3b).

Procedure B: Positive Control Swabs


(1a) Add 13 drops of Reagent

(2a) Introduce the CerTest *S. pneumoniae* Positive Control swab and mix 1 minute. Extract the liquid

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

Procedure C: Negative Control Swab.


- Hold Reagent vertically. Add slowly 13 free falling drops of Reagent into the testing tube (1).
- Immediately remove Negative control swab from the pouch and put the swab into the testing tube with reagents, mix 1 minute and extract as much liquid possible from the swab, squeezing the sides of the tube as the swab is withdrawn. Discard the swab (2).
- Remove the CerTest *S. pneumoniae*+*Legionella* device from its sealed bag just before using.
- Use a separate pipette and device for each sample or control. Dispense 3 drops from the testing tube, into the circular window marked with the letter A (3) and other 3 drops from the same testing tube into the circular window marked with the letter B (4).
- 10. Read the result at 15 minutes.** Do not read the test result later than 15 minutes.

Procedure C: Negative Control Swab


Positive and negative controls should be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.

INTERPRETATION OF THE RESULTS (please refer to the illustration below)

INVALID RESULTS

	A (<i>S. pneumoniae</i>)	B (<i>Legionella</i>)	Interpretation of the results
1.	-	-	There is no <i>S. pneumoniae</i> or <i>L. pneumophila</i> presence. No infection caused by <i>S. pneumoniae</i> or <i>L. pneumophila</i> serogroup 1. Negative control result.
2.	+	+	There is <i>S. pneumoniae</i> and <i>L. pneumophila</i> presence. Infection caused by <i>S. pneumoniae</i> and <i>L. pneumophila</i> serogroup 1. Positive controls result.
3.	+	-	There is <i>S. pneumoniae</i> presence. Infection caused by <i>S. pneumoniae</i> . <i>S. pneumoniae</i> Positive control result.
4.	-	+	There is <i>L. pneumophila</i> presence. Infection caused by <i>L. pneumophila</i> serogroup 1. <i>Legionella</i> Positive control result.
5.	Any other result	Any other result	Invalid result either A or B, we recommend repeating the assay using the same sample with another test.

INVALID: Total absence of any control coloured line (BLUE) regardless the appearance or not of the test lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. The blue lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

External Positive and Negative Controls are included in the kit (see procedure B and C). The use of positive and negative controls is recommended to assure functionality of reagents and proper performance of assay procedure.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. CerTest *S. pneumoniae*+*Legionella* should be used only with human urine samples. The use of other samples (e. g. plasma, serum or other body fluids or environmental samples (water)) has not been established.
3. Positive results determine the presence of *S. pneumoniae* and/or *Legionella* serogroup 1 in urine samples; nevertheless, a positive result should be followed up with additional laboratory techniques to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.

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4. A negative result is not meaningful because of it is possible the antigens concentration in the urine sample is lower than the detection limit values. If the symptoms or situation still persist, a *S. pneumoniae* and/or *L. pneumophila* serogroup 1 determination should be carried out from a culture or other methods.
5. Excretion of *S. pneumoniae* and *Legionella* antigen in urine may vary depending on the individual patient. For *legionella* antigen excretion may begin as early as 3 days after onset of symptoms and persist for up to 1 year afterwards. A positive CerTest *S. pneumoniae*+*Legionella* result can occur due to current or past infection and therefore is not definitive for infection without other supporting evidence.
6. It has been proven that although the test is compatible with boric acid, its presence increases the reactivity and can lead to false positives, so it is not recommended as a preservative.

EXPECTED VALUES

Streptococcus pneumoniae is the main cause of community acquired pneumonia, but also responsible cause of other diseases: meningitis, bloodstream infection, otitis. In United States of America, some analysis estimated that this pathogen was responsible of more than 4 million of illness episodes, 400000 hospitalization and 20000 deaths per year. The performance of the CerTest *Streptococcus pneumoniae* test has been evaluated by using positive and negative urine samples from European patients, as well as spiked urines.

More than 50 *Legionella* species have been characterized, and 25 species are known to cause human disease. Most human infections are caused by *Legionella pneumophila*, which is responsible for approximately 90% of the identified clinical cases, and the predominant serogroup is serogroup 1.

In recent years, hot springs have been recognized as a major source of outbreaks of legionellosis. Several outbreaks and a number of cases of *Legionella* infection that were associated with hot springs have been reported in Japan, Spain and France.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit)

Detection limit value of CerTest *S. pneumoniae*+*Legionella* is: 0.25ng/mL of CWPS for strip A (*S. pneumoniae*) and 12.5ng/mL (pool of several serovars of *Legionella pneumophila*) for strip B (*Legionella*).

Clinical sensitivity and specificity

An evaluation was performed, with urine samples, comparing the results obtained by the immunocromatographic test (CerTest *S. pneumoniae*+*Legionella*, CerTest) and other commercial rapid test (BinaxNOW® *Streptococcus pneumoniae* Antigen Card, Alere) for strip A. The results were as follows:

		IC test: BinaxNOW® <i>Streptococcus pneumoniae</i> Antigen Card		
		+	-	Total
IC test: CerTest <i>S. pneumoniae</i> + <i>Legionella</i> (<i>S. pneumoniae</i>)	+	33	1	34
	-	4	73	77
	Total	37	74	111

CerTest <i>S. pneumoniae</i> + <i>Legionella</i> (<i>S. pneumoniae</i>) vs Binax NOW® <i>Streptococcus pneumoniae</i> Antigen Card		
	Mean Value	95% confidence interval
Sensitivity	89.20%	74.6 – 97.0%
Specificity	98.60%	92.7 – 100.0%
PPV	97.10%	84.7 – 99.9%
NPV	94.80%	87.2 – 98.6%

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Another evaluation was performed, with urine samples, comparing the results obtained by an immunochromatographic test (CerTest *S. pneumoniae*+*Legionella*, CerTest) and other commercial rapid test (Binax NOW® *Legionella* Urinary Antigen Card, Alere) for strip B. The results were as follows:

		IC test: Binax NOW® <i>Legionella</i> Urinary Antigen Card		
		+	-	Total
IC test: CerTest <i>S. pneumoniae</i> + <i>Legionella</i> (<i>Legionella</i>)	+	31	1	32
	-	0	117	117
		31	118	149

CerTest <i>S. pneumoniae</i> + <i>Legionella</i> (<i>Legionella</i>) vs Binax NOW® <i>Legionella</i> Urinary Antigen Card		
	Mean Value	95% confidence interval
Sensitivity	100.00%	88.8 – 100.0%
Specificity	99.20%	95.4 – 100.0%
PPV	96.90%	83.8 - 99.9%
NPV	100.00%	96.9 – 100.0%

The results showed a high sensitivity and specificity to detect *S. pneumoniae* and *Legionella pneumophila* serogroup 1 using CerTest *S. pneumoniae*+*Legionella*.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *S. pneumoniae*+*Legionella*; no cross reactivity against other pathogens occasionally present in urine:


Legionella (Strip A) *Streptococcus pneumoniae* (Strip B)

Reproducibility Study

Evaluation studies were performed to determine reproducibility of the CerTest *S. pneumoniae*+*Legionella*, including inter-day, inter-laboratory, inter and intra lot, showing high reproducibility in all cases.

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
1. Montagnani F, Fanetti A, Stolzuoli L, Croci L, Arena F, Zanchi A, Cellesi C. J. "Pneumococcal disease in a paediatric population in a hospital of central Italy: a clinical and microbiological case series from 1992 to 2006", *Infect.* 2008 Mar;56 (3):179-84. doi: 10.1016/j.jinf.2007.12.002. Epub 2008 Jan 25.
2. Liu Y, Wang H, Chen M, Sun Z, Zhao R, Zhang L, Wang H, Zhang H, Wang L, Chu Y, Liu Y, Ni Y. "Serotype distribution and antimicrobial resistance patterns of *Streptococcus pneumoniae* isolated from children in China younger than 5 years". *Diagn Microbiol Infect Dis.* 2008 Jul;61(3):256-63. doi: 10.1016/j.diagmicrobio.2008.02.004. Epub 2008 Mar 20.
3. VPD Surveillance Manual, 5th Edition, 2012, 11 Pneumococcal Disease: Chapter 11-1.
4. ROIG, J. et. al. Comparative study of *Legionella pneumophila* and other nosocomial-acquired pneumoniae. *Chest.* 1991;99:344-50.
5. BERDAL, B.P., et al. Detection of *Legionella pneumophila* antigen in urine by enzyme-linked immunospecific assay. *J. Clin. Microbiol.* 1979;9:575-578.
6. WHITE A., et al. Rapid diagnosis of Legionnaires' disease. *Trans Am Clin. Climatol. Assoc.* 1982;93:50-62.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

IVD	<i>In vitro</i> diagnostic device Producto para diagnóstico <i>in vitro</i>	Keep dry Almacenar en lugar seco	Use by Fecha de caducidad	Manufacturer Fabricante	LOT	Batch code Número de lote	
i	Consult instructions for use Consultar las instrucciones de uso	Temperature limitation Limitación de temperatura	Contains sufficient for <n> test Contiene <n> test	DIL	Sample diluent Diluyente de muestra	REF	Catalogue number Número de referencia
Control I	Positive Control. Control Positivo	Control II	Negative Control. Control Negativo				

For information purposes only

For information purposes only

CerTest
BIOTEC S.L.

Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840, San Mateo de Gállego, Zaragoza (SPAIN)
www.certest.es