# REAL TIME PCR DETECTION KIT

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## ENGLISH

## 1. Intended of use

VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit is designed for specific identification and differentiation of human *Chlamydomphila pneumoniae*, *Mycoplasma pneumoniae* and *Legionella pneumophila* in clinical samples from patients with signs and symptoms of respiratory infection. This test is intended for use as an aid in the diagnosis of *Chlamydomphila pneumoniae*, *Mycoplasma pneumoniae* and/or *Legionella pneumophila* in combination with clinical and epidemiological risk factors. DNA is extracted from clinical specimens, multiplied using Real Time amplification and detected using specific primers and a fluorescent reporter dye probe for *Chlamydomphila pneumoniae*, *Mycoplasma pneumoniae* and/or *Legionella pneumophila*.

## 2. Summary and Explanation

Community-acquired pneumonia (CAP) is a major respiratory disease with a high prevalence in the general population, clinical heterogeneity and variable severity. Pneumonia usually causes symptoms for 3–4 weeks, and daily activities may be impaired for a further 3 weeks on average. *Chlamydomphila pneumoniae*, *Mycoplasma pneumoniae* and *Legionella pneumophila*, are some of the causes of community-acquired pneumonia.

*Legionella pneumophila*, the bacterium responsible for Legionnaires' disease, was identified in 1976 after a large outbreak at a hotel in Philadelphia, USA. The most common form of transmission of *Legionella* is inhalation of contaminated aerosols produced in conjunction with water sprays. Infection can also occur by aspiration of contaminated water or ice, particularly in susceptible hospital patients. Legionnaires' disease has an incubation period of 2 to 10 days. Untreated Legionnaires' disease usually worsens during the first week. Combining urine antigen testing with culture or molecular assays currently provides the best algorithm for diagnosis of *Legionella* disease.

*Chlamydomphila pneumoniae* cause illness by damaging the lining of the respiratory tract (throat, windpipe, and lungs). *C. pneumoniae* respiratory infection occurs worldwide and in all age groups. The seroprevalence to *Chlamydomphila pneumoniae* is low in infants but it can be higher than 50% in adults. Seroepidemiological studies show that 50 to 75% of adults have antibodies against *Chlamydomphila pneumoniae*. Most people are infected and reinfected throughout their life. However, not everyone who is exposed to *Chlamydomphila pneumoniae* develops pneumonia. *Chlamydomphila pneumoniae* has been associated with the establishment of atherosclerotic disease and heart attacks.

*Mycoplasma pneumoniae* infection is a mild illness that is most common in young adults and school-aged children. Outbreaks of *Mycoplasma pneumoniae* occur mostly in crowded environments, when small droplets of water that contain the bacteria get into the air by coughing and sneezing while in close contact with others. The incubation period is usually between 1 to 4 weeks.

The diagnosis of CAP caused by *Mycoplasma pneumoniae* and *Chlamydomphila pneumoniae* is traditionally based on cultures and serology, which have special requirements and are time-consuming.

### 3. Principle of the procedure

VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit is designed for the diagnosis of *Chlamydomphila pneumoniae*, *Mycoplasma pneumoniae* and *Legionella pneumophila* in clinical samples. After DNA isolation, the identification of *Chlamydomphila pneumoniae*, *Mycoplasma pneumoniae* and/or *Legionella pneumophila* is performed by the amplification of a conserved region of the *argR* gene for *Chlamydomphila pneumoniae*, *CARDS* gene for *Mycoplasma pneumoniae* and *mip* gene for *Legionella pneumophila*, using specific primers and a fluorescent-labelled probe.

VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence could be measured on real time PCR platforms.

VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase, retrotranscriptase) in an stabilized format, as well as an internal control to monitor PCR inhibition. *Legionella pneumophila* DNA targets are amplified and detected in channel FAM, *Chlamydomphila pneumoniae* DNA in channel ROX, *Mycoplasma pneumoniae* DNA in channel Cy5 and the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

### 4. Reagents provided

VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table 1 and Table 2:

Reference	Reagent/Material	Description	Color	Amount
VS-CML1SL/ VS-CML1SH	<i>C. pneumoniae</i> , <i>M. pneumoniae</i> & <i>L. pneumophila</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	6/12 X 8-well strip
VS-RB02	Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
VS-CML1C	<i>C. pneumoniae</i> , <i>M. pneumoniae</i> & <i>L. pneumophila</i> Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
VS-NC1	Negative control	Non template control	Violet	1 vial x 1 mL
VS-H2O	Water RNase/DNase free	Water RNase/DNase free	White	1 vial x 1 mL
VS-OCS	Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit with Ref. VS-CML106L, VS-CML106H, VS-CML112L and VS-CML112H.

Reference	Reagent/Material	Description	Color	Amount
VS-CML1PL/ VS-CML1PH	<i>C. pneumoniae</i> , <i>M. pneumoniae</i> & <i>L. pneumophila</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1 plate
VS-RB02	Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
VS-CML1C	<i>C. pneumoniae</i> , <i>M. pneumoniae</i> & <i>L. pneumophila</i> Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
VS-NC1	Negative control	Non template control	Violet	1 vial x 1 mL
VS-H2O	Water RNase/DNase free	Water RNase/DNase free	White	1 vial x 1 mL
VS-OCS	Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 X 8-cap strip

Table 2. Reagents and materials provided VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit with Ref VS-CML113L and VS-CML113H.

## 5. Reagents and equipment to be supplied by the user

The following list includes materials that are required for using but not included in the VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler) (to check compatibility see Annex I).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes.
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposal gloves.

## 6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until expiration date stated in the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate in aliquots to minimize freeze and thaw cycles.
- Keep components away from sunlight.

## 7. Precautions for users

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.

- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious as well as all reagents and materials that have been exposed to the samples and must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

## 8. Test procedure

### 8.1. DNA EXTRACTION

Perform the sample preparation according to the recommendations appearing in the instructions for use of extraction kit used.

For DNA extraction from clinical samples you can use your manually or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions for use. We have validated the following extraction kits/systems:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- RIDA® Xtract (r-Biopharm).
- Maxwell® 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® 16 instrument (Promega).

### 8.2. LYOPHILIZED POSITIVE CONTROL

*C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Positive Control contains high copies template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Positive Control (red vial) adding 100 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate in aliquots to minimize freeze and thaw cycles.

### 8.3. PCR PROTOCOL

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run. Peel off protective aluminum seal from plates or strips.

1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Positive Control (red vial) or Negative Control (violet vial) in different wells and close the wells with the caps provided. Centrifuge briefly.

Load the plate or the strips in the thermocycler.

3) Set up your thermocycler.

Program your thermocycler following the conditions below and start the run:

Cycles	Step	Time	Temperature
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 3. PCR protocol

Fluorogenic data should be collected during the extension step (\*) through the FAM (*Legionella pneumophila*), ROX (*Chlamydomphila pneumoniae*), Cy5 (*M. pneumoniae*) and HEX, JOE or VIC channels (Internal Control (IC)). Depending on the equipment used select the proper detection channel (see Annex 2). In Applied Biosystems 7500 Fast Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none.

## 9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in negative control well and the presence of signal for *Legionella pneumophila*, *Chlamydomphila pneumoniae* and *Mycoplasma pneumoniae* positive control well. Check Internal Control signal to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software itself of the used real time PCR equipment according to manufacturer’s instructions.

Using the following table read and analyze the results:

<i>L. pneumophila</i> (FAM)	<i>C. pneumoniae</i> (ROX)	<i>M. pneumoniae</i> (Cy5)	Internal control (HEX)	Negative Control	Positive Control	Interpretation
+	+	+	+/-	-	+	<i>L. pneumophila</i> , <i>C. pneumoniae</i> and <i>M. pneumoniae</i> Positive
-	-	-	+	-	+	<i>L. pneumophila</i> , <i>C. pneumoniae</i> and <i>M. pneumoniae</i> Negative
+	-	-	+/-	-	+	<i>L. pneumophila</i> Positive, <i>C. pneumoniae</i> and <i>M. pneumoniae</i> Negative
+	+	-	+/-	-	+	<i>L. pneumophila</i> and <i>C. pneumoniae</i> Positive, and <i>M. pneumoniae</i> Negative
+	-	+	+/-	-	+	<i>L. pneumophila</i> and <i>M. pneumoniae</i> Positive, and <i>C. pneumoniae</i> Negative
-	+	-	+/-	-	+	<i>C. pneumoniae</i> Positive, <i>L. pneumophila</i> and <i>M. pneumoniae</i> Negative
-	+	+	+/-	-	+	<i>C. pneumoniae</i> and <i>M. pneumoniae</i> Positive, <i>L. pneumophila</i> Negative
-	-	+	+/-	-	+	<i>M. pneumoniae</i> Positive, <i>L. pneumophila</i> and <i>C. pneumoniae</i> Negative
-	-	-	-	-	+	Experiment fail
+	+	+	+	+	-	Experiment fail

Table 4. Sample interpretation


+: Amplification curve

-: No amplification curve

A sample is considered positive if the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

## 10. Limitations of the test

- The result of the test should be evaluated by a health care professional and evaluated in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated with urine, sputum and bronchoalveolar lavage (BAL) specimens.
- The quality of the test depends on the quality of the sample; proper DNA from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Legionella pneumophila*, *Chlamydomphila pneumoniae* and/or *M. pneumoniae*, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

## 11. Quality control

VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

## 12. Performance characteristics

### 12.1. CLINICAL SENSITIVITY AND SPECIFICITY

VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit was evaluated with 20 clinical specimens (urine, sputum and bronchoalveolar lavage) from symptomatic patients. *Legionella pneumophila* could be detected in 9 samples. Eleven specimens resulted negative for *Chlamydomphila pneumoniae*, *Mycoplasma pneumoniae* and *Legionella pneumophila*. The results were compared with those obtained by a commercial Real Time PCR Kit (RIDA®GENE Legionella (R-biopharm)), showing a total concordance.

The results show a high sensitivity and specificity to detect *Chlamydomphila pneumoniae*, *Mycoplasma pneumoniae* and *Legionella pneumophila* using VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit.

### 12.2. ANALYTICAL SENSITIVITY

VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit has a detection limit of  $\geq 10$  DNA copies per reaction for *Chlamydomphila pneumoniae*, *Mycoplasma pneumoniae* and *Legionella pneumophila* (Figure 2, 3 and 4).

Figure 2. Dilution series of *Legionella pneumophila* ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (FAM channel).


Figure 3. Dilution series of *Chlamydomphila pneumoniae* ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (ROX channel).



Figure 4. Dilution series of *Mycoplasma pneumoniae* ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Cy5 channel).

### 12.3. ANALYTICAL SPECIFICITY

The specificity of the *C. pneumoniae*, *M. pneumoniae* and *L. pneumophila* assay was confirmed by testing a panel consisting of different microorganisms representing the most common respiratory pathogens.

No cross-reactivity of *Chlamydomphila pneumoniae* was detected against any of the following microorganisms tested:

Cross-reactivity testing					
<i>Bordetella pertussis</i>	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-	Respiratory syncytial virus (RSV)	-
<i>Legionella bozemanii</i>	-	Methicillin-resistant <i>Staphylococcus aureus</i>	-	Human adenovirus	-
<i>Legionella micdadei</i>	-	<i>Haemophilus influenzae</i>	-	Human coronavirus 229E	-
<i>Legionella dumoffii</i>	-	<i>Moraxella catarrhalis</i>	-	Human rhinovirus	-
<i>Legionella longbeachae</i>	-	Influenza B/Brisbane/60/2008 virus	-	MERS Coronavirus	-
<i>Legionella pneumophila</i>	-	Influenza B/Florida/04/06 virus	-	Influenza A/California/7/2009(H1N1) virus	-
<i>Mycoplasma pneumoniae</i>	-	Influenza B/Phuket/3073/2013	-	Influenza A/Perth/16/2009(H3N2) virus	-
<i>Streptococcus pneumoniae</i>	-	Human parainfluenza 1, 2, 3 and 4 viruses	-	Influenza A/New Caledonia/20/99(H1N1) virus	-
Influenza A/Switzerland/9715293/2013	-	Influenza A/Turkey/Germany R2485+86/2014	-	Human bocavirus	-
Human Metapneumovirus A and B	-				

Table 5. Reference pathogenic microorganisms used in this study

No cross-reactivity of *Mycoplasma pneumoniae* was detected against any of the following microorganisms tested:

Cross-reactivity testing					
<i>Bordetella pertussis</i>	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-	Respiratory syncytial virus (RSV)	-
<i>Legionella bozemanii</i>	-	Methicillin-resistant <i>Staphylococcus aureus</i>	-	Human Metapneumovirus A and B	-
<i>Legionella micdadei</i>	-	<i>Haemophilus influenzae</i>	-	Human coronavirus 229E	-
<i>Legionella dumoffii</i>	-	<i>Moraxella catarrhalis</i>	-	Human rhinovirus	-
<i>Legionella longbeachae</i>	-	Influenza B/Brisbane/60/2008 virus	-	MERS Coronavirus	-
<i>Legionella pneumophila</i>	-	Influenza B/Florida/04/06 virus	-	Influenza A/California/7/2009(H1N1) virus	-
<i>Chlamydomphila pneumoniae</i>	-	Influenza B/Phuket/3073/2013	-	Influenza A/Perth/16/2009(H3N2) virus	-
<i>Streptococcus pneumoniae</i>	-	Human parainfluenza 1, 2, 3 and 4 viruses	-	Influenza A/New Caledonia/20/99(H1N1) virus	-
Influenza A/Switzerland/9715293/2013	-	Influenza A/Turkey/Germany R2485+86/2014	-	Human adenovirus	-
Human bocavirus	-				

Table 6. Reference pathogenic microorganisms used in this study.

No cross-reactivity of *Legionella pneumophila* was detected against any of the following microorganisms tested:

Cross-reactivity testing					
<i>Bordetella pertussis</i>	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-	Respiratory syncytial virus (RSV)	-
<i>Legionella bozemanii</i>	-	Methicillin-resistant <i>Staphylococcus aureus</i>	-	Human Metapneumovirus A and B	-
<i>Legionella micdadei</i>	-	<i>Haemophilus influenzae</i>	-	Human coronavirus 229E	-
<i>Legionella dumoffii</i>	-	<i>Moraxella catarrhalis</i>	-	Human rhinovirus	-
<i>Legionella longbeachae</i>	-	Influenza B/Brisbane/60/2008 virus	-	MERS Coronavirus	-
<i>Chlamydomphila pneumoniae</i>	-	Influenza B/Florida/04/06 virus	-	Influenza A/California/7/2009(H1N1) virus	-
<i>Mycoplasma pneumoniae</i>	-	Influenza B/Phuket/3073/2013	-	Influenza A/Perth/16/2009(H3N2) virus	-
<i>Streptococcus pneumoniae</i>	-	Human parainfluenza 1, 2, 3 and 4 viruses	-	Influenza A/New Caledonia/20/99(H1N1) virus	-
Influenza A/Switzerland/9715293/2013	-	Influenza A/Turkey/Germany R2485+86/2014	-	Human bocavirus	-
Human adenovirus	-				

Table 7. Reference pathogenic microorganisms used in this study.

#### 12.4. ANALYTICAL REACTIVITY

The reactivity of VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit for *Chlamydomphila pneumoniae* was evaluated against *Chlamydomphila pneumoniae* showing positive results.

The reactivity of VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit for *Mycoplasma pneumoniae* was evaluated Mycoplasma pneumoniae showing positive results.

The reactivity of VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit for *L. pneumophila* was evaluated against *Legionella pneumophila* sg1, sg3 and sg6, showing positive results.

**ANNEX 1:****COMPATIBILITY OF THE MOST COMMON REAL TIME PCR EQUIPMENT**

Low profile strips can be used in all PCR thermocyclers equipped with low profile block, like systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with high or regular profile block, like systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Applied Biosystems	7500 Fast Real-Time PCR System
Applied Biosystems	7500 Fast Dx Real-Time PCR System
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
Bio-Rad	CFX96 Touch™ Real-Time PCR Detection System
Roche	LightCycler® 480 Real-Time PCR System
Roche	LightCycler® 96 Real-Time PCR System
Agilent Technologies	AriaMx Real-Time PCR System
DNA-Technology	DTlite Real-Time PCR System
DNA-Technology	DTprime Real-time Detection Thermal Cycler
Qiagen	Rotor-Gene® Q*
Cepheid	SmartCycler®

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Bio-Rad	CFX96 Touch™ Deep Well Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ 5 Real-Time PCR Detection System
Eppendorf	Mastercycler™ep <i>realplex</i>
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
Abbott	Abbott m2000 RealTime System
BIONEER	Exicycler™ 96
DNA-Technology	DTlite Real-Time PCR System
DNA-Technology	DTprime Real-time Detection Thermal Cycler
Qiagen	Rotor-Gene® Q*
Cepheid	SmartCycler®

Table A1/A2. Compatible low and high profile Real Time PCR systems.

\* The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into specific Rotor-Gene® Q or SmartCycler® tubes.

ANNEX 2:

**DETECTION CHANNELS OF MOST COMMON REAL TIME PCR EQUIPMENT**

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

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### 13. Bibliography/Bibliografía

1. D.A. Wilson *et al.* Detection of Legionella pneumophila by Real-Time PCR for the mip Gene. *Journal of Clinical Microbiology* 2003; 41(7):3327-3330.
2. K.A. Thurman *et al.* Detection of Mycoplasma pneumoniae, Chlamydomphila pneumoniae, and Legionella spp. in clinical specimens using a single-tube multiplex real-time PCR assay. *Diagnostic Microbiology and Infectious Disease* 2011; 70(1):1-9.
3. T.P. Atkinson *et al.* Epidemiology, clinical manifestations, pathogenesis and laboratory detection of Mycoplasma pneumoniae infections. *FEMS Microbiology Reviews* 2008; 32(6):956-973.
4. M.H. diaz *et al.* Detection of Mycoplasma pneumoniae and Chlamydomphila pneumoniae directly from respiratory clinical specimens using a rapid real-time polymerase chain reaction assay. *Diagnostic Microbiology and Infectious Disease* 2012; 73(3):278-280.
5. M. Dunne *et al.* Laboratory Tests for Legionnaire’s Disease. *Infectious Disease Clinics in North America* 2016; 31(1):167-178.
6. Centers for Disease Control and Prevention (<https://www.cdc.gov/>)
7. World Health Organization (<http://www.who.int/en/>)

### 14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96 Touch™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, Rotor-Gene® Q (Qiagen) and SmartCycler® (Cepheid). When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

---

VIASURE *C. pneumoniae*, *M. pneumoniae* & *L. pneumophila* Real Time PCR Detection Kit ha sido validado en los siguientes equipos: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96 Touch™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, Rotor-Gene® Q (Qiagen) y SmartCycler® (Cepheid). Cuando se utiliza el equipo Applied Biosystems 7500 Fast con tiras, se recomienda colocar el soporte adecuado para reducir el riesgo de aplastar el tubo (Ref. PN 4388506).

- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™ is registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

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