REAL TIME PCR DETECTION KIT

Respiratory Viral Panel I

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VIASURE <i>Respiratory Viral Panel I</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-RPA106L
VIASURE <i>Respiratory Viral Panel I</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-RPA106H
VIASURE <i>Respiratory Viral Panel I</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-RPA112L
VIASURE <i>Respiratory Viral Panel I</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-RPA112H

VIASURE

ENGLISH

1. Intended of use

VIASURE *Respiratory Viral Panel I* Real Time PCR Detection Kit is designed for simultaneous detection of Influenza A, Influenza B, and Respiratory Syncytial (RSV) viruses and subtyping of Influenza A (H1N1)pdm09, H3N2, H5N1, and H7N9 in respiratory samples from patients with signs and symptoms of respiratory viral infection. This test is intended for use as an aid in the diagnosis of Flu A, Flu B and/or RSV in combination with clinical and epidemiological risk factors and evaluation of infections with Influenza A subtypes: (H1N1)pdm09, H3N2, H5N1, and H7N9. RNA is extracted from specimens, amplified using RT-PCR and detected using fluorescent reporter dye probes specific for Flu A (as well as, (H1N1)pdm09, H3N2, H5N1, and H7N9 subtypes), Flu B and RSV.

2. Summary and Explanation

Influenza virus is an enveloped, single stranded RNA virus that belongs to the *Orthomyxoviridae* family and causes the majority of viral lower respiratory tract infections. Influenza A and B are a significant cause of morbidity and mortality worldwide, considering that elderly and compromised individuals are especially at risk of developing severe illness and complications such as pneumonia. People feel some or all of these symptoms: fever or feeling feverish/chills, cough, sore throat, nasal stuffiness and discharge, myalgia, headaches, and anorexia. The influenza viruses can be spread from person to person in two different ways: through the air (large droplets and aerosols from sneezing and coughing), and by direct or indirect contact.

Influenza A viruses has been classified into subtypes based on the combinations of the envelope proteins hemagglutinin (HA) and neuraminidase (NA). Currently, influenza A(H1N1) and A(H3N2) are the circulating seasonal influenza A virus subtypes. This seasonal A(H1N1) virus is the same virus that caused the 2009 influenza pandemic. Influenza A(H3N2v) viruses were first detected in people in 2011, being its infections associated with prolonged exposure to pigs at agricultural fairs. Besides, humans can also be infected with avian influenza virus subtypes A(H5N1) and A(H7N9) since their emergence in China in 2003 and 2013, respectively. Influenza A(H5N1) virus has spread from Asia to Europe and Africa, and has become entrenched in poultry populations in some countries. Outbreaks have resulted in millions of poultry infections and several hundred human cases, who often display severe pneumonia and greater than 50% mortality. Avian Influenza A(H7N9) virus spreads faster than H5N1 and commonly resulted in severe respiratory illness, although its mortality rate (20%) is lower than that attributed to H5N1 virus. Whereas, Influenza B is only divided into 2 antigenically and genetically distinct lineages, Victoria and Yamagata.

Human respiratory syncytial viruses (RSV) belong to the *Paramyxoviridae* family and are the most important viral agents of acute respiratory infections. RSV is an enveloped, nonsegmented, negative, single stranded linear RNA genome virus. Respiratory syncytial virus is a common contributor of respiratory infections causing bronchitis, pneumonia, and chronic obstructive pulmonary infections in people of all ages. People often feel some or all of these symptoms: rhinorrhea, low-grade fever, cough, sore throat, headache, and wheezing. RSV is transmitted

via large nasopharyngeal secretion droplets from infected individuals, close contact, or self-inoculation after touching contaminated surfaces.

Diagnosis can be problematic, as a wide range of pathogens can cause acute respiratory infections presenting with similar clinical syndromes. Real-time PCR assays have been shown to be a sensitive and specific diagnostic tool for the detection of Influenza A, Influenza B and RSV viruses, as well as, of Influenza A subtyping.

3. Principle of the procedure

VIASURE *Respiratory Viral Panel I* Real Time PCR Detection Kit is designed for the diagnosis of Influenza A, Influenza B, and RSV viruses and subtyping of Influenza A (H1N1)pdm09, H3N2, H5N1, and/or H7N9 in clinical samples. The detection is done in one step real time RT format where the reverse transcription and the subsequent amplification of specific target sequence occur in the same reaction well. The isolated RNA target is transcribed generating complementary DNA by reverse transcriptase which is followed by the amplification of a conserved region of the *M1* gene for Flu A and Flu B and the *N* gene for RSV (Flu A, Flu B & RSV, strip 1), as well as, the *hemagglutinin* gene for subtyping of Influenza A ((H1N1)pdm09, H3N2, H5N1, and/or H7N9) (Flu Typing II, strip 2), using specific primers and a fluorescent-labelled probes.

VIASURE *Respiratory Viral Panel I* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE *Respiratory Viral Panel I* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase, retrotranscriptase) in an stabilized format, as well as an internal control to monitor PCR inhibition. Each kit includes two kind of strips and each one corresponds to one different assay. The first strip contains the multiplex reaction mix for the detection of Flu A, Flu B and RSV (VS-ABR1SL/V S-ABR1SH *Flu A, Flu B & RSV* 8-well strips). Influenza A RNA targets are amplified and detected in channel FAM, Influenza B RNA in channel ROX, and RSV RNA in channel Cy5. The second strip contains the multiplex reaction mix for the subtyping of Influenza A ((H1N1)pdm09, H3N2, H5N1, and/or H7N9) (VS-HXN1SL/ VS-HXN1SH *Flu Typing II* 8-well strips). (H1N1)pdm09 RNA targets are amplified and detected in channel FAM, H5N1 RNA in HEX, VIC or JOE channel (depending on the equipment), H3N2 RNA in channel ROX, and H7N9 RNA in channel Cy5. The first strip (VS-ABR1SL/V S-ABR1SH *Flu A, Flu B & RSV* 8-well strips) contains an internal control (IC) which is amplified and detected in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

4. Reagents provided

VIASURE *Respiratory Viral Panel I* Real Time PCR Detection Kit includes the following materials and reagents detailed in Table 1:

Reference	Reagent/Material	Description	Colour	Amount
VS-ABR1SL/ VS-ABR1SH	<i>Flu A, Flu B & RSV</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	3/6 x 8-well strip
VS-HXN1SL/ VS-HXN1SH	<i>Flu Typing II</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	3/6 x 8-well strip
VS-RB02	Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
VS-RP01C	Respiratory Viral Panel Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
VS-NC1	Negative control	Non template control	Violet	1 vial x 1 mL
VS-H2O	Water RNase/DNase free	Water RNase/DNase free	White	1 vial x 1 mL
VS-OCS	Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE *Respiratory Viral Panel I* Real Time PCR Detection Kit with Ref. VS-RPA106L, VS-RPA106H, VS-RPA112L and VS-RPA112H.

5. Reagents and equipment to be supplied by the user

The following list includes materials that are required for using but not included in the VIASURE *Respiratory Viral Panel I* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler) (to check compatibility see Annex I).
- RNA extraction kit.
- Centrifuge for 1.5 mL tubes.
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposal gloves.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until expiration date stated in the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate in aliquots to minimize freeze and thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- Make sure to use a well for the detection of Flu A, Flu B and RSV and another well for the subtyping assay. Be careful not to mix them throughout the process.

- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious as well as all reagents and materials that have been exposed to the samples and must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

8. Test procedure

8.1. RNA EXTRACTION

Perform the sample preparation according to the recommendations appearing in the instructions for use of extraction kit used.

For RNA extraction from respiratory samples you can use your manually or automatic routine optimized system. Also, you can use any commercially available RNA extraction kit and follow the manufacturer's instructions for use. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- RIDA® Xtract (r-Biopharm).
- Maxwell® 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® 16 instrument (Promega).
- Total Nucleic Acid Isolation (TNAI) Kit, using COBAS® AmpliPrep (ROCHE).

8.2. LYOPHILIZED POSITIVE CONTROL

Respiratory Viral Panel Positive Control contains high copies template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Respiratory Viral Panel* Positive Control (red vial) adding 100 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate in aliquots to minimize freeze and thaw cycles.

8.3. PCR PROTOCOL

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

1) **Reconstitute the number of wells you need.**

Add 15 µL of Rehydration Buffer (blue vial) into each well.

2) Adding samples and controls.

Add 5 µL of RNA sample, reconstituted *Respiratory Viral Panel* Positive Control (red vial) or Negative Control (violet vial) in different wells and close the wells with the caps provided. Centrifuge briefly.

Load the plate or the strips in the thermocycler.

3) Set up your thermocycler.

Program your thermocycler following the conditions below and start the run:

Cycles	Step	Time	Temperature
1	Reverse transcription	15 min	45°C
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 2. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (Influenza A and H1N1), HEX, JOE or VIC (Internal Control (IC) and H5N1), ROX (Influenza B and H3N2), and Cy5 channels (RSV and H7N9). Depending on the equipment used select the proper detection channel (see Annex 2). In Applied Biosystems 7500 Fast Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none.

9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in negative control well and the presence of signal for *Respiratory Viral Panel* positive control well. The analysis of the samples is done by the software itself of the used real time PCR equipment according to manufacturer’s instructions.

Interpretation of results for VS-ABR1SL/VS-ABR1SH Flu A, Flu B & RSV 8-well strips:

- A sample is considered positive for **Influenza A** if there is an amplification signal in FAM channel, the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.
- A sample is considered positive for **Influenza B** if there is an amplification signal in ROX channel, the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.
- A sample is considered positive for **RSV** if there is an amplification signal in CY5 channel, the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.

- A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive (HEX channel). An inhibition of the PCR reaction can be excluded by the amplification of internal control.
- The experiment is considered failed if there is an amplification signal in the Negative Control well and/or there is not amplification signal in the Positive Control well. We recommend to repeat the assay again.
- In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

Interpretation of results for VS-HXN1SL/ VS-HXN1SH Flu Typing II 8-well strips:

- A sample is considered positive for **(H1N1)pdm09** if there is an amplification signal in FAM channel and the Ct value obtained is less than 40.
- A sample is considered positive for **H5N1** if there is an amplification signal in HEX/VIC/JOE channel and the Ct value obtained is less than 40.
- A sample is considered positive for **H3N2** if there is an amplification signal in ROX channel and the Ct value obtained is less than 40.
- A sample is considered positive for **H7N9** if there is an amplification signal in Cy5 channel and the Ct value obtained is less than 40.
- A sample is considered negative, if the sample shows no amplification signal in the detection system.
- The experiment is considered failed if there is an amplification signal in the Negative Control well and/or there is not amplification signal in the Positive Control well. We recommend to repeat the assay again.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Flu A, Flu B & RSV).



Figure 2. Correct run of negative and positive control run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Flu Typing II).

10. Limitations of the test

- The result of the test should be evaluated by a health care professional and evaluated in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated with throat samples dissolved in transport medium.
- The quality of the test depends on the quality of the sample; proper RNA from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Respiratory Viral Panel*, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

11. Quality control

VIASURE *Respiratory Viral Panel I* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well included in the strip 1 (VS-ABR1SL/Vs-ABR1SH Flu A, Flu B & RSV 8-well strip) confirms the correct performance of the technique.

12. Performance characteristics

12.1. CLINICAL SENSITIVITY AND SPECIFICITY

The clinical performance of VIASURE *Respiratory Viral Panel I* Real Time PCR Detection Kit was tested using 256 or 100 respiratory specimens (throat swabs) from symptomatic patients. These results were compared with those obtained by molecular detection method (CLART® PneumoVir DNA array” assay (Genomica)).

The results were as follows:

VIASURE Respiratory Viral Panel I Real Time PCR Detection Kit	CLART® PneumoVir DNA array assay (Genomica)			
		+	-	Total
	+	105	0	105
	-	2*	149	151
	Total	107	149	256

Table 3. Comparative results for Flu A.

*The low amount of template RNA in these respiratory samples is below the detection limit of the method used.

VIASURE Respiratory Viral Panel I Real Time PCR Detection Kit	CLART® PneumoVir DNA array assay (Genomica)			
		+	-	Total
	+	51	0	51
	-	1*	204	205
	Total	52	204	256

Table 4. Comparative results for Flu B.

*The low amount of template RNA in this respiratory sample is below the detection limit of the method used.

VIASURE Respiratory Viral Panel I Real Time PCR Detection Kit	CLART® PneumoVir DNA array assay (Genomica)			
		+	-	Total
	+	41	0	41
	-	0	215	215
	Total	41	215	256

Table 5. Comparative results for RSV.

VIASURE Respiratory Viral Panel I Real Time PCR Detection Kit	CLART® PneumoVir DNA array assay (Genomica)			
		+	-	Total
	+	54	1*	55
	-	5*	40	45
	Total	59	41	100

Table 6. Comparative results for (H1N1)pdm09.

*The low amount of template RNA in this respiratory sample is below the detection limit of the method used.

VIASURE Respiratory Viral Panel I Real Time PCR Detection Kit	CLART® PneumoVir DNA array assay (Genomica)			
		+	-	Total
	+	38	0	38
	-	3*	59	62
	Total	41	59	100

Table 7. Comparative results for H3N2.

*The low amount of template RNA in this respiratory sample is below the detection limit of the method used.

The results show a high sensitivity and specificity to detect Flu A, Flu B and RSV, as well as, Influenza A ((H1N1)pdm09, H3N2) subtypes using VIASURE Respiratory Viral Panel I Real Time PCR Detection Kit.

12.2. ANALYTICAL SENSITIVITY

VIASURE Respiratory Viral Panel I Real Time PCR Detection Kit has a detection limit of ≥ 50 RNA copies per reaction for H5N1 (Figure 7) and ≥ 10 RNA copies per reaction for Flu A, Flu B, RSV, (H1N1)pdm09, H3N2 and H7N9 (Figure 3, 4, 5, 6, 8 and 9)

Figure 3. Dilution series of Influenza A (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Multiplex reaction mix Flu A, Flu B & RSV, channel FAM).


Figure 4. Dilution series of Influenza B (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Multiplex reaction mix *Flu A*, *Flu B* & *RSV*, channel ROX).


Figure 5. Dilution series of RSV (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Multiplex reaction mix *Flu A*, *Flu B* & *RSV*, channel Cy5).


Figure 6. Dilution series of Influenza A(H1N1)pdm09 (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Multiplex reaction mix *Flu Typing II*, channel FAM).


Figure 7. Dilution series of Influenza A(H5N1) (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Multiplex reaction mix *Flu Typing II*, channel HEX).


Figure 8. Dilution series of Influenza A(H3N2) (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Multiplex reaction mix *Flu Typing II*, channel ROX).



Figure 9. Dilution series of Influenza A(H7N9) (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96 Touch™ Real-Time PCR Detection System (Multiplex reaction mix *Flu Typing II*, channel Cy5).

12.3. ANALYTICAL SPECIFICITY

The specificity of the Respiratory Viral Panel I assay was confirmed by testing a panel consisting of 27 microorganisms representing the most common respiratory pathogens. No cross-reactivity was detected against any of the following microorganisms tested, except the targeted pathogens of each assay:

Cross-reactivity testing					
<i>Bordetella pertussis</i>	-	Methicillin-resistant <i>Staphylococcus aureus</i>	-	Influenza B/Florida/04/06 virus	-/+
<i>Legionella bozemanii</i>	-	<i>Haemophilus influenzae</i>	-	Influenza B/Phuket/3073/2013 virus	-/+
<i>Legionella micdadei</i>	-	<i>Moraxella catarrhalis</i>	-	Human parainfluenza 1, 2, 3 and 4 viruses	-
<i>Legionella dumoffii</i>	-	Influenza A/California/7/2009(H1N1) virus	-/+	Human metapneumovirus A and B	-
<i>Legionella longbeachae</i>	-	Influenza A/Perth/16/2009(H3N2) virus	-/+	Human coronavirus 229E	-
<i>Legionella pneumophila</i>	-	Influenza A/New Caledonia/20/99(H1N1) virus	-/+	MERS Coronavirus	-
<i>Mycoplasma pneumoniae</i>	-	Influenza A/Switzerland/9715293/2013 (H3N2) virus	-/+	Human rhinovirus	-
<i>Streptococcus pneumoniae</i>	-	Influenza A/Turkey/Germany R2485+86/2014 virus	-/+	Human Adenovirus	-
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-	Influenza B/Brisbane/60/2008-like virus	-/+	Respiratory Syncytial virus RSV	-/+

Table 8. Reference pathogenic microorganisms used in this study.

12.4. ANALYTICAL REACTIVITY

The reactivity of VIASURE *Respiratory Viral Panel I* Real Time PCR Detection Kit for Flu A was evaluated against strains: A/California/7/2009(H1N1)pdm09-like virus, A/Perth/16/2009(H3N2)-like virus, A/New Caledonia/20/99(H1N1), A/Switzerland/9715293/2013 and A/Turkey/Germany R2485+86/2014 showing positive results.

The reactivity of VIASURE *Respiratory Viral Panel I* Real Time PCR Detection Kit for Flu B was evaluated against strains: B/Brisbane/60/2008-like virus (B/Victoria lineage), B/Florida/04/06 and B/Phuket/3073/2013 (B/Yamagata lineage), showing positive results.

The reactivity of VIASURE *Respiratory Viral Panel I* Real Time PCR Detection Kit for RSV was evaluated against strains: Human Respiratory Syncytial Virus (RSV) A and B, showing positive results.

The reactivity of VIASURE *Respiratory Viral Panel I* Real Time PCR Detection Kit for Influenza A(H1N1)pdm09 was evaluated against strains: A/California/7/2009(H1N1)pdm09-like virus, showing positive results.

The reactivity of VIASURE *Respiratory Viral Panel I* Real Time PCR Detection Kit for Influenza A(H3N2) was evaluated against strains: A/Perth/16/2009(H3N2)-like virus and A/Switzerland/9715293/2013 virus, showing positive results.

ANNEX 1:

COMPATIBILITY OF THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with low profile block, like systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with high or regular profile block, like systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Applied Biosystems	7500 Fast Real-Time PCR System
Applied Biosystems	7500 Fast Dx Real-Time PCR System
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
Bio-Rad	CFX96 Touch™ Real-Time PCR Detection System
Roche	LightCycler® 480 Real-Time PCR System
Roche	LightCycler® 96 Real-Time PCR System
Agilent Technologies	AriaMx Real-Time PCR System
DNA-Technology	DTlite Real-Time PCR System
DNA-Technology	DTprime Real-time Detection Thermal Cycler
Qiagen	Rotor-Gene® Q*
Cepheid	SmartCycler®*

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Bio-Rad	CFX96 Touch™ Deep Well Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ 5 Real-Time PCR Detection System
Eppendorf	Mastercycler™ ep <i>realplex</i>
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
Abbott	Abbott m2000 RealTime System
BIONEER	Exicycler™ 96
DNA-Technology	DTlite Real-Time PCR System
DNA-Technology	DTprime Real-time Detection Thermal Cycler
Qiagen	Rotor-Gene® Q*
Cepheid	SmartCycler®*

Table A1/A2. Compatible low and high profile Real Time PCR systems.

* The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into specific Rotor-Gene® Q or SmartCycler® tubes.

ANNEX 2

DETECTION CHANNELS OF MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

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
13. Bibliography/Bibliografía

1. G. Neumann *et al.* Transmission of Influenza A viruses. *Virology* 2015; 234-246.
2. W.P. Glezen *et al.* The burden of Influenza B: A structured literature review. *American Journal of Public Health* 2013; 103(3): 43-51.
3. Y. Yang *et al.* Simultaneous typing and HA/NA subtyping of influenza A and B viruses including the pandemic influenza A/H1N1 2009 by multiplex real-time RT-PCR. *Journal of Virological Methods* 2010; 167(1): 37-44.
4. R.L. Kuo *et al.* Influenza A/B virus detection and influenza A virus subtyping with emphasis on the novel H7N9 virus by using multiplex real-time RT-PCR. *Journal of Virological Methods* 2014; 208:41-46.
5. K. Liderot *et al.* Secondary bacterial infections in patients with seasonal influenza A and pandemic H1N1. *BioMed Research International* 2013; 2013:376219.
6. Fan *et al.* Detection of a novel avian influenza A (H7N9) virus in humans by multiplex one-step real-time RT-PCR assay. *BMC Infectious Diseases* 2014; 14: 541.
7. World Health Organization. WHO information for molecular diagnosis of influenza virus—update. Available: http://www.who.int/influenza/gisrs_laboratory/molecular_diagnosis/en/. Accessed 2015 Dec 30.
8. S. Subhash Bawage *et al.* Recent Advances in Diagnosis, Prevention, and Treatment of Human Respiratory Syncytial Virus. *Advances in Virology* 2013.
9. C. E. French *et al.* Risk of nosocomial respiratory syncytial virus infection and effectiveness of control measures to prevent transmission events: a systematic review. *Influenza and Other Respiratory Viruses* 2016;10(4): 268-290.
10. F. de-Paris *et al.* Optimization of one-step duplex real-time RT-PCR for detection of influenza and respiratory syncytial virus in nasopharyngeal aspirates. *Journal of Virological Methods* 2012; 186(1-2): 189- 192.
11. A. Hu *et al.* Simultaneous detection, subgrouping, and quantitation of respiratory syncytial virus A and B by real-time PCR. *Journal of Clinical Microbiology* 2003; 41(1): 149-154.

14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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VIASURE *Respiratory Viral Panel I* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System, Bio-Rad CFX96 Touch™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, Rotor-Gene® Q (Qiagen) and SmartCycler® (Cepheid). When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

VIASURE *Respiratory Viral Panel I* Real Time PCR Detection Kit ha sido validado en los siguientes equipos: Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System, Bio-Rad CFX96 Touch™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, Rotor-Gene® Q (Qiagen) y SmartCycler® (Cepheid). Cuando se utiliza el equipo Applied Biosystems 7500 Fast con tiras, se recomienda colocar el soporte adecuado para reducir el riesgo de aplastar el tubo (Ref. PN 4388506).

- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

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