CALPROTECTIN TURBILATEX Combo

Calprotectin
latex turbidimetric assay

Cat.TL-022CP100ED, TL-022CP200ED, TL-022CP400ED

CERTEST BIOTEC S.L.

For information purposes only

ENGLISH

INTENDED USE

Calprotectin Turbilatex is a latex turbidimetric assay **only for the quantitative detection of calprotectin in human solid stool samples** (not to be used for body fluid as blood, serum, plasma, urine, cerebrospinal fluid, oral fluid, synovial fluid or empyema fluid).

This assay is simple and widely applicable. This product is optimized for several automated analyser.

Test results should **exclusively be used to differentiate IBD patients with inflammation from IBD patients without inflammation and from irritable bowel syndrome (IBS).**

For professional *in vitro* diagnostic use only.

SYNTHESIS

Calprotectin is a neutrophil cytosolic protein with antimicrobial properties, which is present at increased concentration in stool samples during bowel inflammation. The stability of the protein to degradation keeps it stable in faeces for up to 7 days at room temperature, making it an ideal analyte. Calprotectin is released by activation of leukocytes, giving increased levels in plasma, cerebral spinal fluid, synovial fluid, urine or stools as a consequence of disease in the relevant organ(s). Calprotectin inhibits zinc-dependent enzyme systems, as a result kills microbes and induces apoptosis in normal and cancer cells. In the presence of calcium, calprotectin is remarkably resistant to proteolytic degradation and so is stable in stools kept at room temperature for 7 days.

PRINCIPLE

Calprotectin Turbilatex is a latex turbidimetric assay **only for the quantitative detection of Calprotectin in human solid stool samples.** This assay is simple, widely applicable and test results are obtained in a very short time.

The intended use of the test is **exclusively to differentiate IBD patients with inflammation from IBD patients without inflammation and from irritable bowel syndrome (IBS).**

Calprotectin latex turbidimetric assay is based on agglutination reactions. These involve *in vitro* aggregation of microscopic latex particles. This aggregation consists in the specific reaction between antigen and antibodies, antigen contained in the sample and the antibodies anti-antigen coated on polystyrene latex particles. The sample is mixed with a suspension containing antibodies against the antigen bound to latex particles. If antigen is present in the sample it will react with the antibodies and form an aggregate. If no antigen is present in the sample the mixture will keep its appearance as a smooth suspension. Such turbidity is measured as an increase in absorbance at the determinate wave and is proportional to the quantity of antigen contained in the sample.

PRECAUTIONS

- For professional *in vitro* diagnostic use only.
- A trained person in Turbidimetric technique and autoanalyzer use is required.
- The analyser should be ready before performing any assay.
- Read and follow the instructions for use provided with the kit.
- If the result exceeds the measurement range, use the sample diluent to dilute the sample and repeat the assay again.
- Do not use after expiration date.
- Do not use the reagents if pack is damaged or opened.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The solutions should be discarded in a proper container after testing following local regulations.

STORAGE AND STABILITY

Kit components must be stored at refrigerated temperature (2-8°C). Do not freeze.

Reagents are stable up to the expiration date printed on the label, always considering that reagent containers must be properly closed to avoid any contamination, kept away from the sunlight and conserved at 2-8°C.

REAGENTS

Materials provided**Cat.: TL-022CP100ED**

- TL-022CP01 Calprotectin Turbilatex Reagent 1, (1 x 27mL)
- TL-022CP02 Calprotectin Turbilatex Reagent 2, (1 x 5mL)
- TL-022CP70 Calprotectin Turbilatex Calibrator 0, (1 x 1mL)
- TL-022CP71 Calprotectin Turbilatex Calibrator 1, (1 x 1mL)
- TL-022CP72 Calprotectin Turbilatex Calibrator 2, (1 x 1mL)
- TL-022CP73 Calprotectin Turbilatex Calibrator 3, (1 x 1mL)
- TL-022CP74 Calprotectin Turbilatex Calibrator 4, (1 x 1mL)
- TL-022CP75 Calprotectin Turbilatex Calibrator 5, (1 x 1mL)
- TL-022CP08 Calprotectin Turbilatex Control 1, (2 x 1mL)
- TL-022CP09 Calprotectin Turbilatex Control 2, (2 x 1mL)

Cat.: TL-022CP200ED

- TL-022CP01 Calprotectin Turbilatex Reagent 1, (2 x 27mL)
- TL-022CP02 Calprotectin Turbilatex Reagent 2, (1 x 8mL)
- TL-022CP70 Calprotectin Turbilatex Calibrator 0, (1 x 1mL)
- TL-022CP71 Calprotectin Turbilatex Calibrator 1, (1 x 1mL)
- TL-022CP72 Calprotectin Turbilatex Calibrator 2, (1 x 1mL)
- TL-022CP73 Calprotectin Turbilatex Calibrator 3, (1 x 1mL)
- TL-022CP74 Calprotectin Turbilatex Calibrator 4, (1 x 1mL)
- TL-022CP75 Calprotectin Turbilatex Calibrator 5, (1 x 1mL)
- TL-022CP08 Calprotectin Turbilatex Control 1, (2 x 1mL)
- TL-022CP09 Calprotectin Turbilatex Control 2, (2 x 1mL)

Cat.: TL-022CP400ED

- TL-022CP01 Calprotectin Turbilatex Reagent 1, (4 x 27mL)
- TL-022CP02 Calprotectin Turbilatex Reagent 2, (2 x 8mL)
- TL-022CP70 Calprotectin Turbilatex Calibrator 0, (1 x 1mL)
- TL-022CP71 Calprotectin Turbilatex Calibrator 1, (1 x 1mL)
- TL-022CP72 Calprotectin Turbilatex Calibrator 2, (1 x 1mL)
- TL-022CP73 Calprotectin Turbilatex Calibrator 3, (1 x 1mL)
- TL-022CP74 Calprotectin Turbilatex Calibrator 4, (1 x 1mL)
- TL-022CP75 Calprotectin Turbilatex Calibrator 5, (1 x 1mL)
- TL-022CP08 Calprotectin Turbilatex Control 1, (2 x 1mL)
- TL-022CP09 Calprotectin Turbilatex Control 2, (2 x 1mL)

Materials and instruments that could be required but no provided

- Specimen collection container.
- TL-022CP03E Calprotectin Turbilatex Sample Diluent.
- TL-0221213E Microtubes.
- Disposable gloves and laboratory equipment.
- Automated analyser.
- MST-0006MC Sample Dilution Vials MC.
- MST-0008C Sample Dilution Vials C.

SPECIMEN COLLECTION AND PREPARATION

Collect sufficient quantity of human solid stool samples. These samples should be collected in clean and dry containers (no preservatives or transport media). The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. Homogenize stool samples as thoroughly as possible prior to preparation.

To process the collected stool samples:

Procedure A: using catalogue reference TL-022CP03E.

Use a centrifugation tube for each sample to be tested. Label centrifugation tube with name or number of patient.

1. Homogenize the sample. Add 20mg of sample into centrifugation tube.
2. Add 2mL of sample diluent (TL-022CP03E).
3. Shake vigorously the tube in order to assure good sample dispersion (vortex) until sample is completely dissolved. Rest tubes on the bench for 10 minutes to get a proper calprotectin extraction.
4. Centrifuge for 15 minutes at 10000g or 10 minutes at 15000g.
5. Take the supernatant to automated analyser vial.

Procedure B: using catalogue reference MST-0006MC/MST-0008C.

Follow the indication of the instruction for use of the product.

ANALYTICAL PROCEDURE

Allow reagents and stool samples to reach room temperature (15-30°C) prior to testing.

R1 and R2 are ready to use.

Preparation of the calibration curve

For calibration use only the following materials:


TL-022CP70, TL-022CP71, TL-022CP72, TL-022CP73, TL-022CP74 and TL-022CP75, Calprotectin Turbilatex Calibrator vials: liquid calibrator containing recombinant human calprotectin at different concentrations.

Concentration is indicated on the label of each vial.

Calibrate the system at least once a month is extremely recommended. Recalibrate the system when reagent lot is changed or when the controls are out of the assigned range given on the control label or certificate of analysis provided with the kit.

Allow Calprotectin Turbilatex Calibrator to reach room temperature (15-30°C) prior to testing.

Calprotectin Turbilatex Calibrator vials (calibrator 0 to 5) are ready to use.

Picture 1. Prepared calibrators for performing the calibration curve (5 points and blank)

Quality control

For quality control only use the following materials:


TL-022CP08 Calprotectin Turbilatex Control 1 and TL-022CP09 Calprotectin Turbilatex Control 2: liquid controls at two different concentrations of recombinant human calprotectin. Concentration is indicated on the label of the vial. The use of control materials at two different concentrations is recommended in order to verify test precision.

Controls should be assayed every day before running patient faecal samples extract to validate the calibration curve.

If the obtained results are out of the tolerance range, the equipment, the reagents or the technique must be reviewed.

Allow Calprotectin Turbilatex C1 & C2 to reach room temperature (15-30°C) prior to testing.

Calprotectin Turbilatex Control vials (Control 1 and Control 2) are ready to use.

Picture 2. Prepared controls C1 and C2 (check control concentration on the label).

Analytical procedure, see annex information.

INTERPRETATION OF RESULTS

Positive results: higher or equal than the cut-off fixed by the clinical lab.

Recommended: 50 µg of hCp/g of stool for diagnostic procedures and 200 µg of hCp/g of stool for screening procedures.

Positive results determine the abnormal presence of human Calprotectin (hCp) in stool samples.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

Limit of Detection (LOD): **see annex information.**

Calprotectin Turbilatex cut-off value:

Calprotectin concentration values lower than 50 µg of hCp/g of stool are considered normal values and that is not indicative of an inflammation of gastrointestinal tract.

Calprotectin concentration values higher than 200 µg of hCp/g of stool are indicative of an inflammation of gastrointestinal tract.

Prozone

Studies have been made up to a concentration of 8000 µg of hCp/g of stool and no false negative results have been observed. Studies using higher concentrations have not been carried out.

Clinical sensitivity and specificity

An evaluation was performed comparing a turbidimetric assay (Calprotectin Turbilatex, CerTest) and another commercial immunoassay (Calprest®, Eurospital). The results were as follows:

	Sensitivity	Specificity
Calprotectin Turbilatex vs Calprest®	94%	>99%

The results showed a high sensitivity and specificity to detect human calprotectin (hCp) using Calprotectin Turbilatex.

Interferences

An evaluation was performed to determine the interferences of Calprotectin Turbilatex. No interferences were founded against other substances tested: BSA, Hemin, Ascorbic acid.

Cross reactivity

An evaluation was performed to determine the cross reactivity; no cross reactivity was founded against other faecal markers occasionally present in faeces, such as: bovine and pig haemoglobin, bovine and pig transferrin, bovine lactoferrin and human haemoglobin, transferrin and lactoferrin.

LIMITATIONS

1. Calprotectin Turbilatex should **only be used in human solid stool samples (not to be used for body fluid such as blood, serum, plasma, urine, cerebrospinal fluid, oral fluid, synovial fluid or empyema fluid)**. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
2. Positive results determine the presence of hCp in faecal samples. A positive result should be followed up with additional diagnostic invasive procedures, a colonoscopy and a biopsy in order to confirm the diagnosis and to establish the inflammation extent.
3. If symptoms or situation still persist, calprotectin determination should be carried out invasive techniques. Negative results do not exclude IBD with inflammation, some diseases such as celiac sprue and microscopic colitis polyps that mainly involve mononuclear inflammation.
4. Stool samples from patients with non-steroidal anti-inflammatory drug (NSAID) treatment could show positive result.
5. Neonatal faecal calprotectin levels have been reported higher than those in normal children with a mean of 167µg/g (range 22-860 µg/g).

EXPECTED VALUES

Some studies established equal or higher 50 µg hCp/g faeces as cut-off value to allow detect adult patients with gastrointestinal inflammatory pathology that will require to diagnosis additional diagnostic invasive procedures.

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1. Angriman I. et al. Enzymes in feces: Useful markers of chronic inflammatory bowel disease. *Clinica Chimica Acta* 381 Feb 2007, p. 63-68.
2. Quail, M.A. et al. Fecal Calprotectin Complements Routine Laboratory Investigations in Diagnosing Childhood Inflammatory Bowel Disease. *Inflamm Bowel Dis*, Vol 15 No 5; May 2009, p. 756-759.
3. Gaya D.R., et al. Faecal calprotectin in the assessment of Crohn's disease activity. *Q J Med* 2005, Vol 98, May 2005, p. 435-441.
4. Langhorst, M.D. et al. Noninvasive Markers in the Assessment of Intestinal Inflammation in Inflammatory Bowel Diseases: Performance of Fecal Lactoferrin, Calprotectin and PMN-Elastase, CRP, and Clinical Indices. *Am. J. Gastroenterol.* 2008; Vol 103, p. 162-169.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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The application notes for use our turbidimetric products with the different equipments could be consulted in our website:

www.certest.es

For information purposes only

For information purposes only

CerTest
BIOTEC S.L.

Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840, San Mateo de Gállego, Zaragoza (SPAIN)
www.certest.es

F-381 rev00