CERTEST

Entamoeba

ONE STEP
Entamoeba CARD TEST

CERTEST BIOTEC S.L.

CERTEST *Entamoeba*

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ENGLISH

INTENDED USE

CerTest *Entamoeba* one step card test is a coloured chromatographic immunoassay for the qualitative detection of *Entamoeba histolytica* and *Entamoeba dispar* in stool samples.

CerTest *Entamoeba* card test offers a simple and a highly sensitive screening assay to make a presumptive diagnosis of amoebiasis.

INTRODUCTION

Entamoeba histolytica, *E. dispar*, and *E. moshkovskii* are morphologically identical but biochemically and genetically different. These parasites colonize the human gut, but only *E. histolytica* is thought to be capable of causing disease.

Amoebiasis is the infection of the human gastrointestinal tract by *Entamoeba histolytica*, a protozoan parasite that is capable of invading the intestinal mucosa and may spread to other organs, mainly the liver. The acceptance of *E. dispar* as a distinct but closely related protozoan species has had profound implications for the epidemiology of amoebiasis, since most asymptomatic infections found worldwide are now attributed to this noninvasive amoeba. Roughly, 90% of people infected with *Entamoeba* infections have *E. dispar*.

The disease may manifest itself as an acute, chronic or as an asymptomatic infection. The leading symptom is dysentery whose main consequence one galling of the colon, diarrhoea and abdominal pain. Complications of the disease may be hepatic abscesses, pulmonary abscesses or even cerebral abscesses which, if untreated, usually end in death.

Infection occurs by water pollution, vegetables, raw fruits and other food or poorly cooked or washed with infective cysts from contaminated faeces. It is possible that flies and cockroaches carrying cysts from faeces to food

Amoebic dysentery occurs frequently in tropical countries but are also cases in temperate and cold areas. In Africa, tropical Asia and Latin America, more than two thirds of the population has these intestinal parasites, although most can be virtually asymptomatic infections.

TEST PRINCIPLE

CerTest *Entamoeba* is based on the principle of a qualitative immunochromatographic assay for the determination of *Entamoeba histolytica* and *Entamoeba dispar* in stool samples.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Entamoeba* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Entamoeba*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*Entamoeba* monoclonal antibodies-red polystyrene microspheres), which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*Entamoeba* antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no *Entamoeba* antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-*Entamoeba* antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appears. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the CerTest *Entamoeba* test. Do not use any other commercial kit component.


- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool samples as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

1. Take out the cap of the stool collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2), to collect faecal sample (approx. 125mg) and add it to the stool collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

VERY IMPORTANT: Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.

MATERIALS

MATERIALS PROVIDED

- CerTest *Entamoeba* card tests
- Instructions for use
- Stool collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE


Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove CerTest *Entamoeba* card test from its sealed bag just before using it.
3. Take the stool collection tube, cut the end of the cap (4) and dispense 4 drops in the circular window marked with the letter S (5). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS (please refer to the illustration below)

	<i>Entamoeba</i>	Interpretation of the results
1.	- GREEN	There is no <i>Entamoeba histolytica</i> or <i>E. dispar</i> presence. No infection caused by <i>Entamoeba histolytica</i> or <i>E. dispar</i> .
2.	+ GREEN-RED	There is <i>Entamoeba histolytica</i> or <i>E. dispar</i> presence. Infection caused by <i>Entamoeba histolytica</i> or <i>E. dispar</i> .
3.	Any other result	Invalid result, we recommend repeating the assay using the same sample with another test.

NEGATIVE: Only the GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.
3. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
4. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
5. CerTest *Entamoeba* should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
6. After one week of infection, the number of parasites in faeces is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms.
7. Positive results determine the presence of *Entamoeba histolytica* and/or *E. dispar* in faecal samples; nevertheless, a positive result should be followed up with additional laboratory techniques (biochemical methods or microscopy). A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
8. A negative result is not meaningful because of it is possible the antigens concentration in the stool samples is lower than the detection limit value. If the symptoms or situation still persist, an *Entamoeba histolytica* and/or *E. dispar* determination should be carried out with another technique (for example, microscopy).

EXPECTED VALUES

Entamoeba histolytica is the third main cause of morbidity and mortality due to parasitic disease in humans after malaria and schistosomiasis. It is estimated that *E. histolytica* and *E. dispar* parasitize approximately 10% of the world population, of which 90% are asymptomatic infections. About between 50000 and 100000 deaths every year are due to this disease.

In developed countries, amoebiasis tends to be more common in older and immunodepressed individuals. However, in tropical regions, the epidemiology of amoebiasis is completely different and is more common among the general population and particularly among patients attending health care centres with diarrhoea.

Out of 10% of the world's population infected by *E. histolytica*, only 1% develops invasive form of the disease. Only approximately 10% (pathogenic infections) need to be treated. In developed countries, *E. histolytica* infection is not endemic and occurs mostly after traveling to areas of endemicity.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity

An evaluation was performed, using faecal samples, comparing the results obtained by an immunochromatographic test (CerTest *Entamoeba*, CerTest) and other commercial qPCR tests (VIASURE *Entamoeba histolytica* Real Time Detection Kit and VIASURE *Entamoeba dispar* Real Time Detection Kit, CerTest). The results were as follows:

		qPCR test: VIASURE <i>Entamoeba histolytica</i> Real Time Detection Kit		
		+	-	Total
IC test: CerTest <i>Entamoeba</i>	+	4	1	5
	-	1	107	108
	Total	5	108	113

	Mean Value <i>E. histolytica</i>	95% confidence interval
Sensitivity	80%	28.4 - 99.5%
Specificity	99%	94.9 - 100%
PPV	80%	28.4 - 99.5%
NPV	97%	94.9 - 100%

		qPCR test: VIASURE <i>Entamoeba dispar</i> Real Time Detection Kit		
		+	-	Total
IC test: CerTest <i>Entamoeba</i>	+	15	1	16
	-	6	107	113
	Total	21	108	129

	Mean Value <i>E. dispar</i>	95% confidence interval
Sensitivity	71%	47.8 - 88.7%
Specificity	99%	94.9 - 100%
PPV	94%	69.8 - 99.8%
NPV	95%	88.8 - 98.0%

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	Mean Value <i>E. histolytica</i> + <i>E. dispar</i>	95% confidence interval
Sensitivity	73%	52.2 - 88.4%
Specificity	99%	96.7 - 100%
PPV	90%	69.6 - 98.8%
NPV	97%	93.6 - 98.7%

The results showed a high sensitivity and specificity to detect *Entamoeba histolytica* and *E. dispar* using CerTest *Entamoeba*.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *Entamoeba*; no cross reactivity against gastrointestinal pathogens occasionally present in faeces:

<i>Campylobacter jejuni</i>	<i>Helicobacter pylori</i>	<i>Shigella boydii</i>
<i>Campylobacter coli</i>	<i>Listeria monocytogenes</i>	<i>Shigella dysenteriae</i>
<i>Clostridium difficile</i>	<i>Salmonella enteritidis</i>	<i>Shigella flexneri</i>
<i>Escherichia coli</i> O157:H7	<i>Salmonella paratyphi</i>	<i>Shigella sonnei</i>
<i>Cryptosporidium parvum</i>	<i>Salmonella typhi</i>	<i>Staphylococcus aureus</i>
<i>Giardia lamblia</i>	<i>Salmonella typhimurium</i>	

The logo for CerTest BIOTEC is centered on the page. It consists of a solid blue circle. Inside the circle, the word "CerTest" is written in a white, sans-serif font. The "C" is a standard uppercase letter, while the "T" is a stylized, taller uppercase letter. Below "CerTest", the word "BIOTEC" is written in a smaller, white, all-caps, sans-serif font.

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
1. BLESSMANN J., et al. "Longitudinal Study of Intestinal *Entamoeba histolytica* Infections in Asymptomatic Adult Carriers". *Journal of clinical microbiology*, Oct. 2003, p. 4745–4750 Vol. 41, No. 10.
2. GONIN P., et al. "Detection and Differentiation of *Entamoeba histolytica* and *Entamoeba dispar* Isolates in Clinical Samples by PCR and Enzyme-Linked Immunosorbent Assay". *Journal of clinical microbiology*, Jan. 2003, p. 237–241 Vol. 41, No. 1.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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