CERTEST

Legionella


ONE STEP *Legionella pneumophila*
CARD TEST

CERTEST BIOTEC S.L.

For information purposes only

For information purposes only

CERTEST *Legionella*

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ENGLISH

INTENDED USE

CerTest *Legionella* one step card test is a coloured chromatographic immunoassay for the qualitative detection of *Legionella pneumophila* serogroup 1 in human urine samples.

CerTest *Legionella* card test offers a highly sensitive and simple screening assay to make a presumptive diagnosis of legionellosis in infected humans.

INTRODUCTION

Legionnaires' Disease, named after the outbreak in 1976 at the American Legion convention in Philadelphia, is caused by *Legionella pneumophila* and is characterized as an acute febrile respiratory illness ranging in severity from mild illness to fatal pneumonia. The disease occurs in both epidemic and endemic forms and sporadic cases are not easily differentiated from other respiratory infections by clinical symptoms. An estimated 25,000 to 100,000 cases of *Legionella* infection occur in the United States annually. The resulting mortality rate, ranging from 25% to 40%, can be lowered if the disease is diagnosed rapidly and appropriate antimicrobial therapy is instituted early. Known risk factors include immunosuppression, cigarette smoking, alcohol consumption and concomitant pulmonary disease. The young and the elderly are particularly susceptible. *Legionella pneumophila* is responsible for 80-90% of reported cases of Legionella infection with serogroup 1 accounting for greater than 70% of all legionellosis. Current methods for the laboratory detection of pneumonia caused by *Legionella pneumophila* require a respiratory specimen (e.g. expectorated sputum, bronchial washing, transtracheal aspirate, lung biopsy) or paired sera (acute and convalescent) for an accurate diagnosis.

The CerTest *Legionella* allows for early diagnosis of *Legionella pneumophila* serogroup 1 infection through detection of a specific soluble antigen present in the urine of patients with Legionnaires' Disease. *Legionella pneumophila* serogroup 1 antigen has been detected in urine as early as three days after the onset of symptoms. The test is rapid, giving a result within 15 minutes, and utilizes a urine specimen which is convenient for collection, transport, and subsequent detection of early, as well as later, stages of disease.

TEST PRINCIPLE

CerTest *Legionella* is based on the principle of a qualitative immunochromatographic assay for the determination of *Legionella pneumophila* serogroup 1 in urine samples.

The strip consists of a nitrocellulose membrane pre-coated with polyclonal antibodies on the test line (T), in the results window, against *Legionella pneumophila* (*L. pneumophila*) and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (polyclonal antibodies anti-*L. pneumophila*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to blue polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*L. pneumophila* polyclonal antibodies-red polystyrene microspheres), which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*L. pneumophila* antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no *L. pneumophila* antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-*L. pneumophila* antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control blue-conjugate complex and the control line will always appear. The presence of this blue line serves as 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) as an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the CerTest *Legionella* test. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

SPECIMEN COLLECTION AND PREPARATION

Urine specimens should be collected in standard containers. The samples can be stored at room temperature (15-30°C) if assayed within 24 hours of collection. Alternatively, specimens may be stored at 2-8°C for up to 14 days or at -10°C to -20°C for longer periods before testing. When necessary, urine specimens should be shipped in leak-proof containers at 2-8°C or frozen. Allow all specimens to equilibrate to room temperature before testing.

MATERIALS

MATERIALS PROVIDED

- CerTest *Legionella* card tests
- Plastic pipettes
- Testing tubes
- Reagent (sample and controls diluent)
- Positive Control Swab: Inactivated *L. pneumophila* swab + testing tube+ pipette
- Negative Control Swab: *L. pneumophila* negative swab + testing tube+ pipette
- Instructions for use

MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE

Allow tests, urine samples, reagent and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

Procedure A: Patient samples

1. Using a separate testing tube or vial for each sample. Add 3 drops of urine sample (1).
2. Add 1 drop of Reagent into the testing tube or vial and mix (2). Homogenize the sample.
3. Remove CerTest *Legionella* card test from its sealed bag just before using it.
4. Use a separate pipette and device for each sample or control. Dispense 3 drops from the testing tube in the circular window marked with the letter S (3).
5. **Read the results at 15 minutes.** Do not read the test result later than 15 minutes.

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
Procedure B: Positive and Negative Control Swabs.


1. Hold Reagent vertically. Add slowly 13 free falling drops of Reagent into the testing tube (1).
2. Immediately remove Positive control swab from the pouch and put the swab into the testing tube with reagents, mix 1 minute and extract as much liquid possible from the swab, squeezing the sides of the tube as the swab is withdrawn. Discard the swab (2).
3. Remove the CerTest *Legionella* device from its sealed bag just before using.
4. Use a separate pipette and device for each sample or control. Dispense 3 drops from the testing tube, into the circular window marked with the letter S (3).
5. **Read the result at 15 minutes.** Do not read the test result later than 15 minutes.

Repeat the procedure for Negative Control Swab using the Reagent, same used for sample dilution and for Positive Control Swab.

Positive and negative controls should be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.

INTERPRETATION OF RESULTS (please refer to the illustration below)

	<i>Legionella</i>	Interpretation of the results
1.	- BLUE	There is no <i>Legionella pneumophila</i> presence. No infection caused by <i>L. pneumophila</i> serogroup 1. Negative control result.
2.	+ BLUE-RED	There is <i>Legionella pneumophila</i> presence. Infection caused by <i>L. pneumophila</i> serogroup 1. Positive control result.
3.	Any other result	Invalid result, we recommend repeating the assay using the same sample with another test.

NEGATIVE: Only a BLUE line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the BLUE control line, another RED line also appears in the site marked with the letter T (test line) in the results window.

INVALID: Total absence of the control coloured line (BLUE) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

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QUALITY CONTROL

Internal procedural control is included in the test. A blue line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

External Positive and Negative Controls are included in the kit (see procedure B). The use of positive and negative controls is recommended to assure functionality of reagents and proper performance of assay procedure.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. CerTest *Legionella* should be used only with human urine samples. The use of other samples (e.g. plasma, serum or other body fluids or environmental samples (water)) has not been established.
3. Positive results determine the presence of *L. pneumophila* serogroup 1 in urine samples; nevertheless, a positive result should be followed up with additional laboratory techniques (culture or serology) to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
4. A negative result is not meaningful because of it is possible the antigens concentration in the urine sample is lower than the detection limit value. If symptoms or situation still persist, a *L. pneumophila* serogroup 1 determination should be carried out from a culture.
5. Excretion of *Legionella* antigen in urine may vary depending on the individual patient. Antigen excretion may begin as early as 3 days after onset of symptoms and persist for up to 1 year afterwards. A positive CerTest *Legionella* result can occur due to current or past infection and therefore is not definitive for infection without other supporting evidence.
6. It has been proven that although the test is compatible with boric acid, its presence increases the reactivity and can lead to false positives, so it is not recommended as a preservative.

EXPECTED VALUES

More than 50 *Legionella* species have been characterized, and 25 species are known to cause human disease. Most human infections are caused by *Legionella pneumophila*, which is responsible for approximately 90% of the identified clinical cases, and the predominant serogroup is serogroup 1.

In recent years, hot springs have been recognized as a major source of outbreaks of legionellosis. Several outbreaks and a number of cases of *Legionella* infection that were associated with hot springs have been reported in Japan, Spain and France.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit)

Detection limit value of CerTest *Legionella* is: 12.5ng/mL (pool of several serovars of *Legionella pneumophila*).

Clinical sensitivity and specificity

An evaluation was performed, with urine samples, comparing the results obtained by an immunochromatographic test (CerTest *Legionella*, CerTest) and other commercial rapid test (Binax NOW® *Legionella* Urinary Antigen Card, Alere). The results were as follows:

		IC test: Binax NOW® <i>Legionella</i> Urinary Antigen Card		
		+	-	Total
IC test: CerTest <i>Legionella</i>	+	31	1	32
	-	0	117	117
		31	118	149

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CerTest <i>Legionella</i> vs Binax NOW® <i>Legionella</i> Urinary Antigen Card		
	Mean Value	95% confidence interval
Sensitivity	100.00%	88.8 – 100.0%
Specificity	99.20%	95.4 – 100.0%
PPV	96.90%	83.8 - 99.9%
NPV	100.00%	96.9 – 100.0%

The results showed a high sensitivity and specificity to detect *L. pneumophila* serogroup 1 using CerTest *Legionella*.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *Legionella*; no cross reactivity against other pathogens occasionally present in urine:

Streptococcus pneumoniae

Reproducibility Study

Evaluation studies were performed to determine reproducibility of the CerTest *Legionella*, including inter-day, inter-laboratory, inter and intra lot, showing high reproducibility in all cases.

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