---

# FOB TURBILATEX Combo

---

FOB  
latex turbidimetric assay

Cat.TL-022FB100ED, TL-022FB200ED, TL-022FB400ED

---

CERTEST BIOTEC S.L.

---

**CerTest**  
BIOTEC

For information purposes only

For information purposes only

ENGLISH

INTENDED USE

FOB Turbilatex is a latex turbidimetric assay for the quantitative detection of haemoglobin (faecal occult blood) in human stool samples. This assay is simple and widely applicable. This product is optimized for several automated analyser. For professional in vitro diagnostic use only.

SYNTHESIS

Colorectal cancer is the second leading cause of illness and death in Western world. The screening with faecal occult blood tests is based on the concept that important target colonic neoplasm, such as early-stage cancer and large adenomatous polyps. Colorectal cancer is also associated with local acute inflammatory reaction being visualized, in some cases, by white cell neutrophil scanning. Haemoglobin is the iron-containing oxygen-transport protein in the red blood cells of all vertebrates that may be leaked into gastrointestinal tract and then discharged with the faeces in gastrointestinal bleeding diseases. When gastrointestinal blood is lost, the stool will contain a combination of intact or nearly intact haemoglobin, intact heme and heme-derived porphyrins in amounts that depend on the site and amount of bleeding and the transit time through the gut. FOB Turbilatex tests detect intact or nearly intact human haemoglobin, being a very specific technique for detecting loss of blood from the lower intestine.

PRINCIPLE

FOB Turbilatex is a latex turbidimetric assay for the quantitative detection of haemoglobin in human stool samples. This assay is simple, widely applicable and test results are obtained in a very short time. The intended use of the test is to make a presumptive diagnosis of faecal occult blood (gastrointestinal bleeding). FOB latex turbidimetric assay is based on agglutination reactions. These involve in vitro aggregation of microscopic latex particles. This aggregation consists in the specific reaction between antigen and antibodies, antigen contained in the sample and the antibodies anti-antigen coated on polystyrene latex particles. The sample is mixed with a suspension containing antibodies against the antigen bound to latex particles. If antigen is present in the sample it will react with the antibodies and form an aggregate. If no antigen is present in the sample the mixture will keep its appearance as a smooth suspension. Such turbidity is measured as an increase in absorbance at the determinate wave and is proportional to the quantity of antigen contained in the sample.

PRECAUTIONS

- For professional in vitro diagnostic use only.
- The analyser should be ready before performing any assay.
- Read and follow the instructions for use provided with the kit.
- If the result exceeds the measurement range, use the sample diluent to dilute the sample and repeat the assay again.
- Do not use after expiration date.
- Do not use the reagents if pack is damaged or opened.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The solutions should be discarded in a proper container after testing following local regulations.

STORAGE AND STABILITY

Kit components must be stored at refrigerated temperature (2-8°C). Do not freeze. Reagents are stable up to the expiration date printed on the label, always considering that reagent containers must be properly closed to avoid any contamination, must be kept away from the sunlight and conserved at 2-8°C.

REAGENTS

Materials provided

Cat.: TL-022FB100ED

Cat.: TL-022FB200ED

Cat.: TL-022FB400ED

- TL-022FB01 FOB Turbilatex Reagent 1, (1 x 20mL)
- TL-022FB02 FOB Turbilatex Reagent 2, (1 x 6mL)
- TL-022FB70 FOB Turbilatex Calibrator 0, (1 x 0.3mL)
- TL-022FB71 FOB Turbilatex Calibrator 1, (1 x 0.3mL)
- TL-022FB72 FOB Turbilatex Calibrator 2, (1 x 0.3mL)
- TL-022FB73 FOB Turbilatex Calibrator 3, (1 x 0.3mL)
- TL-022FB74 FOB Turbilatex Calibrator 4, (1 x 0.3mL)
- TL-022FB75 FOB Turbilatex Calibrator 5, (1 x 0.3mL)
- TL-022FB08 FOB Turbilatex Control 1, (2 x 0.5mL)
- TL-022FB09 FOB Turbilatex Control2, (2 x 0.5mL)
- TL-022FB01 FOB Turbilatex Reagent 1, (2 x 20mL)
- TL-022FB02 FOB Turbilatex Reagent 2, (1 x 11mL)
- TL-022FB70 FOB Turbilatex Calibrator 0, (1 x 0.3mL)
- TL-022FB71 FOB Turbilatex Calibrator 1, (1 x 0.3mL)
- TL-022FB72 FOB Turbilatex Calibrator 2, (1 x 0.3mL)
- TL-022FB73 FOB Turbilatex Calibrator 3, (1 x 0.3mL)
- TL-022FB74 FOB Turbilatex Calibrator 4, (1 x 0.3mL)
- TL-022FB75 FOB Turbilatex Calibrator 5, (1 x 0.3mL)
- TL-022FB08 FOB Turbilatex Control 1, (2 x 0.5mL)
- TL-022FB09 FOB Turbilatex Control2, (2 x 0.5mL)
- TL-022FB01 FOB Turbilatex Reagent 1, (4 x 20mL)
- TL-022FB02 FOB Turbilatex Reagent 2, (2 x 11mL)
- TL-022FB70 FOB Turbilatex Calibrator 0, (1 x 0.3mL)
- TL-022FB71 FOB Turbilatex Calibrator 1, (1 x 0.3mL)
- TL-022FB72 FOB Turbilatex Calibrator 2, (1 x 0.3mL)
- TL-022FB73 FOB Turbilatex Calibrator 3, (1 x 0.3mL)
- TL-022FB74 FOB Turbilatex Calibrator 4, (1 x 0.3mL)
- TL-022FB75 FOB Turbilatex Calibrator 5, (1 x 0.3mL)
- TL-022FB08 FOB Turbilatex Control 1, (2 x 0.5mL)
- TL-022FB09 FOB Turbilatex Control2, (2 x 0.5mL)

Materials and instruments that could be required but no provided

- Specimen collection container.
- TL-022FB03E FOB Turbilatex Sample Diluent.
- TL-0221213E Microtubes.
- Disposable gloves and laboratory equipment.
- Automated analyser.

SPECIMEN COLLECTION AND PREPARATION

Collect sufficient quantity of human stool samples. These samples should be collected in clean and dry containers (no preservatives or transport media). The samples can be stored completely dissolved in sample diluent in the refrigerator (2-8°C) for 7 days prior to testing. Homogenize stool samples as thoroughly as possible prior to preparation.

To process the collected stool samples:

Use a centrifugation tube for each sample to be tested. Label centrifugation tube with name or number of patient.

- 1. Homogenize the sample. Add 20 mg of sample into centrifugation tube.

2. Add 2 mL of sample diluent (TL-022FB03E).
3. Shake vigorously the tube in order to assure good sample dispersion (vortex) until sample is completely dissolved.
4. Centrifuge for 15 minutes at 10000 g or 10 minutes at 15000 g.
5. Take the supernatant to automated analyser vial.

**ANALYTICAL PROCEDURE**

**Allow reagents and stool samples to reach room temperature (15-30°C) prior to testing.**

R1 and R2 are ready to use.

**Preparation of the calibration curve**

For calibration use only the following materials:


TL-022FB70, TL-022FB71, TL-022FB72, TL-022FB73, TL-022FB74 and TL-022FB75, FOB Turbilatex Calibrator vials: liquid calibrator containing human haemoglobin at different concentrations.

Concentration is indicated on the label of each vial.

Calibrate the system every week is extremely recommended. Recalibrate the system when reagent lot is changed or when the controls are out of the assigned range given on the control label or certificate of analysis provided with the kit.

**Allow FOB Turbilatex Calibrator to reach room temperature (15-30°C) prior to testing.**

FOB Turbilatex Calibrator vials (calibrator 0 to 5) are ready to use.

**Picture 1.** Prepared calibrators for performing the calibration curve (5 points and blank)

**Quality control**

For quality control only use the following materials:


TL-022FB08 FOB Turbilatex Control 1 and TL-022FB09 FOB Turbilatex Control 2; liquid controls at two different concentrations of human haemoglobin. Concentration is indicated on the label of the vial. The use of control materials at two different concentrations is recommended in order to verify test precision.

Controls should be assayed every day before running patient faecal samples extract to validate the calibration curve.

If the obtained results are out of the tolerance range, the equipment, the reagents or the technique must be reviewed.

**Allow FOB Turbilatex C1 & C2 to reach room temperature (15-30°C) prior to testing.**

FOB Turbilatex Control vials (Control 1 and Control 2) are ready to use.

**Picture 2.** Prepared controls C1 and C2 (check control concentration on the label).

**Analytical procedure, see annex information.**

**INTERPRETATION OF RESULTS**

Positive results: higher or equal than the cut-off fixed by the clinical lab.

Recommended: 5 µg of hHb/g of stool for diagnostic procedures and 20 µg of hHb/g of stool for screening procedures.

Positive results determine the abnormal presence of human Haemoglobin (hHb) in stool samples.

**PERFORMANCE CHARACTERISTICS**

**Analytical sensitivity**

Limit of Detection (LOD): **see annex information.**

FOB Turbilatex cut-off value:

Haemoglobin concentration values low than 50 ng of hHb/mL are considered normal values and that is not indicative of an inflammation of gastrointestinal tract.

Haemoglobin concentration values higher than 200 ng of hHb/mL are indicative of a bleeding in the gastrointestinal tract.

FOB Turbilatex cut-off value: 50-200 ng hHb/mL (X100 dilution factor, 5-20 µg hHb/g of stool).

**Prozone**

Studies have been made up to a concentration of 10 µg of hHb/mL and no false negative results have been observed. Studies using higher concentrations have not been carried out.

**Clinical sensitivity and specificity**

An evaluation was performed comparing a turbidimetric assay (FOB Turbilatex, CerTest) and another commercial immunoassay (EIKEN FOB Latex, EIKEN Chemicals, CO. LTD). The results were as follows:

	Sensitivity	Specificity
FOB Turbilatex vs EIKEN FOB Latex	96%	>99%

The results showed a high sensitivity and specificity to detect human haemoglobin (hHb) using FOB Turbilatex.

**Interferences and cross reactivity**

An evaluation was performed to determine the cross reactivity; no cross reactivity was founded against other faecal markers occasionally present in the faeces, such as: bovine lactoferrin, human calprotectin, lactoferrin and transferrin.

The cross reactions with bovine and pig haemoglobins are not present at concentrations below 1 µg/mL and 2.5 µg/mL respectively.

**LIMITATIONS**

1. FOB Turbilatex should only be used in human stool samples. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
2. Positive results determine the presence of human haemoglobin in faecal samples. A positive result should be followed up with additional diagnostic invasive procedures to determine the exact cause and source of the blood in the stool.
3. If symptoms or situation still persist, haemoglobin determination should be carried out invasive techniques. Negative results do not exclude bleeding, as some polyps and colorectal cancers may bleed intermittently or not during certain stages of the disease. Moreover, blood may not be uniformly distributed in stool samples.
4. This test may be less sensitive for detecting upper gastrointestinal bleeding because blood degrades as it passes through the gastrointestinal tract.
5. Patients should not collect samples during their menstrual period, if they have bleeding haemorrhoids, blood in urine or if they have strained during bowel movement.

**EXPECTED VALUES**

Colorectal cancer affects both men and women of all racial and ethnic groups, being most often found in people aged 50 years or older. Regarding men, colorectal cancer is the third most common cancer after prostate and lung cancers. Regarding women, colorectal cancer is the third most common cancer after breast and lung cancers.

Some studies established equal or higher 5.1 µg hHb/g faeces as cut-off values to allow detect adult patients with gastrointestinal bleeding that will require to diagnosis additional diagnostic invasive procedures. In FOB Turbilatex test, this value is equivalent to 50 ng hHb/mL of diluent.

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1. Towler BP, Irwig L, Glasziou P, Weller D, Kewenter J. Screening for colorectal cancer using the faecal occult blood test, Hemoccult. Cochrane Database Syst Rev. 2000;(2): CD001216.
2. Ransohoff DF and Lang CA. Screening for colorectal cancer with the Fecal Occult Blood Test: a background paper. Ann Intern Med. 1997; 126: 811-822.
3. Ransohoff DF and Lang CA. Suggested technique for Faecal Occult Blood testing and interpretation in colorectal cancer screening. Ann Intern Med. 1997; 126: 808-810.

**SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO**

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
The application notes for use our turbidimetric products with the different equipments could be consulted in our website:

[www.certest.es](http://www.certest.es)

For information purposes only

For information purposes only

For information purposes only

**CerTest**  
BIOTEC S.L.

Pol. Industrial Río Gállego II, Calle J, Nº 1,  
50840, San Mateo de Gállego, Zaragoza (SPAIN)  
[www.cerTest.es](http://www.cerTest.es)

F-381 rev00