---

# H. pylori TURBILATEX Combo

---

H. pylori  
latex turbidimetric assay

Cat.TL-022HP100ED, TL-022HP200ED


---

CERTEST BIOTEC S.L.

---

For information purposes only

**ENGLISH**

**INTENDED USE**

*H. pylori* Turbilateral is a diagnostic assay **for the quantitative detection of *Helicobacter pylori* in human stool samples.**

This assay is simple and widely applicable. This product is optimized for several automated analyser.

For professional *in vitro* diagnostic use only.

**SYNTHESIS**

*H. pylori* is a spiral-shaped bacterium, found in the gastric mucous layer or adherent to the epithelial lining of the stomach. *H. pylori* causes more than 90% of duodenal ulcers and up to 80% of gastric ulcers.

The importance of *Helicobacter pylori* testing has increased greatly since the strong correlation between the presence of bacteria and confirmed gastrointestinal diseases (stomach and duodenum) like gastritis, peptic ulcer disease and gastric carcinoma.

**PRINCIPLE**

*H. pylori* Turbilateral is a latex turbidimetric assay **for the quantitative detection of *H. pylori* in human stool samples.** This assay is simple, widely applicable and test results are obtained in a very short time.

The intended use for this assay is the detection of specific *H. pylori* antigen at very low concentrations.

*H. pylori* latex turbidimetric assay is based on agglutination reactions. These involve *in vitro* aggregation of microscopic latex particles. This aggregation consists in the specific reaction between antigen and antibodies, antigen contained in the sample and the antibodies anti-antigen coated on polystyrene latex particles. The sample is mixed with a suspension containing antibodies against the antigen bound to latex particles. If antigen is present in the sample it will react with the antibodies and form an aggregate. If no antigen is present in the sample the mixture will keep its appearance as a smooth suspension. Such turbidity is measured as an increase in absorbance at the determinate wave and is proportional to the quantity of antigen contained in the sample.

**PRECAUTIONS**

- For professional *in vitro* diagnostic use only.
- The analyser should be ready before performing any assay.
- Read and follow the instructions for use provided with the kit.
- If the result exceeds the measurement range, use the sample diluent to dilute the sample and repeat the assay again.
- Do not use after expiration date.
- Do not use the reagents if pack is damaged or opened.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The solutions should be discarded in a proper container after testing following local regulations.

**STORAGE AND STABILITY**

Kit components must be stored at refrigerated temperature (2-8°C). Do not freeze.

Reagents are stable up to the expiration date printed on the label, always considering that reagent containers must be properly closed to avoid any contamination, must be kept away from the sunlight and conserved at 2-8°C.

**REAGENTS**

**Materials provided**

**Cat.: TL-022HP100ED**

- TL-022HP01 *H. pylori* Turbilateral Reagent 1, (1 x 31mL)
- TL-022HP02 *H. pylori* Turbilateral Reagent 2, (1 x 2.5mL)
- TL-022HP70 *H. pylori* Turbilateral Calibrator 0, (1 x 0.3mL)
- TL-022HP71 *H. pylori* Turbilateral Calibrator 1, (1 x 0.3mL)
- TL-022HP72 *H. pylori* Turbilateral Calibrator 2, (1 x 0.3mL)
- TL-022HP73 *H. pylori* Turbilateral Calibrator 3, (1 x 0.3mL)
- TL-022HP74 *H. pylori* Turbilateral Calibrator 4, (1 x 0.3mL)
- TL-022HP75 *H. pylori* Turbilateral Calibrator 5, (1 x 0.3mL)
- TL-022HP08 *H. pylori* Turbilateral Control 1, (2 x 0.5mL)
- TL-022HP09 *H. pylori* Turbilateral Control 2, (2 x 0.5mL)

**Cat.: TL-022HP200ED**

- TL-022HP01 *H. pylori* Turbilateral Reagent 1, (2 x 31mL)
- TL-022HP02 *H. pylori* Turbilateral Reagent 2, (1 x 5mL)
- TL-022HP70 *H. pylori* Turbilateral Calibrator 0, (1 x 0.3mL)
- TL-022HP71 *H. pylori* Turbilateral Calibrator 1, (1 x 0.3mL)
- TL-022HP72 *H. pylori* Turbilateral Calibrator 2, (1 x 0.3mL)
- TL-022HP73 *H. pylori* Turbilateral Calibrator 3, (1 x 0.3mL)
- TL-022HP74 *H. pylori* Turbilateral Calibrator 4, (1 x 0.3mL)
- TL-022HP75 *H. pylori* Turbilateral Calibrator 5, (1 x 0.3mL)
- TL-022HP08 *H. pylori* Turbilateral Control 1, (2 x 0.5mL)
- TL-022HP09 *H. pylori* Turbilateral Control 2, (2 x 0.5mL)

**Materials and instruments that could be required but no provided**

- Specimen collection container.
- TL-022HP03E *H. pylori* Turbilateral Sample Diluent.
- TL-0221213E Microtubes.
- Disposable gloves and laboratory equipment.
- Automated analyser.

**SPECIMEN COLLECTION AND PREPARATION**

Collect sufficient quantity of human stool samples. These samples should be collected in clean and dry containers (no preservatives or transport media). The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. Homogenize stool samples as thoroughly as possible prior to preparation.

**To process the collected stool samples:**

Use a microtube for each sample to be tested. Label microtube with name or number of patient.

1. Homogenize the sample. Add 20 mg of sample into a microtube.
2. Add 2 mL of sample diluent (TL-022HP03E).
3. Shake vigorously the microtube in order to assure good sample dispersion (vortex) until sample is completely dissolved.
4. Rest the microtube with sample up for 20 minutes to get a proper extraction.
5. Take the supernatant to automated analyser vial.


**ANALYTICAL PROCEDURE**

**Allow reagents and stool samples to reach room temperature (15-30°C) prior to testing.**

R1 and R2 are ready to use.

**Preparation the calibration curve.**


For calibration use only the following materials:  
 TL-022HP70, TL-022HP71, TL-022HP72, TL-022HP73, TL-022HP74 and TL-022HP75, *H. pylori* Turbilatex Calibrator vials: liquid calibrator containing recombinant *H. pylori* antigen at different concentrations.  
 Concentration is indicated on the label of each vial.  
 Calibrate the system every week is extremely recommended. Recalibrate the system when reagent lot is changed or when the controls are out of the assigned range given on the control label or certificate of analysis provided with the kit.  
**Allow *H. pylori* Turbilatex Calibrator to reach room temperature (15-30°C) prior to testing.**  
*H. pylori* Turbilatex Calibrator vials (calibrators 0 to 5) are ready to use.

**Picture 1.** Prepared calibrators for performing the calibration curve (5 points and blank)

**Quality control**

For quality control only use the following materials:  
 TL-022HP08 *H. pylori* Turbilatex Control 1 and TL-022HP09 *H. pylori* Turbilatex Control 2: liquid controls at two different concentrations of recombinant *H. pylori* antigen.  
 Concentration is indicated on the label of the vial. The use of control materials at two different concentrations is recommended in order to verify test precision.  
 Controls should be assayed every day before running patient faecal samples extract to validate the calibration curve.  
 If the obtained results are out of the tolerance range, the equipment, the reagents or the technique must be reviewed.  
**Allow *H. pylori* Turbilatex C1 & C2 to reach room temperature (15-30°C) prior to testing.**  
*H. pylori* Turbilatex Control vials (Control 1 and C2) are ready to use.

**Picture 2.** Prepared controls C1 and C2 (check control concentration on the label).

**Analytical procedure, see annex information**

| INTERPRETATION OF RESULTS                                                                                                                                                                                                                                                                                                                           |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Negative results: Lower to the limit of detection of 0.5 ng <i>H. pylori</i> antigen/mL.<br>Negative results determine the absence of <i>H. pylori</i> antigen in higher concentrations than the detection limit.                                                                                                                                   |
| Grey area results: Between 0.5 and 1 ng <i>H. pylori</i> antigen/mL.<br>Results between 0.5 and 1 ng/mL are between the detection and quantification limit. A second sampling is recommended to clarify the diagnosis. In case, the second result also shows a value in this range, the recommendation is to follow up the patient some time later. |
| Positive results: Higher to the limit of quantification of 1 ng <i>H. pylori</i> antigen/mL.<br>Positive results determine the presence of <i>H. pylori</i> in human stool samples.                                                                                                                                                                 |

**PERFORMANCE CHARACTERISTICS**

**Analytical sensitivity**

Limit of detection (LOD): 0.5 ng of anti-*H. pylori* antigen/mL (50 ng/g of stool)., determined x10 replicates with stabilizing diluent and reported as LOB +2xSD.

**Prozone**  
Studies have been made up to a concentration of 2 µg of anti-*H. pylori* antigen/mL and no false negative results have been observed. Studies using higher concentrations have not been carried out.

**Clinical sensitivity and specificity**

An evaluation was performed comparing a turbidimetric assay (*H. pylori* Turbilatex, CerTest) and an immunochromatographic test (CerTest *H. pylori*, CerTest). The results were as follows:

|                                                         | Sensitivity | Specificity |
|---------------------------------------------------------|-------------|-------------|
| <i>H. pylori</i> Turbilatex vs CerTest <i>H. pylori</i> | 88%         | >98%        |

The results showed a high sensitivity and specificity to detect *H. pylori* using *H. pylori* Turbilatex.

**Interferences and cross reactivity**

An evaluation was performed to determine the cross reactivity of the specific antibody anti- *H. pylori*. No cross reactivity was founded against other bacteria present on faecal samples, such as: *Campylobacter coli*, *Salmonella enteritidis*, *Shigella dysenteriae*, *Campylobacter jejuni*, *Salmonella paratyphi*, *Shigella flexneri*, *Clostridium difficile*, *Salmonella typhi*, *Shigella sonnei*, *Escherichia coli*, O157:H7, *Salmonella typhimurium*, *Staphylococcus aureus*, *Listeria monocytogenes*, *Shigella boydii*, *Yersinia enterocolitica*.

**LIMITATIONS**

1. *H. pylori* Turbilatex should only be used in human stool samples. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
2. Positive results determine the presence of *H. pylori* in faecal samples. A positive result should be followed up with additional invasive techniques (endoscopy) to determine the presence of these bacteria. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
3. Positive results can be obtained even when the symptoms had disappeared.
4. Negatives results could be not meaningful because the antigen concentration could be lower than the detection limit value. If the symptoms or situation still persist, a *H. pylori* determination should be carried out, on a sample from an enrichment culture or using an invasive technique.

**EXPECTED VALUES**

Studies have found that more than 90% of patients with duodenal ulcer and 80% of patients with gastric ulcer are infected with *Helicobacter pylori*.

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