---

# CERTEST GDH+ Lactoferrin

---

ONE STEP  
GDH+Lactoferrin COMBO CARD  
TEST

---

CERTEST BIOTEC S.L.

---

For Information Purposes Only

**ENGLISH****INTENDED USE**

CerTest GDH+Lactoferrin one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection of *Clostridium difficile* Glutamate Dehydrogenase (GDH) and human lactoferrin (hLf) in stool samples.

CerTest GDH+Lactoferrin combo card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of *Clostridium difficile* infection and intestinal inflammation presence.

**INTRODUCTION**

*Clostridium difficile* (*C. difficile*), a Gram-positive spore bearing anaerobic bacterium is the major aetiological agent of diarrhoea and colitis associated with antibiotics. *C. difficile* is the most common cause of health care-associated diarrhoea in developed countries and is a major source of nosocomial morbidity and mortality worldwide.

Disease due to *C. difficile* develops when the organism is allowed to proliferate in the colon, most commonly after antibiotic use has eliminated competing flora. *C. difficile* can release two high-molecular-weight toxins, toxin A and toxin B, which are responsible for the clinical manifestations, which range from mild, self-limited watery diarrhoea to fulminant pseudomembranous colitis, toxic megacolon and death.

*Clostridium difficile* Glutamate Dehydrogenase (GDH) is an enzyme produced in large quantities by all toxigenic and non-toxigenic strains, making it an excellent marker for the organism.

The toxigenic culture (TC) is used as the gold standard technique to determine *Clostridium difficile* infection. This method consists in culture and isolation of *C. difficile* from faeces, followed by toxin testing of the isolate, a labour-intensive assay to obtain a result.

Lactoferrin is a glycoprotein component of neutrophil secondary granules, a primary component of the acute inflammatory response and is released from fecal leukocytes. This protein is resistant to proteolysis in the faeces and may serve as a marker of inflammation in the intestine. The major cause of faecal neutrophils in patients with chronic diarrhoea is chronic inflammatory bowel disease of the colon (i.e., Crohn's Disease and ulcerative colitis). Lactoferrin has been also studied as a predictor of infection with invasive enteropathogens in children with diarrhoea. Bacterial inflammatory diarrhea may be caused by *Shigella*, *Salmonella*, *Campylobacter* and *Clostridium difficile*.

**TEST PRINCIPLE**

CerTest GDH+Lactoferrin is based on the principle of a qualitative immunochromatographic assay for the determination *Clostridium difficile* Glutamate Dehydrogenase (GDH) and human lactoferrin (hLf) in stool samples.

**Strip A** consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against GDH and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-GDH) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

**Strip B** consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against lactoferrin and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-human lactoferrin) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is GDH positive, the antigen of the diluted sample reacts with the red-coloured conjugates complex (anti-GDH monoclonal antibodies-red polystyrene microspheres) in the strip A, and if the sample is lactoferrin positive, the antigens of the diluted sample react with the red-coloured conjugated complex (anti-human lactoferrin monoclonal antibodies-red polystyrene microspheres) in strip B, which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-GDH antibodies present on the membrane of strip A (test line) and the anti-human lactoferrin antibodies present on the membrane of strip B (test line) capture the coloured conjugate and the red line will be visible in both strips. These bands are used to interpret the result.

If the sample is negative, there is no GDH and lactoferrin antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-GDH and the anti-human lactoferrin antibodies present on the membranes (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, in both strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on both membranes will capture control green-conjugate complex and both control lines will always appear. The presence of these green lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

## CERTEST GDH+Lactoferrin

One Step test to detect Clostridium difficile antigen GDH and human lactoferrin in combo card format/  
Prueba combo de un solo paso para detección del antígeno GDH de Clostridium difficile y lactoferrina humana en formato cassette

### STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

### PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest GDH+Lactoferrin combo card test. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

### SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. For longer storage, maximum 6 months, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Ensure only the amount needed is thawed because of freezing and defrosting cycles are not recommended for GDH antigen. Homogenise stool sample as thoroughly as possible prior to preparation.


**Please take note that there are two different stool collection tubes.**

**a. "GDH" stool collection tubes (without reducer piece).**

**Specimen preparation (see illustration):**


1. Take out the cap of the stool collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2), to collect faecal sample (approx. 125mg) and add it to the collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

**b. "Lactoferrin" stool collection tubes (with a reducer piece).** 

**Specimen preparation (see illustration):**


1. Take out the cap of the stool collection tube and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (1), to collect faecal sample and add it to the stool collection tube. For liquid samples, add approx. 15µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (2). The stool collection tube with dilute sample can be stored for 7 days in the refrigerator (2-8°C) prior to testing.

**MATERIALS**

**MATERIALS PROVIDED**

- CerTest GDH+Lactoferrin combo card tests
- Instructions for use
- "GDH" stool collection tubes with diluent
- "Lactoferrin" stool collection tubes with diluent

**MATERIALS REQUIRED BUT NO PROVIDED**


- Specimen collection container
- Disposable gloves
- Timer

**TEST PROCEDURE**

**Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.**

1. Proceed to shake both stool collection tubes in order to assure good sample dispersion.
2. Remove CerTest GDH+Lactoferrin combo card test from its sealed bag just before using it.
3. Take the "GDH" stool collection tube, cut the end of the cap (4) and dispense 4 drops in the circular window marked with the letter A (5). Avoid adding solid particles with the liquid.
4. Take the "Lactoferrin" stool collection tube, cut the end of the cap (6) and dispense 4 drops in the circular window marked with the letter B (7). Avoid adding solid particles with the liquid.
5. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.


If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of extraction diluent until seeing the liquid running through the reaction zone.

One Step test to detect *Clostridium difficile* antigen GDH and human lactoferrin in combo card format/  
 Prueba combo de un solo paso para detección del antígeno GDH de *Clostridium difficile* y lactoferrina humana en formato cassette


### GDH Strip A-procedure

Add 4 drops in A-circular window


### Lactoferrin Strip B-procedure


Add 4 drops in B-circular window

### INTERPRETATION OF THE RESULTS (please refer to the illustration below)

### INVALID RESULTS

	A (GDH)	B (Lactoferrin)	Interpretation of the results
1.	- GREEN	- GREEN	There are no GDH and human lactoferrin presence which might mean no <i>Clostridium difficile</i> infection.
2.	+ GREEN-RED	+ GREEN-RED	There are GDH and human lactoferrin presence which might mean <i>Clostridium difficile</i> infection (clinical severe disease), intestinal inflammation associated with other infections ( <i>Campylobacter</i> , <i>Salmonella</i> and <i>Shigella spp.</i> ) or non-infectious causes of intestinal inflammation (ulcerative colitis, gastrointestinal bleeding, duodenal ulcer, etc.).
3.	+ GREEN-RED	- GREEN	Infection caused by <i>Clostridium difficile</i> which might mean milder diarrhoea or possible carrier of disease.
4.	- GREEN	+ GREEN-RED	There is human lactoferrin presence which might mean an intestinal inflammation associated with other infections ( <i>Campylobacter</i> , <i>Salmonella</i> and <i>Shigella spp.</i> ) or non-infectious causes of intestinal inflammation (ulcerative colitis, gastrointestinal bleeding, duodenal ulcer, etc.).
5.	Any other result	Any other result	Invalid result, either A or B, we recommend repeating the assay using the same sample with another test.

## CERTEST GDH+Lactoferrin

One Step test to detect *Clostridium difficile* antigen GDH and human lactoferrin in combo card format/  
Prueba combo de un solo paso para detección del antígeno GDH de *Clostridium difficile* y lactoferrina humana en formato cassette

**INVALID:** Total absence of any control coloured line (GREEN) regardless the appearance or not of the test lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

### NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration antigen present in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

### QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

### LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
3. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
4. CerTest GDH+Lactoferrin should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
5. Stool from patients with active inflammatory bowel diseases that usually involve significant neutrophilic inflammation of the intestine, such as Crohn's disease and ulcerative colitis, would be positive for faecal lactoferrin.
6. Lactoferrin is a component of breast milk, the test will be positive in breast fed children and should not be used to evaluate neonates receiving breast milk.
7. Positive results determine the presence of *Clostridium difficile* and human lactoferrin in faecal samples and corroborate an intestinal inflammation associated with advanced *Clostridium difficile* infection (severe disease); nevertheless, it can be due to toxigenic or non-toxigenic strains of *Clostridium difficile*. A positive result should be followed up with additional laboratory techniques (toxigenic culture) to determine the strain. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
8. A negative result is not meaningful because of it is possible the antigens concentration in the stool samples is lower than the detection limit value. If the symptoms or situation still persist, a *Clostridium difficile* determination should be carried out on a sample from an enrichment culture. Negative results do not exclude inflammation, some diseases such as celiac sprue and microscopic colitis polyps that involve mainly monocuclear inflammation.

### EXPECTED VALUES

*Clostridium difficile* is associated with 95-100% cases of pseudomembranous colitis, 60-75% cases of antibiotic-associated colitis and 35% of cases of antibiotic-associated diarrhoea cases.

Some studies established equal or higher 10µg hLf/g faeces as cut-off values to allow detect adult patients with intestinal inflammatory pathology that will require to diagnosis additional diagnostic invasive procedures.

### PERFORMANCE CHARACTERISTICS

#### Analytical sensitivity

Detection limit value of GDH is 0.8ng/mL for *Clostridium difficile* antigen GDH and cut-off value of Lactoferrin is 100ng/mL (10µg hLf/g faeces) for human lactoferrin.

#### Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest GDH+Lactoferrin, CerTest) and other available commercial tests (*C. DIFF QUIK CHEK Complete*® and IBD EZ VUE®, TechLab). The results were as follows:

# CERTEST GDH+Lactoferrin

One Step test to detect *Clostridium difficile* antigen GDH and human lactoferrin in combo card format/  
Prueba combo de un solo paso para detección del antígeno GDH de *Clostridium difficile* y lactoferrina humana en formato cassette

		IC test: <i>C. DIFF QUIK CHEK Complete</i> ®		
		+	-	Total
IC test: CerTest GDH+Lactoferrin (Strip A: GDH)	+	26	0	26
	-	0	48	48
Total		26	48	74

	Sensitivity	Specificity	PPV	NPV
CerTest GDH+Lactoferrin (Strip A) vs <i>C. DIFF QUIK CHEK Complete</i> ®	>99%	>99%	>99%	>99%

		IBD EZ VUE®		
		+	-	Total
IC test: CerTest GDH+Lactoferrin (Strip B: Lactoferrin)	+	20	0	20
	-	0	44	44
Total		20	44	64

	Sensitivity	Specificity	PPV	NPV
CerTest GDH+Lactoferrin (Strip B) vs IBD EZ VUE®	>99%	>99%	>99%	>99%

The results showed a high sensitivity and specificity to detect *Clostridium difficile* glutamate dehydrogenasa (GDH) and human lactoferrin using CerTest GDH+Lactoferrin.

### Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest GDH+Lactoferrin; no cross reactivity against gastrointestinal pathogens and other faecal markers occasionally present in faeces:

- |                                   |                               |                                |
|-----------------------------------|-------------------------------|--------------------------------|
| <i>Bovine and pig haemoglobin</i> | <i>Human transferrin</i>      | <i>Shigella boydii</i>         |
| <i>Bovine and pig transferrin</i> | <i>Human calprotectin</i>     | <i>Shigella dysenteriae</i>    |
| <i>Bovine and pig lactoferrin</i> | <i>Listeria monocytogenes</i> | <i>Shigella flexneri</i>       |
| <i>Campylobacter jejuni</i>       | <i>Salmonella enteritidis</i> | <i>Shigella sonnei</i>         |
| <i>Campylobacter coli</i>         | <i>Salmonella paratyphi</i>   | <i>Staphylococcus aureus</i>   |
| <i>Escherichia coli</i> O157:H7   | <i>Salmonella typhi</i>       | <i>Yersinia enterocolitica</i> |
| <i>Helicobacter pylori</i>        | <i>Salmonella typhimurium</i> |                                |

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		IC test: C. DIFF QUIK CHEK Complete®		
		+	-	Total
IC test: CerTest GDH+Lactoferrin (Tira A: GDH)	+	26	0	26
	-	0	48	48
	Total	26	48	74

	Sensibilidad	Especificidad	VPP	VPN
CerTest GDH+Lactoferrin (Tira A) vs C. DIFF QUIK CHEK Complete®	>99%	>99%	>99%	>99%

		IBD EZ VUE®		
		+	-	Total
IC test: CerTest GDH+Lactoferrin (Tira B: Lactoferrin)	+	20	0	20
	-	0	44	44
	Total	20	44	64

	Sensibilidad	Especificidad	VPP	VPN
CerTest GDH+Lactoferrin (Tira B) vs IBD EZ VUE®	>99%	>99%	>99%	>99%

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
1. WREN, M.W.D, et al. "Laboratory diagnosis of *Clostridium difficile* infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory". *British Journal of Biomedical Science*, 66 (1), 2009.
2. VAISHNAVI, Ch., "Clinical spectrum & pathogenesis of *Clostridium difficile* associated diseases". *Indican J. Med. Res.* 131, April 2010, pp 487-499.
3. Angriman I. et al. Enzymes in feces: Useful markers of chronic inflammatory bowel disease. *Clinica Chimica Acta* 381 Feb 2007, p. 63-68.
4. Guerrant R. et al. Measurement of Fecal Lactoferrin as a Marker of Fecal Leukocytes. *Journal of Clinical Microbiology*, Vol 30 No 5; May 1992, p. 1238-1242.
5. Langhorst, M.D. et al. Noninvasive Markers in the Assessment of Intestinal Inflammation in Inflammatory Bowel Diseases: Performance of Fecal Lactoferrin, Calprotectin and PMN-Elastase, CRP, and Clinical Indices. *Am. J. Gastroenterol.* 2008; Vol 103, p. 162-169.

### SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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