---

**CERTEST**  
*Rotavirus+*  
*Norovirus GI+GII*

---

ONE STEP  
*Rotavirus, Norovirus GI and*  
*Norovirus GII*  
COMBO CARD TEST

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## ENGLISH

### INTENDED USE

CerTest *Rotavirus+Norovirus* GI+GII one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection of *Rotavirus*, *Norovirus* genogroup I and II (GI and GII) in stool samples.

CerTest *Rotavirus+Norovirus* GI+GII combo card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of *Rotavirus* and/or *Norovirus* infection.

### INTRODUCTION

*Rotavirus*, *Adenovirus* and *Astrovirus* are most common and major causes of severe gastroenteritis in infants and young children. Pattern also observed in adults. They are transmitted by faecal-oral contact. The main symptoms of viral gastroenteritis are watery diarrhoea and vomiting. The affected person may also have headache, fever and abdominal cramps ("stomach ache"). In general, the symptoms begin 1 to 2 days following infection with a virus that causes gastroenteritis and may last for 1 to 10 days, depending on which virus causes the illness (*Rotavirus* 3 days, *Adenovirus* 5-8 days and *Astrovirus* 3 days).

*Noroviruses* are the leading cause of epidemic gastroenteritis, including foodborne outbreaks, in the United States. *Noroviruses* belong to the family *Caliciviridae* and can be grouped into five genogroups (GI through GV), which are further divided into at least 34 genotypes. Genotypes GI, GII and GIV infect humans, causing gastroenteritis, whereas GIII and GV typically infect animals.

The symptoms of *Norovirus* related diseases are those typical of gastroenteritis, that is, vomiting, watery diarrhea and abdominal cramps. Vomiting is a characteristic symptom in the majority of *Norovirus* infections (64% of adults and 81% of children). Other symptoms such as general malaise, low grade fever, nausea and fatigue are also present in over 90% of cases. The incubation period of the disease is generally between 12 and 48 hours, while infection lasts between 12 and 60 hours. Infection may also be asymptomatic, and thus contribute to the spread of the virus in the community. As a rule, the disease does not have serious consequences, and most patients recover within 1-2 days without complications. Debilitated patients and persons with weaker immune systems such as children, elderly or chronic patients may be affected by more serious forms of the disease. Specially, dehydration may represent a serious complication for children, the elderly and persons with a precarious metabolic balance or cardio circulatory instability.

### TEST PRINCIPLE

CerTest *Rotavirus+Norovirus* GI+GII is based on the principle of a qualitative immunochromatographic assay for the determination of *Rotavirus*, *Norovirus* GI and *Norovirus* GII in stool samples.

**Strip A** consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Rotavirus* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Rotavirus*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

**Strip B** consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Norovirus* GI and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Norovirus* GI) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

**Strip C** consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Norovirus* GII and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Norovirus* GII) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

If the sample is *Rotavirus* positive, the antigen of the diluted sample reacts with the red-coloured conjugates complex (anti-*Rotavirus* monoclonal antibodies-red polystyrene microspheres) in the strip A, if the sample is *Norovirus* GI positive, the antigens of the diluted sample react with the red-coloured conjugates complex (anti-*Norovirus* GI monoclonal antibodies-red polystyrene microspheres) in the strip B and if the sample is *Norovirus* GII positive, the antigens of the diluted sample react with the red-coloured conjugated complex (anti-*Norovirus* GII monoclonal antibodies-red polystyrene microspheres) in the strip C, which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*Rotavirus* antibodies present on the membrane of strip A (test line), the anti-*Norovirus* GI antibodies present on the membrane of strip B (test line) and the anti-*Norovirus* GII antibodies present on the membrane of strip C (test line) capture the coloured conjugate and the red lines will be visible in the three strips. These bands are used to interpret the result.

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If the sample is negative, there is no *Rotavirus*, *Norovirus GI* and *Norovirus GII* presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-*Rotavirus*, the anti-*Norovirus GI* and the anti-*Norovirus GII* antibodies present on the membranes (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in the three strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on three membranes will capture control green-conjugate complex and three control lines will always appear. The presence of these green lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

### STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

### PRECAUTIONS


- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest *Rotavirus+Norovirus GI+GII* combo card test. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

### SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

#### Specimen preparation (see illustration):

1. Take out the cap of the stool collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2), to collect faecal sample (approx. 125mg) and add it to the stool collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

### MATERIALS

#### MATERIALS PROVIDED

- CerTest *Rotavirus+Norovirus GI+GII* combo card tests
- Instructions for use
- Stool collection tubes with diluent

#### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer


## TEST PROCEDURE

Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.


1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove CerTest *Rotavirus+Norovirus* GI+GII combo card test from its sealed bag just before using it.
3. Take the stool collection tube, cut the end of the cap (4) and dispense 4 drops in the circular window marked with the letter A (5), 4 drops, using the same tube, in the circular window marked with the letter B (6), and 4 drops, using the same tube, in the circular window marked with the letter C (7). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of extraction diluent until seeing the liquid running through the reaction zone.


### Rotavirus Strip A-procedure

Add 4 drops in A-circular window


### Norovirus GI Strip B-procedure

Add 4 drops in B-circular window


### Norovirus GII Strip C-procedure


Add 4 drops in C-circular window


## INTERPRETATION OF THE RESULTS (please refer to the illustration below)


A: Green→Negative Rotavirus  
B: Green→Negative Norovirus GI  
C: Green→Negative Norovirus GII


A: Green/Red→Positive Rotavirus  
B: Green/Red→Positive Norovirus GI  
C: Green/Red→Positive Norovirus GII


A: Green/Red→Positive Rotavirus  
B: Green/Red→Positive Norovirus GI  
C: Green→Negative Norovirus GII


A: Green/Red→Positive Rotavirus  
B: Green→Negative Norovirus GI  
C: Green/Red→Positive Norovirus GII


A: Green/Red→Positive Rotavirus  
B: Green→Negative Norovirus GI  
C: Green→Negative Norovirus GII

A: Green→Negative Rotavirus  
B: Green/Red→Positive Norovirus GI  
C: Green/Red→Positive Norovirus GII


A: Green→Negative Rotavirus  
B: Green/Red→Positive Norovirus GI  
C: Green→Negative Norovirus GII

A: Green→Negative Rotavirus  
B: Green→Negative Norovirus GI  
C: Green/Red→Positive Norovirus GII

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	<b>A</b> <b>(Rotavirus)</b>	<b>B</b> <b>(Norovirus GI)</b>	<b>C</b> <b>(Norovirus GII)</b>	<b>Interpretation of the results</b>
<b>1.</b>	- GREEN	- GREEN	- GREEN	There are no Rotavirus, <i>Norovirus</i> GI and <i>Norovirus</i> GII. No infection caused by <i>Rotavirus</i> , <i>Norovirus</i> GI or <i>Norovirus</i> GII.
<b>2.</b>	+ GREEN-RED	+ GREEN-RED	+ GREEN-RED	There are <i>Rotavirus</i> , <i>Norovirus</i> GI and <i>Norovirus</i> GII presence. Infection caused by <i>Rotavirus</i> , <i>Norovirus</i> GI and <i>Norovirus</i> GII.
<b>3.</b>	+ GREEN-RED	+ GREEN-RED	- GREEN	There is <i>Rotavirus</i> and <i>Norovirus</i> GI presence. Infection caused by <i>Rotavirus</i> and <i>Norovirus</i> GI.
<b>4.</b>	+ GREEN-RED	- GREEN	+ GREEN-RED	There is <i>Rotavirus</i> and <i>Norovirus</i> GII presence. Infection caused by <i>Rotavirus</i> and <i>Norovirus</i> GII.
<b>5.</b>	+ GREEN-RED	- GREEN	- GREEN	There is <i>Rotavirus</i> presence. Infection caused by <i>Rotavirus</i> .
<b>6.</b>	- GREEN	+ GREEN-RED	+ GREEN-RED	There is <i>Norovirus</i> GI and <i>Norovirus</i> GII presence. Infection caused by <i>Norovirus</i> GI and <i>Norovirus</i> GII.
<b>7.</b>	- GREEN	+ GREEN-RED	- GREEN	There is <i>Norovirus</i> GI presence. Infection caused by <i>Norovirus</i> GI.
<b>8.</b>	- GREEN	- GREEN	+ GREEN-RED	There is <i>Norovirus</i> GII presence. Infection caused by <i>Norovirus</i> GII.
<b>9.</b>	Any other result			Invalid result either A, B or C, we recommend repeating the assay using the same sample with another test.

**INVALID:** Total absence of any control coloured line (GREEN) regardless the appearance or not of the test lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

## NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

## QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

## LIMITATIONS

- The test must be carried out within 2 hours after opening the sealed bag.
- An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
- The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
- After one week of infection, the number of viruses in faeces is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms.
- CerTest *Rotavirus+Norovirus* GI+GII should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
- Positive results determine the presence of *Rotavirus*, *Norovirus* GI and *Norovirus* GII in faecal samples, nevertheless, a positive result should be followed up with additional laboratory techniques (biochemical methods, PCR, microscopy) to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
- A negative result is not meaningful because of it is possible the antigens concentration in the stool samples is lower than the detection limit value. If the symptoms or situation still persist, a *Rotavirus* and a *Norovirus* determination should be carried out with another technique (for example microscopy).

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## EXPECTED VALUES

Acute gastroenteritis is a common disorder in young children; moreover, the associated dehydration is a leading cause of admission to hospital in industrialized countries. Acute diarrhoea is a major health problem throughout the world and main source of mortality in developing countries.

Enteric viruses have been recognized as the most significant etiological agents of the disease and yet, four categories of viruses are being considered clinically relevant: group A *Rotavirus* (family *Reoviridae*), *Norovirus* (family *Caliciviridae*), *Adenovirus* and *Astrovirus*. Several studies proved co-infections in infants 46% with acute watery diarrhoea.

In recent years, global public health concern caused by *Norovirus* was raised referring to significant morbidity and occasional mortality, which frequently occurred in schools, hospitals, cruises and other semi-closed institutions. The pathogen is highly contagious with a low infectious dose and spreads by ingestion of contaminated food and/or water or inhaling the agent particle by contact with patients. It has been deduced that approximately 90% of the outbreaks of acute nonbacterial gastroenteritis in the United States were caused by *Norovirus* and 23 million cases were affected annually. In addition to human cost, *Norovirus* infections cause economic losses that approach or exceed \$650,000 in supplies and staff time off. Although previous studies showed the disease was mild and self-limiting, recent studies have showed its ability in causing more severe complications than previously expected.

## PERFORMANCE CHARACTERISTICS

### Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed using an immunochromatographic test (CerTest *Rotavirus+Norovirus* GI+GII, CerTest) and these results were confirmed with commercial available tests (Ridascreen®*Rotavirus* ELISA Test, r-Biopharm) for *Rotavirus* and (Simple Norovirus, Operon) for *Norovirus* GI+GII and confirmed by PCR. The results were as follows:

		Ridascreen® <i>Rotavirus</i> ELISA Test		
		+	-	Total
IC Test: CerTest <i>Rotavirus+Norovirus</i> GI+GII ( <i>Rotavirus</i> )	+	18	1	19
	-	0	43	43
Total		18	44	62

	Sensitivity	Specificity	PPV	NPV
CerTest <i>Rotavirus+Norovirus</i> GI+GII ( <i>Rotavirus</i> ) vs Ridascreen® <i>Rotavirus</i> ELISA Test	>99%	98%	>94%	>99%

		IC test: Simple Norovirus		
		+	-	Total
IC test: CerTest <i>Rotavirus+Norovirus</i> GI+GII ( <i>Norovirus</i> GI)	+	2	0	2
	-	0	48	48
Total		2	48	50

	Sensitivity	Specificity	PPV	NPV
CerTest <i>Rotavirus+Norovirus</i> GI+GII ( <i>Norovirus</i> GI) vs Simple Norovirus	>99%	>99%	>99%	>99%

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The results showed a high sensitivity and specificity to detect *Rotavirus*, genogroup I and II *Noroviruses* using CerTest *Rotavirus+Norovirus* GI+GII.

### Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *Rotavirus+Norovirus* GI+GII; no cross reactivity against gastrointestinal pathogens occasionally present in faeces:

<i>Adenovirus</i>	<i>Hepatitis A</i>	<i>Salmonella paratyphi</i>
<i>Astrovirus</i>	<i>Giardia lamblia</i>	<i>Salmonella typhimurium</i>
<i>Campylobacter coli</i>	<i>Helicobacter pylori</i>	<i>Salmonella typhi</i>
<i>Campylobacter jejuni</i>	<i>Listeria monocytogenes</i>	<i>Shigella boydii</i>
<i>Clostridium difficile</i>	<i>Norovirus GI (Strip A,C)</i>	<i>Shigella dysenteriae</i>
<i>Cryptosporidium parvum</i>	<i>Norovirus GII (Strip A,B)</i>	<i>Shigella flexneri</i>
<i>Entamoeba histolytica</i>	<i>Rotavirus (Strip B,C)</i>	<i>Shigella sonnei</i>
<i>Enterovirus</i>	<i>RSV</i>	<i>Staphylococcus aureus</i>
<i>Escherichia coli</i> O157:H7	<i>Salmonella enteritidis</i>	<i>Yersinia enterocolitica</i>
<i>Escherichia coli</i> O111		

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1. CUKOR G., and BLACKLOW N. R., "Human Viral Gastroenteritis", Microbiological Reviews, Vol. 48 No 2, June 1984, pp. 157-179.
2. ESTES, M. K. and COHEN, J.; "Rotavirus Gene Structure and Function", Microbiological Reviews, Vol. 53 No 4, Dec. 1989, pp. 410-449.
3. PAI C. H., SHAHRABADI M. S., and INCE B., "Rapid Diagnosis of Rotavirus Gastroenteritis by a Commercial Latex Agglutination Test", Journal of Clinical Microbiology, Vol. 22 No 5, Nov. 1985, pp. 846-850.
4. CUKOR, G.,PERRON, D.M., and BLACKLOW, N. R.: "Detection of Rotavirus in Human Stools by Using Monoclonal Antibody", Journal of Clinical Microbiology, Vol. 19, 888-892.
5. NEEL K. KRISHNA, B.A, "Identification of Structural Domains Involved in Astrovirus Capsid Biology", Viral Immunol. 2005; 18(1): 17-26.
6. BON, F. et al. "Prevalence of group A Rotavirus, human calicivirus, Astrovirus type 40 and 41 infections among children with acute gastroenteritis in Dijon, France." J. Clin. Microbiol. 37 No 9 3055-3058 (1999).
7. LI-JUAN LIU, et al. "Identification of Norovirus as the Top Enteric Viruses Detected in Adult Cases with Acute Gastroenteritis". *Am. J. Trop. Med, Hyg*, 82(4), 2010, p. 717-722.
8. GIUSEPPINA LA ROSA, et al. "Recreational and drinking water as a source of norovirus gastroenteritis outbreaks: a review and update". *Environmental biotechnology* 4 (1) 2008, p. 15-24.
9. THOMAS R. FRIEDEN, et al. "Emergence of new norovirus strain GII.4 Sydney-United States, 2012". *MMWR Morbidity and Mortality Weekly Rep.* 2013 Jan 25;62 (3):55.

## SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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