CERTEST

Influenza A

ONE STEP *Influenza A*
CARD TEST

CERTEST BIOTEC S.L.

For Information Only

Purposes

For Information Purposes Only

ENGLISH**INTENDED USE**

CerTest *Influenza A* one step card test is a coloured chromatographic immunoassay for the qualitative detection of *Influenza* type A from nasal swab, nasopharyngeal wash or aspirate specimens.

CerTest *Influenza A* card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of *Influenza* type A respiratory infection.

INTRODUCTION

Influenza is caused by a virus that attacks mainly the upper respiratory tract - the nose, throat and bronchi and rarely also the lungs. The infection usually lasts for about a week. It is characterized by sudden onset of high fever, myalgia, headache and severe malaise, non-productive cough, sore throat and rhinitis. Most people recover within one to two weeks without requiring any medical treatment. In the very young, the elderly and people suffering from medical conditions such as lung diseases, diabetes, cancer, kidney or heart problems, *Influenza* poses a serious risk. In these people, the infection may lead to severe complications of underlying diseases, pneumonia and death.

The currently circulating *Influenza* viruses that cause human disease are divided into two groups: A and B. *Influenza A* has 3 subtypes which are important for humans: A (H3N2), A (H1N1) and A (H5N1), of which the former is currently associated with most deaths. *Influenza* viruses are defined by 2 different protein components, known as antigens, on the surface of the virus. They are spike-like features called haemagglutinin (H) and neuraminidase (N) components.

TEST PRINCIPLE

CerTest *Influenza A* is based on the principle of a qualitative immunochromatographic assay for the determination of *Influenza* type A from nasal swab, nasopharyngeal wash or aspirate specimens.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Influenza A* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Influenza A*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*Influenza A* monoclonal antibodies-red polystyrene microspheres), which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*Influenza A* antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no *Influenza A* antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-*Influenza A* antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appear. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contains preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the CerTest *Influenza A* test. Do not use any other commercial kit component.

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- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

SPECIMEN COLLECTION AND PREPARATION

Samples should be process as soon as possible after collection. If this is not possible, the samples can be stored in the refrigerator (2-8 °C) for 8 hours prior to testing.

Specimen preparation (see illustration):

Nasal swab method:

1. Remove the swab from its packing.
2. Collect specimen with a sterile swab from one nostril.
3. Insert the swab into the nostril to the nasopharynx, rotating against the nasal wall (to ensure swab contains cells as well as mucus).
4. Repeat procedure using other nostril.
5. Process the swab as soon as possible after collecting the specimen.

Specimen collection

Nasopharyngeal aspirate method (suction apparatus, sterile suction catheter):

For adult:

1. Place the irrigator up to the nose.
2. Let the sterile saline water run into the nose (2.5 mL). It will run out the opposite side. Tilt and twist the irrigator side to side and up and down directing the water flow into all portions of the nasal cavity.
3. Collect the wash in a clean specimen container, tilt the head forward and allow the water with mucus to run out of the nostril into the specimen container. Repeat the mucus collection for the other nostril and collect it into the same container.

Specimen collection

For children:

1. Use an aspiration bulb or bulb syringe to instill the saline water into one nostril leaning the children head.
2. Aspirate the mix of mucus-saline water into the bulb and transfer it into a clean container.
3. Repeat for the other nostril and transfer the fluid into the same specimen container.

MATERIALS

MATERIALS PROVIDED

- CerTest *Influenza A* card tests
- Reagent B (sample diluent)
- Swabs
- Disposable pipettes
- Testing tubes
- Instructions for use
- CerTest *Influenza A* Positive Control swab + instructions for use

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer
- Vortex

TEST PROCEDURE

Allow tests, samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.


- Procedure A using nasal swab samples:

1. Add 15 drops (1) Reagent B and immediately put the swab into the tube.
2. Mix the solution by rotating the swab forcefully against the side of the tube at least 1 minute. Best results are obtained when the specimen is vigorously extracted in the solution (2). Extract as much liquid as possible from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
3. Remove CerTest *Influenza A* card test from its sealed bag just before using it.
4. Use a separate pipette and test for each sample or control. Dispense exactly 4 drops from the testing tube, into the circular window marked with the letter S (3).

5. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.


If the test does not run due to the type of sample, stir the sample added in the sample window (S) with the pipette. If it doesn't work, dispense a drop of Reagent B until seeing the liquid running through the reaction zone.


- Procedure B using nasopharyngeal wash or aspirate samples:

1. Add 6 drops (1) of the nasopharyngeal wash or aspirate samples with a pipette and 9 drops (2) of Reagent B in a testing tube. Mixer with vortex at least 1 minute to homogenize. Best results are obtained when the specimen is vigorously extracted in the solution (3).
2. Remove CerTest *Influenza A* card test from its sealed bag just before using it.
3. Dispense exactly 4 drops from the testing tube, in the circular window marked with the letter S (4).
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to the type of sample, stir the sample added in the sample window (S) with the pipette. If it doesn't work, dispense a drop of Reagent B until seeing the liquid running through the reaction zone.

INTERPRETATION OF THE RESULTS (please refer to the illustration below)

| | Influenza A | Interpretation of the results |
|-----------|--------------------|---|
| 1. | - | There is no <i>Influenza</i> type A presence. No infection caused by <i>Influenza</i> type A. |
| | GREEN | |
| 2. | + | There is <i>Influenza</i> type A presence. Infection caused by <i>Influenza</i> type A (FLU A). |
| | GREEN-RED | |
| 3. | Any other result | Invalid result, we recommend repeating the assay using the same sample with another test. |

NEGATIVE: Only the GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window

INVALID: Total absence of any control coloured line (GREEN) regardless the appearance or not of the test lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. CerTest *Influenza A* should be used only with nasal swab, nasopharyngeal wash and aspirate samples. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper nasal specimens must be obtained.
3. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
4. Positive results determine the presence of *Influenza* type A. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
5. A negative result is not meaningful because of it is possible the antigens concentration in the nasal sample is lower than the detection limit value. If the symptoms or situation still persist, it is recommended that all negative results undergo confirmatory testing using other method and/or virus identification by cell culture or PCR.

EXPECTED VALUES

Yearly *Influenza* epidemics can seriously affect all age groups, but the highest risk of complications occur among children younger than age two, adults age 65 or older and people of any age with certain medical conditions, such as chronic heart, lung, kidney, liver, blood or metabolic diseases (such as diabetes) or weakened immune systems.

Influenza affects between 20% and 30% of the United States population annually. Three fourths of the infected individuals develop an acute respiratory illness and one third of these seek medical attention. In a typical year, more than 20,000 to 40,000 Americans die and more than 100,000 are hospitalized because of complications related to influenza.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity

An evaluation, with nasal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest *Influenza A*, CerTest) and another commercial available immunochromatographic test (BinaxNOW® Influenza A&B). The results were as follows:

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The results showed a high sensitivity and specificity to detect *Influenza* type A using CerTest *Influenza A*.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *Influenza A*; no cross reactivity against organisms that cause other respiratory infections:

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Virus Respiratorio Sincitial

Adenovirus

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1. Daniel E. Noyola, et al. Comparison of a New Neuraminidase Detection Assay with an Enzyme Immunoassay, Immunofluorescence, and Culture for Rapid Detection of Influenza A and B Viruses in Nasal Wash Specimens. *Journal of Clinical Microbiology*, Mar. 2000, p. 1161–1165 Vol. 38, No. 3
2. Adriana Weinberg and Miranda L. Walker. Evaluation of Three Immunoassay Kits for Rapid Detection of Influenza Virus A and B. *Clinical And Diagnostic Laboratory Immunology*, Mar. 2005, p. 367–370 Vol. 12, No. 3
3. Andreea C. Cazacu, Gail J. Demmler et al. Comparison of a New Lateral-Flow Chromatographic Membrane Immunoassay to Viral Culture for Rapid Detection and Differentiation of Influenza A and B Viruses in Respiratory Specimens. *Journal of Clinical Microbiology*, Aug. 2004, p. 3661–3664 Vol. 42, No. 8.
4. Kenneth E. Irmen and James J. Kelleher. Use of Monoclonal Antibodies for Rapid Diagnosis of Respiratory Viruses in a Community Hospital. *Clinical and Diagnostic Laboratory Immunology*, May 2000, p. 396–403 Vol. 7, No. 3.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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