# VIASURE

## Real Time PCR Detection Kits

by CerTest  
BIOTEC

### *Sexual Health Panel I*

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VIASURE Sexual Health Panel I Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-SP0112L
VIASURE Sexual Health Panel I Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-SP0112H
VIASURE Sexual Health Panel I Real Time PCR Detection Kit 96-well plate, low profile	VS-SP0113L
VIASURE Sexual Health Panel I Real Time PCR Detection Kit 96-well plate, high profile	VS-SP0113H

For Research Use Only (RUO)

*This product has no declared clinical intended purpose and is not for clinical diagnostic use. No claim or representation is intended to provide information for the diagnosis, prevention, or treatment of a disease.*

## ENGLISH

### 1. Intended use

VIASURE Sexual Health Panel I Real Time PCR Detection Kit is designed for the specific identification and differentiation of *Neisseria gonorrhoeae*, *Chlamydia trachomatis*, *Mycoplasma genitalium*, *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and/or *Mycoplasma hominis*; Herpes virus 1, Herpes virus 2 and/or *Treponema pallidum*; *Candida albicans*, *Gardnerella vaginalis* and/or *Trichomonas vaginalis* in urogenital specimens, urine samples and rectal samples from patients with signs and symptoms of sexually transmitted diseases (STDs). This test is intended to be used for research purposes, without any medical objective is not regarded as devices for performance evaluation. DNA is extracted from clinical specimens, multiplied using Real Time amplification and detected using specific primers and a fluorescent reporter dye probe for *N. gonorrhoeae*, *C. trachomatis*, *M. genitalium*, *T. vaginalis*, *U. urealyticum*, *U. parvum* and *M. hominis*; HSV 1, HSV 2 and *T. pallidum*; *C. albicans*, *G. vaginalis* and *T. vaginalis*.

### 2. Summary and Explanation

Sexually transmitted infections (STIs) represent a group of diseases that affect the sexual and reproductive health of millions of people, being a public problem of interest. Etiological agents responsible for STIs include fungi, bacteria, parasites and viruses. Some of these microorganisms are eliminated after a period of time, while others are recurrent and some remain in the body asymptotically, allowing the progress of the disease and generating consequences such as inflammations of the genito-urinary tract, infertility and even the development of cancer.

*Trichomonas vaginalis* infection is one of the most common sexually transmitted diseases (STDs) in the world. *Trichomonas vaginalis* is a flagellated pathogen protozoan belonging to the order *Trichomonadida* located in both, male and female urogenital tract but it has also been isolated from the respiratory tract of infants and adults. In women it can be found in the vagina and in the urethra, while in men it can be found in the urethra, the prostate and the epididymis. *Trichomonas vaginalis* infection has been associated with vaginitis, cervicitis and urethritis, premature rupture of membranes and premature delivery in pregnant women. *Trichomonas vaginalis* infection has also been associated with an increased risk of HIV acquisition and transmission in women.

*Mycoplasma hominis* colonizes the lower urogenital tract and is associated with urogenital infections, particularly bacterial vaginosis and non-gonococcal urethritis. It is also involved in extra genital infections, such as postpartum or post-abortion fever, in post-cesarean wound infections or after a hysterectomy. In neonates, it can cause meningitis, brain abscesses and eye infections. In adults, bacteremia, septic arthritis, osteitis, endocarditis, mediastinitis, brain abscesses and respiratory infections have been described. Most patients have predisposing factors, including immunosuppression, trauma, respiratory problems or post-manipulation and/or surgery of the genitourinary tract. It is considered a bacterium present in the normal micro flora in 20% of the male population and in 40% of the female population.

*Mycoplasma genitalium* is a facultative anaerobic organism and a recognized cause of nongonococcal urethritis in men. In women, *M. genitalium* has been associated with cervicitis, endometritis, pelvic inflammatory disease

(PID), infertility, susceptibility to human immunodeficiency virus (HIV), and adverse birth outcomes, indicating a consistent relationship with female genital tract pathology. The long-term reproductive consequences of *M. genitalium* infection in asymptomatic individuals need to be investigated further.

*Ureaplasma* spp. are often isolated from human genital mucosa of individuals with a lack of symptoms. In humans, two major species, namely *Ureaplasma parvum* and *Ureaplasma urealyticum* (family *Mycoplasmataceae*, genus *Ureaplasma*) are part of the genital flora of men and women and are present in almost 70% of sexually active population. These bacteria cause inflammation and lead to chorioamnionitis, preterm deliveries, and premature rupture of membranes. Recently, their pathogenic roles on spermatogenesis and subsequent male infertility have been also suggested. Although both male and female carriers are asymptomatic, *Ureaplasma* spp. is occasionally isolated from the neonatal lesions of bronchopulmonary dysplasia, intraventricular haemorrhage, and necrotizing enterocolitis.

*Neisseria gonorrhoeae* is an obligate human pathogen and is the etiological agent of gonorrhea. Syndromes include cervicitis in women, and urethritis, pharyngitis and proctitis in both sexes. If untreated, women may experience severe sequelae of pelvic inflammatory disease, chronic pelvic pain, ectopic pregnancy and tubal infertility, while men may develop epididymitis, prostatitis and urethral stricture.

*Chlamydia trachomatis* is a small intracellular bacterium that requires living cells to multiply. There are 18 serotypes; D-K are the ones that cause sexually transmitted infections, as well as neonatal infections. The majority of genital chlamydial infections in both males and females are asymptomatic. When symptoms do occur, lower urogenital tract infection can manifest as cervicitis in females and urethritis in males and females. Whether symptomatic or asymptomatic, untreated chlamydia can ascend to the upper genital tract. In males, this can cause epididymitis, which is not thought to be an important cause of long-term sequelae. However, in females, upper tract infection can result in pelvic inflammatory disease (PID), a spectrum of clinical disorders involving infection and inflammation of the uterus, fallopian tubes, ovaries, or adjacent peritoneum.

Genital herpes can be caused by two very similar viruses, herpes simplex virus HSV-1 or HSV-2, being the infections caused by first one more often. Herpes simplex virus (HSV)-1 and -2 are large, double-stranded DNA viruses that cause lifelong persistent infections characterized by periods of quiescence and recurrent disease. These two HSV types cannot be distinguished clinically. In fact, they share a high degree of genetic homology, but they also have specific regions with small nucleotide variations which may allow discrimination. HSV-1 and HSV-2 infection occurs via inoculation of virus particles into susceptible mucosal surfaces. Afterwards, these neurotropic viruses can become latent in the local sensory ganglion, periodically reactivating to cause symptomatic lesions, or undergo asymptomatic viral release, with the potential for disease transmission and infection. The same treatment is used for both HSV-1 and HSV-2 infections, the location of the lesions and the chronicity of the infection (primary or recurrent) determine dosage and frequency.

The spirochete *Treponema pallidum*, the etiologic agent of syphilis, causes a multistage sexually transmitted infection (STI). Pathogenic treponemes cause venereal syphilis, yaws, endemic syphilis, and pinta—multistage, infections that, although similar, can be differentiated based on clinical, epidemiologic, and geographic criteria. Only venereal syphilis is transmitted by sexual activity. The pathogenic treponemes are uncultivable, slow-growing microorganisms with identical flat-wave morphologies. They poorly tolerate desiccation, elevated

temperature, and ambient oxygen tension, traits that explain why efficient transmission requires close personal contact.

*Candida albicans*, is a common opportunistic fungal pathogen that dwells on human mucosal surfaces, cancause fungal infections, especially in immunocompromised and high-risk surgical patients. Although *Candida albicans* is still the major species isolated from clinical samples in the majority of individuals, it is well known that some other non-albicans *Candida* spp. such as *Candida glabrata*, *Candida krusei*, *Candida parapsilosis*, and *Candida tropicalis* infections are significantly widespread.

*Gardnerella vaginalis* is a facultative Gram-variable anaerobic bacterium belonging to the family *Bifidobacteriaceae* and is the dominant species involved in bacterial vaginosis (BV). Although BV is a polymicrobial condition with no single causative agent, 100 % of women with BV are colonized with *G. vaginalis*. *G. vaginalis* is also a normal commensal of the vaginal micro flora in healthy women.

Conventional diagnostic assays, lack sensitivity, require viable organisms and thus special shipment conditions and, sometimes, invasive sampling. As nucleic acid amplification tests, as real time PCR, allow us to overcome some of these limitations. Several molecular diagnostic assays have recently been commercialized to assist the syndromic diagnosis of STIs. In addition, their implicit multiplexing capacity allows for the detection of multiple pathogens in a single sample.

### 3. Principle of the procedure

VIASURE Sexual Health Panel / Real Time PCR Detection Kit is designed for the identification and differentiation of *N. gonorrhoeae*, *C. trachomatis*, *M. genitalium*, *T. vaginalis*, *U. urealyticum*, *U. parvum* and/or *M. hominis*; HSV 1, HSV 2 and/or *T. pallidum*; *C. albicans*, *G. vaginalis* and/or *T. vaginalis* in urogenital specimens, urine samples and rectal samples. After DNA isolation, the identification of *N. gonorrhoeae*, *C. trachomatis*, *M. genitalium*, *T. vaginalis*, *U. urealyticum*, *U. parvum* and/or *M. hominis*; HSV 1, HSV 2 and/or *T. pallidum*; *C. albicans*, *G. vaginalis* and/or *T. vaginalis* is performed by the amplification of a conserved region of the specific genes (Table 3, Pathogens detected in each reaction tube and target genes)using specific primers and fluorescent-labelled probes.

VIASURE Sexual Health Panel / Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE Sexual Health Panel / Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase) in a stabilized format, as well as an internal control to monitor PCR inhibition (if present). Each DNA targets are amplified and detected in specific channels (FAM, HEX, ROX, and/or Cy5) and the internal control (IC) in HEX and/or Cy5 channels (if present) (depending on the equipment used select the proper detection channel, see Annex 2) (Table 5).

## 4. Reagents provided

VIASURE Sexual Health Panel / Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables 1 and 2:

Reagent/Material	Description	Colour	Amount
<i>Sexual Health Panel I</i> 8-well strips (2 Hyperplex reaction)	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control (if present) in stabilized format	White	6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Sexual Health Panel I</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNAse free	RNAse/DNAse free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE Sexual Health Panel / Real Time PCR Detection Kit with Ref. VS-SP0112L and VS-SP0112H.

Reagent/Material	Description	Colour	Amount
<i>Sexual Health Panel I</i> 96-well plate (2 Hyperplex reaction)	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control (if present) in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Sexual Health Panel I</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNAse free	RNAse/DNAse free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	12 X 8-cap strip

Table 2. Reagents and materials provided in VIASURE Sexual Health Panel / Real Time PCR Detection Kit with Ref. VS-SP0113L and VS-SP0113H.

*Sexual Health Panel I* 8-well strip (2 Hyperplex) contains in each reaction well the following reaction mixes for the detection of the specific targets (Table 3):

	Code	Controls	Reaction mix placed into each well Pathogens and target genes
1	NMT	IC	<i>Neisseria gonorrhoeae, Chlamydia trachomatis &amp; Mycoplasma genitalium</i>
			<i>Neisseria gonorrhoeae</i> ( <i>porA</i> and <i>Opa</i> genes), <i>Chlamydia trachomatis</i> (a region within ORF2 of the chlamydial plasmid) and <i>Mycoplasma genitalium</i> ( <i>MgPa adhesin gene</i> )
2	URE	-	<i>Trichomonas vaginalis, Ureaplasma urealyticum, Ureaplasma parvum &amp; Mycoplasma hominis</i>
			<i>Trichomonas vaginalis</i> ( <i>T. vaginalis</i> -specific 2-kb repeated sequence), <i>Ureaplasma urealyticum</i> and <i>Ureaplasma parvum</i> ( <i>ureasa gene</i> ), <i>Mycoplasma hominis</i> ( <i>yidC gene</i> )
3	HHT	IC	<i>Herpes virus 1, Herpes virus 2 &amp; Treponema pallidum</i>
			<i>Herpes virus 1</i> ( <i>US4gene</i> ), <i>Herpes virus 2</i> ( <i>US6gene</i> ) and <i>Treponema pallidum</i> (16S rRNA gene)
4	CGT	IC	<i>C. albicans, G. vaginalis &amp; T. vaginalis</i>
			<i>Trichomonas vaginalis</i> ( <i>T. vaginalis</i> -specific 2-kb repeated sequence), <i>Candida albicans</i> (5.8S rRNA gene) and <i>Gardnerella vaginalis</i> (16S rRNA gene)

Table 3. *Sexual Health Panel I* 8-well strips provided in VIASURE *Sexual Health Panel I* Real Time PCR Detection Kit. Reaction mix placed into each well, pathogens detected and target genes. IC: Internal control. Note that the first reaction mix is placed into the ends of the strip.

## 5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Sexual Health Panel I* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortexer.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.

VIASURE *Sexual Health Panel I* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTlite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

## 6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

## 7. Precautions for users

- This VIASURE Sexual Health Panel / Real Time PCR Detection kit is for Research Use Only.
- The product is intended for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-SP0113L and VS-SP0113H). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- **Note that the first reaction mix is placed into the ends of the strip.**
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.

## 8. Test procedure

### 8.1. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For DNA extraction from clinical samples you can use your manual or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- Maxwell® 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® 16 instrument (Promega).
- ZP02006 MagPurix Bacterial DNA Extraction Kit, using the MagPurix 12A instrument (Zinexts Life Science Corp.).
- US121 NX-48 Urine/Swab DNA Kit using the Nextractor® NX-48 (Genolution Inc.).
- MagDEA® Dx SV, using magLEAD® 6gC (Precision System Science)).
- NucleoSpin ® RNA Virus (MACHEREY-NAGEL).
- NucleoMag® Pathogen (MACHEREY-NAGEL).
- Invisorb® Spin Universal Kit (Stratec).

### 8.2. Lyophilized positive control

*Sexual Health Panel / Positive Control* contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Sexual Health Panel / Positive Control* (red vial) by adding 200 µL of the supplied Water RNase/DNAse free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

### 8.3. PCR protocol

Determine and separate the number of required strips, including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plate or strips. Note that the first reaction mix is placed into the ends of the strip (Table 3) and each strip includes two Hyperplex reactions, wells 1 to 4: NMT, URE, HHT and CGT reaction mixes.

- 1) Reconstitute the number of strips (Hyperplex reactions) you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted Sexual Health Panel / Positive Control (red vial) or Negative Control (violet vial) in the 8 wells of each strip and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips or 96-well plate.

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 4. PCR protocol

Fluorogenic data should be collected during the extension step (\*) through the FAM, ROX, Cy5 and/or HEX (JOE or VIC) channels following Table 5. Depending on the equipment used select the proper detection channel (see Annex 2).

The reaction mixes places in each well of Hyperplex strip allow the detection of the specific target pathogens in the following channels (Table 5).

Code	Channels			
	FAM	HEX	ROX	Cy5
1 NMT	<i>Chlamydia trachomatis</i>	<i>Mycoplasma genitalium</i>	<i>Neisseria gonorrhoeae</i>	Internal Control
2 URE	<i>Trichomonas vaginalis</i>	<i>Ureaplasma urealyticum</i>	<i>Ureaplasma parvum</i>	<i>Mycoplasma hominis</i>
3 HHT	Herpes virus 1	Internal Control	Herpes virus 2	<i>Treponema pallidum</i>
4 CGT	<i>Trichomonas vaginalis</i>	Internal Control	<i>Candida albicans</i>	<i>Gardnerella vaginalis</i>
4 CGT	<i>Trichomonas vaginalis</i>	Internal Control	<i>Candida albicans</i>	<i>Gardnerella vaginalis</i>
3 HHT	Herpes virus 1	Internal Control	Herpes virus 2	<i>Treponema pallidum</i>
2 URE	<i>Trichomonas vaginalis</i>	<i>Ureaplasma urealyticum</i>	<i>Ureaplasma parvum</i>	<i>Mycoplasma hominis</i>
1 NMT	<i>Chlamydia trachomatis</i>	<i>Mycoplasma genitalium</i>	<i>Neisseria gonorrhoeae</i>	Internal Control

Table 5. Pathogens identified in each detection channel.

## 9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control of each well from Hyperplex reaction (4 wells) and the presence of signal for *Sexual Health Panel I* positive control of each well from Hyperplex reaction (4 wells). Check Internal Control signal (with the exception of URE reaction mix) to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

A sample is considered positive if the Ct value obtained is less than 40 and the internal control (if available) shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive (if available). An inhibition of the PCR reaction can be excluded by the amplification of internal control (if available). Using the following table read and analyze the results:

Reaction Mix	Pathogens	Channels			
		FAM	HEX	ROX	Cy5
1	NMT	<i>Chlamydia trachomatis</i>	Positive		
		<i>Mycoplasma genitalium</i>		Positive	
		<i>Neisseria gonorrhoeae</i>			Positive
		Internal Control			Positive/Negative
2	URE	<i>Trichomonas vaginalis</i> *	Positive		
		<i>Ureaplasma urealyticum</i>		Positive	
		<i>Ureaplasma parvum</i>			Positive
		<i>Mycoplasma hominis</i>			Positive
3	HHT	Herpes virus 1	Positive		
		Internal Control		Positive/Negative	
		Herpes virus 2			Positive
		<i>Treponema pallidum</i>			Positive
4	CGT	<i>Trichomonas vaginalis</i> *	Positive		
		Internal Control		Positive/Negative	
		<i>Candida albicans</i>			Positive
		<i>Gardnerella vaginalis</i>			Positive

Table 6. Sample interpretation. Positive: Amplification curve. Empty: No amplification curve.

\* Any patient sample displaying an amplification curve for *Trichomonas vaginalis* with URE and/or CGT reaction mixes should be considered as positive for *Trichomonas vaginalis*.


Table 7. Correct run of negative and positive controls run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

The result is considered invalid if there is signal of amplification in negative controls or absence of signal in the positive wells. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells (with the exception of URE reaction mix), we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

## 10. Limitations of the test

- The test must be used for RUO purposes (Exclusive Use in Research), without any medical objective.
- Although this assay can be used with other types of samples, it has been validated only with STD swabs (urogenital male or genital female origin), endocervical specimens, urethral samples, urogenital specimens,

rectal-vaginal specimens, rectal sample, vaginal specimen, urine samples, serum sample, semen and throat samples (for URE and NMT); samples in transport medium, endocervical swabs, urine, plasma, cell lysates and synthetic CSF (for HHT); and endocervical swab, vaginal swab, and vagino-rectal swabs, urine, plasma, brochoalveolar lavage (BAL) saline, sputum and samples in transport medium (for CGT).

- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *N. gonorrhoeae*, *C. trachomatis*, *M. genitalium*, *T. vaginalis*, *U. urealyticum*, *U. parvum* and/or *M. hominis*; HSV 1, HSV 2 and/or *T. pallidum*; *C. albicans*, *G. vaginalis* and/or *T. vaginalis*, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

## 11. Quality control

VIASURE Sexual Health Panel / Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique (with the exception of URE reaction mix).

## 12. Preliminary results of the test characteristics

### 12.1. Clinical sensitivity and specificity

The clinical performance of ***Sexually Transmitted Diseases Reaction mix*** (*Neisseria gonorrhoeae*, *Chlamydia trachomatis* & *Mycoplasma genitalium* and *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* & *Mycoplasma hominis*) was tested using 130 different specimens (endocervical specimens, urogenital specimens, urethral samples, vaginal specimens, rectal-vaginal specimens, urine samples, serum samples, semen and throat samples) from symptomatic patients.

These results of both Multiplex reactions were compared with those obtained with a molecular detection method (FTD Urethritis plus (Fast Track Diagnostics)).

The results were as follows:

<i>Sexually Transmitted Diseases Reaction mix</i>	FTD Urethritis plus (Fast Track Diagnostics)			
		+	-	Total
+	12	2*‡	14	
-	4*#	112	116	
Total	16	114	130	

Table 8. Comparative results for *Neisseria gonorrhoeae*.

\*The low amount of template DNA in these genitourinary samples is below the detection limit of the method used. †1/2 sample was confirmed by Seegene test as positive. # Neither identified by routine diagnosis.

Sexually Transmitted Diseases Reaction mix	FTD Urethritis plus (Fast Track Diagnostics)			
		+	-	Total
	+	37	4*‡	41
	-	3#	86	89
	Total	40	90	130

Table 9. Comparative results for *Chlamydia trachomatis*.

\*The low amount of template DNA in these genitourinary samples is below the detection limit of the method used. ‡ 4/4 samples were confirmed by routine diagnosis. # 2/3 identified by routine diagnosis.

Sexually Transmitted Diseases Reaction mix	FTD Urethritis plus (Fast Track Diagnostics)			
		+	-	Total
	+	7	0	7
	-	1*	122	123
	Total	8	122	130

Table 10. Comparative results for *Mycoplasma genitalium*.

\*The low amount of template DNA in this genitourinary sample is below the detection limit of the method used.

Sexually Transmitted Diseases Reaction mix	FTD Urethritis plus (Fast Track Diagnostics)			
		+	-	Total
	+	39	0	39
	-	3*#	88	91
	Total	42	88	130

Table 11. Comparative results for *Trichomonas vaginalis*.

\*The low amount of template DNA in these genitourinary samples is below the detection limit of the method used. # 2/3 identified by routine diagnosis.

Sexually Transmitted Diseases Reaction mix	FTD Urethritis plus (Fast Track Diagnostics)			
		+	-	Total
	+	26	4*‡	30
	-	3*#	97	100
	Total	29	101	130

Table 12. Comparative results for *Ureaplasma urealyticum*.

\*The low amount of template DNA in these genitourinary samples is below the detection limit of the method used. ‡ 2/4 samples were confirmed by routine diagnosis as positive. # 1/3 identified positive by routine diagnosis.

Sexually Transmitted Diseases Reaction mix	FTD Urethritis plus (Fast Track Diagnostics)			
		+	-	Total
	+	53	16*	69
	-	0	61	61
	Total	53	77	130

Table 13. Comparative results for *Ureaplasma parvum*.

\*The low amount of template DNA in these genitourinary samples is below the detection limit of the method used. These Real Time assays might be show different analytical sensitivity.

<b>Sexually Transmitted Diseases Reaction mix</b>	FTD Urethritis plus (Fast Track Diagnostics)			
		+	-	Total
	+	51	0	51
	-	1*	78	79
	Total	52	78	130

Table 14. Comparative results for *Mycoplasma hominis*.

\*The low amount of template DNA in these genitourinary samples is below the detection limit of the method used.

The clinical performance of *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Reaction mix was tested using 8 endocervical swabs from symptomatic patients. The comparison was carried out using the results obtained from the routine analysis performed in the hospital with the molecular detection method "CLART® ENTHERPEX" (Genómica). Two samples were HSV-1 positive and 6 were detected as HSV-2 positive. The concordance of both kits was 100%. The clinical performance of *Herpes virus 1*, *Herpes virus 2* & *Treponema pallidum* Reaction mix was tested using 89 samples from QCMD, UK NEQAS and INSTAND EQA programs. These samples included urine, plasma, cell lysates, synthetic CSF, genital, anal or oropharyngeal ulcer-causing pathogens. These results were compared with those obtained with a molecular detection method (FTD Genital ulcer (Fast Track Diagnostics)) and/or with the corresponding reports. All samples were detected correctly. 23/81 were HSV-1 positive, 21/81 were HSV-2 positive and 15/81 were *Treponema pallidum* positive.

The clinical performance of *C. albicans*, *G. vaginalis* & *T. vaginalis* Reaction mix was tested using QCMD EQA programs, and the results were compared with the corresponding reports. The clinical performance for *Candida albicans* was tested with the QCMD 2017 *Candida* spp. EQA Programme and QCMD 2017 Immunocompromised EQA Pilot Study. All *C. albicans* samples (7 samples positives) were properly detected. Other species of *Candida*, including *C. glabrata*, *C. parapsilosis*, and *C. Krusei*, as well as *Candida* negative samples were reported as negative, and other opportunistic microorganisms but pathogens for immunocompromised patients as *Pneumocystis jirovecii*, *Aspergillus fumigatus*, and JC virus were reported as negative. The clinical performance for *Gardnerella vaginalis* and *Trichomonas vaginalis* was tested with the QCMD 2017 Sexually Transmitted Infections I EQA Pilot Study and the QCMD 2018 Trichomonas vaginalis EQA Pilot Study. All *G. vaginalis* (2 samples positives) and *T. vaginalis* (9 samples positives) samples were properly detected. Other sexually transmitted pathogens including *Mycoplasma hominis*, *Mycoplasma genitalium* and *Ureaplasma urealyticum*, as well as negative samples were reported as negative.

The results show a high sensitivity and specificity to detect *N. gonorrhoeae*, *C. trachomatis*, *M. genitalium*, *T. vaginalis*, *U. urealyticum*, *U. parvum* and *M. hominis*; HSV 1, HSV 2 and *T. pallidum*; *C. albicans*, *G. vaginalis* and *T. vaginalis*; using VIASURE Sexual Health Panel / Real Time PCR Detection Kit.

## 12.2. Analytical sensitivity

VIASURE Sexual Health Panel / Real Time PCR Detection Kit has a detection limit of  $\geq 10$  DNA copies per reaction for *N. gonorrhoeae*, *C. trachomatis*, *M. genitalium*, *T. vaginalis*, *U. urealyticum*, *U. parvum* and *M. hominis*; HSV 1, HSV 2 and *T. pallidum*; *C. albicans*, *G. vaginalis* and *T. vaginalis*.

### 12.3. Analytical specificity

The specificity of the *Sexual Health Panel I* assay was confirmed by testing a panel consisting of different microorganisms representing the most common enteric and genitourinary pathogens or flora present in the gastrointestinal tract or urogenital system. No cross-reactivity was detected between almost any of the following microorganisms tested, except the targeted pathogens of each assay.

Cross-reactivity testing					
<i>Candida parapsilosis</i>	-	<i>Enterococcus faecium</i>	-	<i>Pseudomonas aeruginosa</i>	-
<i>Candida tropicalis</i>	-	<i>Escherichia coli</i> 0.1285; O18:H7:K1	-	<i>Serratia marcescens</i>	-
<i>Candida glabrata</i>	-	<i>Gardnerella vaginalis</i>	-/+	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-
<i>Candida krusei</i>	-	<i>Haemophilus influenzae</i> MinnA	-	<i>Stenotrophomonas maltophilia</i>	-
<i>Candida dubliniensis</i>	-	<i>Haemophilus ducreyi</i> class 1	-	<i>Streptococcus agalactiae</i>	-
<i>Candida albicans</i>	-/+	<i>Klebsiella oxytoca</i>	-	<i>Streptococcus pneumoniae</i>	-
<i>Aspergillus fumigatus</i>	-	<i>Klebsiella pneumoniae</i>	-	<i>Treponema pallidum</i>	-/+
<i>Bacteroides fragilis</i>	-	<i>Listeria ivanovii</i>	-	<i>Trichomonas vaginalis</i>	-/+
<i>Acinetobacter baumannii</i>	-	<i>Listeria monocytogenes</i>	-	<i>Ureaplasma parvum</i>	-/+
<i>Chlamydia trachomatis</i> (LGV)	-/+	<i>Listeria innocua</i>	-	<i>Ureaplasma urealyticum</i>	-/+
<i>Chlamydia trachomatis</i> (SW)	-/+	<i>Mycoplasma genitalium</i>	-/+	Cytomegalovirus	-
<i>Chlamydia trachomatis</i> (Genovar F)	-/+	<i>Mycoplasma hominis</i>	-/+	Herpes simplex virus 1	-/+
<i>Enterobacter cloacae</i>	-	<i>Neisseria gonorrhoeae</i>	-/+	Herpes simplex virus 2	-/+
<i>Enterobacter aerogenes</i>	-	<i>Neisseria meningitidis</i>	-	Human papillomavirus 16	-
<i>Enterococcus faecalis</i>	-	<i>Proteus mirabilis</i>	-	Human papillomavirus 18	-

Table 27. Reference pathogenic microorganisms used in this study

### 12.4. Analytical reactivity

The reactivity of VIASURE *Sexual Health Panel I* Real Time PCR Detection Kit for ***Sexually Transmitted Diseases Reaction mix*** (*Neisseria gonorrhoeae*, *Chlamydia trachomatis* & *Mycoplasma genitalium* and *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* & *Mycoplasma hominis*) was evaluated against *Neisseria gonorrhoeae* (*Neisseria gonorrhoeae* strains St 49226 and Lvl Ng PorA), *Chlamydia trachomatis* (*Chlamydia trachomatis* (SW and LGV), Swedish strain, and Genovar F, and Serovars D, E, I, K and J), *Mycoplasma genitalium*, *Trichomonas vaginalis*, *Ureaplasma urealyticum*, *Ureaplasma parvum* and *Mycoplasma hominis*, showing positive result.

The reactivity of VIASURE ***Herpes virus 1, Herpes virus 2 & Treponema pallidum Reaction mix***, was evaluated against *Herpes virus 1*(HSV 1 MacIntyre strain), *Herpes virus 2* (HSV 2 MS strain) and *Treponema pallidum*, showing positive result.

The reactivity of VIASURE ***C. albicans, G. vaginalis & T. vaginalis Reaction mix***, was evaluated against *Candida albicans*, *Gardnerella vaginalis* and *Trichomonas vaginalis*, showing positive result.

## ANNEX 1

## COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast Real-Time PCR System <sup>(1)</sup>
Applied Biosystems	7500 Fast Dx Real-Time PCR System <sup>(1)</sup>
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	StepOne Plus™ Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	StepOne™ Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	Mini Opticon™ Real-Time PCR Detection System <sup>(6)</sup>
Cepheid	SmartCycler® <sup>(3)</sup>
Qiagen	Rotor-Gene® Q <sup>(3)</sup>
Roche	LightCycler® 480 Real-Time PCR System <sup>(4)</sup>
Roche	LightCycler® 96 Real-Time PCR System <sup>(4)</sup>
Roche	Cobas z480 Analyzer <sup>(4)</sup>

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	ABI PRISM 7000 <sup>(6)</sup>
Applied Biosystems	ABI PRISM 7700 <sup>(5)</sup>
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System <sup>(6)</sup>
Bio-Rad	MyiQ™2 Real-Time PCR Detection System <sup>(6)</sup>
Cepheid	SmartCycler® <sup>(3)</sup>
DNA-Technology	DTprime Real-time Detection Thermal Cycler <sup>(2)</sup>
DNA-Technology	DTlite Real-Time PCR System <sup>(2)</sup>
Eppendorf	Mastercycler™ ep realplex
Qiagen	Rotor-Gene® Q <sup>(3)</sup>
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System <sup>(2)</sup>
VIASURE	VIASURE 96 Real Time PCR System <sup>(2)</sup>

(1)Select Ramp Speed "Standard".

(2)See Annex 3 to check optical measurement exposure setting.

(3)The product should be reconstituted following the appropriate procedure (see Test Procedure, section 8.3) and transferred into the specific tubes designed to perform on Rotor-Gene® Q or SmartCycler® instruments.

(4)Shell Frame grid plate which fits in these Roche qPCR System is necessary.

(5)No detection in Cy5 channel.

(6)Detection in FAM and HEX channels only.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

## ANNEX 2

## DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline. Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation is required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	RQX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 Fl for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 µl) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

## ANNEX 3

### OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". Set exposition values as follow:

Universal exposition values are as follow:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -1000, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel – 500.

For "VIASURE Sexual Health Panel / Real Time PCR Detection Kit", it is recommend set the exposition values as follow:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -150, HEX channel – 3000, ROX channel – 2000 and Cy5 channel - 1500.
- DTlite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 150, HEX channel - 500, ROX channel – 500 and Cy5 channel – 100.

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<i>Sexually Transmitted Diseases Reaction mix</i>	FTD Urethritis plus (Fast Track Diagnostics)			
		+	-	Total
	+	53	16*	69
	-	0	61	61
	Total	53	77	130

Table 13. Comparativa de resultados para *Ureaplasma parvum*.

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<i>Sexually Transmitted Diseases Reaction mix</i>	FTD Urethritis plus (Fast Track Diagnostics)			
		+	-	Total
	+	51	0	51
	-	1*	78	79
	Total	52	78	130

Table 14. Comparativa de resultados para *Mycoplasma hominis*.

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## 13. Bibliography/Bibliografía

1. D. León *et al.* Molecular detection of sexually transmitted agents in a symptomatic group of men and its relationship with sexual behaviour. *Revista Chilena de Infectología* 2016; 33(5):505-512.
2. A. K. Wendel *et al.* Trichomonas vaginalis Polymerase Chain Reaction Compared with Standard Diagnostic. *Clinical Infectious Diseases* 2020; 35:576-80.
3. G. Kusdian *et al.* The biology of Trichomonas vaginalis in the light of urogenital tract infection. Molecular and Biochemical Parasitology Journal (2015).
4. T. Posse *et al.* Bacteriemia por Mycoplasma hominis: un agente etiológico subestimado. *Revista Argentina Microbiología*. 2017;02.009.
5. F. Kazumasa *et al.* Rapid and simple detection of Ureaplasma species from vaginal swab samples using a loop-mediated isothermal amplification method. *American Journal of Reproductive Immunology* 2017;e12771.
6. S. Ona *et al.* *Mycoplasma genitalium*: An Overlooked Sexually Transmitted Pathogen in Women? *Infectious Diseases in Obstetrics and Gynecology* 2016; 2016: 4513089.
7. N. Lai-King *et al.* The laboratory diagnosis of *Neisseria gonorrhoeae*. *Canadian Journal of Infectious Diseases and Medical Microbiology*. 2005 Jan-Feb; 16(1): 15-25.
8. CDC Grand Rounds: Chlamydia prevention: challenges and strategies for reducing disease burden and sequelae. *Centers for Disease Control and Prevention (CDC)*.
9. G. Fernández *et al.* Usefulness of a novel multiplex real-time PCR assay for the diagnosis of sexually-transmitted infections. *Enfermedades Infecciosas y Microbiología Clínica* 2016.
10. S. Alexander *et al.* Evaluation of strategies for confirming *Neisseria gonorrhoeae* nucleic acid amplification tests. *Journal of Medical Microbiology* 2011; 60: 909-912.
11. J. Schirm *et al.* *Trichomonas vaginalis* detection using real-time TaqMan PCR. *Journal of Microbiological Methods* 2007; 68(2): 243-247.
12. J. Yi *et al.* Detection and biovar discrimination of *Ureaplasma urealyticum* by real-time PCR. *Molecular and Cellular Probes* 2005; 19(4): 255-260.
13. J.S. Jensen *et al.* Use of TaqMan 5' nuclease real-time PCR for quantitative detection of *Mycoplasma genitalium* DNA in males with and without urethritis who were attendees at a sexually transmitted disease clinic. *Journal of Clinical Microbiology* 2004; 42(2): 683-692.
14. T.A. Halse *et al.* A multiplexed real-time PCR assay for rapid detection of *Chlamydia trachomatis* and identification of serovar L-2, the major cause of Lymphogranuloma venereum in New York. *Molecular and Cellular Probes* 2006 Oct;20(5):290-7. Epub 2006 Mar 6.

15. D. Férandon *et al.* Development of a real-time PCR targeting the *yidC* gene for the detection of *Mycoplasma hominis* and comparison with quantitative culture. *Clinical Microbiology and Infection* 2011; 17(2): 155-159.
16. J. D. Radolf *et al.* *Treponema pallidum*, the syphilis spirochete: making a living as a stealth pathogen. *Nature Reviews Microbiology* 2016; 14(12): 744-759.
17. M. Lieveld *et al.* A high resolution melting (HRM) technology-based assay for cost-efficient clinical detection and genotyping of herpes simplex virus (HSV)-1 and HSV-2. *Journal of Virological Methods* 2017; Volume 248, 181-186.
18. M. Costa-Silva *et al.* Cross-sectional study of *Treponema pallidum* PCR in diagnosis of primary and secondary syphilis. *International Journal of Dermatology* 2017; doi: 10.1111/ijd.13823.
19. M.A. Minaya *et al.* Molecular evolution of herpes simplex virus 2 complete genomes: Comparison between primary and recurrent infections. *American Society for Microbiology* 2017; doi:10.1128/JVI.00942-17.
20. C.O. Onyango *et al.* Evaluation of a TaqMan Array Card for Detection of Central Nervous System Infections. *Journal of Clinical Microbiology* 2017 Jul;55(7): 2035-2044.
21. A. K. Wendel *et al.* Trichomonas vaginalis Polymerase Chain Reaction Compared with Standard Diagnostic. *Clinical Infectious Diseases* 2020; 35:576-80.
22. G. Kusdian *et al.* The biology of Trichomonas vaginalis in the light of urogenital tract infection. *Molecular and Biochemical Parasitology Journal* (2015).
23. J. Schirm *et al.* *Trichomonas vaginalis* detection using real-time TaqMan PCR. *Journal of Microbiological Methods* 2007; 243-247.
24. H. Fukumoto *et al.* Development of a new real-time PCR system for simultaneous detection of bacteria and fungi in pathological samples. *International Journal of Clinical and Experimental Pathology* 2015; 8(11):15479-15488.
25. C. Cox *et al.* New assay for *Gardnerella vaginalis* loads correlates with Nugent scores and has potential in the diagnosis of bacterial vaginosis. *Journal of Medical Microbiology* 2015; 64, 978-984.
26. J.G. Kusters *et al.* A multiplex real-time PCR assay for routine diagnosis of bacterial vaginosis. *European Journal of Clinical Microbiology & Infectious Diseases* 2015; 34:1779-1785.
27. Y. Gong *et al.* *Candida albicans* Heat Shock Proteins and Hsps-Associated Signaling Pathways as Potential Antifungal Targets. *Frontiers in Cellular and Infection Microbiology* 2017; 7:520.
28. M. Shirani *et al.* Chemical composition and antifungal effect of hydroalcoholic extract of *Allium tripedale* (Tautv.) against *Candida* species. *Current Medical Mycology* 2017; (1): 6-12.

## 14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: January 2019

For information purposes only

CerTest Biotec, S.L.

Pol. Industrial Rio Gállego II · Calle J, Nº1  
50840, San Mateo de Gállego, Zaragoza (Spain)  
[www.certest.es](http://www.certest.es)

VIASURE online

F-362 rev01

VIASURE

Real Time PCR Detection Kit

CerTest  
BIOTEC

For information purposes only