VIASURE

Real Time PCR Detection Kits

by **CerTest**
BIOTEC

Shigella/EIEC

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| | |
|---|------------|
| VIASURE <i>Shigella/EIEC</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile | VS-SHY106L |
| VIASURE <i>Shigella/EIEC</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile | VS-SHY106H |
| VIASURE <i>Shigella/EIEC</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile | VS-SHY112L |
| VIASURE <i>Shigella/EIEC</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile | VS-SHY112H |
| VIASURE <i>Shigella/EIEC</i> Real Time PCR Detection Kit 96-well plate, low profile | VS-SHY113L |
| VIASURE <i>Shigella/EIEC</i> Real Time PCR Detection Kit 96-well plate, high profile | VS-SHY113H |
| VIASURE <i>Shigella/EIEC</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene® | VS-SHY136 |
| VIASURE <i>Shigella/EIEC</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene® | VS-SHY172 |

ENGLISH

1. Intended use

VIASURE *Shigella*/EIEC Real Time PCR Detection Kit is designed for specific identification of *Shigella* and enteroinvasive *Escherichia coli* (EIEC) in human stool samples from patients with signs and symptoms of *Shigella*/EIEC infection. This test is intended for use as an aid in the diagnosis of *Shigella*/EIEC infection in combination with clinical and epidemiological risk factors. DNA is extracted from stool specimens, multiplied using Real Time amplification and detected using specific primers and a fluorescent reporter dye probe for *Shigella*/EIEC.

2. Summary and Explanation

Shigella species are gram-negative organisms which annually cause an estimated 165 million cases of shigellosis worldwide, resulting in 1 million deaths. Shigellosis is a dysentery characterized by severely bloody and mucus-containing diarrhea. The disease is caused by any of the four species of *Shigella* (*S. dysenteriae*, *S. flexneri*, *S. boydii* and *S. sonnei*) and by enteroinvasive *Escherichia coli* (EIEC). The differentiation of *Shigella* and enteroinvasive *E. coli* is complicated due to the ability to cause dysentery using of the same method of invasion.

Shigella is a significant agent of foodborne illness, especially with foods that require hand processing and/or are prepared from raw or previously cooked products without re-heating. The low infectious dose (10 cells), allows the disease to be spread effectively by infected food or water, and also by person-to person contact.

Several virulence factors encoded on the chromosome and/or plasmids have been reported, among which the invasion plasmid antigen H encoded by *ipaH* gene helps diffusion of the pathogen from cell to cell. This gene has a multiple-copy nature and has been located on both the chromosome and the invasive plasmid in all *Shigella* species and in EIEC. For this reason, Real-Time PCR assays for *Shigella*/EIEC detection has been developed using the *ipaH* gene as target.

3. Principle of the procedure

VIASURE *Shigella*/EIEC Real Time PCR Detection Kit is designed for the diagnosis of shigellosis caused by *Shigella*/EIEC in human stool samples. After DNA isolation, the identification of *Shigella*/EIEC is performed by the amplification of a conserved region of the *ipaH* gene using specific primers and a fluorescent-labeled probe.

VIASURE *Shigella*/EIEC Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence could be measured on real time PCR platforms.

VIASURE *Shigella*/EIEC Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase) in a stabilized format, as well as an internal control to monitor PCR inhibition. *Shigella*/EIEC DNA targets are amplified and detected in FAM channel and the

internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

4. Reagents provided

VIASURE *Shigella*/EIEC Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables 1, 2 and 3. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table 1 includes materials and reagents to be used with 8-well strips compatible devices (See Annex 1). Table 2 includes materials and reagents to be used with 96-well plate compatible devices (See Annex 1). Table 3 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips.

| Reagent/Material | Description | Colour | Amount |
|--|--|-------------|---------------------|
| <i>Shigella</i> /EIEC 8-well strips | A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format | White | 6/12 x 8-well strip |
| Rehydration Buffer | Solution to reconstitute the stabilized product | Blue | 1 vial x 1.8 mL |
| <i>Shigella</i> /EIEC Positive Control | Non-infectious synthetic lyophilized DNA | Red | 1 vial |
| Negative control | Non template control | Violet | 1 vial x 1 mL |
| Water RNase/DNase free | RNase/DNase free water | White | 1 vial x 1 mL |
| Tear-off 8-cap strips | Optical caps for sealing wells during thermal cycling | Transparent | 6/12 X 8-cap strip |

Table 1. Reagents and materials provided in VIASURE *Shigella*/EIEC Real Time PCR Detection Kit with Ref. VS-SHY106L, VS-SHY106H, VS-SHY112L and VS-SHY112H.

| Reagent/Material | Description | Color | Amount |
|--|--|-------------|------------------|
| <i>Shigella</i> /EIEC 96-well plate | A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format | White | 1 plate |
| Rehydration Buffer | Solution to reconstitute the stabilized product | Blue | 1 vial x 1.8 mL |
| <i>Shigella</i> /EIEC Positive Control | Non-infectious synthetic lyophilized DNA | Red | 1 vial |
| Negative control | Non template control | Violet | 1 vial x 1 mL |
| Water RNase/DNase free | RNase/DNase free water | White | 1 vial x 1 mL |
| Tear-off 8-cap strips | Optical caps for sealing plate during thermal cycling | Transparent | 12 X 8-cap strip |

Table 2. Reagents and materials provided in VIASURE *Shigella*/EIEC Real Time PCR Detection Kit with Ref VS-SHY113L and VS-SHY113H.

| Reagent/Material | Description | Colour | Amount |
|---------------------------------------|--|-------------|---------------------|
| <i>Shigella/EIEC</i> 4-well strips | A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format | Transparent | 9/18 x 4-well strip |
| Rehydration Buffer | Solution to reconstitute the stabilized product | Blue | 1 vial x 1.8 mL |
| <i>Shigella/EIEC</i> Positive Control | Non-infectious synthetic lyophilized DNA | Red | 1 vial |
| Negative control | Non template control | Violet | 1 vial x 1 mL |
| Water RNase/DNase free | RNase/DNase free water | White | 1 vial x 1 mL |
| 4-cap strips | Optical caps for sealing wells during thermal cycling | Transparent | 9/18 X 4-cap strip |

Table 3. Reagents and materials provided in VIASURE *Shigella/EIEC* Real Time PCR Detection Kit with Ref. VS-SHY136 and VS-SHY172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Shigella/EIEC* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortexer.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments).

VIASURE *Shigella/EIEC* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne™ Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTLite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.

- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- The product is intended for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-SHY113L, VS-SHY113H, VS-SHY136 and VS-SHY172). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-SHY136 and VS-SHY172 (compatible with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

8.1. Sample preparation

Stool samples should be collected in clean containers and processed as soon as possible to guarantee the quality of the test. We recommend to use fresh samples.

For longer storage, the samples must be frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation. Freezing and thawing cycles are not recommended.

Perform the sample preparation according to the recommendations appearing in the instructions for use of extraction kit used.

8.2. DNA extraction

For DNA extraction from human stool samples you can use your manually or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions for use. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- QIAamp DNA Mini kit (QIAGEN).
- QIAamp DNA Stool Mini Kit (QIAGEN).
- Invisorb® Spin Universal Kit (Strattec).
- UltraClean® Tissue & Cells DNA Isolation Kit (Mobio).
- RIDA® Xtract (R-biopharm).
- Maxwell®RSC Blood DNA Kit, using the Maxwell® 16 instrument (Promega).

8.3. Lyophilized positive control

Shigella/EIEC Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Shigella/EIEC* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

8.4. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA extracted from each sample, reconstituted *Shigella/EIEC* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® kits).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

| Cycles | Step | Time | Temperature |
|--------|--|--------|-------------|
| 1 | Polymerase activation | 2 min | 95°C |
| 45 | Denaturation | 10 seg | 95°C |
| | Annealing/Extension (Data collection*) | 50 seg | 60°C |

Table 4. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*Shigella*/EIEC) and HEX, JOE or VIC channels (IC). Depending on the equipment used select the proper detection channel (see Annex 2). In the Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne Plus™ Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that the passive reference option for ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *Shigella*/EIEC in the positive control well. Check Internal Control signal to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer’s instructions.

Using the following table read and analyze the results:

| <i>Shigella</i> /EIEC (FAM) | Internal control (HEX) | Negative Control | Positive Control | Interpretation |
|-----------------------------|------------------------|------------------|------------------|--------------------------------------|
| + | +/- | - | + | <i>Shigella</i>/EIEC Positive |
| - | + | - | + | <i>Shigella</i>/EIEC Negative |
| - | - | - | + | Experiment fail |
| + | + | + | - | Experiment fail |

Table 5. Sample interpretation

+: Amplification curve


-: No amplification curve

A sample is considered positive if the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples, it has been validated only with DNA extracted from human faecal samples.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.

- There is a possibility of false positive results due to cross-contamination by *Shigella*/EIEC either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

11. Quality control

VIASURE *Shigella*/EIEC Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

12. Performance characteristics


12.1. Clinical sensitivity and specificity

A total amount of 307 faecal samples from symptomatic patients were tested using VIASURE *Shigella*/EIEC Real Time PCR Detection Kit. *Shigella*/EIEC was detected in 1 sample. This result was confirmed with that obtained by a commercial Real Time PCR Kit RIDA®GENE EHEC/EPEC (R-Biopharm), which detects and differentiates EHEC, STEC (EHEC), EPEC, *Shigella dysenteriae* type 1 and *Shigella*/EIEC.

12.2. Analytical sensitivity

VIASURE *Shigella*/EIEC Real Time PCR Detection Kit has a detection limit of ≥ 10 DNA copies per reaction. (Figure 2).

Figure 2. Dilution series of *Shigella*/EIEC (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel FAM).

12.3. Analytical specificity

The specificity of the *Shigella*/EIEC assay was confirmed by testing a panel consisting of different microorganisms representing the most common enteric pathogens or flora present in the intestine. No cross-reactivity was detected between almost any of the following microorganisms tested.

| Cross-reactivity testing | | | | | |
|--|---|---|---|------------------------------------|---|
| <i>Helicobacter pylori</i> | - | <i>Campylobacter coli</i> | - | <i>Candida albicans</i> | - |
| <i>Helicobacter hepaticus</i> | - | <i>Campylobacter jejuni</i> subsp. <i>jejuni</i> | - | <i>Arcobacter butzleri</i> | - |
| <i>Helicobacter cinaedi</i> | - | <i>Campylobacter upsaliensis</i> | - | <i>Pseudomonas aeruginosa</i> | - |
| <i>Helicobacter heilmannii</i> | - | <i>Proteus vulgaris</i> | - | <i>Enterococcus faecalis</i> | - |
| <i>Salmonella typhi</i> | - | <i>Aeromonas hydrophila</i> subsp. <i>hydrophila</i> | - | <i>Yersinia enterocolitica</i> O:3 | - |
| <i>Salmonella paratyphi A</i> | - | <i>Citrobacter freundii</i> | - | <i>Yersinia enterocolitica</i> O:9 | - |
| <i>Salmonella paratyphi B</i> | - | <i>Staphylococcus aureus</i> subsp. <i>aureus</i> | - | <i>Bacteroides fragilis</i> | - |
| <i>Salmonella typhimurium</i> | - | <i>Serratia liquefaciens</i> | - | <i>Cryptosporidium parvum</i> | - |
| <i>Salmonella bongori</i> | - | <i>Vibrio parahaemolyticus</i> | - | <i>Giardia intestinalis</i> | - |
| <i>Salmonella enteritidis</i> | - | <i>Clostridium difficile</i> | - | <i>Entamoeba histolytica</i> | - |
| <i>Salmonella enterica</i> subsp. <i>entérica</i> | - | <i>Clostridium perfringens</i> | - | Adenovirus serotypes 40 | - |
| <i>Salmonella pullorum</i> | - | Enterotoxigenic <i>Escherichia coli</i> | - | Adenovirus serotypes 41 | - |
| <i>Salmonella gallinarum</i> | - | Enteropathogenic <i>Escherichia coli</i> | - | Rotavirus A | - |
| <i>Campylobacter lari</i> | - | <i>Klebsiella oxytoca</i> | - | Norovirus Genotypes I and II | - |
| <i>Campylobacter fetus</i> | - | <i>Listeria monocytogenes</i> | - | Astrovirus Genotype I-VIII | - |

Table 6. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE *Shigella*/EIEC Real Time PCR Detection Kit was evaluated against *Shigella dysenteriae*, *Shigella flexneri*, *Shigella boydii* and *Shigella sonnei* showing positive results.

ANNEX 1

COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

| Table A.1 LOW PROFILE BLOCK THERMOCYCLERS | |
|---|---|
| Manufacturer | Model |
| Agilent Technologies | AriaMx/AriaDx Real-Time PCR System |
| Applied Biosystems | 7500 Fast Real-Time PCR System ⁽¹⁾ |
| Applied Biosystems | 7500 Fast Dx Real-Time PCR System ⁽¹⁾ |
| Applied Biosystems | QuantStudio™ 12K Flex 96-well Fast |
| Applied Biosystems | QuantStudio™ 6 Flex 96-well Fast |
| Applied Biosystems | QuantStudio™ 7 Flex 96-well Fast |
| Applied Biosystems | QuantStudio™ 3 Fast Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | QuantStudio™ 5 Fast Real-Time PCR System |
| Applied Biosystems | StepOne Plus™ Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | StepOne™ Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | ViiA™ 7 Fast Real-Time PCR System |
| BIONEER | Exicycler™ 96 |
| Bio-Rad | CFX96™ / CFX96™ IVD Real-Time PCR Detection System |
| Bio-Rad | Mini Opticon™ Real-Time PCR Detection System ⁽⁶⁾ |
| Cepheid | SmartCycler® ⁽³⁾ |
| Qiagen | Rotor-Gene® Q ⁽³⁾ |
| Roche | LightCycler® 480 Real-Time PCR System ⁽⁴⁾ |
| Roche | LightCycler® 96 Real-Time PCR System ⁽⁴⁾ |
| Roche | Cobas z480 Analyzer ⁽⁴⁾ |

| Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS | |
|--|--|
| Manufacturer | Model |
| Abbott | Abbott m2000 RealTime System |
| Applied Biosystems | 7300 Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | 7500 Real-Time PCR System |
| Applied Biosystems | 7900 HT Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | ABI PRISM 7000 ⁽⁶⁾ |
| Applied Biosystems | ABI PRISM 7700 ⁽⁵⁾ |
| Applied Biosystems | QuantStudio™ 12K Flex 96-well |
| Applied Biosystems | QuantStudio™ 6 Flex 96-well |
| Applied Biosystems | QuantStudio™ 7 Flex 96-well |
| Applied Biosystems | QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | QuantStudio™ 5 Real-Time PCR System |
| Applied Biosystems | ViiA™ 7 Real-Time PCR System |
| Analytik Jena Biometra | Optical |
| Analytik Jena Biometra | qTOWER 2.0 |
| BIONEER | Exicycler™ 96 |
| Bio-Rad | CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System |
| Bio-Rad | iCycler iQ™ Real-Time PCR Detection System |
| Bio-Rad | iCycler iQ™5 Real-Time PCR Detection System |
| Bio-Rad | MyiQ™ Real-Time PCR Detection System ⁽⁶⁾ |
| Bio-Rad | MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾ |
| Cepheid | SmartCycler® ⁽³⁾ |
| DNA-Technology | DTprime Real-time Detection Thermal Cycler ⁽²⁾ |
| DNA-Technology | DTlite Real-Time PCR System ⁽²⁾ |
| Eppendorf | Mastercycler™ep realplex |
| Qiagen | Rotor-Gene® Q ⁽³⁾ |
| Stratagene / Agilent Technologies | Mx3000P™ Real Time PCR System |
| Stratagene / Agilent Technologies | Mx3005P™ Real Time PCR System |
| VIASURE | VIASURE 48 Real Time PCR System ⁽²⁾ |
| VIASURE | VIASURE 96 Real Time PCR System ⁽²⁾ |

(1) Select Ramp Speed "Standard".
 (2) See Annex 3 to check optical measurement exposure setting.
 (3) The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Rotor-Gene® Q or SmartCycler® tubes.
 (4) Shell Frame grid plate which fits in these Roche qPCR System is necessary.
 (5) No detection in Cy5 channel.
 (6) Detection in FAM and HEX channels only.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

ANNEX 2

DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

| REAL-TIME PCR THERMOCYCLER | VIASURE CHANNEL | DETECTION CHANNEL | OBSERVATIONS |
|-------------------------------|-----------------|-------------------|---|
| Bio-Rad CFX96™ | FAM | FAM | Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it. |
| | HEX | HEX | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| ABI 7500 Applied Biosystems | FAM | FAM | Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected. |
| | HEX | VIC | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| Roche Lightcycler®480II | FAM | 465/510 | Colour Compensation is required |
| | HEX | 533/580 | |
| | ROX | 533/610 | |
| | Cy5 | 618/660 | |
| Smartcycler® Cepheid | FAM | Channel 1 | |
| | HEX | Channel 2 | |
| | ROX | Channel 3 | |
| | Cy5 | Channel 4 | |
| Abbott m2000rt | FAM | FAM | |
| | HEX | VIC | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| Mx3000P™ Mx 3005P™ Stratagene | FAM | FAM | Passive reference option for ROX must be none |
| | HEX | VIC | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| AriaMx Agilent | FAM | FAM | |
| | HEX | HEX | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| Rotor-Gene®Q Qiagen | FAM | Green | In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition". |
| | HEX | Yellow | |
| | ROX | Orange | |
| | Cy5 | Red | |
| Mic Real Time PCR Cycler bms | FAM | Green | In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10) |
| | HEX | Yellow | |
| | ROX | Orange | |
| | Cy5 | Red | |
| Exicycler™ 96 BIONEER | FAM | FAM | |
| | HEX | JOE | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |

Table A3: Detection fluorescence channels of different Real Time PCR systems.

ANNEX 3

OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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

13. Bibliography/Bibliografía

1. A. Ub-Din *et al.* Relationship among *Shigella* spp. and enteroinvasive *Escherichia coli* (EIEC) and their differentiation. *Brazilian Journal of Microbiology* 2014; 45(4): 1131-1138.
 2. V. D. Thiem *et al.* Detection of *Shigella* by a PCR assay targeting the *ipaH* gene suggests increased prevalence of shigellosis in Nha Trang, Vietnam. *Journal of Clinical Microbiology* 2004; 42(5): 2031-2035.
 3. W. S. Lin *et al.* A quantitative PCR assay for rapid detection of *Shigella* species in fresh produce. *Journal of Food Protection* 2010; 73(2): 221-233.
 4. W. Mokhtari *et al.* Real-Time PCR using SYBR Green for the detection of *Shigella* spp. in food and stool samples. *Molecular and Cellular Probes* 2013, 27; 53-59.
- S. Ghosh *et al.* Genetic characterization of *Shigella* spp. isolated from diarrhoeal and asymptomatic children. *Journal of Medical Microbiology* 2014, 63; 903-910.

14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico in vitro

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| Tabla A.2 TERMOCICLADORES CON BLOQUE DE PERFIL ALTO | |
|---|---|
| Fabricante | Modelo |
| Abbott | Abbott m2000 RealTime System |
| Applied Biosystems | 7300 Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | 7500 Real-Time PCR System |
| Applied Biosystems | 7900 HT Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | ABI PRISM 7000 ⁽⁶⁾ |
| Applied Biosystems | ABI PRISM 7700 ⁽⁵⁾ |
| Applied Biosystems | QuantStudio™ 12K Flex 96-well |
| Applied Biosystems | QuantStudio™ 6 Flex 96-well |
| Applied Biosystems | QuantStudio™ 7 Flex 96-well |
| Applied Biosystems | QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | QuantStudio™ 5 Real-Time PCR System |
| Applied Biosystems | ViiA™ 7 Real-Time PCR System |
| Analytik Jena Biometra | Optical |
| Analytik Jena Biometra | qTOWER 2.0 |
| BIONEER | Exicycler™ 96 |
| Bio-Rad | CFX96™ / CFX96™ IVD Real-Time PCR Detection System |
| Bio-Rad | iCycler iQ™ Real-Time PCR Detection System |
| Bio-Rad | iCycler iQ™5 Real-Time PCR Detection System |
| Bio-Rad | MyiQ™ Real-Time PCR Detection System ⁽⁶⁾ |
| Bio-Rad | MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾ |
| Cepheid | SmartCycler® ⁽³⁾ |
| DNA-Technology | DTprime Real-time Detection Thermal Cycler ⁽²⁾ |
| DNA-Technology | DTlite Real-Time PCR System ⁽²⁾ |
| Eppendorf | Mastercycler™ep realplex |
| Qiagen | Rotor-Gene® Q ⁽³⁾ |
| Stratagene / Agilent Technologies | Mx3000P™ Real Time PCR System |
| Stratagene / Agilent Technologies | Mx3005P™ Real Time PCR System |
| VIASURE | VIASURE 48 Real Time PCR System ⁽²⁾ |
| VIASURE | VIASURE 96 Real Time PCR System ⁽²⁾ |

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™, StepOnePlus™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: May 2019

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
www.certest.es

VIASURE online

F-362 rev01

VIASURE

CE

IVD

Real Time PCR Detection Kits

CerTest
BIOTEC