CERTEST

H. pylori+

Transferrin

ONE STEP *H. pylori*+Transferrin
COMBO CARD TEST

CERTEST BIOTEC S.L.

ENGLISH**INTENDED USE**

CerTest *H. pylori*+Transferrin one step combo card test is a colored chromatographic immunoassay for the simultaneous qualitative detection of *H. pylori* and transferrin in stool samples.

CerTest *H. pylori*+Transferrin combo card test offers a simple, a highly sensitive and non-invasive screening assay to determine, on the one hand, the presence of *H. pylori* that causes gastric and duodenal ulcers and on the other, the presence of gastrointestinal bleeding and moreover, to monitor patients treatment (*H. pylori* eradication).

INTRODUCTION

Helicobacter pylori (*H. pylori*) is a spiral-shaped bacterium that is found in the gastric mucous layer or adherent to the epithelial lining of the stomach. *H. pylori* causes more than 90% of duodenal ulcers and up to 80% of gastric ulcers. The importance of *Helicobacter pylori* testing has increased due to the strong correlation between the presence of bacteria and confirmed gastrointestinal diseases (stomach and duodenum) like gastritis, peptic ulcer disease and gastric carcinoma. Patients suffering from gastro-duodenal ulcer may develop complications such as upper gastrointestinal bleeding or perforation. A precise diagnosis of *H. pylori* after eradication reduces significantly the risk of ulcer reactivation, and thus prevents new episodes of upper gastrointestinal bleeding. The invasive and non-invasive methods are used for diagnosis of infection with *H. pylori* in patients with symptoms of gastrointestinal disease.

Transferrin is a blood-derived component. In gastrointestinal bleeding diseases, it may be leaked into gastrointestinal tract and then discharged with the faeces. Transferrin is stable in faeces and a good marker to detect gastrointestinal bleeding.

This immunochromatographic assay detects simultaneously *H. pylori* and human transferrin in stool samples, obtaining more accurate testing results regarding the *H. pylori* infection.

TEST PRINCIPLE

CerTest *H. pylori*+Transferrin is based on the principle of a qualitative immunochromatographic assay for the determination of *H. pylori* and human transferrin in stool samples.

Strip A consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *H. pylori* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*H. pylori*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

Strip B consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against human transferrin and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-human transferrin) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is *H. pylori* positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*H. pylori* monoclonal antibodies-red polystyrene microspheres) in the strip A, and if the sample is human transferrin positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-human transferrin monoclonal antibodies-red polystyrene microspheres) in the strip B, which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*H. pylori* antibodies present on the membrane of strip A (test line) and the anti-human transferrin antibodies present on the membrane of strip B (test line) capture the coloured conjugate and the red line will be visible in both strips. These bands are used to interpret the result.

If the sample is negative, there is no *H. pylori* and human transferrin presence and yet, the *H. pylori* and human transferrin may be present in a concentration lower than the detection limit and the cut-off value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-*H. pylori* and the anti-human transferrin antibodies present on the membranes (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in both strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on both membranes will capture control green-conjugate complex and both control lines will always appear. The presence of these green lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the CerTest *H. pylori*+Transferrin test. Do not use any other commercial kit component.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

SPECIMEN COLLECTION AND PREPARATION


Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8 °C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

Please take note that there are two different stool collection tubes.

- a. “*H. pylori*” stool collection tubes (without reducer piece). 

Specimen preparation (see illustration):


1. Take out the cap of the collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2) making sure that, at each insertion, only the stick's screw is covered with the sample to take the appropriate amount of faecal sample (approx. 50mg) and add it to the collection tube. Any additional sample that exceeds the stick's screw could cause wrong results. For liquid samples, add approx. 125µL in the collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

- b. “Transferrin” stool collection tube (with a reducer piece). 

Specimen preparation (see illustration):

1. Take out the cap of the stool collection tube and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (1), to collect faecal sample and add it to the stool collection tube. For liquid samples, add approx. 15µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (2). The stool collection tube with dilute sample can be stored until 2 days at room temperature or until 5 days in the refrigerator (2-8°C) prior to testing.


Precautions: Patients suffering from menstrual period, bleeding hemorrhoids, blood in urine or strain during bowel movement should not collect samples.

MATERIALS

MATERIALS PROVIDED

- CerTest *H. pylori*+Transferrin combo card tests
- Instructions for use
- "*H. pylori*" stool collection tubes with diluent
- "Transferrin" stool collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE


Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay. Please take note that there are two procedures by using different stool collection tubes:

1. Proceed to shake each stool collection tube in order to assure good sample dispersion.
2. Remove the CerTest *H. pylori*+Transferrin combo card test from its sealed bag just before using it.
3. Take the "*H. pylori*" stool collection tube, cut the end of the cap (4) and dispense 3 drops in the circular window marked with the letter A (5). Avoid adding solid particles with the liquid.
4. Take the "Transferrin" stool collection tube, cut the end of the cap (6) and dispense 3 drops in the circular window marked with the letter B (7). Avoid adding solid particles with the liquid.
5. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.


If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

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



INTERPRETATION OF RESULTS (please refer to the illustration below)

1. 


A: Green→Negative *H. pylori*
B: Green→Negative Transferrin
2. 

A: Green-Red→Positive *H. pylori*
B: Green→Negative Transferrin
3. 

A: Green-Red→Positive *H. pylori*
B: Green-Red→Positive Transferrin
4. 

A: Green→Negative *H. pylori*
B: Green-Red→Positive Transferrin

INVALID RESULTS

5. 

	A (<i>H. pylori</i>)	B (Transferrin)	Interpretation of the results
1.	- GREEN	- GREEN	No infection caused by <i>H. pylori</i> and no gastrointestinal bleeding problem.
2.	+ GREEN-RED	- GREEN	Infection caused by <i>H. pylori</i> , without gastrointestinal bleeding.
3.	+ GREEN-RED	+ GREEN-RED	Infection caused by <i>H. pylori</i> , with gastrointestinal bleeding, which might mean gastrointestinal diseases (stomach and duodenum) like gastritis, peptic ulcer disease and gastric carcinoma.
4.	- GREEN	+ GREEN-RED	No infection caused by <i>H. pylori</i> but gastrointestinal bleeding problem.
5.	Any other result	Any other result	Invalid result either A or B, we recommend repeating the assay using the same sample with another test.

INVALID: Total absence of any control coloured line (GREEN) regardless the appearance or not of the test lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of *H. pylori* and human transferrin present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
3. The intensity of test line may vary from very strong at high *H. pylori*/human transferrin concentration to faint when the *H. pylori*/human transferrin concentration is close to the detection limit and cut-off value of the test.

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- CerTest *H. pylori*+Transferrin should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
- Positive results determine the presence of *H. pylori* and/or human transferrin in faecal samples; nevertheless, a positive result should be followed up with additional invasive techniques (endoscopy) to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations. Transferrin positive result can be also due to several causes besides colorectal bleeding, such as hemorrhoids, blood in urine or stomach irritations.
- A negative result is not meaningful because of it is possible the *H. pylori* and/or human transferrin concentration in the stool sample is lower than the detection limit or cut-off value. If the symptoms or situation still persist, *H. pylori* determination should be carried out on a sample from an enrichment culture or other any technique. Transferrin negative results do not exclude bleeding, as some polyps and colorectal cancers may bleed intermittently or not during certain stages of the disease. Moreover, blood may not be uniformly distributed in stool samples.
- Patients suffering from menstrual period, bleeding hemorrhoids, blood in urine or strain during bowel movement should not collect samples to detect transferrin (false positive result).
- Mucous stool samples may cause non-specific reactions in the test- Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
- Faecal samples whose antigen concentration is near to the cut-off value of the CerTest *H. pylori*+Transferrin test (only for Transferrin strip) could be lost if the sample has been diluted and preserved until 2 days at room temperature or until 5 days at 2-8°C until to be tested. For faecal samples near to the cut-off value is better to dilute the sample and make immediately the test.

EXPECTED VALUES

Studies have found that more than 90% of patients with duodenal ulcer and 80% of patients with gastric ulcer are infected with *Helicobacter pylori*. CerTest *H. pylori* has been compared with different methods: cultures, Urea Breath Test and Urease Test, demonstrating an overall accuracy of >92%.

Peptic ulcer is mostly the common cause of acute intestinal hemorrhage. *H. pylori* is a pathogenetic agent in chronic active gastritis and is strongly associated with ulcer disease.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

Detection limit range of CerTest *H. pylori*+Transferrin (*H. pylori* strip) is: 0.78ng/mL – 0.09ng/mL of *H. pylori* recombinant outer membrane protein.

Cut-off value of CerTest *H. pylori*+Transferrin (Transferrin strip) is 4ng/mL (0.4µg hTf/g of faeces) for human transferrin.

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest *H. pylori*+Transferrin, CerTest), and qPCR technique (VIASURE *Helicobacter pylori* Real Time Detection kit, CerTest) and another commercial available immunochromatographic test (Human Hexagon, OBTI). The samples were taken from patients with the same as *Helicobacter pylori* infection symptoms arrived in Gastroenterology area in a Spanish hospital. The results were as follows:

		qPCR test: Viasure <i>Helicobacter pylori</i> Real Time Detection kit		
		+	-	Total
IC test: <i>H. pylori</i> +Transferrin (<i>H. pylori</i>)	+	54	1	55
	-	1	60	61
Total		55	61	116

CerTest <i>H. pylori</i> +Transferrin (<i>H. pylori</i>) vs Viasure <i>Helicobacter pylori</i> Real Time Detection Kit		
	Mean Value	95% CI (confidence interval)
Sensitivity	98.2%	90.3 – 100.0%
Specificity	98.4%	91.2 – 100.0%
PPV	98.2%	90.3 - 100.0%
NPV	98.4%	91.2 –100.0%

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CerTest <i>H. pylori</i> +Transferrin (Transferrin) vs Human Hexagon		
	Mean Value	95% CI (confidence interval)
Sensitivity	100.0%	69.2 – 100.0%
Specificity	100.0%	76.8 – 100.0%
PPV	100.0%	69.2 - 100.0%
NPV	100.0%	76.8 –100.0%

The results showed a high sensitivity and specificity to detect *H. pylori* and human transferrin using CerTest *H. pylori*+Transferrin.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *H. pylori*+Transferrin; no cross reactivity against gastrointestinal pathogens and other faecal marker transferrin occasionally present in faeces:

Calprotectin	Pig Haemoglobin	<i>Shigella boydii</i>
<i>Campylobacter coli</i>	Lactoferrin	<i>Shigella dysenteriae</i>
<i>Campylobacter jejuni</i>	<i>Listeria monocytogenes</i>	<i>Shigella flexneri</i>
<i>Clostridium difficile</i>	<i>Salmonella enteritidis</i>	<i>Shigella sonnei</i>
<i>Escherichia coli</i>	<i>Salmonella paratyphi</i>	<i>Staphylococcus aureus</i>
Human Haemoglobin	<i>Salmonella typhi</i>	<i>Yersinia enterocolitica</i>
Bovine Haemoglobin	<i>Salmonella typhimurium</i>	

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