# VIASURE

## Real Time PCR Detection Kits

by CerTest  
BIOTEC

### Gastrointestinal Panel IV

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VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit 12 x 8-well strips, low profile

VS-GP0412L

VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit 12 x 8-well strips, high profile

VS-GP0412H

VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit 96-well plate, low profile

VS-GP0413L

VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit 96-well plate, high profile

VS-GP0413H

## ENGLISH

### 1. Intended use

VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit is designed for the specific identification and differentiation of *Clostridium difficile* Toxin A and B; *Salmonella*, *Campylobacter* and/or *Yersinia enterocolitica*; *Campylobacter coli*, *Campylobacter lari* and/or *Campylobacter jejuni*; and/or Enterohemorrhagic *Escherichia coli* (EHEC), Shiga toxin-producing *Escherichia coli* (STEC), Enteropathogenic *Escherichia coli* (EPEC) and/or Enteroinvasive *Escherichia coli* (EIEC)/*Shigella* in human stool samples from patients with signs and symptoms of gastrointestinal infection. This test is intended for use as an aid in the diagnosis of the *Clostridium difficile* Toxin A and B; *Salmonella*, *Campylobacter* and/or *Y. enterocolitica*; *C. coli*, *C. lari* and/or *C. jejuni*; and/or STEC, EPEC and/or EIEC/*Shigella* in combination with clinical and epidemiological risk factors. DNA is extracted from stool specimens, multiplied using Real Time amplification and detected using specific primers and a fluorescent reporter dye probe for *Clostridium difficile* Toxin A and B; *Salmonella*, *Campylobacter* and *Yersinia enterocolitica*; *C. coli*, *C. lari* and *C. jejuni*; and STEC, EPEC and EIEC/*Shigella*.

### 2. Summary and Explanation

Virus, bacteria and/or parasites can cause gastroenteritis, an inflammation of the gastrointestinal tract involving both the stomach and the small intestine. Gastrointestinal infections are self-limited and resolve within a few days. However, in a healthcare setting and in specific populations (newborns/infants, immunocompromized patients or elderly populations), they are potentially serious. Rapid diagnosis, appropriate treatment and infection control measures are therefore particularly important in these contexts.

*Clostridium difficile* is a gram-positive, sporogenic, anaerobic bacillus that belongs to the Clostridiaceae family. It has been identified as a causal agent of antibiotic-associated diarrhea (AAD) and its infection is associated with high morbidity and mortality in the elderly. The major risk factors for *C. difficile* infection are broad spectrum antibiotics exposure, hospitalization, and advanced age. The severity of its infection ranges from mild diarrhea and pseudomembranous colitis to toxic megacolon, perforations of the colon and occasionally, sepsis and even death. The main routes of transmission are the fecal-oral or aerosols. In fact, infected persons with acute diarrhea can lead to heavy contamination of the environment with spores, which can persist in dust or on surfaces for months and be transmitted to other hospitalized patients or to healthcare workers once again. Toxigenic strains of *C. difficile* can colonize the gut, replicate and produce enterotoxin A and cytotoxin B, encoded by *tcdA* and *tcdB* genes which are included in a 19,6 kb pathogenicity locus. Toxins A and B share significant sequence homology and have similar domains, except the carboxyl terminal that differs significantly between the two toxins and it is the receptor binding portion.

The *Salmonella* genus is divided taxonomically into six *Salmonella enterica* subspecies and *Salmonella bongori*. *Salmonella* causes two types of diseases: enteric (typhoid) fever and acute gastroenteritis commonly referred to as salmonellosis. Transmission of *Salmonella* occurs through contaminated food (poultry, poultry products, beef, pork, eggs, milk, and seafood), water or contact with infected animals. Patients infected with *Salmonella* frequently suffer nausea, vomiting, abdominal cramps/pains, achiness, fever, intense diarrhea and/or headache.

Campylobacter species are gram-negative, nonspore forming, spiral, or curved-shaped bacteria. Among the more than 26 species currently classified in the Campylobacter genus, most human diseases are attributed to three major food-borne species: *Campylobacter jejuni*, *Campylobacter coli*, and *Campylobacter lari*. These three species can be isolated from poultry and are of greatest concern to the poultry industry.

*Campylobacter* is considered one of the most common causes of diarrheal illness worldwide. Poultry is a major reservoir and source of transmission of *Campylobacter* to humans. In particular *C. jejuni* is the common species found in beef, *C. coli* is often isolated from pork, and *C. lari* is predominant in shorebirds. Other risk factors include consumption of animal products and water, contact with animals, and even person-to-person transmission (fecal-oral or via fomites).

Infection with *Campylobacter* causes gastroenteritis characterised by fever, vomiting, headaches, and abdominal pain with watery or bloody diarrhea, for a median duration of 6 days. Besides gastroenteritis, these three species can cause periodontitis, septicemia, and second trimester intrauterine growth restriction. Furthermore, *C. jejuni* infection may lead to autoimmune conditions such as Guillain-Barré syndrome (GBS) and Miller Fisher syndrome (MFS).

The genus *Yersinia* has three well known human and animal pathogens: *Yersinia enterocolitica*, *pestis* and *pseudotuberculosis*. There are six biotypes of *Y. enterocolitica*; five of which are considered pathogenic in humans (biotypes 1B, 2, 3, 4 and 5). In addition, there are 60-70 serotypes, among which O:3, O:9, O:8, O:5,27 are mainly associated with human disease. *Yersinia enterocolitica* is a foodborne pathogen and its clinical manifestations typically include nausea, vomiting, abdominal pain, diarrhea and fever. There is a strong evidence that the food of animal origin especially pork and dairy products are responsible for human infections.

*Escherichia coli* is (*E. coli*) a gram-negative microorganism that can be an innocuous resident of the gastrointestinal tract, but it also has the pathogenic capacity to cause enteric disease, and extraintestinal diseases, as urinary tract infections (UTIs) and sepsis/meningitis. Pathogenic variants of *E. coli* (pathovars or pathotypes) cause much morbidity and mortality worldwide, due to they have low infectious doses and are transmitted through ubiquitous mediums, including food and water. Of the strains that cause diarrhoeal diseases, six pathotypes are now recognized: Enterohaemorrhagic *E. coli* (EHEC), Enterotoxigenic *E. coli* (ETEC), Enteroinvasive *E. coli* (EIEC), Enteropathogenic *E. coli* (EPEC), Enteroaggregative *E. coli* (EAEC), and Diffusely adherent *E. coli* (DAEC). Furthermore, different *E. coli* strains may belong to more than one pathotype group owing to the expression of different virulence factors.

Enterohaemorrhagic *E. coli* (EHEC) is a subset of Shiga toxin-producing *E. coli* (STEC), also called verotoxin-producing *E. coli*. STEC are a diverse group of food-borne pathogens which cause a wide spectrum of human diseases, ranging from mild diarrhoea to severe human diseases, including hemorrhagic colitis (HC) and a life-threatening complication hemolytic uremic syndrome (HUS). Their virulence is related in part to their capacity to produce Stx1 and/or Stx2, potent cytotoxins that inhibit host cell protein synthesis. In addition, typical EHEC are often characterised by the production of an outer membrane protein called intimin, which is encoded by the eae gene. This protein mediates both tight attachment of bacteria to enterocytes, as well as lesions (both attaching and effacing (A/E lesions)) in the colon. STEC and EHEC strains can be transmitted to humans through person-to-person contact; consumption of raw or undercooked meat, raw milk and other dairy products;

ingestion of other food or drinking water contaminated with animal faeces; direct contact with domestic cattle and other ruminants recognised as a major reservoir, and contaminated bathing/ recreational water. The more clinically relevant STEC strains belong to the serotype O157:H7, followed by O26:H11, O103:H2, O111:H8, and O145:H28.

Enteropathogenic *E. coli* (EPEC) also contains eae as EHEC, but without shiga-like toxin. EPEC adhere to small bowel enterocytes, and destroy the normal microvillar architecture, inducing the characteristic attaching and effacing lesion. EPEC is an important cause of potentially fatal infant diarrhoea in developing countries that is often accompanied by fever, vomiting, and dehydration in children under 2 years of age. It is transmitted via the fecal-oral route through contaminated surfaces, weaning fluids, and human carriers.

Enteroinvasive *E. coli* (EIEC) are biochemically, genetically and pathogenically closely related to *Shigella* spp. Both bacteria express the invasion plasmid antigen H (*ipaH*) gene which is related to invasion. There are four *Shigella* species responsible for human disease (*S. dysenteriae*, *S. flexneri*, *S. sonnei* and *S. boydii*), which cause varying degrees of dysentery. This infection is characterized by fever, abdominal cramps and diarrhoea containing blood and mucous. Severe shigellosis complications are often associated with the Shiga toxin-producing serotype *S. dysenteriae* 1 and can range from local intestinal disorders to systemic manifestations. Instead, EIEC might cause an invasive inflammatory colitis, and occasionally dysentery, but in most cases EIEC elicits watery diarrhoea that is indistinguishable from that due to infection by other *E. coli* pathogens. Conventional transmission of EIEC and *Shigella* is mediated via the fecal-oral route mainly through contaminated food or water or direct person-to person spread.

Culture, antigen detection, microscopy has been considered to be the gold standard methods for diagnosis of these pathogens, however there are time-consuming and not very sensitive. Real-time PCR assay is less labor-intensive and has higher sensitivity and specificity, making it an attractive alternative.

### 3. Principle of the procedure

VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit is designed for the diagnosis of the *Clostridium difficile* Toxin A and B; *Salmonella*, *Campylobacter* and/or *Yersinia enterocolitica*; *Shigella/EIEC*; and/or STEC, EPEC and/or EIEC/*Shigella* in human stool samples. After DNA isolation, the identification of *C. difficile* Toxin A and B; *Salmonella*, *Campylobacter* and/or *Yersinia enterocolitica*; *C. coli*, *C. lari* and/or *C. jejuni*; and/or STEC, EPEC and/or EIEC/*Shigella* is performed by the amplification of a conserved region of the specific genes (Table 1, Pathogens detected in each reaction tube and target genes) using specific primers and fluorescent-labelled probes.

Gastrointestinal Panel IV 8-well strip contains in each reaction well the following reaction mixes for the detection of the specific targets (Table 1):

	<b>Code</b>	<b>Controls</b>	<b>Reaction mix placed into each well</b>
			Pathogens and target genes
 1 2 3 4 4 3 2 1 •	<b>CDA</b>	IC	<b>Clostridium difficile toxins A+B</b>
			<i>Clostridium difficile</i> Toxin A ( <i>tcdA</i> gene) and Toxin B ( <i>tcdB</i> gene)
2 3 4 4 3 2 1 •	<b>SCY</b>	IC	<b>Salmonella, Campylobacter &amp; Y. enterocolitica</b>
			<i>Salmonella</i> ( <i>invA</i> gene), <i>Campylobacter</i> (16S rRNA gene) and <i>Yersinia enterocolitica</i> ( <i>ail</i> gene)
2 3 4 4 3 2 1 •	<b>CLJ</b>	IC	<b>Campylobacter coli, C. lari &amp; C. jejuni</b>
			<i>C. coli</i> ( <i>CeuE</i> gene), <i>C. lari</i> ( <i>Gyrasa A</i> gene) and <i>C. jejuni</i> ( <i>hip0</i> gene)
2 3 4 4 3 2 1 •	<b>EEE</b>	IC	<b>E. coli EHEC, EPEC &amp; EIEC</b>
			<i>EHEC</i> , <i>EPEC</i> and/or <i>EIEC/Shigella</i> ( <i>stx1</i> , <i>stx2</i> , <i>eae</i> and/or <i>IpaH</i> genes)

Table 1. Gastrointestinal Panel IV 8-well strips provided in VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit. Reaction mix placed into each well, pathogens detected and target genes. IC: Internal control. Note that the first reaction mix is placed into the ends of the strip.

VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase) in a stabilized format, as well as an internal control to monitor PCR inhibition. Each DNA targets are amplified and detected in specific channels (FAM, HEX, ROX, and/or Cy5) and the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2)(Table 5).

#### 4. Reagents provided

VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit includes the following materials and reagents detailed in Table 2 and Table 3:

Reagent/Material	Description	Colour	Amount
Gastrointestinal Panel IV 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Gastrointestinal Panel IV Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNAse free	RNAse/DNAse free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 2. Reagents and materials provided in VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit with Ref. VS-GP0412L and VS-GP0412H.

Reagent/Material	Description	Color	Amount
Gastrointestinal Panel IV 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
Gastrointestinal Panel IV Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNAse free	RNAse/DNAse free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 X 8-cap strip

Table 3. Reagents and materials provided in VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit with Ref. VS-GP0413L and VS-GP0413H.

## 5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortexer.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.

VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTlite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When

using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

## 6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

## 7. Precautions for users

- The product is intended for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-GP0413L and VS-GP0413H). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- **Note that the first reaction mix is placed into the ends of the strip.**
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.

## 8. Test procedure

### 8.1. Sample preparation

Stool samples should be collected in clean containers and processed as soon as possible to guarantee the quality of the test. We recommend to use fresh samples.

For longer storage, the samples must be frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation. Freezing and thawing cycles are not recommended.

Perform the sample preparation according to the recommendations appearing in the instructions for use of extraction kit used.

#### 8.1.1. DNA extraction

For DNA extraction from human stool samples you can use your manually or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions for use. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- QIAamp DNA Stool Mini Kit (QIAGEN).
- Invisorb® Spin Universal Kit (Stratec).
- RIDA® Xtract (r-Biopharm).
- NucleoMag® Pathogen (Macherey Nagel).
- NucleoSpin® RNA Virus (Macherey Nagel).
- Maxwell® RSC Blood DNA Kit, using the Maxwell® 16 instrument (Promega).
- MagDEA Dx SV kit, using the magLEAD® 6gC instrument (Precision System Science Co.)
- ZP02011 MagPurix Viral/Pathogen Nucleic Acids Extraction Kit A, using the MagPurix 12A instrument (Zinexts Life Science Corp.) (Device intended to isolate viral nucleic acids).
- ZP02012 MagPurix Viral/Pathogen Nucleic Acids Extraction Kit B, using the MagPurix 12A instrument (Zinexts Life Science Corp.) (Device intended to isolate viral nucleic acids).

### 8.2. Lyophilized positive control

*Gastrointestinal Panel IV Positive Control* contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Gastrointestinal Panel IV Positive Control* (red vial) by adding 200 µL of the supplied Water RNase/DNAse free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

### 8.3. PCR protocol

Determine and separate the number of required strips, including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plate or strips. **Note that the first reaction mix is placed into the ends of the strip (Table 1) and each strip includes two Gastrointestinal Panel IV reactions, wells 1 to 4: CDA, SCY, CLJ and EEE reaction mixes.**

- 1) Reconstitute the number of strips you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *Gastrointestinal Panel IV* Positive Control (red vial) or Negative Control (violet vial) in the 8 wells of each strip and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips or 96-well plate.

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 4. PCR protocol

Fluorogenic data should be collected during the extension step (\*) through the FAM, ROX, Cy5 and/or HEX (JOE or VIC) channels following Table 5. Depending on the equipment used select the proper detection channel (see Annex 2).

The reaction mixes places in each well of *Gastrointestinal Panel IV* strip allow the detection of the specific target pathogens in the following channels (Table 5).

Code	Channels			
	FAM	HEX	ROX	Cy5
1 CDA	Toxina B	Internal Control	-	Toxina A
2 SCY	Salmonella	Internal Control	<i>Yersinia enterocolitica</i>	Campylobacter
3 CLJ	<i>Campylobacter jejuni</i>	Internal Control	<i>Campylobacter lari</i>	<i>Campylobacter coli</i>
4 EEE	<i>stx1</i> and <i>stx2</i> genes	Internal Control	<i>IpaH</i> gene	eae gene
4 EEE	<i>stx1</i> and <i>stx2</i> genes	Internal Control	<i>IpaH</i> gene	eae gene
3 CLJ	<i>Campylobacter jejuni</i>	Internal Control	<i>Campylobacter lari</i>	<i>Campylobacter coli</i>
2 SCY	Salmonella	Internal Control	<i>Yersinia enterocolitica</i>	Campylobacter
1 CDA	Toxina B	Internal Control	-	Toxina A

Table 5. Pathogens identified in each detection channel.

## 9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative controls of each well from *Gastrointestinal Panel IV* reaction (4 wells) and the presence of signal for *Gastrointestinal Panel IV* positive control of each well from *Gastrointestinal Panel IV* reaction (4 wells). Check Internal Control signal (CLJ, CDA, SCY and EEE reaction mixes) to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.


A sample is considered positive if the Ct value obtained is less than 40 and the Internal Control shows or not an amplification signal. Sometimes, the detection of internal is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.

A sample is considered negative, if the sample shows no amplification signal in the detection system but the Internal Control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control (if available). Using the following table read and analyze the results:

Reaction Mix	Pathogens	Channels			
		FAM	HEX (IC/EC)	ROX	Cy5
1	CDA	Toxina B	Positive		
		Internal Control		Positive/Negative	
		Toxina A			Positive
2	SCY	Salmonella	Positive		
		Internal Control		Positive/Negative	
		<i>Yersinia enterocolitica</i>			Positive
		Campylobacter			Positive
3	CLJ	Campylobacter <i>jejuni</i>	Positive		
		Internal Control		Positive/Negative	
		<i>Campylobacter lari</i>			Positive
		<i>Campylobacter coli</i>			Positive
4	EEE	EHEC Detected	Positive	Positive/Negative	Positive
		STEC (EHEC) Detected		Positive/Negative	Positive
		EPEC Detected		Positive/Negative	Positive
		EIEC/ <i>Shigella</i> Detected		Positive/Negative	Positive
		<i>Shigella dysenteriae</i> type 1 Detected	Positive	Positive/Negative	Positive

Table 6. Sample interpretation. Positive : Amplification curve. Empty: No amplification curve.


Table 7. Correct run of negative and positive controls run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

The result is considered invalid if there is signal of amplification in negative controls or absence of signal in the positive wells. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells, we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

## 10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated with DNA extracted from human faecal samples.

- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Clostridium difficile* toxins A+B; *Salmonella*, *Campylobacter* and/or *Yersinia enterocolitica*; *C. coli*, *C. lari* and/or *C. jejuni*; and/or STEC, EPEC and/or EIEC/*Shigella* either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

## 11. Quality control

VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

## 12. Performance characteristics

### 12.1. Clinical sensitivity and specificity

***Clostridium difficile* toxins A+B Reaction mix** was tested using human stool samples from symptomatic patients. These results were compared with those obtained with a molecular detection method (RealStar *C. difficile* PCR Kit (Altona)).

The results were as follows:

	RealStar <i>C. difficile</i> PCR Kit (Altona)			
		+	-	Total
<b><i>Clostridium difficile</i> toxins A+B Reaction mix</b>	+	62	19*	81
	-	0	34	34
	Total	62	53	115

Table 8. Comparative results for Toxin A

\* The low amount of template DNA detected in these faecal samples is below the detection limit of the method used.

	RealStar <i>C. difficile</i> PCR Kit (Altona)			
		+	-	Total
<b><i>Clostridium difficile</i> toxins A+B Reaction mix</b>	+	77	5*	82
	-	2*	31	33
	Total	79	36	115

Table 9. Comparative results for Toxin B

\* The low amount of template DNA detected in these faecal samples is below the detection limit of the method used.

The clinical performance of ***Salmonella*, *Campylobacter* & *Y. enterocolitica* Reaction mix** was tested using 400 faecal samples from symptomatic patients. The results were compared with those obtained by commercial Real Time RT-PCR kit ("RIDA®GENE Bacterial Stool Panel" (r-Biopharm) for *Salmonella* and *Campylobacter*). The results were as follows:

<b>Salmonella, Campylobacter &amp; Y. enterocolitica Reaction mix</b>	RIDA®GENE Bacterial Stool Panel (r-Biopharm)			
		+	-	Total
	+	20	1*#	21
	-	0	379	379
	<b>Total</b>	<b>20</b>	<b>380</b>	<b>400</b>
	Table 10. Comparative results for <i>Salmonella</i>			

Table 10. Comparative results for *Salmonella*

\* The low amount of template DNA detected in these faecal samples is below the detection limit of the method used and so it may produce random positive results.

# Besides, this sample was culture positive for *Salmonella* G.

<b>Salmonella, Campylobacter &amp; Y. enterocolitica Reaction mix</b>	RIDA®GENE Bacterial Stool Panel (r-Biopharm)			
		+	-	Total
	+	58	13*#	71
	-	0	329	329
	<b>Total</b>	<b>58</b>	<b>342</b>	<b>400</b>
	Table 11. Comparative results for <i>Campylobacter</i>			

Table 11. Comparative results for *Campylobacter*

\* The low amount of template DNA detected in these faecal samples is below the detection limit of the method used and so it may produce random positive results.

# Besides, these 13 positive samples have been evaluated by an additional commercial Real Time PCR Kit (Mericon-Campylobacter spp, QIAGEN), confirming our results.

<b>Salmonella, Campylobacter &amp; Y. enterocolitica Reaction mix</b>	RIDA®GENE Bacterial Stool Panel (r-Biopharm)			
		+	-	Total
	+	4	0	4
	-	0	396	396
	<b>Total</b>	<b>4</b>	<b>396</b>	<b>400</b>
	Table 12. Comparative results for <i>Yersinia enterocolitica</i>			

Table 12. Comparative results for *Yersinia enterocolitica*

**Campylobacter coli, C. lari & C. jejuni Reaction mix** was tested using 62 faecal samples (for *C. lari* are spiked clinical samples) from symptomatic patients. These results were compared with those obtained by Routine Diagnostic procedures (Microbiological culture and MALDI-TOF).

The results were as follows:

<b>Campylobacter coli, C. lari &amp; C. jejuni Reaction mix</b>	Culture + MALDI-TOF			
		+	-	Total
	+	10	0	10
	-	2	50	52
	<b>Total</b>	<b>12</b>	<b>50</b>	<b>62</b>
	Table 13. Comparative results for <i>C. coli</i>			

Table 13. Comparative results for *C. coli*

<b>Campylobacter coli, C. lari &amp; C. jejuni Reaction mix</b>	Culture + MALDI-TOF			
		+	-	Total
	+	30	0	30
	-	3	29	32
	<b>Total</b>	<b>33</b>	<b>29</b>	<b>62</b>
	Table 14. Comparative results for <i>C. jejuni</i>			

Table 14. Comparative results for *C. jejuni*

<b>Campylobacter coli, C. lari &amp; C. jejuni Reaction mix</b>	Culture + MALDI-TOF			
		+	-	Total
	+	5	0	5
	-	0	57	57
	<b>Total</b>	<b>5</b>	<b>57</b>	<b>62</b>
	Table 15. Comparative results for <i>C. lari</i>			

Table 15. Comparative results for *C. lari*

A small cohort of 25 faecal samples from symptomatic patients and with clinical suspicion of *E. coli* infection were tested using ***E. coli* EHEC, EPEC & EIEC Reaction mix**. The results were confirmed with that obtained by a commercial Real Time PCR Kit (RIDA®GENE EHEC/EPEC real-time multiplex PCR (R-biopharm)). We found 6 positive samples for STEC, 1 positive sample for *Shigella dysenteriae* serotype 1, 2 specimens for EHEC and 4 samples for EPEC using R-biopharm and VIASURE kits.

The results show a high sensitivity and specificity to detect *Clostridium difficile* toxins A+B; *Salmonella*, *Campylobacter* and/or *Yersinia enterocolitica*; *Campylobacter coli*, *C. lari* & *C. jejuni*; and/or STEC, EPEC and/or EIEC/*Shigella* using VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit.

## 12.2. Analytical sensitivity

VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit has a detection limit of  $\geq 10$  DNA copies per reaction for *Clostridium difficile* toxins A+B; *Salmonella*, *Campylobacter* and/or *Yersinia enterocolitica*; *C. coli*, *C. lari* and/or *C. jejuni*; and/or *stx1*, *stx2*, *eae* and/or *IpaH* genes.

## 12.3. Analytical specificity

The specificity of the Gastrointestinal Panel IV assay was confirmed by testing a panel consisting of different microorganisms representing the most common enteric pathogens or flora present in the intestine. No cross-reactivity was detected between almost any of the following microorganisms tested, except the targeted pathogens of each assay.

Cross-reactivity testing					
<i>Helicobacter pylori</i>	-	<i>Campylobacter coli</i>	-/+	<i>Candida albicans</i>	-
<i>Helicobacter hepaticus</i>	-	<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	-/+	<i>Arcobacter butzleri</i>	-
<i>Helicobacter cinaedi</i>	-	<i>Campylobacter upsaliensis</i>	-/+	<i>Pseudomonas aeruginosa</i>	-
<i>Helicobacter heilmannii</i>	-	<i>Proteus vulgaris</i>	-	<i>Enterococcus faecalis</i>	-
<i>Shigella flexneri</i>	-/+	<i>Citrobacter freundii</i>	-	<i>Yersinia enterocolitica</i> O:3	-/+
<i>Shigella dysenteriae</i>	-/+	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-	<i>Yersinia enterocolitica</i> O:9	-/+
<i>Salmonella typhi</i>	-/+	<i>Serratia liquefaciens</i>	-	<i>Bacteroides fragilis</i>	-
<i>Salmonella paratyphi A</i>	-/+	<i>Vibrio parahaemolyticus</i>	-	Adenovirus serotypes 40/41	-
<i>Salmonella paratyphi B</i>	-/+	<i>Clostridium difficile</i>	-/+	Rotavirus A	-
<i>Salmonella typhimurium</i>	-/+	<i>Clostridium perfringens</i>	-	Norovirus Genotypes I and II	-
<i>Salmonella bongori</i>	-/+	Enterotoxigenic <i>Escherichia coli</i>	-	Astrovirus Genotype I-VIII	-
<i>Salmonella enteritidis</i>	-/+	<i>Blastocystis hominis</i>	-	Sapovirus	-
<i>Salmonella enterica</i> subsp. <i>enterica</i>	-/+	Enterohemorrhagic <i>Escherichia coli</i>	-/+	<i>Entamoeba histolytica</i>	-
<i>Salmonella pullorum</i>	-/+	<i>Klebsiella oxytoca</i>	-	<i>Entamoeba dispar</i>	-
<i>Salmonella gallinarum</i>	-/+	<i>Listeria monocytogenes</i>	-	<i>Cryptosporidium parvum/hominis</i>	-
<i>Campylobacter lari</i>	-/+	<i>Aeromonas caviae</i>	-	<i>Giardia intestinalis</i>	-
<i>Campylobacter fetus</i>	-/+	<i>Aeromonas hydrophila</i> subsp. <i>Hydrophila</i>	-	<i>Dientamoba fragilis</i>	-

Table 16. Reference pathogenic microorganisms used in this study.

## 12.4. Analytical reactivity

The reactivity of VIASURE Gastrointestinal Panel IV Real Time PCR Detection Kit for ***Clostridium difficile* Reaction mix**, was evaluated against Toxins A (*Clostridium difficile* CECT 531 and *Clostridium difficile* ribotypes 001, 002, 005, 014, 027, 057, 106 and 078/126) and Toxins B (*Clostridium difficile* CECT 531 and *Clostridium difficile* ribotypes 001, 002, 005, 014, 017, 027, 057, 106, and 078/126), showing positive results.

The reactivity of ***Salmonella, Campylobacter & Y. enterocolitica* Reaction mix**, was evaluated against Salmonella (*Salmonella typhi*, *Salmonella paratyphi A*, *Salmonella paratyphi B*, *Salmonella typhimurium*, *Salmonella bongori*, *Salmonella enteritidis*, *Salmonella enterica* subsp. *enterica*, *Salmonella pullorum*, *Salmonella gallinarum*, *Salmonella mbandaka*, *Salmonella braenderup* and *Salmonella Tennessee*); Campylobacter (*Campylobacter jejuni* subsp. *jejuni*, *Campylobacter coli*, *Campylobacter lari*, *Campylobacter upsaliensis*, *Campylobacter fetus*, *Campylobacter concisus*, *Campylobacter hyoilestinalis*, *Campylobacter gracilis*, *Campylobacter helveticus*, *Campylobacter curvus* and *Campylobacter rectus*) and *Yersinia enterocolitica* (*Yersinia enterocolitica* serotypes O:3, O:8 and O:9), showing positive result.

The reactivity of ***Campylobacter coli, C. lari & C. jejuni* Reaction mix**, was evaluated *Campylobacter coli*, *Campylobacter lari* and *Campylobacter jejuni*, showing positive result.

The reactivity of ***E. coli EHEC, EPEC & EIEC* Reaction mix** was evaluated against *E. coli* strains of each pathotype and *Shigella* species which contain the following virulence factors: *stx1* and/or *stx2* and *eae* genes (EHEC), *stx1* and/or *stx2* genes (STEC) and *eae* gene (EPEC) and *ipaH* (EIEC, *S. dysenteriae*, *S. flexneri*, *S. sonnei* and *S. boydii*), showing positive results.

## ANNEX 1

**COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT**

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast Real-Time PCR System <sup>(1)</sup>
Applied Biosystems	7500 Fast Dx Real-Time PCR System <sup>(1)</sup>
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 3 Fast Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	QuantStudio™ 5 Fast Real-Time PCR System
Applied Biosystems	StepOne Plus™ Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	StepOne™ Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	Mini Opticon™ Real-Time PCR Detection System <sup>(6)</sup>
Cepheid	SmartCycler® <sup>(3)</sup>
Qiagen	Rotor-Gene® Q <sup>(3)</sup>
Roche	LightCycler ®480 Real-Time PCR System <sup>(4)</sup>
Roche	LightCycler ®96 Real-Time PCR System <sup>(4)</sup>
Roche	Cobas z480 Analyzer <sup>(4)</sup>

(1)Select Ramp Speed "**Standard**".

(2)See Annex 3 to check optical measurement exposure setting.

(3)The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Rotor-Gene® Q or SmartCycler® tubes.

(4)Shell Frame grid plate which fits in these Roche qPCR System is necessary.

(5)No detection in Cy5 channel.

(6)Detection in FAM and HEX channels only.

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	ABI PRISM 7000 <sup>(6)</sup>
Applied Biosystems	ABI PRISM 7700 <sup>(5)</sup>
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System <sup>(6)</sup>
Bio-Rad	MyiQ™2 Real-Time PCR Detection System <sup>(6)</sup>
Cepheid	SmartCycler® <sup>(3)</sup>
DNA-Technology	DTprime Real-time Detection Thermal Cycler <sup>(2)</sup>
DNA-Technology	DTlite Real-Time PCR System <sup>(2)</sup>
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q <sup>(3)</sup>
Stratagene / Agilent Technologies	Mx3000PTM Real Time PCR System
Stratagene / Agilent Technologies	Mx3005PTM Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System <sup>(2)</sup>
VIASURE	VIASURE 96 Real Time PCR System <sup>(2)</sup>

Table A1/A2. Compatible low and high profile Real Time PCR systems.

## ANNEX 2

**DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT**

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation is required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 µl) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

## ANNEX 3

### OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -500\*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

\*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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
## 13. Bibliography/Bibliografía

1. S. D. Bélanger et al. Rapid detection of *Clostridium difficile* in feces by Real-Time PCR. *Journal of Clinical Microbiology* 2003; 41(2): 730-734.
2. J. G. Bartlett et al. Clinical recognition and diagnosis of *Clostridium difficile* infection. *Clinical Infectious Diseases* 2008; 46: S12-S18.
3. R. Mutters et al. Quantitative detection of *Clostridium difficile* in hospital environmental samples by Real-Time polymerase chain reaction. *Journal of Hospital Infection* 2009; 71: 43-48.
4. R. A. Luna et al. Rapid stool based diagnosis of *Clostridium difficile* infection by Real-Time PCR in a children's hospital. *Journal of Clinical Microbiology* 2011; 49(3): 851-857.
5. E. de Jong et al. Clinical and laboratory evaluation of a real-time PCR for *Clostridium difficile* toxin A and B genes. *European Journal of Clinical Microbiology and Infectious Diseases* 2012; 31: 2219-2225.
6. J. Hoofar, et al. Automated 5' Nuclease PCR Assay for Identification of *Salmonella enterica*. *Journal of Clinical Microbiology* 2000; 38: 3429-3435.
7. D. Rodríguez-Lázaro et al. A rapid and direct real time PCR-based method for identification of *Salmonella* spp. *Journal of Microbiological Methods* 2003; 54: 381- 390.
8. F. Barletta et al. Multiplex real-time PCR for detection of *Campylobacter*, *Salmonella*, and *Shigella*. *Journal of Clinical Microbiology* 2013; 51: 2822-2829
9. N.O. Kaakoush et al. Global Epidemiology of *Campylobacter* Infection. *Clinical Microbiology Reviews* 2015; 28(3): 687-720.
10. S.M. Man. The clinical importance of emerging *Campylobacter* species. *Nature Reviews Gastroenterology & Hepatology* 2011; 8(12): 669-685.
11. R.F. de Boer et al. Detection of *Campylobacter* species and *Arcobacter butzleri* in stool samples by use of real-time multiplex PCR. *Journal of Clinical Microbiology* 2013; 51(1): 253-259.
12. N. Botteldoorn et al. Quantification of *Campylobacter* spp. in chicken carcass rinse by real-time PCR. *Journal of Applied Microbiology* 2008; 105(6):1909-1918.
13. M. Arrausi-Subiza, et al. Evaluation of different enrichment methods for pathogenic *Yersinia* species detection by real time PCR. *BMC Veterinary Research* 2014; 10(1): 192.
14. A.D. Jourdan et al. Development of a fluorogenic 5' nuclease PCR assay for detection of the *ail* gene of pathogenic *Yersinia enterocolitica*. *Applied and Environmental Microbiology* 2000; 66(9):3750-3755.
15. S.T. Lambertz et al. Real-time PCR method for detection of pathogenic *Yersinia enterocolitica* in food. *Applied and Environmental Microbiology* 2008;74(19):6060-6067.
16. J. Liang et al. Two novel *ail*-positive biotype 1A strains of *Yersinia enterocolitica* isolated in China with unequal adhesion and invasion properties. *Infection, Genetics and Evolution* 2014;27:83-88.
17. A. Ub-Din et al. Relationship among *Shigella* spp. and enteroinvasive *Escherichia coli* (EIEC) and their differentiation. *Brazilian Journal of Microbiology* 2014; 45(4): 1131-1138.
18. N.O. Kaakoush et al. Global Epidemiology of *Campylobacter* Infection. *Clinical Microbiology Reviews* 2015; 28(3): 687-720.

19. K. C. Liu et al. Simultaneous Identification of *Campylobacter jejuni*, *Campylobacter coli*, and *Campylobacter lari* with SmartCycler-Based Multiplex Quantitative Polymerase Chain Reaction. *Foodborne pathogens and disease*, 2017 Volume 14, Number 7.
20. M. E. Patrick et al. Features of illnesses caused by five species of *Campylobacter*, Foodborne Diseases Active Surveillance Network (FoodNet) – 2010–2015. *Epidemiology & Infection* 2017. doi:10.1017/S0950268817002370.
21. C. A. Whitehouse et al. Antimicrobial Resistance in *Campylobacter* Species: Mechanisms and Genomic Epidemiology. *Advances in Applied Microbiology*, 2018, Volume 103.
22. AM. Mayr et al. Rapid detection and differentiation of *Campylobacter jejuni*, *Campylobacter coli*, and *Campylobacter lari* in food, using multiplex real-time PCR. *Journal of Food Protection*, 2010;73(2):241-50.
23. E.M. Nielsen et al. Detection and characterization of verocytotoxin-producing *Escherichia coli* by automated 5' nuclease PCR assay. *Journal of Clinical Microbiology* 2003; 41(7): 2884-2893.
24. J. Antikainen et al. A quantitative polymerase chain reaction assay for rapid detection of 9 pathogens directly from stools of travelers with diarrhea. *Clinical gastroenterology and hepatology: the official clinical practice journal of the American Gastroenterological Association* 2013; 11(10):1300-1307.e3.
25. L. Ercoli, et al. Prevalence and characteristics of verotoxigenic *Escherichia coli* strains isolated from pigs and pork products in Umbria and Marche regions of Italy. *International Journal of Food Microbiology* 2016; 232: 7-14.
26. T.E. Grys et al. Rapid and sensitive detection of Shiga toxin-producing *Escherichia coli* from nonenriched stool specimens by real-time PCR in comparison to enzyme immunoassay and culture. *Journal of Clinical Microbiology* 2009; 47(7): 2008-12.
27. J. Madic et al. Detection of Shiga toxin-producing *Escherichia coli* serotypes O26:H11, O103:H2, O111:H8, O145:H28, and O157:H7 in raw-milk cheeses by using multiplex real-time PCR. *Applied Environmental Microbiology* 2011; 77(6): 2035-41.
28. K. Verstraete et al. A qPCR assay to detect and quantify Shiga toxin-producing *E. coli* (STEC) in cattle and on farms: a potential predictive tool for STEC culture-positive farms. *Toxins (Basel)* 2014; 6(4): 1201-21.
29. Q. Chen et al. Rapid genetic typing of diarrheagenic *Escherichia coli* using a two-tube modified molecular beacon based multiplex real-time PCR assay and its clinical application. *Annals of Clinical Microbiology and Antimicrobials* 2014;13:30.
30. J.B. Kape. Pathogenic *Escherichia coli*. *Nature Reviews Microbiology* 2004; 2(2): 123-40.
31. M.A. Croxen et al. Recent advances in understanding enteric pathogenic *Escherichia coli*. *Clinical Microbiology Reviews* 2013; 26(4): 822-80.
32. E.J. Won et al. Multiplex Real-Time PCR Assay Targeting Eight Parasites Customized to the Korean Population: Potential Use for Detection in Diarrheal Stool Samples from Gastroenteritis Patients. *PLoS One* 2016;11(11): e0166957.
33. L. Chui et al. Comparison of Shiga toxin-producing *Escherichia coli* detection methods using clinical stool samples. *The Journal of Molecular Diagnostics* 2010; 12(4): 469-475.

## 14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico in vitro

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

- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: November 2019

**CerTest Biotec, S.L.**

Pol. Industrial Río Gállego II · Calle J, Nº1  
50840, San Mateo de Gállego, Zaragoza (Spain)  
[www.certest.es](http://www.certest.es)

VIASURE online

F-362 rev01

**VIASURE**

Real Time PCR Detection Kits

**CerTest**  
BIOTEC