VIASURE


Real Time PCR Detection Kits

by CerTest
BIOTEC

Malaria differentiation

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VIASURE <i>Malaria differentiation</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-MAD106L
VIASURE <i>Malaria differentiation</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-MAD106H
VIASURE <i>Malaria differentiation</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-MAD112L
VIASURE <i>Malaria differentiation</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-MAD112H
VIASURE <i>Malaria differentiation</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-MAD136

ENGLISH

1. Intended use

VIASURE *Malaria differentiation* Real Time PCR Detection Kit is designed for the specific detection and differentiation of DNA of *Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium ovale*, *Plasmodium malariae* and/or *Plasmodium knowlesi* in whole blood samples from patients with signs and symptoms of malaria infection. This test is intended for use as an aid in the diagnosis of these five *Plasmodium* species responsible for malaria in combination with clinical and epidemiological risk factors. DNA is extracted from specimens, multiplied using Real Time amplification and detected using specific primers and a fluorescent reporter dye probe for *P. falciparum*, *P. vivax*, *P. ovale*, *P. malariae* and *P. knowlesi*.

2. Summary and Explanation

Malaria is a major threat to public health worldwide in communicable diseases, since its prevalence and mortality exceeds those of any other.

Malaria is a treatable infectious disease transmitted by vectors and caused by the asexual form of eukaryotic parasites of the genus *Plasmodium*. It is an apicomplex of the *Hematozoa* class, family *Plasmodiidae*. The parasites are transmitted to people by the bites of infected mosquitoes of the genus *Anopheles*.

It is known that there are five species that are capable of infecting humans: *Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium malariae*, *Plasmodium ovale* and *Plasmodium knowlesi*, the latter infected only primates, but has been responsible for cases of infection in humans in Malaysia and Borneo. The two most frequent species are *P. vivax* and *P. falciparum*. *P. falciparum* infection is the most common cause of severe malaria, is responsible for the majority of deaths from malaria and can cause other complications such as anemia and in children of pregnant women with malaria, can cause low birth weight. The classic symptomatology of malaria is a triad: fever, chills and sweating, and sometimes with other nonspecific symptoms such as headache and pain in other parts of the body, without associated complications. In a few cases severe malaria is present with neurological complications or life-threatening complications.

The most severe cases (malaria tropica), usually caused by the most aggressive species, *P. falciparum*, usually include high fever, chills, diarrhea, headache and in a few hours can evolve into a severe condition with hepatic, renal, and coagulation disorders, pulmonary and cerebral edema, encephalopathy, coma and death. Even mild cases can quickly evolve into a deadly form, so early diagnosis and treatment are essential. The incubation time takes 7 days or more after the first exposure.

The genetic diversity gives *Plasmodium* the ability to evade the immune response of the host and produce variants resistant to drugs and vaccines, being this largely responsible for the successful survival of the parasite in the evolutionary history, as well as the failure of the measures used with the aim of eradicating it.

Since it was first described in 1880, the diagnosis of this disease has been made by observing the different forms of the parasite in the microscopic examination of peripheral blood smears stained with various dyes. Today, this technique is still the reference method. However, the laboriousness that requires the training of a good

microscopist and the difficulty of observing low parasitemias has led to the development of new, simpler but more sensitive and specific techniques. Diagnosing malaria in time can be vital for the patient, since the appearance of complications is closely related to the delay in the initiation of treatment.

In addition, molecular diagnostic tests can detect amounts of *Plasmodium* when the parasitaemia is below the minimum level of detection in the examination of the blood sample.

3. Principle of the procedure

VIASURE *Malaria differentiation* Real Time PCR Detection Kit is designed for the diagnosis of *P. falciparum*, *P. vivax*, *P. ovale*, *P. malariae* and *P. knowlesi* in whole blood samples. After DNA isolation, the identification of *P. falciparum*, *P. vivax*, *P. ovale*, *P. malariae* and *P. knowlesi* is performed by the amplification of a conserved region of the 18S rRNA gene using specific primers and a fluorescent-labeled probe.

VIASURE *Malaria differentiation* Real Time PCR Detection Kit is based on the 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bounded to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence can be measured on Real Time PCR platforms.

VIASURE *Malaria differentiation* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase) in a stabilized format, as well as an internal control to monitor PCR inhibition. Each kit includes two kind of strips and each one corresponds to one different assay. The first strip contains the multiplex reaction mix for the detection of *P. malariae*, *P. knowlesi* and *P. ovale* (*P. malariae*, *P. knowlesi* & *P. ovale* 4/8-well strips). *P. malariae* is amplified and detected in FAM channel, *P. knowlesi* is amplified and detected in ROX channel, *P. ovale* is amplified and detected in Cy5 channel and the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2). The second strip contains the multiplex reaction mix for the detection of *P. falciparum* and *P. vivax* (*P. falciparum* + *P. vivax* 4/8-well strips). *P. falciparum* is amplified and detected in FAM channel, *P. vivax* is amplified and detected in ROX channel and the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

4. Reagents provided

VIASURE *Malaria differentiation* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables 1 and 2. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table 1 includes materials and reagents to be used with 8-well strips compatible devices (See Annex 1). Table 2 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips.

Reagent/Material	Description	Colour	Amount
<i>P. malariae</i> , <i>P. knowlesi</i> & <i>P. ovale</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	3/6 x 8-well strip
<i>P. falciparum</i> + <i>P. vivax</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	3/6 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Malaria differentiation</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE *Malaria differentiation* Real Time PCR Detection Kit with Ref. VS-MAD106L, VS-MAD106H, VS-MAD112L and VS-MAD112H.

Reagent/Material	Description	Colour	Amount
<i>P. malariae</i> , <i>P. knowlesi</i> & <i>P. ovale</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	Transparent	9x 4-well strip
<i>P. falciparum</i> + <i>P. vivax</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	Transparent	9x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Malaria differentiation</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	18 x 4-cap strip

Table 2. Reagents and materials provided in VIASURE *Malaria differentiation* Real Time PCR Detection Kit with Ref. VS-MAD136. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Malaria differentiation* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips (if available).
- Vortexer.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments).

VIASURE *Malaria differentiation* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTlite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- The product is indented for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-MAD136). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For reference VS-MAD136 (compatible with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.

- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Make sure to use a well for the detection of *P. malariae*, *P. knowlesi* and *P. ovale* and another well for *P. falciparum* and *P. vivax* assay. Be careful not to mix them throughout the process.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

8.1. DNA extraction

Perform the sample preparation following the manual system protocol or according to the recommendations appearing in the instructions for use of extraction kit used.

For DNA extraction from whole blood samples you can use your manual or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- Invisorb® Spin Universal Kit (Strattec).
- MagDEA Dx SV kit, using the magLEAD® 6gC instrument (Precision System Science Co.)
- EZ1 DNA Blood Kit, using the BioRobot EZ1 DSP Workstation (QIAGEN).

8.2. Lyophilized positive control

Malaria differentiation Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Malaria differentiation* Positive Control (red vial) by adding 200 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

8.3. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

2) Adding samples and controls.

Add 5 µL of DNA extracted from each sample, reconstituted *Malaria differentiation* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® kits).

Load the plate or the strips in the thermocycler.

3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Polymerase activation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 3. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*P. malariae* and *P. falciparum*), ROX (*P. knowlesi* and *P. vivax*), Cy5 (*P. ovale*) and HEX, JOE or VIC channels (IC). Depending on the equipment used select the proper detection channel (see Annex 2). In the Applied Biosystems 7500 Fast Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that the passive reference option for ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *Malaria differentiation* in the positive control well. Check Internal Control signal to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

Interpretation of results for *P. malariae*, *P. knowlesi* & *P. ovale* 4/8-well strips:

<i>P. malariae</i> (FAM)	<i>P. knowlesi</i> (ROX)	<i>P. ovale</i> (Cy5)	Internal control (HEX)	Negative Control	Positive Control	Interpretation
+	+	+	+/-	-	+	<i>P. malariae</i> , <i>P. knowlesi</i> and <i>P. ovale</i> Positive
-	-	-	+	-	+	<i>P. malariae</i> , <i>P. knowlesi</i> and <i>P. ovale</i> Negative
+	-	-	+/-	-	+	<i>P. malariae</i> Positive, <i>P. knowlesi</i> and <i>P. ovale</i> Negative
+	+	-	+/-	-	+	<i>P. malariae</i> and <i>P. knowlesi</i> Positive, and <i>P. ovale</i> Negative
+	-	+	+/-	-	+	<i>P. malariae</i> and <i>P. ovale</i> Positive, and <i>P. knowlesi</i> Negative
-	+	-	+/-	-	+	<i>P. knowlesi</i> Positive, <i>P. malariae</i> and <i>P. ovale</i> Negative
-	+	+	+/-	-	+	<i>P. knowlesi</i> and <i>P. ovale</i> Positive, <i>P. malariae</i> Negative
-	-	+	+/-	-	+	<i>P. ovale</i> Positive, <i>P. malariae</i> and <i>P. knowlesi</i> Negative
-	-	-	-	-	+	Experiment fail
+	+	+	+	+	-	Experiment fail

Table 4. Sample interpretation

+: Amplification curve

-: No amplification curve

Interpretation of results for *P. falciparum* + *P. vivax* 4/8-well strips:

<i>P. falciparum</i> (FAM)	<i>P. vivax</i> (ROX)	Internal control (HEX)	Negative Control	Positive Control	Interpretation
+	+	+/-	-	+	<i>P. falciparum</i> and <i>P. vivax</i> Positives
-	-	+	-	+	<i>P. falciparum</i> and <i>P. vivax</i> Negative
+	-	+/-	-	+	<i>P. falciparum</i> Positive and <i>P. vivax</i> Negative
-	+	+/-	-	+	<i>P. vivax</i> Positive and <i>P. falciparum</i> Negative
-	-	-	-	+	Experiment fail
+	+	+	+	-	Experiment fail

Table 5. Sample interpretation

+: Amplification curve

-: No amplification curve

A sample is considered positive if the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System. (*P. malariae*, *P. knowlesi* & *P. ovale*).



Figure 2. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System. (*P. falciparum* + *P. vivax*).

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is

recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated only with DNA extracted from whole blood samples.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Malaria differentiation* Positive Control, either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

11. Quality control

VIASURE *Malaria differentiation* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *Malaria differentiation* Real Time PCR Detection Kit was tested using 226 blood samples from symptomatic patients. These results were compared with those obtained with a molecular detection method (FTD *Malaria differentiation* (Fast Track Diagnostics)).

The results were as follows:

VIASURE <i>Malaria differentiation</i> Real Time PCR Detection Kit	FTD <i>Malaria differentiation</i> (Fast Track Diagnostics)		
		+	-
	+	2	0
	-	0	224
	Total	2	224

Table 6. Comparative results for *P. malariae*.

VIASURE Malaria differentiation Real Time PCR Detection Kit	FTD Malaria differentiation (Fast Track Diagnostics)			
		+	-	Total
	+	2	0	2
	-	0	224	224
	Total	2	224	226

Table 7. Comparative results for *P. ovale*.

VIASURE Malaria differentiation Real Time PCR Detection Kit	FTD Malaria differentiation (Fast Track Diagnostics)			
		+	-	Total
	+	5	0	5
	-	0	221	221
	Total	5	221	226

Table 8. Comparative results for *P. vivax*.

VIASURE Malaria differentiation Real Time PCR Detection Kit	FTD Malaria differentiation (Fast Track Diagnostics)			
		+	-	Total
	+	55	0	55
	-	0	171	171
	Total	55	171	226

Table 9. Comparative results for *P. falciparum*.

Sensitivity, specificity, PPV and NPV values for VIASURE Malaria differentiation Real Time PCR detection kit (CerTest Biotec) compared with FTD Malaria differentiation (Fast Track Diagnostics) are shown in next table:

	<i>P. malariae</i>	<i>P. ovale</i>	<i>P. vivax</i>	<i>P. falciparum</i>
SE (%)	>99	>99	>99	>99
SP (%)	>99	>99	>99	>99
PPV (%)	>99	>99	>99	>99
NPV (%)	>99	>99	>99	>99

Table 10. Sensitivity (SE), Specificity (SP), Positive Predictive Value (PPV) and Negative Predictive Value (NPV) calculated for VIASURE assay and FTD Malaria differentiation (Fast Track Diagnostics) as reference. (CI = Confidence Interval 95%).

The results show a high sensitivity and specificity to detect *P. falciparum*, *P. vivax*, *P. ovale* and *P. malariae* using VIASURE Malaria differentiation Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE *Malaria differentiation* Real Time PCR Detection Kit has a detection limit of ≥ 10 DNA copies per reaction for *P. falciparum*, *P. vivax*, *P. ovale*, *P. malariae* and *P. knowlesi* (Figure 3, 4, 5, 6 and 7).

Figure 3. Dilution series of *P. malariae* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *P. malariae*, *P. knowlesi* & *P. ovale*, channel FAM).


Figure 4. Dilution series of *P. knowlesi* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *P. malariae*, *P. knowlesi* & *P. ovale*, channel ROX).


Figure 5. Dilution series of *P. ovale* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *P. malariae*, *P. knowlesi* & *P. ovale*, channel Cy5).


Figure 6. Dilution series of *P. falciparum* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *P. falciparum* + *P. vivax*, channel FAM).



Figure 7. Dilution series of *P. vivax* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Multiplex reaction mix *P. falciparum* + *P. vivax*, channel ROX).

12.3. Analytical specificity

The specificity of the *Malaria differentiation* assay was confirmed by testing a panel consisting of different microorganisms representing the most common vector-borne pathogens, and other 18S-possessing homologous eukaryotes. No cross-reactivity was detected between almost any of the following microorganisms tested, except the targeted pathogens of each assay.

Cross-reactivity testing					
<i>Anaplasma marginale</i>	-	Chikungunya virus strain S27 Petersfield	-	<i>Theileria annulata</i>	-
<i>Aspergillus fumigatus</i>	-	Chikungunya virus Martinique isolate	-	Tick Borne Encephalitis Virus (TBEV) strain Neudorfl	-
<i>Bartonella henselae</i> strain Houston-1	-	Chikungunya virus strain F24	-	<i>Toxoplasma gondii</i> (Type II)	-
<i>Borrelia hermsii</i>	-	Dengue 1 virus strain Hawaii	-	<i>Treponema phagedenis</i>	-
<i>Borrelia lusitanae</i>	-	Dengue 2 virus strain New Guinea C	-	<i>Trypanosoma cruzi</i>	-
<i>Borrelia valaisiana</i>	-	Dengue 3 virus strain H87	-	Usutu Virus	-
<i>Borrelia azfeli</i> strain P-Ko/1984	-	Dengue 4 virus strain H241	-	West Nile virus strain NY99	-
<i>Borrelia bavariensis</i>	-	Japanese encephalitis	-	West Nile virus Heja	-
<i>Borrelia bisetti</i>	-	Japanese Encephalitis virus Nakayama	-	West Nile virus Ug37	-
<i>Borrelia burgdorferi</i> sensu stricto strain IRS	-	<i>Leptospira</i>	-	Yellow Fever virus strain 17D	-
<i>Borrelia burgdorferi</i> sensu stricto strain B31	-	<i>Pneumocystis jirovecii</i>	-	Yellow Fever Virus French Neurotropic	-
<i>Borrelia garinii</i>	-	<i>Plasmodium falciparum</i> 3D7	-/+	Zika Virus strain 11474/16 (French Polynesia)	-
<i>Borrelia japonica</i>	-	<i>Rickettsia conorii</i> strain Moroccan	-	Zika Virus strain 11468/16(French Polynesia)	-
<i>Borrelia miyamotoi</i>	-	Rift Valley Fever Virus AR21229	-	Zika Virus (African)	-
<i>Borrelia spielmanii</i>	-	Rift Valley Fever Virus MP12	-	Zika virus strain PF13/251013-18 (Asian)	-
<i>Coxiella burnetii</i> strain Nine Mile Q	-	St Louis Encephalitis virus	-	Zika virus strain FB-GWUH-2016	-
Chikungunya virus WHO IS (R91064)	-				

Table 11. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE *Malaria differentiation* Real Time PCR Detection Kit was evaluated against *Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium ovale*, *Plasmodium malariae* and *Plasmodium knowlesi*, showing positive results.

ANNEX 1

COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast Real-Time PCR System ⁽¹⁾
Applied Biosystems	7500 Fast Dx Real-Time PCR System ⁽¹⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 3 Fast Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Fast Real-Time PCR System
Applied Biosystems	StepOne Plus™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	StepOne™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	Mini Opticon™ Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
Qiagen	Rotor-Gene® Q ⁽³⁾
Roche	LightCycler® 480 Real-Time PCR System ⁽⁴⁾
Roche	LightCycler® 96 Real-Time PCR System ⁽⁴⁾
Roche	Cobas z480 Analyzer ⁽⁴⁾

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

(1) Select Ramp Speed "Standard".

(2) See Annex 3 to check optical measurement exposure setting.

(3) The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Rotor-Gene® Q or SmartCycler® tubes.

(4) Shell Frame grid plate which fits in these Roche qPCR System is necessary.

(5) No detection in Cy5 channel.

(6) Detection in FAM and HEX channels only.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

ANNEX 2

DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation is required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

ANNEX 3

OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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VIASURE <i>Malaria differentiation</i> Real Time PCR Detection Kit	FTD <i>Malaria differentiation</i> (Fast Track Diagnostics)			
		+	-	Total
	+	2	0	2
	-	0	224	224
	Total	2	224	226

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VIASURE <i>Malaria differentiation</i> Real Time PCR Detection Kit	FTD <i>Malaria differentiation</i> (Fast Track Diagnostics)			
		+	-	Total
	+	5	0	5
	-	0	221	221
	Total	5	221	226

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VIASURE <i>Malaria differentiation</i> Real Time PCR Detection Kit	FTD <i>Malaria differentiation</i> (Fast Track Diagnostics)			
		+	-	Total
	+	55	0	55
	-	0	171	171
	Total	55	171	226

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13. Bibliography/Bibliografía

1. S. E. Shokoples *et al.* Multiplexed Real-Time PCR Assay for Discrimination of *Plasmodium* Species with Improved Sensitivity for Mixed Infections. *Journal of clinical microbiology*, 2009, p. 975–980.
2. M. Rougemont *et al.* Detection of Four *Plasmodium* Species in Blood from Humans by 18S rRNA Gene Subunit-Based and Species-Specific Real-Time PCR Assays. *Journal of Clinical Microbiology*, 2004, p. 5636–5643.
3. P. CS. Divis *et al.* A TaqMan real-time PCR assay for the detection and quantitation of *Plasmodium knowlesi*. *Malaria Journal*, 2010, 9:344.
4. A. Moody *et al.* Rapid Diagnostic Tests for Malaria Parasites. *Clinical Microbiology Reviews*, 2002. 15(1): 66–78.
5. PF Mens *et al.* Detection and identification of human *Plasmodium* species with real-time quantitative nucleic acid sequence-based amplification. *Malaria Journal*, 2006. 3:5:80.

6. I. Safeukui *et al.* Evaluation of FRET real-time PCR assay for rapid detection and differentiation of *Plasmodium* species in returning travellers and migrants. *Malaria Journal*, 2008. 7: 70.
7. N. Tangpukdee *et al.* Malaria Diagnosis: A Brief Review. *The Korean Journal of Parasitology*, 2009. Vol. 47, No. 2: 93-102.
8. N. Steenkeste *et al.* Towards high-throughput molecular detection of *Plasmodium*: new approaches and molecular markers., *Malaria Journal*, 2009. 29;8:86.
9. J. Maltha *et al.* Malaria rapid diagnostic tests: *Plasmodium falciparum* infections with high parasite densities may generate false positive *Plasmodium vivax* pLDH lines. *Malaria Journal*, 2010.10;9:198.
10. N. Bourgeois *et al.* Comparison of three real-time PCR methods with blood smears and rapid diagnostic test in *Plasmodium* sp. Infection. *Clinical Microbiology and Infection*, 2010 Volume 16, Issue 8, Pages 1305–1311.
11. E. Kamau *et al.* Multiplex qPCR for Detection and Absolute Quantification of Malaria. *Ploze One*, 2013. 29.
12. A. K *et al.* Rapid diagnostic tests for diagnosing uncomplicated non-falciparum or *Plasmodium vivax* malaria in endemic countries. *Cochrane Database Syst Rev*, 2014. 18;(12):CD011431.
13. S. A. Bickersmith *et al.* A sensitive, specific and reproducible real-time polymerase chain reaction method for detection of *Plasmodium vivax* and *Plasmodium falciparum* infection in field-collected anophelines. *Memórias do Instituto Oswaldo Cruz*, 2015. vol.110 no.4.
14. F. D. Krampa *et al.* Recent Progress in the Development of Diagnostic Tests for Malaria., *Diagnostics (Basel)*, 2017; 7(3): 54.
15. N. Siwal *et al.* Malaria diagnosis by PCR revealed differential distribution of mono and mixed species infections by *Plasmodium falciparum* and *P. vivax* in India. *Ploze One*, 2018. 22.
16. J. M. Roth *et al.* *Plasmodium* Detection and Differentiation by Direct-on-Blood PCR Nucleic Acid Lateral Flow Immunoassay. *The Journal of Molecular Diagnostics*, 2018. Vol. 20, nº 1.
17. <https://www.who.int/malaria/areas/diagnosis/rapid-diagnostic-tests/en/>
18. <https://www.malariasite.com/rdt/>
19. https://www.cdc.gov/malaria/diagnosis_treatment/diagnostic_tools.html#tabs-2-2
20. <https://www.who.int/tdr/publications/documents/rdt3.pdf?ua=1>

14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Tabla A.2 TERMOCICLADORES CON BLOQUE DE PERFIL ALTO	
Fabricante	Modelo
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: November 2019

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
www.cer-test.es

VIASURE online

F-362 rev01

VIASURE

CE

IVD

Real Time PCR Detection Kits

CerTest
BIOTEC