CERTEST Calprotectin 50+200

ONE STEP Calprotectin 50+200
COMBO CARD TEST

CERTEST BIOTEC S.L.

CERTEST Calprotectin 50+200

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ENGLISH

INTENDED USE

CerTest Calprotectin one step combo card test is a coloured chromatographic immunoassay (cut-off qualitative test) for the simultaneous qualitative determination of human calprotectin (hCp) in stool samples, to the following cut-off values: A: 50µg/g faeces and B: 200µg/g faeces.

CerTest Calprotectin 50+200 combo card test offers a simple, a highly sensitive and non-invasive screening assay to determine inflammatory disease, to monitor treatment response in these patients, to predict risk of relapse and to monitor small bowel transplantation or graft rejection.

INTRODUCTION

Calprotectin is a neutrophil cytosolic protein with antimicrobial properties, which is present at increased concentration in stool samples during bowel inflammation. The stability of the protein to degradation keeps it stable in faeces for up to 7 days at room temperature, making it an ideal analyte. Calprotectin is released by activation of leukocytes, giving increased levels in plasma, cerebral spinal fluid, synovial fluid, urine or stools as a consequence of disease in the relevant organ(s). Calprotectin inhibits zinc-dependent enzyme systems, as a result kills microbes and induces apoptosis in normal and cancer cells. In the presence of calcium, calprotectin is a remarkably resistant to proteolytic degradation, and so is stable in stools kept at room temperature for 7 days.

This is a non-invasive marker of intestinal inflammation (for example in Ulcerative Colitis (UC) and Crohn's Disease (CD)).

PRINCIPLE OF THE TEST

CerTest Calprotectin 50+200 is based on the principle of a cut-off qualitative immunochromatographic assay for the determination of human calprotectin in stool samples.

Strip A consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against human calprotectin and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-human calprotectin) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

Strip B consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against human calprotectin and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-human calprotectin) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is human calprotectin positive (concentration $\geq 50\mu\text{g/g}$), the antigens of the diluted sample react with the red-coloured conjugates complex (anti-human calprotectin monoclonal antibodies-red polystyrene microspheres) in the strip A, and if the sample is human calprotectin positive (concentration $\geq 200\mu\text{g/g}$), the antigens of the diluted sample react with the red-coloured conjugated complex (anti-human calprotectin monoclonal antibodies-red polystyrene microspheres) in strip B, which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane, by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-human calprotectin antibodies present on the membrane of strip A (test line) and the anti-human calprotectin antibodies present on the membrane of strip B (test line) capture the coloured conjugate and the red line will be visible in both strips. These bands are used to interpret the result.

If the sample is negative, there is no human calprotectin presence and yet, the antigens may be present in a concentration lower than the cut off value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-human calprotectin antibodies present on the membranes (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in both strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on both membranes will capture control green-conjugate complex and both control lines will always appear. The presence of these green lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.


- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest Calprotectin 50+200 test. Do not use any other commercial kit component.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. For longer storage, maximum 6 months, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

1. Take out the cap 1 of the stool collection tube and use the stick to pick up sufficient sample quantity. Introduce the stick once into 4 different parts of the stool sample (1); making sure that, in each insertion, the stick's screw is covered with the sample to take the appropriate amount of stool. Then, introduce the sample in the stool collection tube (2). For liquid samples, add approx. 20µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample (2). Shake the tube in order to assure good sample dispersion (3). The stool collection tube with dilute sample can be stored for 7 days in the refrigerator (2-8°C) prior to testing.

MATERIALS

MATERIALS PROVIDED

- CerTest Calprotectin 50+200 combo card tests
- Instructions for use
- Stool collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE


Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove the CerTest Calprotectin 50+200 combo card test from its sealed bag just before using it.
3. Take the stool collection tube, open the cap 2 (4) and dispense 3 drops in the circular window marked with the letter A (5) and 3 drops, using the same tube, in the circular window marked with the letter B (6). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of extraction buffer until seeing the liquid running through the reaction zone.

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



Calprotectin 200 Strip B-procedures


Add 3 drops in B-circular window


INTERPRETATION OF THE RESULTS (please refer to the illustration below)

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












	A (Calprotectin 50 (concentration ≥ 50µg/g))	B (Calprotectin 200 (concentration ≥ 200µg/g))	Interpretation of the results
1.	-	-	Calprotectin marker is not present in patient sample (<50µg/g), which might mean neither active gastrointestinal inflammation, nor risk of relapse (CD or UC relapse), nor graft or transplantation rejection.
	GREEN	GREEN	
2.	+	+	High quantity of calprotectin marker (≥200µg/g) is present in patient sample, which might mean active gastrointestinal inflammation, risk of relapse in clinical remission, graft rejection or transplantation rejection.
	GREEN-RED	GREEN-RED	
3.	+	-	Calprotectin marker is present in patient sample (≥50µg/g and <200µg/g), which might mean active gastrointestinal inflammation.
	GREEN-RED	GREEN	
4.	-	+	Invalid result. If B is positive, A has to be also positive.
	GREEN	GREEN-RED	
5.	Any other result	Any other result	Invalid result, either A or B, we recommend repeating the assay using the same sample with another test.

INVALID: Total absence of any control coloured line (GREEN) regardless the appearance or not of the tests lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay using a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of human calprotectin present in the specimen.

QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours after opening of the sealed bag.
2. An excess of stool sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
3. The intensity of test line may vary from very strong at high human calprotectin concentration to faint when the human calprotectin concentration is close to cut-off value of the test.
4. CerTest Calprotectin 50+200 should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
5. Positive results determine the presence of human calprotectin in faecal samples; nevertheless, it can be due to several causes (Inflammatory Bowel Disease, colorectal cancer and some other enteropathies). A positive result should be followed up with additional diagnostic invasive procedures, colonoscopy and biopsy in order to confirm the diagnosis and to establish the inflammation extent.
6. A negative result is not meaningful because of it is possible the human calprotectin concentration in the stool sample is lower than the cut-off values. If symptoms or situation still persist, calprotectin determination should be carried out invasive techniques. Negative results do not exclude inflammation, some diseases such as celiac sprue and microscopic colitis polyps that mainly involve mononuclear inflammation.
7. Stool samples from patients with non-steroidal anti-inflammatory drug (NSAID) treatment might show positive result.
8. Neonatal faecal calprotectin levels have been reported higher than those in normal children with a mean of 167µg/g (range 22-860µg/g).
9. Stool from patients with active inflammatory bowel disease that usually involve significant neutrophilic inflammation of the intestine, such as Crohn's disease and ulcerative colitis, would be positive for faecal calprotectin. CerTest Calprotectin 50+200 could be sensitive for this diagnosis in patients with chronic diarrhoea.
10. Bloody stool samples can contain components that may cause non-specific reactions in the test. Every bloody stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
11. Mucous stool samples may cause non-specific reactions in the test. Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.

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EXPECTED VALUES

Some studies established equal or higher 50µg hCp/g faeces as cut-off value to allow detect adult patients with gastrointestinal inflammatory pathology that will require to diagnosis additional diagnostic invasive procedures.

The clinical course of inflammatory bowel disease is characterized by a succession of relapses and remissions. Some studies established equal or higher 200µg hCp/g faeces as cut-off value to allow predict clinical relapse of disease activity in patients with Ulcerative Colitis and in Crohn's Disease and detect some neonatal patients (higher levels than normal children) with gastrointestinal inflammatory pathology that will require to diagnosis additional diagnostic invasive procedures.

Faecal calprotectin level evolution in individuals showed that faecal calprotectin level can predict rejection days before histopathological diagnosis.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

Cut-off values of CerTest Calprotectin 50+200 are 500ng/mL (50µg hCp/g faeces) for strip A and 2000ng/mL (200µg hCp/g faeces) for strip B.

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest Calprotectin 50+200, CerTest) and another commercial immunoassay (Calprest®, Eurospital). The results were as follows:

		ELISA test: Calprest®		
		+	-	Total
IC test: CerTest Calprotectin 50+200	+	34	2	36
	-	2	26	28
	Total	36	28	64

	Sensitivity	Specificity	PPV	NPV
CerTest Calprotectin 50+200 vs ELISA test: Calprest®	>94%	93%	>94%	93%

The results showed a high sensitivity and specificity to detect human calprotectin using CerTest Calprotectin 50+200.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest Calprotectin 50+200; no cross reactivity against other gastrointestinal pathogens, other organism, substances and/or faecal markers occasionally present in faeces:

Adenovirus	Clostridium perfringens	Human Lactoferrin	Salmonella typhi
Astrovirus	Cryptosporidium	Human Transferrin	Salmonella typhimurium
Bovine Haemoglobin	Entamoeba dispar	Legionella pneumophila	Shigella boydii
Bovine Lactoferrin	Entamoeba histolytica	Listeria monocytogenes	Shigella dysenteriae
Bovine Transferrin	Escherichia coli O111	Norovirus GI	Shigella flexneri
Campylobacter coli	Escherichia coli O26	Norovirus GII	Shigella sonnei
Campylobacter jejuni	Escherichia coli O157	Pig Haemoglobin	Streptococcus pneumococcal
Clostridium difficile antigen GDH	Giardia	Rotavirus	Streptococcus pyogenes
Clostridium difficile Toxin A	Helicobacter pylori	Salmonella enteritidis	Yersinia Enterocolitica O:3
Clostridium difficile Toxin B	Human Haemoglobin	Salmonella paratyphi A	Yersinia Enterocolitica O:9

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1. Angriman I. et al. Enzymes in feces: Useful markers of chronic inflammatory bowel disease. *Clinica Chimica Acta* 381 Feb 2007, p. 63-68
2. Quail, M.A. et al. Fecal Calprotectin Complements Routine Laboratory Investigations in Diagnosing Childhood Inflammatory Bowel Disease. *Inflamm Bowel Dis*, Vol 15 No 5; May 2009, p. 756-759
3. Gaya D.R., et al. Faecal calprotectin in the assessment of Crohn's disease activity. *Q J Med* 2005, Vol 98, May 2005, p. 435-441.
4. Langhorst, M.D. et al. Noninvasive Markers in the Assessment of Intestinal Inflammation in Inflammatory Bowel Diseases: Performance of Fecal Lactoferrin, Calprotectin and PMN-Elastase, CRP, and Clinical Indices. *Am. J. Gastroenterol.* 2008; Vol 103, p. 162-169.
5. Tibble JA., et al., "Surrogate markers of intestinal inflammation are predictive of relapse in patients with inflammatory bowel disease". *Gastroenterology*, Vol 119 (1), Jul 2000, p. 15-22.
6. Gisbert J.P. et al. "Utilidad de la calprotectina y lactoferrina fecal en la predicción de la recidiva de la enfermedad inflamatoria intestinal (EII)". *Gastroenterología y Hepatología*, Vol. 32, Issue 3, March 2009, p. 226-227.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/

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