# Transferrin TURBILATEX<sup>®</sup> Combo

**Transferrin latex turbidimetric assay**

Cat. TL-022TF100ED, TL-022TF200ED

IU-TL-022TF100ED-200ED en-es v.04

## CONTENT/CONTENIDO

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## ENGLISH

## INTENDED USE

Transferrin Turbilatex® is a latex turbidimetric assay **for the quantitative detection of transferrin in human stool samples.**

This assay is simple and widely applicable. This product is optimized for several automated analyser.

For professional *in vitro* diagnostic use only.

## SYNTHESIS

Transferrin, an iron-binding glycoprotein, is essential for transporting iron through the bloodstream, impacting both iron metabolism and overall homeostasis (1). Clinically, it serves as a valuable marker in gastrointestinal (GI) pathologies, particularly those associated with chronic bleeding. Conditions like colorectal cancer, inflammatory bowel disease (IBD), and peptic ulcers often involve subtle GI bleeding that can lead to iron deficiency anemia when persistent. The presence of transferrin in feces, particularly due to its stability against enzymatic degradation in the GI tract, indicates blood loss within the GI system that may not be detected otherwise (2).

In the context of colorectal cancer, transferrin detection is significant because early-stage cancers often present with microscopic bleeding. Similarly, IBD patients may experience chronic GI bleeding, which can lead to complications like severe anemia and impact their quality of life (3). By quantifying fecal transferrin, clinicians can gain insights into the presence and extent of bleeding, enabling a proactive approach to diagnosis and treatment. Thus, transferrin testing is vital in managing GI disorders, providing key information on bleeding sources and helping tailor treatments for underlying pathologies.

## PRINCIPLE

**Transferrin Turbilatex® Combo** is a latex turbidimetric assay **for the quantitative detection of transferrin in human stool samples.** This assay is simple, widely applicable and test results are obtained in a very short time.

The intended use of the test is the detection of specific transferrin antigen at very low concentrations.

Transferrin latex turbidimetric assay is based on agglutination reactions. These involve *in vitro* aggregation of microscopic latex particles. This aggregation consists in the specific reaction between antigen and antibodies, antigen contained in the sample and the antibodies anti-antigen coated on polystyrene latex particles. The sample is mixed with a suspension

containing antibodies against the antigen bound to latex particles. If antigen is present in the sample, it will react with the antibodies and form an aggregate. If no antigen is present in the sample the mixture will keep its appearance as a smooth suspension. Such turbidity is measured as an increase in absorbance at determinate wavelength and is proportional to the quantity of antigen contained in the sample.

## PRECAUTIONS

---

- For professional *in vitro* diagnostic use only.
- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in Turbidimetric technique, autoanalyzer use and *in vitro* diagnostic procedures.
- The analyser should be ready before performing any assay.
- Read the instructions for use carefully before use the kit.
- Do not use after expiration date.
- Do not use the reagents if pack is damaged or opened.
- Do not freeze to maintain the integrity and performance of the test.
- Do not mix screw caps from different reagent vials.
- Components provided in the kit are approved for use in **Transferrin Turbilatex® Combo**. Do not mix or use with components from other product kit or lots. Do not use with other commercial kit components (from another supplier). The test's satisfactory performance is ensured only when components from the same lot of **Transferrin Turbilatex® Combo** are used.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- All the specimens should be considered potentially hazardous as well as all the materials that have been exposed to the samples and they must be handled according to the local or national safety regulations. Must be handled in the same manner as an infectious agent. Use proper infection control practices. These practices should include, but are not limited to, personal protective equipment (PPE), such as laboratory coat, surgical or appropriate mask, or face shield, disposable gloves, and eye protection. Take necessary precautions during the collection, transport, storage, handling, and disposal of the samples.
- **Transferrin Turbilatex® Combo** do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which are in concentrations higher than the value established in the mentioned regulation for

their declaration. Material Safety Data Sheet are not included with this device, nevertheless it could be consulted upon request to Certest Biotec.

- **Transferrin Turbilatex® Combo** has been validated with the analyser Biolis from Tokyo Boeki. For different analysers, transferrin should be measured according to the specific application parameters for each specific chemistry analyzer. More information about the application notes available can be found on our website.
- Homogenise the reagents before use, make sure that the reagents are free of bubbles.
- Evaporation of calibrators and controls can lead to incorrect results. Follow the storage recommendations detailed in the following section.
- If the values obtained for the controls are not valid, repeat the quality control measurement with new controls. If the control values are still invalid, recalibrate the equipment.
- If the result exceeds the measurement range, use the sample diluent to dilute the sample and repeat the assay again.
- The solutions should be discarded in a proper container after testing following local regulations.

## STORAGE AND STABILITY

Kit components must be stored at refrigerated temperature (2-8°C). **Do not freeze.**

Reagents are stable up to the expiration date printed on the label, always considering that reagent containers must be properly closed to avoid any contamination, kept away from the sunlight and conserved at 2-8°C.

## REAGENTS

### Materials provided

<u>References</u>	<u>Cat.:</u> <u>TL-022TF100ED</u>	<u>Cat.:</u> <u>TL-022TF200ED</u>
- TL-022TF01 Transferrin Turbilatex® Reagent 1	1x22mL	2 x 22mL
- TL-022TF02 Transferrin Turbilatex® Reagent 2	1 x 8mL	1 x 13mL
- TL-022TF70/ TL-022TF71/ TL- 022TF72- TL-022TF73/ TL- 022TF74/ TL-022TF75 Transferrin Turbilatex® Calibrators	6 x 1mL	
- TL-022TF08/ TL-022TF09 Transferrin Turbilatex® Control 1&2	4 x 1mL	

### Materials and instruments that could be required but no provided

- Specimen collection container.
- Disposable gloves and laboratory equipment.

- Automated analyser.
- MST-0019U Universal Turbilatex® Sample Collection Vial.

## SPECIMEN COLLECTION AND PREPARATION

Collect enough quantity of human solid stool samples. These samples should be collected in clean and dry containers (no preservatives or transport media). The samples can be stored in the refrigerator (2-8°C) for 3 days prior to testing. If not immediately tested, freeze the stored samples at -20 °C maximum 6 months. In this case, the sample will be totally thawed, and brought to room temperature (15-30°C) before testing. Homogenize stool samples as thoroughly as possible prior to preparation.

**To process the collected stool samples using catalogue reference MST-0019U follow the indication of the instruction for use of the product.**

## ANALYTICAL PROCEDURE

R1 and R2 are ready to use.

### Preparation of the calibration curve


For calibration use only the following materials:

TL-022TF70, TL-022TF71, TL-022TF72, TL-022TF73, TL-022TF74 and TL-022TF75, Transferrin Turbilatex® Calibrator vials: liquid calibrator containing recombinant transferrin antigen at different concentrations.

Concentration is indicated on the label of each vial.

Calibrate the system at least once a month is extremely recommended. Recalibrate the system when reagent lot is changed or when the controls are out of the assigned range given on the control label or certificate of analysis provided with the kit.

Transferrin Turbilatex® Calibrator vials (calibrator 0 to 5) are ready to use.

**Picture 1.** Prepared calibrators for performing the calibration curve (5 points and blank)

## Quality control


For quality control only use the following materials:

TL-022TF08 Transferrin Turbilatex® Control 1 and TL-022TF09 Transferrin Turbilatex® Control 2: liquid controls at two different concentrations of recombinant transferrin antigen. Concentration is indicated on the label of the vial. The use of control materials at two different concentrations is recommended in order to verify test precision.

Controls should be assayed every day before running patient faecal samples extract to validate the calibration curve.

If the obtained results are out of the tolerance range, the equipment, the reagents or the technique must be reviewed. Transferrin Turbilatex® Control vials (Control 1 and Control 2) are ready to use.

### Control vials

Picture 2. Prepared controls C1 and C2 (check control concentration on the label)

## INTERPRETATION OF RESULTS

Positive results: higher or equal than the cut-off fixed by the clinical lab.

Recommended: 10 ng of transferrin/mL (1 µg of transferrin/g of stool) for diagnostic procedures.

Positive results determine the abnormal presence of transferrin in stool samples.

**Transferrin Turbilatex® Combo** cut-off value:

Transferrin concentration values lower than 10 ng of hTf/mL are considered normal values and that is not indicative of bleeding in the gastrointestinal tract.

Transferrin concentration values equal or higher than 10 ng hTf/mL are considered abnormal values and that is indicative of bleeding in the gastrointestinal tract.

## PERFORMANCE CHARACTERISTICS

### Analytical sensitivity

Limit of detection (LOD): 1.4 ng hTf/mL. The lower limit of detection of **Transferrin Turbilatex® Combo** was determined on 20 samples and 2 sample replicates as the mean value + 2 SD.

Limit of quantification (LOQ): 2 ng hTf/mL. The lower limit of quantification is defined as the lowest actual amount of analysis that can be reliably detected; imprecision is < 20% as CV%.

### Measuring range

**Transferrin Turbilatex® Combo** assay measuring range is 2-1250 ng hTf/mL. Samples higher concentrated than 250 ng hTf/mL must be diluted for proper quantification by the user, using additional sample buffer.

### Prozone

Studies have been made up to a concentration of 10000 ng hTf/mL and no false negative results have been observed. Studies using higher concentrations have not been carried out.

Samples with concentrations up to 1250 ng hTf/mL can be measured without inhibitory prozone effect.

### Clinical sensitivity and specificity

An evaluation was performed comparing a turbidimetric assay (**Transferrin Turbilatex® Combo**, CerTest) and an immunochromatographic test (CerTest Transferrin, CerTest). The results were as follows:

	CerTest Transferrin			Total
		+	-	
Transferrin Turbilatex® Combo	+	36	0	36
	-	2	32	34
	Total	38	32	70

Table 1 **Transferrin Turbilatex® Combo** performance versus CerTest Transferrin.

Transferrin Turbilatex® vs CerTest Transferrin		
	Mean Value	95% confidence interval
Sensitivity	94.7%	82.3-99.4
Specificity	100.0%	90.3-100.0%
PPV	100.0%	89.1-100.0%
NPV	94.1%	80.3-99.3%
LR+	61.77	3.941-968.1
LR-	0.065	0.02-0.216

Table 2 **Transferrin Turbilatex® Combo** clinical results. True positive (TP), false positive (FP), false negative (FN), true negative values (TN), and total samples (N). Sensitivity (Sens), specificity (Spec), Predictive Positive Values (PPV), Predictive Negative Values (NPV), 95% interval of confidence (CI), Likelihood positive ratio (LR+), Likelihood negative ratio (LR-) for **Transferrin Turbilatex® Combo**.

The results showed a high sensitivity and specificity to detect transferrin using **Transferrin Turbilatex® Combo**.

## Interferences

An evaluation was performed to determine the interferences of **Transferrin Turbilatex® Combo**. No interferences were founded against other substances tested: BSA, Hemin, Ascorbic acid.

## Cross reactivity

An evaluation was performed to determine the cross reactivity; no cross reactivity was founded against other faecal markers occasionally present in faeces, such as:

bovine and pig haemoglobin, bovine transferrin, human haemoglobin, calprotectin and lactoferrin.

## Precision

An evaluation was performed to determine the accuracy of **Transferrin Turbilatex® Combo**. **Transferrin Turbilatex® Combo** shows high in its reproducibility results on known concentration of Transferrin.

## Linearity

**Transferrin Turbilatex® Combo** using calibrator kit is linear in the calibration range of 3.7-250 ng hTf/mL.

The analytical performance of **Transferrin Turbilatex® Combo** has been validated using the Tokyo Boeki Biolis 24i analyser. For more information about other compatible equipment, please refer to the corresponding application. The different turbidimetric equipment with which our turbidimetric products are compatible can be found in **Annex I**.

## LIMITATIONS

---

1. **Transferrin Turbilatex® Combo** should **only be used in human stool samples**. The use of other samples has not been established. Proper faecal specimens must be obtained.
2. Positive results determine the presence of transferrin in faecal samples. A positive result should be followed up with additional invasive procedures to determine the exact cause and source of the blood in the stool.

3. If symptoms or situation still persist, transferrin determination should be carried out invasive techniques. Negative results do not exclude bleeding, as some polyps and colorectal cancers may bleed intermittently or not during certain stages of the disease. Moreover, blood may not be uniformly distributed in stool samples.
4. Patients should not collect samples during their menstrual period, if they have bleeding haemorrhoids, blood in urine or if they have strained during bowel movement.

## EXPECTED VALUES

---

Colorectal cancer (CRC) is a leading cause of cancer incidence and mortality globally. According to the World Health Organization, it is the third most common cancer worldwide, accounting for approximately 10% of all cancer cases and is the second leading cause of cancer-related deaths worldwide (4). In Europe, it ranks as the second most common cancer after breast cancer, with over 500,000 new cases reported annually (5). The increasing incidence among individuals under 50 has prompted a shift in screening recommendations, with guidelines now advocating routine screening starting at age 45 (6).

Globally, enhancing early detection through screening strategies is crucial, as early-stage CRC is more treatable, leading to better patient outcomes (7). Implementing standardized screening protocols can significantly reduce CRC-related morbidity and mortality across diverse populations.

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1. Ogun AS, Adeyinka A. Biochemistry, Transferrin. StatPearls. 2022 Nov 16;
2. Jin P, Wu Z tao, Meng M ming, Wang X, Wang X wei, Gong L juan, et al. Combined fecal transferrin test and immuno fecal occult blood test for detecting colorectal cancer and advanced adenoma in asymptomatic and symptomatic populations. J Cancer Sci Ther. 2012;4(8):243–8.
3. Chardalias L, Papaconstantinou I, Gklavas A, Politou M, Theodosopoulos T. Iron Deficiency Anemia in Colorectal Cancer Patients: Is Preoperative Intravenous Iron Infusion Indicated? A Narrative Review of the Literature. Cancer Diagnosis & Prognosis. 2023 Mar 3;3(2):163.
4. World Health Organization. Colorectal cancer. Fact sheet. [Internet]. 2023 [cited 2024 Oct 29]. Available from: <https://www.who.int/news-room/fact-sheets/detail/colorectal-cancer>
5. International Agency for Research on Cancer. Cancer Today. Absolute numbers, Incidence, Both sexes, in 2022 [Internet]. 2022 [cited 2024 Oct 29]. Available from: [https://gco.iarc.who.int/today/en/dataviz/pie?mode=cancer&group\\_populations=1&populations=908](https://gco.iarc.who.int/today/en/dataviz/pie?mode=cancer&group_populations=1&populations=908)
6. American Cancer Society. Colorectal Cancer Early Detection, Diagnosis, and Staging [Internet]. 2024. Available from: [https://seer.cancer.gov/csr/1975\\_2016/](https://seer.cancer.gov/csr/1975_2016/),
7. Young GP, Symonds EL, Allison JE, Cole SR, Fraser CG, Halloran SP, et al. Advances in Fecal Occult Blood Tests: The FIT Revolution. Dig Dis Sci. 2015 Mar 1;60(3):609–22.

## SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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The application notes for use our turbidimetric products with the different equipments could be consulted in our website: [www.certest.es](http://www.certest.es)

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