
CerTest
BIOTEC

CerTest Turbilatex

PANCREATIC ELASTASE TURBILATEX® Combo

Pancreatic Elastase latex turbidimetric assay

Cat.TL-022EL100ED, TL-022EL200ED, TL-022EL400ED, TL-022EL800ED

IU-TL -022EL100ED-200ED-400ED-800ED en-es v.1020

ENGLISH

INTENDED USE

Pancreatic Elastase Turbilatex® is a latex turbidimetric assay **only for the quantitative detection of Pancreatic Elastase E1 in human stool samples.**

This assay is simple and widely applicable. This product is optimized for several automated analyser.

Test results should **exclusively be used to evaluate exocrine pancreatic function in stool samples.**

For professional *in vitro* diagnostic use only.

SYNTHESIS

Exocrine pancreatic insufficiency (IPE) is one of the main complications in pancreas illness. It is due to the inability of the pancreas to synthesize this enzyme in the required amount. Pancreatic elastase E1 is produced in the acinar cells of the pancreas. This enzyme is involved in food digestion, and its value is clinically relevant. Human pancreatic elastase 1 (E1) remains undegraded during intestinal tract, becoming in the perfect marker to determine IPE in stool samples

Pancreatic elastase Turbilatex® is a non-invasive *in vitro* diagnostic method, based on monoclonal antibodies and used to monitor exocrine pancreatic function.

PRINCIPLE

Pancreatic Elastase Turbilatex® is a latex turbidimetric assay **only for the quantitative detection of Pancreatic elastase E1 in human stool samples.** This assay is simple, widely applicable and test results are obtained in a very short time.

The intended use of the test is **exclusively be used to evaluate Pancreatic elastase levels in stool samples associated to Pancreatic insufficiency.**

Pancreatic Elastase latex turbidimetric assay is based on agglutination reactions. These involve *in vitro* aggregation of microscopic latex particles. This aggregation consists in the specific reaction between antigen and monoclonal antibodies, antigen contained in the sample and the antibodies anti-antigen coated on polystyrene latex particles. The sample is mixed with a suspension containing antibodies against the antigen bound to latex particles. If antigen is present in the sample, it will react with the antibodies and form an aggregate. If no antigen is present in the sample the mixture will keep its appearance

as a smooth suspension. Such turbidity is measured as an increase in absorbance at the determined wavelength and is proportional to the quantity of antigen contained in the sample.

PRECAUTIONS

- For professional *in vitro* diagnostic use only.
- A trained person in Turbidimetric technique and autoanalyzer use is required.
- The analyser should be ready before performing any assay.
- Read and follow the instructions for use provided with the kit.
- If the result exceeds the measurement range, use the sample diluent to dilute the sample and repeat the assay again.
- Do not use after expiration date.
- Do not use the reagents if pack is damaged or opened.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The solutions should be discarded in a proper container after testing following local regulations.

STORAGE AND STABILITY

Kit components must be stored at refrigerated temperature (2-8°C). Do not freeze.

Reagents are stable up to the expiration date printed on the label, always considering that reagent containers must be properly closed to avoid any contamination, kept away from the sunlight and conserved at 2-8°C.

REAGENTS

Materials provided

References	Cat.: TL-022EL100ED	Cat.: TL-022EL200ED	Cat.: TL-022EL400ED	Cat.: TL-022EL800ED
- TL-022EL01 Pancreatic Elastase Turbilatex® Reagent 1	1x24mL	2 x 24mL	4 x 24mL	8 x 24mL
- TL-022EL02 Pancreatic Elastase Turbilatex® Reagent 2	1 x 4mL	1 x 6mL	2 x 6mL	4 x 6mL

References	Cat.: <u>TL-022EL100ED</u>	Cat.: <u>TL-022EL200ED</u>	Cat.: <u>TL-022EL400ED</u>	Cat.: <u>TL-022EL800ED</u>
- TL-022EL70/ TL-022EL71/ TL-022EL72- TL-022EL73/ TL- 022EL74/ TL-022EL75 Pancreatic Elastase Turbilatex® Calibrators			6 x 1mL	
- TL-022EL08/ TL-022EL09 Pancreatic Elastase Turbilatex® Control 1&2			4 x 1mL	

Materials and instruments that could be required but no provided

- Specimen collection container.
- TL-022UN03E Universal Turbilatex® Sample Diluent (Required for the Procedure A).
- TL-0221213E Microtubes.
- Disposable gloves and laboratory equipment.
- Automated analyser.
- MST-0019U Universal Turbilatex® Sample Collection Vial (Required for the Procedure B).
- MST-0018MU Universal Turbilatex® Sample Collection Vial (Required for the Procedure B).

SPECIMEN COLLECTION AND PREPARATION

Collect enough quantity of human solid stool samples. These samples should be collected in clean and dry containers (no preservatives or transport media). The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. If not immediately tested, freeze the stored samples at -20 °C maximum 6 months. In this case, the sample will be totally thawed, and brought to room temperature (15-30°C) before testing. Homogenize stool samples as thoroughly as possible prior to preparation.

To process the collected stool samples:

Procedure A: using catalogue reference TL-022UN03E.

Use a microtube for each sample to be tested. Label microtube with name or number of patient.

1. Homogenize the sample. Add 20 mg of sample into microtube.
2. Add 2 mL of sample diluent (TL-022UN03E).

3. Shake vigorously the microtube in order to assure good sample dispersion (vortex) until sample is completely dissolved.
4. Centrifuge for 15 minutes at 10000g or 10 minutes at 15000g.
5. Take the supernatant to automated analyser vial. Diluted sample can be stored for 7 days in the refrigerator (2-8 °C).

Procedure B: using catalogue reference MST-0019U/MST-0018MU.

Follow the indication of the instruction for use of the product.

ANALYTICAL PROCEDURE

R1 and R2 are ready to use.

Preparation of the calibration curve


For calibration use only the following materials:

TL-022EL70, TL-022EL71, TL-022EL72, TL-022EL73, TL-022EL74 and TL-022EL75, Pancreatic Elastase Turbilatex® Calibrator vials: liquid calibrator containing recombinant Pancreatic Elastase E1 at different concentrations.

Concentration is indicated on the label of each vial.

Calibrate the system at least once a month is extremely recommended. Recalibrate the system when reagent lot is changed or when the controls are out of the assigned range given on the control label or certificate of analysis provided with the kit.

Pancreatic Elastase Turbilatex® Calibrator vials (calibrator 0 to 5) are ready to use.

Picture 1. Prepared calibrators for performing the calibration curve (5 points and blank)


Quality control

For quality control only use the following materials:

TL-022EL08 Pancreatic Elastase Turbilatex® Control 1 and TL-022EL09 Pancreatic Elastase Turbilatex® Control 2: liquid controls at two different concentrations of recombinant Pancreatic Elastase. Concentration is indicated on the label of the vial. The use of control materials at two different concentrations is recommended in order to verify test precision. Controls should be assayed every day before running patient faecal samples extract to validate the calibration curve.

If the obtained results are out of the tolerance range, the equipment, the reagents or the technique must be reviewed.

Pancreatic Elastase Turbilatex® Control vials (Control 1 and Control 2) are ready to use.

Picture 2. Prepared controls C1 and C2 (check control concentration on the label)

Analytical procedure, see annex information.

INTERPRETATION OF RESULTS

Positive results: Higher or equal than the cut-off fixed by the clinical lab.

Recommended: 100-200 µg of Pancreatic elastase E1/g of stool for diagnostic procedures.

Positive results determine the presence of normal level of Pancreatic Elastase in faeces.

Pancreatic Elastase Turbilatex® cut-off value:

Pancreatic Elastase E1 values equal or higher than 200 µg of Pancreatic elastase E1/g of stool are indicative of a normal pancreatic exocrine functioning.

Pancreatic Elastase E1 values between 100 and 200 µg of Pancreatic elastase E1/g of stool are indicative of mild to moderate pancreatic exocrine insufficiency.

Pancreatic Elastase E1 values lower than 100 µg of Pancreatic elastase E1/g of stool are indicative of a severe pancreatic exocrine insufficiency.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity

Limit of Detection (LOD): **See annex information.**

Prozone

Studies have been made up to a concentration of 10000 µg of Pancreatic Elastase E1/g of stool and no false negative results have been observed. Studies using higher concentrations have not been carried out.

Data obtained by the analyser Biolis 24i (Tokyo Boeki).

Clinical sensitivity and specificity

An evaluation was performed comparing a turbidimetric assay (Pancreatic Elastase Turbilatex®, CerTest) and another commercial immunoassay (Schebo® Pancreatic Elastase 1, Schebo® Biotech AG). The results were as follows:

	Sensitivity	Specificity
Pancreatic Elastase Turbilatex® vs Schebo® Pancreatic Elastase 1	100 %	93 %

The results showed a high sensitivity and specificity to detect human Pancreatic Elastase using Pancreatic Elastase Turbilatex®.

Interferences

An evaluation was performed to determine the interferences of Pancreatic Elastase Turbilatex®. No interferences were founded against other substances tested: BSA, Hemin, Ascorbic acid.

Cross reactivity

An evaluation was performed to determine the cross reactivity; no cross reactivity was founded against porcine pancreatic elastase.

LIMITATIONS

1. Pancreatic Elastase Turbilatex® should **only be used in human stool samples**. Proper faecal specimens must be obtained.
2. Values in the range between 100-200 µg/g should be consider as mild to moderate pancreatic insufficiency, and they should be reviewed by the specialist.

EXPECTED VALUES

Values of pancreatic elastase over 200 µg/g determine normal activity of the pancreas.

Values of pancreatic elastase between 100-200 µg/g determine mild to moderate pancreatic insufficiency and requires follow up.

Values of pancreatic elastase below 100 µg/g determine severe pancreatic insufficiency and requires additional diagnosis to implement the correct treatment.

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1. Andriulli, A. et al. Exocrine pancreatic insufficiency, as assessed by fecal elastase-1 levels in diabetic patients: an estimate of prevalence in prospective studies. J Diabetes Metab 2014, 5:6.
2. Borowitz, D. et al. Use of fecal elastase-1 to classify pancreatic status in patients with cystic fibrosis. J Pediatr 145 : 322-326. 2004.

3. Cade, A. et al. Evaluation of fecal pancreatic elastase as a measure of pancreatic exocrine function in children with cystic fibrosis. *Pedr Pulmonol* 29 (3): 172-176. 2000.
4. Dominguez-Muñoz, J.E. et al. Potential for screening for pancreatic exocrine insufficiency using the fecal elastase-1 test. *Dig Dis Sci.* 62 (5) : 1119-1130.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO *IN VITRO*

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The application notes for use our turbidimetric products with the different equipments could be consulted in our website: www.certest.es

Note/Nota: Turbilatex® registered brand name/ Turbilatex® marca registrada
Nº: CTM-010179299; community (OAMI)/ Nº: CTM-010179299; comunitaria (OAMI)
Nº: 4,650,387 ; EEUU (UPSTO)

CerTest Turbilatex

CerTest
BIOTEC

One step ahead

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)

Tel. (+34) 976 520 354

Fax (+34) 976 106 268

certest@certest.es | turbilatex@certest.es

www.certest.es