---

CERTEST  
*Clostridium*  
*difficile* Toxin  
A+B

---

ONE STEP *Clostridium difficile*  
antigen Toxin A and Toxin B  
COMBO CARD TEST

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### ENGLISH

#### INTENDED USE

CerTest *Clostridium difficile* Toxin A+B one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection Toxin A and Toxin B of *Clostridium difficile* in stool samples.

CerTest *Clostridium difficile* Toxin A+B combo card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of *Clostridium difficile* infection.

#### INTRODUCTION

*Clostridium difficile* (*C. difficile*), a Gram-positive spore bearing anaerobic bacterium is the major aetiological agent of diarrhoea and colitis associated with antibiotics. *C. difficile* is the most common cause of health care-associated diarrhoea in developed countries and is a major source of nosocomial morbidity and mortality worldwide.

Disease due to *C. difficile* develops when the organism is allowed to proliferate in the colon, most commonly after antibiotic use has eliminated competing flora. *C. difficile* can release two high-molecular-weight toxins, toxin A and toxin B, which are responsible for the clinical manifestations, which range from mild, self-limited watery diarrhoea to fulminant pseudomembranous colitis, toxic megacolon, and death.

The toxigenic culture (TC) is used as the gold standard technique to determine *Clostridium difficile* infection. This method consists in culture and isolation of *C. difficile* from faeces, followed by toxin testing of the isolate, a labour-intensive assay to obtain a result.

#### TEST PRINCIPLE

CerTest *Clostridium difficile* Toxin A+B is based on the principle of a qualitative immunochromatographic assay for the determination of Toxin A and Toxin B of *Clostridium difficile* in stool samples.

**Strip A** consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against Toxin A and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-Toxin A) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

**Strip B** consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against Toxin B and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-Toxin B) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

If the sample is Toxin A positive, the antigens of the diluted sample react with the red-coloured conjugates complex (anti-Toxin A monoclonal antibodies-red polystyrene microspheres) in the strip A, and if the sample is Toxin B positive, the antigens of the diluted sample react with the red-coloured conjugated complex (anti-Toxin B monoclonal antibodies-red polystyrene microspheres) in strip B which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-Toxin A antibodies present on the membrane of strip A (test line) and the anti-Toxin B antibodies present on the membrane of strip B (test line) capture the coloured conjugate and the red line will be visible in both strips. These bands are used to interpret the result.

If the sample is negative, there is no Toxin A and Toxin B presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-Toxin A and the anti-Toxin B antibodies present on the membranes (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in both strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on both membranes will capture control green-conjugate complex and both control lines will always appear. The presence of these green lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

#### STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

### PRECAUTIONS


- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the CerTest *Clostridium difficile* Toxin A+B combo card test. Do not use any other commercial kit component.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

### SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 24 hours prior to testing. For longer storage, the specimen must be kept frozen at -20°C. Freezing and thawing cycles are not recommended. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

#### Specimen preparation (see illustration):

1. Take out the cap of the stool collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2) making sure that, at each insertion, only the stick's screw is covered with the sample to take the appropriate amount of faecal sample and add it to the stool collection tube. Any additional sample that exceeds the stick's screw could cause wrong results. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

### MATERIALS

#### MATERIALS PROVIDED

- CerTest *Clostridium difficile* Toxin A+B combo card tests
- Instructions for use
- Stool collection tubes with diluent

#### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

### TEST PROCEDURE


**Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.**

1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove CerTest *Clostridium difficile* Toxin A+B combo card test from its sealed bag just before using it.
3. Take the stool collection tube, cut the end of the cap (4) and dispense 3 drops in the circular window marked with the letter A (5) and 3 drops, using the same tube, in the circular window marked with the letter B (6). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.


If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of extraction diluent until seeing the liquid running through the reaction zone.

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

## Toxin A Strip A-procedure


## Toxin B Strip B-procedure

## INTERPRETATION OF THE RESULTS (please refer to the illustration below)

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	A (Toxin A)	B (Toxin B)	Interpretation of the results
1.	- GREEN	- GREEN	There are no Toxin A and Toxin B of <i>Clostridium difficile</i> presence. No infection caused by <i>Clostridium difficile</i> .
2.	+ GREEN-RED	+ GREEN-RED	There are Toxin A and Toxin B of <i>Clostridium difficile</i> . Infection caused by <i>Clostridium difficile</i> .
3.	+ GREEN-RED	- GREEN	There is Toxin A of <i>Clostridium difficile</i> . Infection caused by <i>Clostridium difficile</i> .
4.	- GREEN	+ GREEN-RED	There is Toxin B of <i>Clostridium difficile</i> . Infection caused by <i>Clostridium difficile</i> .
5.	Any other result	Any other result	Invalid result either A or B, we recommend repeating the assay using the same sample with another test.

**INVALID:** Total absence of any control coloured line (GREEN) regardless the appearance or not of the test lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

### NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

### QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

### LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
3. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
4. CerTest *Clostridium difficile* Toxin A+B should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
5. Positive results determine the presence of Toxin A and Toxin B of *Clostridium difficile* in faecal samples. A positive result should be followed up with additional laboratory techniques (toxigenic culture) to determine the strain. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
6. A negative result is not meaningful because of it is possible the antigens concentration in the stool samples is lower than the detection limit value. If the symptoms or situation still persist, a *Clostridium difficile* determination should be carried out on a sample from an enrichment culture.
7. Bloody stool samples can contain components that may cause non-specific reactions in the test. Every bloody stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
8. Mucous stool samples may cause non-specific reactions in the test- Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.

### EXPECTED VALUES

*Clostridium difficile* is associated with 95-100% cases of pseudomembranous colitis, 60-75% cases of antibiotic-associated colitis and 35% of cases of antibiotic-associated diarrhoea cases.

### PERFORMANCE CHARACTERISTICS

#### Analytical sensitivity (detection limit)

Detection limit values of CerTest *Clostridium difficile* Toxin A+B (Toxin A strip) is: 2-0.25 ng/mL of *Clostridium difficile* Toxin A

Detection limit values of CerTest *Clostridium difficile* Toxin A+B (Toxin B strip) is: 3.12-0.19 ng/mL of *Clostridium difficile* Toxin B.

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### Clinical sensitivity and specificity (Toxin A+B)

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest *Clostridium difficile* Toxin A+B, CerTest) vs the Evaluation Criteria. Evaluation Criteria consisted on the evaluation in parallel of two rapid (CerTest *Clostridium difficile* Toxin A+B) vs other rapid test from the competitor; the discrepant results were confirmed by qPCR technique  
The results were as follows:

		Evaluation Criteria (Toxin A)		
		+	-	Total
IC test: CerTest <i>Clostridium difficile</i> Toxin A+B (Toxin A)	+	28	0	28
	-	1	97	98
	Total	29	97	126

CerTest <i>Clostridium difficile</i> Toxin A+B (Toxin A) vs Evaluation Criteria		
	Mean Value	95% confidence interval
<b>Sensitivity</b>	96.9%	92.2–99.9%
<b>Specificity</b>	100.0%	96.2–100.0%
<b>PPV</b>	100.0%	87.7–100.0%
<b>NPV</b>	99.0%	94.4–100.0%

		Evaluation Criteria (Toxin B)		
		+	-	Total
IC test: CerTest <i>Clostridium difficile</i> Toxin A+B (Toxin B)	+	36	1	37
	-	0	89	89
	Total	36	90	126

CerTest <i>Clostridium difficile</i> Toxin A+B (Toxin B) vs Evaluation Criteria		
	Mean Value	95% confidence interval
<b>Sensitivity</b>	100.0%	90.3–100.0%
<b>Specificity</b>	98.9%	94.0–100.0%
<b>PPV</b>	97.3%	85.8–99.9%
<b>NPV</b>	100.0%	95.9–100.0%

The results showed a high sensitivity and specificity to detect Toxin A and Toxin B of *Clostridium difficile* using CerTest *Clostridium difficile* Toxin A+B.

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### Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *Clostridium difficile* Toxin A+B; no cross reactivity against gastrointestinal pathogens, other organism, substances and/or faecal occasionally present in faeces:

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1. WREN, M.W.D, et al. "Laboratory diagnosis of *Clostridium difficile* infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory". British Journal of Biomedical Science, 66 (1), 2009.
2. VAISHNAVI, Ch., "Clinical spectrum & pathogenesis of *Clostridium difficile* associated diseases". Indican J. Med. Res. 131, April 2010, pp 487-499.

### SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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