CERTEST

*Giardia*²

ONE STEP
*Giardia*² COMBO CARD TEST

CERTEST BIOTEC S.L.

ENGLISH**INTENDED USE**

CerTest *Giardia*² one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection of *Giardia* (CWP1) and *Giardia* (α 1-*giardin*) in stool samples.

CerTest *Giardia*² combo card test offers a simple and a highly sensitive screening assay to make a presumptive diagnosis of giardiasis.

INTRODUCTION

Giardiasis is a diarrheal illness caused by a very small parasite, *Giardia intestinalis* (also known as *Giardia lamblia* and *Giardia duodenalis*). Once an animal or person is infected with *Giardia*, the parasite lives in the intestine and is passed in the stool. The parasite is protected by an outer shell (cyst) and can survive outside the body and in the environment for a long time.

Giardia, a flagellated protozoan, inhabits the upper part of the small intestine of its host and has a two major states in the life cycle: trophozoites which produces the antigens (α 1-*giardin*) and cyst which produces the antigens (CWP1) After the host ingest the cysts, which are the infective stage, the trophozoites emerge from the cysts in the duodenum and attach to the small intestinal mucosa. They undergo mitotic division in the intracellular lumen, some will encyst to protect themselves and will be eliminated from the host in the feces. The trophozoite is the vegetative form and replicates in the small intestine.

In humans, the clinical effects of *Giardia* infection range from the asymptomatic carrier state to a severe malabsorption syndrome. Possible factors contributing to the variation in clinical manifestations include the virulence of *Giardia* strain, the number of cyst ingested, the age of the host, and the state of the host immune system at the time of infection.

Giardia infection has become one of the most common causes of waterborne disease in humans. The most common symptoms of giardiasis include: diarrhea, loose or watery stool, stomach cramps and upset stomach. These symptoms generally begin 1-2 weeks after infection and may last 2-6 weeks in healthy individuals. Sometimes symptoms last longer and may lead to weight loss and dehydration. Some people will have no symptoms. However, people with weakened immune systems (e.g., persons with HIV/AIDS, cancer patients and transplant patients) or the elderly may have a more serious infection that can lead to severe illness or death.

TEST PRINCIPLE

CerTest *Giardia*² is based on the principle of a qualitative immunochromatographic assay for the determination of *Giardia* (CWP1) antigen strip A and *Giardia* (α 1-*giardin*) antigen strip B in stool samples. Thus, the test can detect trophozoites and cyst.

Strip A consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Giardia* (CWP1) and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Giardia* (CWP1)) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

Strip B consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Giardia* (α 1-*giardin*) and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Giardia* (α 1-*giardin*)) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is *Giardia* (CWP1) positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*Giardia* monoclonal antibodies-red polystyrene microspheres) in the strip A, and if the sample is *Giardia* (α 1-*giardin*) positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*Giardia* monoclonal antibodies-red polystyrene microspheres) in strip B, which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane, by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*Giardia* (CWP1) antibodies present on the membrane of strip A (test line) and the anti-*Giardia* (α 1-*giardin*) antibodies present on the membrane of strip B (test line) capture the coloured conjugate and the red line will be visible in both strips. These bands are used to interpret the result.

If the sample is negative, there is no *Giardia* (CWP1) and *Giardia* (α 1-*giardin*) presence and yet, the antigens may be present in a concentration lower than the detection limit values, for which the reaction will not take place with any red-coloured conjugate complex. The anti-*Giardia* (CWP1) and anti-*Giardia* (α 1-*giardin*) antibodies present on the membranes (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in both strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on both membranes will capture control green-conjugate complex and both control lines will always appear. The presence of these green lines serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS


- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the CerTest *Giardia*² combo card test. Do not use any other commercial kit component.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-3 days prior to testing. For longer storage, maximum 1 month, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

1. Take out the cap of the stool collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2) making sure that, at each insertion, only the stick's screw is covered with the sample to take the appropriate amount of faecal sample and add it to the stool collection tube. Any additional sample that exceeds the stick's screw could cause wrong results. For liquid samples, add approx. 125 μ L in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

VERY IMPORTANT: Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.

MATERIALS

MATERIALS PROVIDED

- CerTest *Giardia*² combo card tests
- Instructions for use
- Stool collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE

Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.



1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove CerTest *Giardia*² combo card test from its sealed bag just before using it.

CERTEST *Giardia*²

One Step test to detect *Giardia* (CWP1 and/or α 1-giardin) in combo card format/
 Prueba combo de un solo paso para detección de *Giardia* (CWP1 y/o α 1-giardin) en formato cassette


- Take the stool collection tube, cut the end of the cap (4) and dispense 3 drops in the circular window marked with the letter A (5) and 3 drops, using the same tube, in the circular window marked with the letter B (6). Avoid adding solid particles with the liquid.
- Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.


Giardia (CWP1) Strip A-procedure

Add 3 drops in A-circular window

Giardia (α 1-giardin) Strip B-procedure

Add 3 drops in B-circular window

INTERPRETATION OF THE RESULTS (please refer to the illustration below)

- A: Green \rightarrow Negative *Giardia* (CWP1)
 B: Green \rightarrow Negative *Giardia* (α 1-giardin)
- A: Green/Red \rightarrow Positive *Giardia* (CWP1)
 B: Green/Red \rightarrow Positive *Giardia* (α 1-giardin)
- A: Green/Red \rightarrow Positive *Giardia* (CWP1)
 B: Green \rightarrow Negative *Giardia* (α 1-giardin)
- A: Green \rightarrow Negative *Giardia* (CWP1)
 B: Green/Red \rightarrow Positive *Giardia* (α 1-giardin)

INVALID RESULTS

-

CERTEST *Giardia*²

One Step test to detect *Giardia* (CWP1 and/or α 1-*giardin*) in combo card format/
Prueba combo de un solo paso para detección de *Giardia* (CWP1 y/o α 1-*giardin*) en formato cassette

	A (<i>Giardia</i> (CWP1))	B (<i>Giardia</i> (α 1- <i>giardin</i>))	Interpretation of the results
1.	- GREEN	- GREEN	There is no <i>Giardia</i> (CWP1, α 1- <i>giardin</i>) presence. No infection caused by <i>Giardia</i> .
2.	+ GREEN-RED	+ GREEN-RED	There is <i>Giardia</i> (CWP1) and <i>Giardia</i> (α 1- <i>giardin</i>) presence. Infection caused by <i>Giardia</i> .
3.	+ GREEN-RED	- GREEN	There is presence (CWP1). Infection caused by <i>Giardia</i> (CWP1).
4.	- GREEN	+ GREEN-RED	There is <i>Giardia</i> (α 1- <i>giardin</i>) presence. Infection caused by <i>Giardia</i> (α 1- <i>giardin</i>).
5.	Any other result	Any other result	Invalid result either A or B, we recommend repeating the assay using the same sample with another test.

INVALID: Total absence of any control coloured line (GREEN) regardless the appearance or not of the test lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.
3. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
4. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
5. After one week of infection, the number of parasites in faeces is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms.
6. CerTest *Giardia*² should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
7. Positive results determine the presence of *Giardia* (CWP1 and/or α 1-*giardin*) in faecal samples; nevertheless, a positive result should be followed up with additional laboratory techniques (biochemical methods or microscopy) to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
8. A negative result is not meaningful because of it is possible the antigens concentration in the stool sample is lower than the detection limit value. If the symptoms or situation still persist, a *Giardia* identification should be carried out with another technique (for example microscopy).
9. Bloody stool samples can contain components that may cause non-specific reactions in the test. Every bloody stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
10. Mucous stool samples may cause non-specific reactions in the test. Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.

EXPECTED VALUES

Giardia infection has become recognized as one of the most common causes of waterborne disease (found in both drinking and recreational water) in humans. *Giardia* is found worldwide.

In one study, it was found that between 60 and 80% of infected children in day care nurseries and their household contact was asymptomatic giardiasis. Asymptomatic individuals are an important reservoir for spread of the infection.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest *Giardia*², CerTest) vs the Evaluation Criteria. Evaluation criteria: three lateral flow test were evaluated (CerTest *Giardia*², CerTest and other two rapid test from the competitors); the discrepant results were confirmed by qPCR technique (Viasure *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time Detection Kits (CerTest)). The results were as follows:

		Evaluation Criteria		
		+	-	Total
IC test: CerTest <i>Giardia</i> ² (CWP1)	+	43	2	45
	-	2	79	81
	Total	45	81	126

CerTest <i>Giardia</i> ² (CWP1) vs Evaluation Criteria		
	Mean Value	95% confidence interval
Sensitivity	95.6%	84.9– 99.5%
Specificity	97.5%	91.4 – 99.7%
PPV	95.6%	84.9 - 99.5%
NPV	97.5%	91.4 – 99.7%

		Evaluation Criteria		
		+	-	Total
IC test: CerTest <i>Giardia</i> ² (α 1-giardin)	+	34	1	35
	-	11	80	91
	Total	45	81	126

CerTest <i>Giardia</i> ² (α 1-giardin) vs Evaluation Criteria		
	Mean Value	95% confidence interval
Sensitivity	75.6%	60.5– 87.1%
Specificity	98.8%	93.3 – 100.0%
PPV	97.1%	85.1 - 99.9%
NPV	87.9%	79.4 – 93.8%

The results showed a high sensitivity and specificity to detect *Giardia* (CWP1) and *Giardia* (α 1-giardin) using CerTest *Giardia*².

CERTEST *Giardia*²

One Step test to detect *Giardia* (CWP1 and/or α 1-giardin) in combo card format/
Prueba combo de un solo paso para detección de *Giardia* (CWP1 y/o α 1-giardin) en formato cassette

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *Giardia*²; no cross reactivity against gastrointestinal pathogens, other organism, substances and/or faecal markers occasionally present in faeces:

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