---

# CERTEST

## *Giardia+*

---

ONE STEP  
*Giardia+* CARD TEST

---

CERTEST BIOTEC S.L.

---

## CERTEST *Giardia*+

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### ENGLISH

#### INTENDED USE

CerTest *Giardia*+ one step card test is a coloured chromatographic immunoassay for the qualitative detection of *Giardia* ( $\alpha$ 1-*giardin* and CWP1) in stool samples.

CerTest *Giardia*+ card test offers a simple and a highly sensitive screening assay to make a presumptive diagnosis of giardiasis.

#### INTRODUCTION

*Giardiasis* is a diarrheal illness caused by a very small parasite, *Giardia intestinalis* (also known as *Giardia lamblia* and *Giardia duodenalis*). Once an animal or person is infected with *Giardia*, the parasite lives in the intestine and is passed in the stool. The parasite is protected by an outer shell (cyst) and can survive outside the body and in the environment for a long time.

*Giardia*, a flagellated protozoan, inhabits the upper part of the small intestine of its host and has a two major states in the life cycle: trophozoites which produces the antigens ( $\alpha$ -1 *giardin*) and cyst with produces the antigens (CWP1) After the host ingest the cysts, which are the infective stage, the trophozoites emerge from the cysts in the duodenum and attach to the small intestinal mucosa. They undergo mitotic division in the intracellular lumen, some will encyst to protect themselves and will be eliminated from the host in the feces. The trophozoite is the vegetative form and replicates in the small intestine.

In humans, the clinical effects of *Giardia* infection range from the asymptomatic carrier state to a severe malabsorption syndrome. Possible factors contributing to the variation in clinical manifestations include the virulence of *Giardia* strain, the number of cyst ingested, the age of the host, and the state of the host immune system at the time of infection.

*Giardia* infection has become one of the most common causes of waterborne disease in humans. The most common symptoms of giardiasis include: diarrhea, loose or watery stool, stomach cramps and upset stomach. These symptoms generally begin 1-2 weeks after infection and may last 2-6 weeks in healthy individuals. Sometimes symptoms last longer and may lead to weight loss and dehydration. Some people will have no symptoms. However, people with weakened immune systems (e.g., persons with HIV/AIDS, cancer patients and transplant patients) or the elderly may have a more serious infection that can lead to severe illness or death.

#### TEST PRINCIPLE

CerTest *Giardia*+ is based on the principle of a qualitative immunochromatographic assay for the determination of *Giardia* ( $\alpha$ 1-*giardin* and CWP1) antigens in stool samples. Thus, the test can detect trophozoites and cyst.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Giardia* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Giardia*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*Giardia* monoclonal antibodies-red polystyrene microspheres), which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*Giardia* antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no *Giardia* antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-*Giardia* antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appear. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

#### STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

#### PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.


- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the CerTest *Giardia*+ test. Do not use any other commercial kit component.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

### SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-3 days prior to testing. For longer storage, maximum 1 month, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

#### Specimen preparation (see illustration):

1. Take out the cap of the stool collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2) making sure that, at each insertion, only the stick's screw is covered with the sample to take the appropriate amount of faecal sample and add it to the stool collection tube. Any additional sample that exceeds the stick's screw could cause wrong results. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

**VERY IMPORTANT: Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.**

### MATERIALS

#### MATERIALS PROVIDED

- CerTest *Giardia*+ card tests
- Instructions for use
- Stool collection tubes with diluent

#### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

### TEST PROCEDURE


**Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.**

1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove CerTest *Giardia*+ card test from its sealed bag just before using it.
3. Take the stool collection tube, cut the end of the cap (4) and dispense 3 drops in the circular window marked with the letter S (5). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.


If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

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### INTERPRETATION OF THE RESULTS (please refer to the illustration below)

	<i>Giardia</i> +	Interpretation of the results
1.	-	There is no <i>Giardia</i> presence. No infection caused by <i>Giardia</i> ( $\alpha$ 1- <i>giardin</i> and/or CWP1).
	GREEN	
2.	+	There is <i>Giardia</i> presence. Infection caused by <i>Giardia</i> ( $\alpha$ 1- <i>giardin</i> and/or CWP1).
	GREEN-RED	
3.	Any other result	Invalid result, we recommend repeating the assay using the same sample with another test.

**NEGATIVE:** Only the GREEN line appears across the results window in the site marked with the letter C (control line).

**POSITIVE:** In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

**INVALID:** Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

### NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

### QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

### LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.
3. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
4. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.

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5. After one week of infection, the number of parasites in faeces is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms.
6. CerTest *Giardia*+ should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
7. Positive results determine the presence of *Giardia* ( $\alpha 1$ -*giardin* and/or CWP1) in faecal samples; nevertheless, a positive result should be followed up with additional laboratory techniques (biochemical methods or microscopy) to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
8. A negative result is not meaningful because of it is possible the antigens concentration in the stool sample is lower than the detection limit value. If the symptoms or situation still persist, a *Giardia* determination should be carried out with another technique (for example microscopy).
9. Bloody stool samples can contain components that may cause non-specific reactions in the test. Every bloody stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
10. Mucous stool samples may cause non-specific reactions in the test- Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.

### EXPECTED VALUES

*Giardia* infection has become recognized as one of the most common causes of waterborne disease (found in both drinking and recreational water) in humans. *Giardia* is found worldwide.

In one study, it was found that between 60 and 80% of infected children in day care nurseries and their household contact was asymptomatic giardiasis. Asymptomatic individuals are an important reservoir for spread of the infection.

### PERFORMANCE CHARACTERISTICS

#### Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest *Giardia*+, CerTest) vs the Evaluation Criteria. Evaluation criteria: three rapid test were evaluated (CerTest *Giardia*+, Certest and other two rapid test from the competitors); the discrepant results were confirmed by qPCR technique (*Viasure Cryptosporidium, Giardia & E. histolytica* Real Time Detection Kits (CerTest)). The results were as follows:

		Evaluation Criteria		
		+	-	Total
IC test: CerTest <i>Giardia</i> +	+	44	2	46
	-	1	79	80
	Total	45	81	126

CerTest <i>Giardia</i> + vs Evaluation Criteria		
	Mean Value	95% confidence interval
Sensitivity	97.8%	88.2–99.9%
Specificity	97.5%	91.4–99.7%
PPV	95.7%	85.2–99.5%
NPV	98.8%	93.2–100.0%

The results showed a high sensitivity and specificity to detect *Giardia* using CerTest *Giardia*+ test.

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### Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *Giardia*+; no cross reactivity against gastrointestinal pathogens, other organism, substances and/or faecal markers occasionally present in faeces:

<i>Adenovirus</i>	<i>Cryptosporidium parvum</i>	<i>Legionella pneumophila</i>	<i>Salmonella typhimurium</i>
<i>Astrovirus</i>	<i>Entamoeba histolytica</i>	<i>Listeria monocytogenes</i>	<i>Shigella boydii</i>
<i>Bovine Haemoglobine</i>	<i>Escherichia coli O157</i>	<i>Norovirus GI</i>	<i>Shigella dysenteriae</i>
<i>Calprotectin</i>	<i>Escherichia coli O:111</i>	<i>Norovirus GII</i>	<i>Shigella flexneri</i>
<i>Campylobacter jejuni</i>	<i>Escherichia coli O:026</i>	<i>Pig Haemoglobine</i>	<i>Shigella sonnei</i>
<i>Clostridium difficile antigen GDH</i>	<i>Helicobacter pylori</i>	<i>Rotavirus</i>	<i>Streptococcus pyogenes</i>
<i>Clostridium difficile Toxin A</i>	<i>Human Haemoglobine</i>	<i>Salmonella enteritidis</i>	<i>Streptococcus pneumococcal</i>
<i>Clostridium difficile Toxin B</i>	<i>Human Lactoferrin</i>	<i>Salmonella paratyphi A</i>	<i>Yersinia Enterocolitica O:3</i>
<i>Clostridium perfringens</i>	<i>Human Transferrin</i>	<i>Salmonella typhi</i>	<i>Yersinia Enterocolitica O:9</i>

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