CERTEST

SARS-CoV-2 + Flu A

ONE STEP SARS-CoV-2 and
Influenza type A Ag COMBO CARD
TEST

CERTEST BIOTEC S.L.

FOR INFORMATION PURPOSES ONLY

CERTEST SARS-CoV-2 + Flu A

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ENGLISH

INTENDED USE

CerTest SARS-CoV-2 + Flu A one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection of nucleoprotein antigen of SARS-CoV-2 and *Influenza* type A from nasopharyngeal swab samples from patients suspected of COVID-19 infection and/or *Influenza A* infection.

CerTest SARS-CoV-2 + Flu A combo card test offers a simple and highly sensitive screening assay as aid in the diagnosis of SARS-CoV-2 and/or *Influenza* type A infection.

INTRODUCTION

An unseen pneumonia spread swiftly in Wuhan, China, in December 2019. Scientists sequenced and recognized a novel β -coronavirus. The disease, caused by the Severe Acute Respiratory Syndrome Corona Virus 2 (SARS-CoV-2), is called as Coronavirus disease 2019 (COVID-19).

Coronaviruses (CoVs) are RNA, enveloped viruses, with a large genome (29.9 kb), amongst which β -CoVs and α -CoVs can infect mammals. Viruses belonging to the Coronaviridae family have the most abundant membrane protein (M) among other proteins, i.e. spike glycoprotein (S), nucleocapsid protein (N) and an envelope protein (E). Spike glycoprotein (S) is one of the targets of T cell response in the immune system. The S protein also promotes the binding of the virus envelope to the ACE2 receptor and the entry of the virus into the cell target. ACE2 receptors are present in cells of the arteries, veins, smooth muscles, small intestine, alveoli of the lungs, hair follicles, cardiac myofibroblasts, skin, brain, and kidney, thus SARS-CoV-2 could potentially infect these tissues.

Clinical forms of this disease comprise from mild to very severe symptoms: pneumonia, fever and respiratory symptoms are the most frequent. Other symptoms of the viral infection include a sore throat, aches, pains, and shortness of breath. In some cases, runny nose, nausea, and diarrhea also occur.

The outbreak of newly discovered coronavirus begins mainly via discharge from the nose or droplets of saliva, once an infected individual sneezes or coughs, and the estimated time for incubation is within 2 weeks. Due to the high infectious rate of SARS-CoV-2, detection of asymptomatic positive patients is probably one of the key points to controlling the outbreak.

Influenza is caused by a virus that attacks mainly the upper respiratory tract - the nose, throat and bronchi and rarely also the lungs. The infection usually lasts for about a week. It is characterized by sudden onset of high fever, myalgia, headache and severe malaise, non-productive cough, sore throat, and rhinitis. Most people recover within one to two weeks without requiring any medical treatment. In the very young, the elderly and people suffering from medical conditions such as lung diseases, diabetes, cancer, kidney or heart problems, *Influenza* poses a serious risk. In these people, the infection may lead to severe complications of underlying diseases, pneumonia and death.

The currently circulating influenza viruses that cause human disease are divided into two groups: A and B. *Influenza A* has 3 subtypes which are important for humans: A (H3N2), A (H1N1) and A (H5N1), of which the former is currently associated with most deaths. *Influenza* viruses are defined by 2 different protein components, known as antigens, on the surface of the virus. They are spike-like features called haemagglutinin (H) and neuraminidase (N) components.

TEST PRINCIPLE

CerTest SARS-CoV-2 + Flu A one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection of nucleoprotein antigen of SARS-CoV-2 and *Influenza* type A from nasopharyngeal swab samples from patients suspected of COVID-19 infection and/or *Influenza A* infection.

Strip A consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against SARS-CoV-2 and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-SARS-CoV-2) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

Strip B consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Influenza* type A and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-type A) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is SARS-CoV-2 positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-SARS-CoV-2 monoclonal antibodies-red polystyrene microspheres) in the strip A, and if the sample is *Influenza A* positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*Influenza* type A monoclonal antibodies-red polystyrene microspheres) in the strip B, which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-SARS-CoV-2 antibodies present on the membrane of strip A

(test line) and the anti-*Influenza* type A antibodies present on the membrane of strip B (test line) capture the coloured conjugate and the one red line will be visible in both strips. These bands are used to interpret the result.

If the sample is negative, there is no SARS-CoV-2 and *Influenza* type A antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-SARS-CoV-2 and anti-*Influenza* type A antibodies present on the membranes (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in both strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on both membranes will capture control green-conjugate complex and both control lines will always appear. The presence of these green lines serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors. Please remember that COVID-19 infection is very dangerous.
- The tests should be discarded in a proper biohazard container after testing.
- Clean up spills thoroughly using an appropriate disinfectant.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest SARS-CoV-2 + Flu A test. Do not use any other commercial kit component.
- Sterile swabs provided in the kits should be only used for taking the nasopharyngeal sample collection. They cannot be reuse.
- Do not touch the head of the sterile swabs provided when opening their primary packaging to avoid contamination.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- All positive results should be processed following local laws and regulations.

SPECIMEN COLLECTION AND PREPARATION

Samples should be process as soon as possible after collection. If this is not possible, the samples can be stored in the refrigerator (2-8 °C) for 8 hours prior to testing.

Specimen preparation (see illustration):

Nasopharyngeal swab method:

1. Remove the swab from its packing.
2. Collect specimen with a sterile swab from one nostril.
3. Insert the swab into the nostril to the nasopharynx, rotating against the wall.
(to ensure swab contains cells as well as mucus).
4. Repeat procedure using other nostril.
5. Process the swab as soon as possible after collecting the specimen.


Follow the test procedure.

Samples previously extracted by transport media:

Even while using direct nasopharyngeal swab is the preferred protocol, the device can be also used with samples from nasopharyngeal swabs previously diluted in transport media such as: VTM, UTM or Saline Buffer.

1. Use the minimum volume of transport media in order to avoid loss in sensitivity. Example: 1.0mL. Using higher volumes will affect sensitivity of the system.
2. Dilute the extracted sample 1:1 in Reagent SARS-CoV-2 + Flu A (sample diluent) provided.

Follow immediately with the point 3 of test procedure

Specimen collection

MATERIALS

MATERIALS PROVIDED

- CerTest SARS-CoV-2 + Flu A combo card tests
- Reagent SARS-CoV-2 + Flu A (sample diluent)
- Swabs
- Disposable pipettes
- Testing tubes
- Instructions for use
- Laboratory rack
- CerTest SARS-CoV-2 Positive Control swab + Instructions for use
- CerTest Influenza A Positive Control swab + Instructions for use

MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE

Allow tests, samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.


1. Add 15 drops (approx. 500µL) (1) Reagent SARS-CoV-2 + Flu A and immediately put the swab into the tube.
2. Mix the solution by rotating the swab forcefully against the side of the tube at least 1 minute. Best results are obtained when the specimen is vigorously extracted in the solution (2). Extract as much liquid as possible from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
3. Remove the CerTest SARS-CoV-2 + Flu A combo card test from its sealed bag just before using it.
4. Use a separate pipette and test for each sample or control. Dispense 3 drops (approx. 120µL) from the testing tube, into the circular window marked with the letter A (3) and 3 drops (approx. 120µL), using the same tube, into the circular window marked with the letter B (4).
5. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to the type of sample, stir the sample added in the sample window (S) with the pipette. If it doesn't work, dispense a drop of Reagent SARS-CoV-2+ Flu A until seeing the liquid running through the reaction zone.

SARS-CoV-2 Strip A-procedure

Add 3 drops in A-circular window


(3)

(1)

Add 15 drops of Reagent SARS-CoV-2 + Flu A


(2)

Put the swab into the tube, rotating 1 minute and extract the liquid

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

INTERPRETATION OF THE RESULTS (please refer to the illustration below)

INVALID RESULTS

	A (SARS-CoV-2)	B (Flu A)	Interpretation of the results
1.	- GREEN	- GREEN	There is no SARS-CoV-2, <i>Influenza</i> type A antigens presence in the sample.
2.	+ GREEN-RED	+ GREEN-RED	There is SARS-CoV-2 and <i>Influenza</i> type A antigens presence in the sample.
3.	+ GREEN-RED	- GREEN	There is SARS-CoV-2 antigen presence in the sample. There is no <i>Influenza</i> type A antigens presence in the sample.
4.	- GREEN	+ GREEN-RED	There is <i>Influenza</i> type A presence in the sample. There is no SARS-CoV-2 antigens presence in the sample.
5.	Any other result	Any other result	Invalid result either A or B, we recommend repeating the assay using the same sample with another test.

INVALID: Total absence of any control coloured line (GREEN) regardless the appearance or not of the test lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

RECOMMENDATIONS

Recommendations World Health Organization: *Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays*- 11 September 2020.

1. To optimize performance, testing should be conducted by trained operators in strict accordance with the test procedure and within the first 5-7 days following the onset of symptoms.

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2. Where possible, all positive samples giving positive results should be transported to laboratories with NAAT (Nucleic acid amplification tests) capability for confirmatory testing.
3. This test could be used to screen at-risk individuals and rapidly isolate positive cases in NAAT-confirmed COVID-19 outbreaks.
4. To monitor trends in disease incidence in communities.
5. For early detection and isolation of positive cases in health facilities, where there is widespread community transmission.
6. A negative result cannot completely exclude an active COVID-19 infection, repeat testing or preferably confirmatory testing should be performed (NAAT) whenever possible, particularly in symptomatic patients.
7. Even CerTest SARS-CoV-2 + Flu A test was not validated using samples from asymptomatic contacts of cases, asymptomatic cases have been demonstrated to have viral loads similar to symptomatic cases, so CerTest SARS-CoV-2 + Flu A could detect as positive.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. CerTest SARS-CoV-2 + Flu A should be used only with nasopharyngeal swab samples. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper nasopharyngeal specimens must be obtained.
3. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
4. Positive results determine the presence of SARS-CoV-2 and/or Influenza type A antigens. A diagnosis of SARS-CoV-2 and/or Influenza type A infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
5. Positive results do not rule out co-infections with other pathogens.
6. Only VTM, UTM and Saline buffer transport media has been validated with CerTest SARS-CoV-2 + Flu A test. When using transport media, the test sensitivity can be reduced due to excessive dilution of sample.
7. A negative result is not meaningful because of it is possible the antigens concentration in the nasopharyngeal samples is lower than the detection limit value. If the symptoms or situation still persist, it is recommended that all negative results undergo confirmatory testing using other method and/or virus identification by cell culture, PCR or other technique.

EXPECTED VALUES

On January 30, 2020, the World Health Organization declared COVID-19 a public health emergency of international concern. Since the onset of the pandemic, SARS-CoV-2 has affected millions of people, causing thousands of victims worldwide.

Recent studies have shown that the primary viral infection found in symptomatic patients was identical to that seen in asymptomatic patients, suggesting that COVID-19 patients with asymptomatic infections are capable of transmitting SARS-CoV-2. During the incubation period, asymptomatic patients are also likely to transmit the virus to healthy people.

Although SARS-CoV-2 infection can affect people of any age, children and the elderly with asthma are in a higher risk category for the disease. Those with underlying health problems such as cancer, chronic respiratory diseases, cardiovascular diseases, and diabetes are also reported to have a higher chance of developing serious illnesses.

Yearly Influenza epidemics can seriously affect all age groups, but the highest risk of complications occur among children younger than age two, adults age 65 or older and people of any age with certain medical conditions, such as chronic heart, lung, kidney, liver, blood or metabolic diseases (such as diabetes) or weakened immune systems.

Influenza affects between 20% and 30% of the United States population annually. Three fourths of the infected individuals develop an acute respiratory illness and one third of these seek medical attention. In a typical year, more than 20,000 to 40,000 Americans die and more than 100,000 are hospitalized because of complications related to influenza.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity

Analytical sensitivity (detection limit)

Detection limit value (typical) of CerTest SARS-CoV-2 + Flu A (strip A: SARS-CoV-2) is 1.00 ng/mL of SARS-CoV-2 recombinant protein or 1×10^3 TCID₅₀/mL of 2019nCoV/USA-WA1/2020.

Detection limit value (typical) of CerTest SARS-CoV-2 + Flu A (strip B: Influenza A) is 6.25 ng/mL of Influenza A recombinant nucleoprotein.

Evaluations for SARS-CoV-2 + Flu A (strip A: SARS-CoV-2)

An evaluation, with nasopharyngeal samples from people suspected of infection by SARS-CoV-2 virus, was performed comparing the results obtained by an immunochromatographic test (CerTest SARS-CoV-2 + Flu A, Certest (SARS-CoV-2 strip) vs qPCR technique (VIASURE). The results were as follows:

IC test: CerTest SARS-CoV-2 + Flu A (SARS-CoV-2)	qPCR technique (SARS-CoV-2)			
		+	-	Total
	+	26	1	27
-	2	233	235	
Total	28	234	262	

CerTest SARS-CoV-2 + Flu A (SARS-CoV-2) vs qPCR technique		
	Mean Value	95% confidence interval
Sensitivity (*)	92.9%	76.5–99.1%
Specificity	99.6%	97.6–100.0%
PPV	96.3%	81.0–99.9%
NPV	99.1%	97.0–99.9%

(*) Taking into account the recommendations for *Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays (11 September 2020)* from the WHO, the sensitivity of the test was calculated with nasopharyngeal samples with high viral load (high viral loads is expected in early symptomatic phases of the illness (with the first 5-7 days of illness) in the range of Ag-RDT test detection).

Evaluations for SARS-CoV-2 + Flu A (strip B: Influenza A)

Results showed in tables below have been obtained from several evaluations developed using an immunochromatographic (CerTest SARS-CoV-2 + Flu A, CerTest (Influenza A)) and others commercial test kit.

1574 respiratory samples have been used for *Influenza A* evaluation (130 positive samples, 11 false positives, 1402 negatives samples and 31 false negatives). The results were as follows:

	Evaluation number	True Positives	False Negatives	True Negatives	False Positives
CerTest SARS-CoV-2 + Flu A (Influenza A) vs XPECT FLU A&B (Remel)	I	7	0	85	1
CerTest SARS-CoV-2 + Flu A (Influenza A) vs Instand EQAS program (March 2016)	II	101	26	955	3
CerTest SARS-CoV-2 + Flu A (Influenza A) vs Binax NOW® Influenza A&B Card (Alere)	III	11	2	87	0
CerTest SARS-CoV-2 + Flu A (Influenza A) vs SD Bioline, Influenza Ag A/B/A (H1N1) (SD Bioline)	IV	1	0	97	0
CerTest SARS-CoV-2 + Flu A (Influenza A) vs MariPOC® respi (ArcDia International)	V	10	3	178	7
	Total	130	31	1402	11

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CerTest SARS-CoV-2 + Flu A (Influenza A) vs other test used in evaluation (I, II, III, IV, V) for Influenza A		
	Mean Value	95% confidence interval
Sensitivity	80.7%	73.8 – 86.5%
Specificity	99.2%	98.6 – 99.6%
PPV	92.2%	86.5 – 96.0%
NPV	97.8%	96.9 – 98.5%

The results showed a high sensitivity and specificity to detect SARS-CoV-2 and *Influenza* type A using CerTest SARS-CoV-2 + Flu A.

Hook Effect

CerTest SARS-CoV-2 + Flu A test does not show inhibitory hook effect at concentration higher than 100,000xLoD for Strip A: SARS-CoV-2 and at concentrations higher than 30,000xLoD for Strip B: Influenza A.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest SARS-CoV-2 + Flu A; no cross reactivity against organisms, pathogens, substances that cause infections:

For Strip A. SARS-CoV-2

<i>Adenovirus</i>	<i>Entamoeba histolytica</i>	<i>Legionella pneumophila</i>	<i>Salmonella typhimurium</i>
<i>Astrovirus</i>	<i>Escherichia coli</i> O157	<i>Listeria monocytogenes</i>	<i>Shigella boydii</i>
<i>Bovine Haemoglobin</i>	<i>Giardia</i>	<i>Norovirus GI</i>	<i>Shigella dysenteriae</i>
<i>Campylobacter jejuni</i>	<i>Helicobacter pylori</i>	<i>Norovirus GII</i>	<i>Shigella flexneri</i>
<i>Clostridium difficile</i> antigen GDH	<i>Human Calprotectin</i>	<i>Pig Haemoglobin</i>	<i>Shigella sonnei</i>
<i>Clostridium difficile</i> Toxin A	<i>Human Haemoglobin</i>	<i>Rotavirus</i>	<i>Streptococcus pneumococcal</i>
<i>Clostridium difficile</i> Toxin B	<i>Human Lactoferrin</i>	<i>Respiratory Syncytial Virus</i>	<i>Streptococcus pyogenes</i>
<i>Coronavirus strain</i> 229E	<i>Human Transferrin</i>	<i>Salmonella enteritidis</i>	<i>Yersinia Enterocolitica</i> O:3
<i>Coronavirus strain</i> NL63	<i>Influenza type A</i>	<i>Salmonella paratyphi A</i>	<i>Yersinia Enterocolitica</i> O:9
<i>Coronavirus strain</i> OC43	<i>Influenza type B</i>	<i>Salmonella typhi</i>	

CerTest SARS-CoV-2 + Flu A (SARS-CoV-2, strip A) showed some cross reaction with SARS and very low with MERS.

For Strip B. Influenza A

<i>Adenovirus</i>	<i>Coronavirus strain</i> OC43	<i>Legionella pneumophila</i>	<i>Salmonella typhi</i>
<i>Astrovirus</i>	<i>Entamoeba histolytica</i>	<i>Listeria monocytogenes</i>	<i>Salmonella typhimurium</i>
<i>Bovine Haemoglobin</i>	<i>Escherichia coli</i> O157	<i>MERS</i>	<i>Shigella boydii</i>
<i>Campylobacter jejuni</i>	<i>Giardia</i>	<i>Norovirus GI</i>	<i>Shigella dysenteriae</i>
<i>Clostridium difficile</i> antigen GDH	<i>Helicobacter pylori</i>	<i>Norovirus GII</i>	<i>Shigella flexneri</i>
<i>Clostridium difficile</i> Toxin A	<i>Human Calprotectin</i>	<i>Pig Haemoglobin</i>	<i>Shigella sonnei</i>
<i>Clostridium difficile</i> Toxin B	<i>Human Haemoglobin</i>	<i>Rotavirus</i>	<i>Streptococcus pneumococcal</i>
<i>Coronavirus SARS-1</i> (SARS)	<i>Human Lactoferrin</i>	<i>Respiratory Syncytial Virus</i>	<i>Streptococcus pyogenes</i>
<i>Coronavirus SARS-2</i> (SARS-CoV-2)	<i>Human Transferrin</i>	<i>Salmonella enteritidis</i>	<i>Yersinia Enterocolitica</i> O:3
<i>Coronavirus strain</i> 229E	<i>Influenza type B</i>	<i>Salmonella paratyphi A</i>	<i>Yersinia Enterocolitica</i> O:9
<i>Coronavirus strain</i> NL63			

Interferences

An evaluation was performed to determine the possible interferences of CerTest SARS-CoV-2 + Flu A.; no interferences against the substances tested were detected:

Potencial interfering substance	Concentration	Result	Potencial interfering substance	Concentration	Result
Metronidazole	3.0 mg/mL	Negative	Loperamide hydrochloride (<i>Fortasec</i>)	1.5 · 10 ⁻² mg/mL	Negative
Ampicillin	3.0 mg/mL	Negative	Heparin (<i>Hibor</i>)	350.0 IU/mL	Negative
Oseltamivir	3.0 · 10 ⁻² mg/mL	Negative	Almagato (<i>Almax</i>)	3.0 mg/mL	Negative
Amantadine	0.3 mg/mL	Negative	Fosfamycin (<i>Monurol</i>)	3.0 · 10 ⁻³ mg/mL	Negative
Ribavirin	3.0 mg/mL	Negative	Acetylcysteine (<i>Fluimucil</i>)	3.0 · 10 ⁻² mg/mL	Negative
Codeine (<i>Toseina</i>)	0.2 mg/mL	Negative	Dexketoprofen trometamol (<i>Enantyum</i>)	0.3 mg/mL	Negative
Benzocaine (<i>Angileptol</i>)	3.0 · 10 ⁻² mg/mL	Negative	Levofloxacin	3.0 · 10 ⁻² mg/mL	Negative
Cloperastine (<i>Flutox</i>)	0.3 mg/mL	Negative	Ciprofloxacin	0.3 mg/mL	Negative
Carbocisteine (<i>Iniston mucolítico</i>)	3.0 · 10 ⁻² mg/mL	Negative	Rifampicin (<i>Rifaldin</i>)	0.3 mg/mL	Negative
Loratadine	0.3 mg/mL	Negative	Phenoxymethylpenicillin potassium	3.0 mg/mL	Negative
Dexchlorpheniramine (<i>Polaramine</i>)	0.3 mg/mL	Negative	Ambroxol hydrochloride (<i>Mucosan</i>)	0.3 mg/mL	Negative
Ebastine (<i>Ebastel</i>)	3.0 mg/mL	Negative	Macrogol 3350 (<i>Movicol</i>)	3.0 mg/mL	Negative
Acetyl Salicylic (<i>Adiro</i>)	0.3 mg/mL	Negative	Lysine Carbocysteinate (<i>Pectox</i>)	3.0 · 10 ⁻² mg/mL	Negative
Ibuprofen (<i>Espidifen</i>)	0.3 mg/mL	Negative	Hydroxyzine dihydrochloride	0.3 mg/mL	Negative
Paracetamol (<i>Dolocatil</i>)	5.0 mg/mL	Negative	Lorazepam	3.0 · 10 ⁻³ mg/mL	Negative
Metamizole (<i>Nolotil</i>)	5.0 mg/mL	Negative	Amoxicillin	3.0 mg/mL	Negative
Prednisone	0.3 mg/mL	Negative	Mercaptopurine	0.3 mg/mL	Negative
Omeprazole	2.0 · 10 ⁻³ mg/mL	Negative			

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