CERTEST

SARS-CoV-2

ONE STEP SARS-CoV-2 Ag CARD
TEST

CERTEST BIOTEC S.L.

For information purposes only

CERTEST SARS-CoV-2

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ENGLISH

INTENDED USE

CerTest SARS-CoV-2 one step card test is a coloured chromatographic immunoassay for the qualitative detection of nucleoprotein antigen of SARS-CoV-2 from nasopharyngeal swab samples from patients suspected of COVID-19 infection.

CerTest SARS-CoV-2 card test offers a simple and a highly sensitive screening assay as an aid in the diagnosis of SARS-CoV-2 infection.

INTRODUCTION

An unseen pneumonia spread swiftly in Wuhan, China, in December 2019. Scientists sequenced and recognized a novel β -coronavirus.

The disease, caused by the Severe Acute Respiratory Syndrome Corona Virus 2 (SARS-CoV-2), is called as Coronavirus disease 2019 (COVID-19).

Coronaviruses (CoVs) are RNA, enveloped viruses, with a large genome (29.9 kb), amongst which β -CoVs and α -CoVs can infect mammals. Viruses belonging to the Coronaviridae family have the most abundant membrane protein (M) among other proteins, i.e. spike glycoprotein (S), nucleocapsid protein (N) and an envelope protein (E). Spike glycoprotein (S) is one of the targets of T cell response in the immune system. The S protein also promotes the binding of the virus envelope to the ACE2 receptor and the entry of the virus into the cell target. ACE2 receptors are present in cells of the arteries, veins, smooth muscles, small intestine, alveoli of the lungs, hair follicles, cardiac myofibroblasts, skin, brain, and kidney, thus SARS-CoV-2 could potentially infect these tissues.

Clinical forms of this disease comprise from mild to very severe symptoms: pneumonia, fever and respiratory symptoms are the most frequent. Other symptoms of the viral infection include a sore throat, aches, pains, and shortness of breath. In some cases, runny nose, nausea, and diarrhea also occur.

The outbreak of newly discovered coronavirus begins mainly via discharge from the nose or droplets of saliva, once an infected individual sneezes or coughs, and the estimated time for incubation is within 2 weeks. Due to the high infectious rate of SARS-CoV-2, detection of asymptomatic positive patients is probably one of the key points to controlling the outbreak.

TEST PRINCIPLE

CerTest SARS-CoV-2 is based on the principle of a qualitative immunochromatographic assay for the qualitative detection of SARS-CoV-2 from nasopharyngeal swab samples from patients suspected of COVID-19 infection.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against SARS-CoV-2 and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-SARS-CoV-2) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-SARS-CoV-2 monoclonal antibodies-red polystyrene microspheres), which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-SARS-CoV-2 antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no SARS-CoV-2 antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-SARS-CoV-2 antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appear. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors. Please remember that COVID-19 infection is very dangerous.

- The tests should be discarded in a proper biohazard container after testing.
- Clean up spills thoroughly using an appropriate disinfectant.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest SARS-CoV-2 test. Do not use any other commercial kit component.
- Sterile swabs provided in the kits should be only used for taking the nasopharyngeal sample collection. They cannot be reuse.
- Do not touch the head of the sterile swabs provided when opening their primary packaging to avoid contamination.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.
- All positive results should be processed following local laws and regulations.

SPECIMEN COLLECTION AND PREPARATION

Samples should be process as soon as possible after collection. If this is not possible, the samples can be stored in the refrigerator (2-8 °C) for 8 hours prior to testing.

Specimen preparation (see illustration):

Nasopharyngeal swab method:

1. Remove the swab from its packing.
2. Collect specimen with a sterile swab from one nostril.
3. Insert the swab into the nostril to the nasopharynx, rotating against the wall.
(to ensure swab contains cells as well as mucus).
4. Repeat procedure using other nostril.
5. Process the swab as soon as possible after collecting the specimen.

Follow the test procedure.

Samples previously extracted by transport media:

Even while using direct nasopharyngeal swab is the preferred protocol, the device can be also used with samples from nasopharyngeal swabs previously diluted in transport media such as: VTM, UTM or Saline Buffer.

1. Use the minimum volume of transport media in order to avoid loss in sensitivity. Example: 1.0 mL. Using higher volumes will affect sensitivity of the system.
2. Dilute the extracted sample 1:1 in Reagent SARS-CoV-2 (sample diluent) provided.

Follow immediately with the point 3 of test procedure

Specimen collection

MATERIALS

MATERIALS PROVIDED

- CerTest SARS-CoV-2 card tests
- Reagent SARS-CoV-2 (sample diluent)
- Swabs
- Disposable pipettes
- Testing tubes
- Instructions for use
- Laboratory rack
- CerTest SARS-CoV-2 Positive Control swab +
Instructions for use

MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer

TEST PROCEDURE


Allow tests, samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Add 15 drops (approx. 500µL) (1) Reagent SARS-CoV-2 and immediately put the swab into the tube.
2. Mix the solution by rotating the swab forcefully against the side of the tube at least 1 minute. Best results are obtained when the specimen is vigorously extracted in the solution (2). Extract as much liquid as possible from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
3. Remove the CerTest SARS-CoV-2 card test from its sealed bag just before using it.
4. Use a separate pipette and test for each sample or control. Dispense 3 drops (approx. 120µL) from the testing tube, into the circular window marked with the letter S (3).
5. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to the type of sample, stir the sample added in the sample window (S) with the pipette. If it doesn't work, dispense a drop of Reagent SARS-CoV-2 until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS (please refer to the illustration below)

| | SARS-CoV-2 | Interpretation of the results |
|----|------------------|---|
| 1. | - GREEN | There is no SARS-CoV-2 antigen presence in the sample. |
| 2. | + GREEN-RED | There is SARS-CoV-2 antigen presence in the sample. |
| 3. | Any other result | Invalid result, we recommend repeating the assay using the same sample with another test. |

NEGATIVE: Only one GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

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NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigen presents in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

RECOMMENDATIONS

Recommendations World Health Organization: Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays- 11 september 2020.

1. To optimize performance, testing should be conducted by trained operators in strict accordance with the test procedure and within the first 5-7 days following the onset of symptoms.
2. Where possible, all positive samples giving positive results should be transported to laboratories with NAAT capability for confirmatory testing.
3. This test could be used to screen at-risk individuals and rapidly isolate positive cases in NAAT- confirmed COVID-19 outbreaks.
4. To monitor trends in disease incidence in communities.
5. For early detection and isolation of positive cases in health facilities, where there is widespread community transmission.
6. A negative result cannot completely exclude an active COVID-19 infection, repeat testing or preferably confirmatory testing should be performed (NAAT) whenever possible, particularly in symptomatic patients.
7. Even CerTest SARS-CoV-2 test was not validated using samples from asymptomatic contacts of cases, asymptomatic cases have been demonstrated to have viral loads similar to symptomatic cases, so CerTest SARS-CoV-2 could detect as positive.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. CerTest SARS-CoV-2 should be used only with nasopharyngeal swab samples. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper nasopharyngeal specimens must be obtained.
3. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
4. Positive results determine the presence of SARS-CoV-2 antigens. A diagnosis of SARS-CoV-2 infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
5. Positive results do not rule out co-infections with other pathogens.
6. Only VTM, UTM and Saline buffer transport media has been validated with CerTest SARS-CoV-2 test. When using transport media, the test sensitivity can be reduced due to excessive dilution of sample.
7. A negative result is not meaningful because of it is possible the antigens concentration in the nasopharyngeal samples is lower than the detection limit value. If the symptoms or situation still persist, it is recommended that all negative results undergo confirmatory testing using other method and/or virus identification by cell culture, PCR or other technique.

EXPECTED VALUES

On January 30, 2020, the World Health Organization declared COVID-19 a public health emergency of international concern. Since the onset of the pandemic, SARS-CoV-2 has affected millions of people, causing thousands of victims worldwide.

Recent studies have shown that the primary viral infection found in symptomatic patients was identical to that seen in asymptomatic patients, suggesting that COVID-19 patients with asymptomatic infections are capable of transmitting SARS-CoV-2. During the incubation period, asymptomatic patients are also likely to transmit the virus to healthy people.

Although SARS-CoV-2 infection can affect people of any age, children and the elderly with asthma are in a higher risk category for the disease. Those with underlying health problems such as cancer, chronic respiratory diseases, cardiovascular diseases, and diabetes are also reported to have a higher chance of developing serious illnesses.

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PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit)

Detection limit value (typical) of CerTest SARS-CoV-2 is 1.0 ng/mL of SARS-CoV-2 recombinant protein or 1×10^3 TCID₅₀/mL of 2019nCoV/USA-WA1/2020.

Clinical sensitivity and specificity

An evaluation, with nasopharyngeal samples from people suspected of infection by SARS-CoV-2 virus, was performed comparing the results obtained by an immunochromatographic test (CerTest SARS-CoV-2, Certest) vs qPCR technique (VIASURE). The results were as follows:

| IC test: CerTest SARS-CoV-2. | qPCR technique | | | Total |
|---------------------------------|----------------|----|-----|-------|
| | | + | - | |
| + | | 26 | 1 | 27 |
| - | | 2 | 233 | 235 |
| Total | | 28 | 234 | 262 |

| CerTest SARS-CoV-2 vs qPCR technique | | |
|--------------------------------------|------------|-------------------------|
| | Mean Value | 95% confidence interval |
| Sensitivity (*) | 92.9% | 76.5–99.1% |
| Specificity | 99.6% | 97.6–100.0% |
| PPV | 96.3% | 81.0–99.9% |
| NPV | 99.1% | 97.0–99.9% |

(*) Taking into account the recommendations for *Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays* (11 September 2020) from the WHO, the sensitivity of the test was calculated with nasopharyngeal samples with high viral load (high viral loads is expected in early symptomatic phases of the illness (with the first 5-7 days of illness) in the range of Ag-RDT test detection).

The results showed a high sensitivity and specificity to detect SARS-CoV-2 using CerTest SARS-CoV-2.

Hook Effect

CerTest SARS-CoV-2 test does not show inhibitory hook effect at concentration higher than 100,000xLoD.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest SARS-CoV-2; no cross reactivity against organism, pathogens, substances that could cause infections:

| | | | |
|--|------------------------------|------------------------------------|------------------------------------|
| <i>Adenovirus</i> | <i>Entamoeba histolytica</i> | <i>Legionella pneumophila</i> | <i>Salmonella typhimurium</i> |
| <i>Astrovirus</i> | <i>Escherichia coli O157</i> | <i>Listeria monocytogenes</i> | <i>Shigella boydii</i> |
| <i>Bovine Haemoglobin</i> | <i>Giardia</i> | <i>Norovirus GI</i> | <i>Shigella dysenteriae</i> |
| <i>Campylobacter jejuni</i> | <i>Helicobacter pylori</i> | <i>Norovirus GII</i> | <i>Shigella flexneri</i> |
| <i>Clostridium difficile antigen GDH</i> | <i>Human Calprotectin</i> | <i>Pig Haemoglobin</i> | <i>Shigella sonnei</i> |
| <i>Clostridium difficile Toxin A</i> | <i>Human Haemoglobin</i> | <i>Rotavirus</i> | <i>Streptococcus pneumococcal</i> |
| <i>Clostridium difficile Toxin B</i> | <i>Human Lactoferrin</i> | <i>Respiratory Syncytial Virus</i> | <i>Streptococcus pyogenes</i> |
| <i>Coronavirus strain 229E</i> | <i>Human Transferrin</i> | <i>Salmonella enteritidis</i> | <i>Yersinia Enterocolitica O:3</i> |
| <i>Coronavirus strain NL63</i> | <i>Influenza type A</i> | <i>Salmonella paratyphi A</i> | <i>Yersinia Enterocolitica O:9</i> |
| <i>Coronavirus strain OC43</i> | <i>Influenza type B</i> | <i>Salmonella typhi</i> | |

CerTest SARS-CoV-2 showed some cross reaction with SARS and very low with MERS.

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Interferences

An evaluation was performed to determine the possible interferences of CerTest SARS-CoV-2; no interferences against the substances tested were detected:

| Potencial interfering substance | Concentration | Result | Potencial interfering substance | Concentration | Result |
|---|---------------------------|----------|--|---------------------------|----------|
| Metronidazole | 3.0 mg/mL | Negative | Loperamide hydrochloride (<i>Fortasec</i>) | 0.15 mg/mL | Negative |
| Ampicillin | 3.0 mg/mL | Negative | Heparin (<i>Hibor</i>) | 350.0 IU/mL | Negative |
| Oseltamivir | $3.0 \cdot 10^{-2}$ mg/mL | Negative | Almagato (<i>Almax</i>) | 3.0 mg/mL | Negative |
| Amantadine | 0.3 mg/mL | Negative | Fosfamycin (<i>Monurol</i>) | $3.0 \cdot 10^{-3}$ mg/mL | Negative |
| Ribavirin | 3.0 mg/mL | Negative | Acetylcysteine (<i>Fluimucil</i>) | $3.0 \cdot 10^{-2}$ mg/mL | Negative |
| Codeine (<i>Toseina</i>) | 0.2 mg/mL | Negative | Dexketoprofen trometamol (<i>Enantyum</i>) | 0.3 mg/mL | Negative |
| Benzocaine (<i>Angileptol</i>) | $3.0 \cdot 10^{-2}$ mg/mL | Negative | Levofloxacin | $3.0 \cdot 10^{-2}$ mg/mL | Negative |
| Cloperastine (<i>Flutox</i>) | 0.3 mg/mL | Negative | Ciprofloxacin | 0.3 mg/mL | Negative |
| Carbocisteine (<i>Iniston mucolítico</i>) | $3.0 \cdot 10^{-2}$ mg/mL | Negative | Rifampicin (<i>Rifaldin</i>) | 0.3 mg/mL | Negative |
| Loratadine | 0.3 mg/mL | Negative | Phenoxymethylpenicillin potassium | 3.0 mg/mL | Negative |
| Dexchlorpheniramine (<i>Polaramine</i>) | 0.3 mg/mL | Negative | Ambroxol hydrochloride (<i>Mucosan</i>) | 0.3 mg/mL | Negative |
| Ebastine (<i>Ebastel</i>) | 3.0 mg/mL | Negative | Macrogol 3350 (<i>Movicol</i>) | 3.0 mg/mL | Negative |
| Acetyl Salicylic (<i>Adiro</i>) | 0.3 mg/mL | Negative | Lysine Carbocysteinate (<i>Pectox</i>) | $3.0 \cdot 10^{-2}$ mg/mL | Negative |
| Ibuprofen (<i>Espidifen</i>) | 0.3 mg/mL | Negative | Hydroxyzine dihydrochloride | 0.3 mg/mL | Negative |
| Paracetamol (<i>Dolocatil</i>) | 5.0 mg/mL | Negative | Lorazepam | $3.0 \cdot 10^{-3}$ mg/mL | Negative |
| Metamizole (<i>Nolotil</i>) | 5.0 mg/mL | Negative | Amoxicillin | 3.0 mg/mL | Negative |
| Prednisone | 0.3 mg/mL | Negative | Mercaptopurine | 0.3 mg/mL | Negative |
| Omeprazole | $2.0 \cdot 10^{-3}$ mg/mL | Negative | | | |

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