# VIASURE

Real Time PCR Detection Kits

by CerTest  
BIOTEC

## **SARS-CoV-2, Flu (A+B) & RSV**

Handbook for the following references/

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## ENGLISH

### 1. Intended use

VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System is an automated real-time RT-PCR test designed for the qualitative detection and differentiation of RNA from the SARS-CoV-2, Influenza A (Flu A), Influenza B (Flu B) and/or Human Respiratory Syncytial Virus A/B (RSV) in respiratory samples from individuals suspected of COVID-19 or other respiratory infection by their healthcare provider. This test is intended to be used as an aid in the identification of the presence of the SARS-CoV-2, Flu A, Flu B and/or RSV viral RNA. The assay uses the BD MAX™ System for automated extraction of RNA and subsequent real-time RT-PCR employing the reagents provided combined with universal reagents and disposables for the BD MAX™ System. RNA is extracted from respiratory specimens, amplified using RT-PCR and detected using fluorescent reporter dye probes specific for SARS-CoV-2, Flu A, Flu B and/or RSV.

### 2. Summary and Explanation

Coronavirus are enveloped non-segmented positive-sense RNA viruses and belong to *Coronaviridae* family. There are six coronavirus species known to cause human diseases. Four viruses (229E, OC43, NL63 and HKU1) cause common cold symptoms and the other two (severe acute respiratory syndrome coronavirus (SARS-CoV) and Middle East respiratory syndrome coronavirus (MERS-CoV)) are zoonotic and producing more severe complications. SARS-CoV and MERS-CoV have caused more than 10,000 cumulative cases in the past two decades, with mortality rates of 34% MERS-CoV and 10% SARS-CoV.

In December 2019, some people that worked at or lived around the Huanan seafood market in Wuhan, Hubei Province, China, have presented pneumonia of unknown cause. Deep sequencing analysis of the respiratory samples indicated a novel coronavirus, which was named firstly 2019 novel coronavirus (2019-nCoV) and lately SARS-CoV-2.

Human-to-human transmission of the SARS-CoV-2 has been confirmed, even in the incubation period without symptoms, and the virus causes severe respiratory illness like those SARS-CoV produced. Although the pneumonia is the principal illness associated, a few patients have developed severe pneumonia, pulmonary edema, acute respiratory distress syndrome, or multiple organ failure and death. Centers of Disease Control and Prevention (CDC) believes that symptoms of SARS-CoV-2 may appear in as few as 2 days or as long as 14 days after exposure, being the most common fever or chills, cough, fatigue, anorexia, myalgia and dyspnea. Less common symptoms are sore throat, nasal congestion, headache, diarrhea, nausea and vomiting. Loss of smell (anosmia) or loss of taste (ageusia) preceding the onset of respiratory symptoms has also been reported. Older adults and people who have severe underlying medical conditions like heart or lung disease or diabetes seem to be at higher risk for developing more serious complications from COVID-19 illness.

CDC recommends upper respiratory tract specimens (nasopharyngeal (NP) and oropharyngeal (OP) swabs, nasal mid-turbinate swab, nasal swab, nasopharyngeal wash/aspirate or nasal wash/aspirate (NW) specimens collected mainly by a healthcare provider) and/or lower respiratory specimens (sputum, endotracheal aspirate, or

bronchoalveolar lavage in patients with more severe respiratory disease) for the identification of SARS-CoV-2 and other respiratory viruses, such as Influenza and RSV.

Influenza viruses belong to the *Orthomyxoviridae* family and cause the majority of viral lower respiratory tract infections. Influenza A and B are a significant cause of morbidity and mortality worldwide, considering that elderly and compromised individuals are especially at risk of developing severe illness and complications such as pneumonia. People feel some or all of these symptoms: fever or feeling feverish/chills, cough, sore throat, nasal stuffiness and discharge, myalgia, headaches, and anorexia. The influenza viruses can be spread from person to person in two different ways: through the air (large droplets and aerosols from sneezing and coughing), and by direct or indirect contact.

Influenza A and B are an enveloped, single stranded RNA viruses that containing eight segmented strands of genome RNA, which typically encodes 11 or 12 viral proteins. The viral envelope, derived from the host plasma membrane, consists of a lipid bilayer containing transmembrane proteins, like hemagglutinin (HA) and neuraminidase (NA), and matrix proteins M1 and M2. Influenza A viruses are further classified into subtypes based on the antigenicity of their "HA" and "NA" molecules, whereas Influenza B is divided into 2 antigenically and genetically distinct lineages, Victoria and Yamagata.

Human respiratory syncytial viruses A and B (RSV) belong to the *Paramyxoviridae* family and are the most important viral agents of acute respiratory infections. RSV is an enveloped, nonsegmented, negative, single stranded linear RNA genome virus. Respiratory syncytial virus is a common contributor of respiratory infections causing bronchitis, pneumonia, and chronic obstructive pulmonary infections in people of all ages. People often feel some or all of these symptoms: rhinorrhea, low-grade fever, cough, sore throat, headache, and wheezing. RSV is transmitted via large nasopharyngeal secretion droplets from infected individuals, close contact, or self-inoculation after touching contaminated surfaces.

Diagnosis can be problematic, as a wide range of pathogens can cause acute respiratory infections presenting with similar clinical syndromes. Real-time PCR assays have been shown to be a sensitive and specific diagnostic tool for the detection of SARS-CoV-2, Flu A, Flu B and RSV viruses.

### 3. Principle of the procedure

VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System is designed for the identification of SARS-CoV-2, Flu A, Flu B and /or RSV in respiratory samples. The detection is done in one step real-time RT-PCR format where the reverse transcription and the subsequent amplification of specific target sequence occur in the same reaction tube. The isolated RNA target is transcribed generating complementary DNA by reverse transcriptase which is followed by the amplification of two conserved regions of N gene (N1 and N2) for SARS-CoV-2, a conserved region of the M1 gene for Flu A and Flu B, and a conserved region of the N gene for RSV using specific primers and fluorescent-labeled probes.

VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the

fluorescent signal which is proportional to the quantity of the target template. This fluorescence is measured on the BD MAX™ System.

VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System is composed of two different reaction tubes. One of the tubes detects and differentiates the RNA from Flu A, Flu B and/or RSV (Transparent Red or 1A foil) and the other tube detects specifically the RNA from SARS-CoV-2 (Transparent Green or 1G foil). Each tube contains all the components necessary for real-time PCR assay (specific primers/probes, dNTPS, buffer, polymerase, reverse-transcriptase) in a stabilized format, as well as an internal control (endogenous in the SARS-CoV-2 reaction tube) to monitor the extraction process and/or inhibition of the polymerase activity. The SARS-CoV-2 assay uses a human housekeeping gene as an endogenous Internal Control (human *RNase P* gene). Human housekeeping genes are involved in basic cell maintenance and, therefore, are expected to be present in all nucleated human cells and maintain relatively constant expression levels. Each RNA targets are amplified and detected in specific channels (475/520, 585/630, and/or 630/665) and the internal control (IC) in channel 530/565. In the Flu A, Flu B and/or RSV assay, Flu A RNA target is amplified and detected in channel 475/520, Influenza B RNA target in channel 585/630, RSV RNA target in channel 630/665 and the internal control (IC) of this assay in channel 530/565. In SARS-CoV-2 assay, N2 target is amplified and detected in channel 475/520, N1 target in channel 630/665 and the endogenous internal control (IC) in channel 530/565.

## 4. Reagents provided

VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System includes the following materials and reagents detailed in Table 1:

| Reference         | Reagent/Material                   | Description                                                                                                       | Color/Barcode                 | Amount                |
|-------------------|------------------------------------|-------------------------------------------------------------------------------------------------------------------|-------------------------------|-----------------------|
| <b>VS-ABR212R</b> | Flu A, Flu B & RSV reaction tube   | A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and internal control in stabilized format            | Transparent Red or 1A foil    | 2 pouches of 12 tubes |
| <b>VS-NCO312</b>  | SARS-CoV-2 (N1 + N2) reaction tube | A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and endogenous internal control in stabilized format | Transparent Green or 1G foil  | 2 pouches of 12 tubes |
| <b>VS-RB09</b>    | Rehydration Buffer tube            | Solution to reconstitute the stabilized product                                                                   | Transparent Orange or 1I foil | 1 pouch of 24 tubes   |

Table 1. Reagents and materials provided in VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System with Cat. N°: VS-FNR124 (444217).

## 5. Reagents and equipment to be supplied by the user

The following list includes the materials and equipment that are required for use but not included in the VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System.

- Real-time PCR instrument: BD MAX™ System.
- BD MAX™ ExK™ TNA-3 (Ref:442827 or 442828)
- BD MAX™ PCR Cartridges (Ref: 437519)

- Vortex.
- Micropipettes (accurate between 2 and 1000 µL).
- Filter tips.
- Powder-free disposable gloves

## 6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- After opening the aluminum pouches which contain the reaction tubes can be used up to 28 days.

## 7. Precautions for users

- The product is indented for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- For *in vitro* diagnostic use.
- Do not use expired reagents and/or materials.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective box is open or broken upon arrival.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches.
- Close protective pouches of reagents promptly with the zip seal after each use. Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and/or kits and/or lots.
- Protect reagents from humidity. Prolonged exposure to humidity may affect product performance.
- Keep components away from light.
- In cases where other PCR tests are conducted in the same general area of the laboratory, care must be taken to ensure that the VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System, BD MAX™ ExK™ TNA-3 extraction kit, any additional reagents required for testing, and the BD MAX™ System are not contaminated. Avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents at all times. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents and cartridges.
- Make sure to use a tube to determine RNA from Influenza A, Influenza B and RSV in Snap-In 2 (green position) and another tube to determine RNA from SARS-CoV-2 in Snap-In 4 (blue position). Be careful not to mix them throughout the entire process.
- To avoid contamination of the environment by amplicons, do not break apart the BD MAX™ PCR Cartridge after use. The seals of the BD MAX™ PCR Cartridge are designed to prevent contamination.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.

- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Wash your hands after finishing the test.
- Samples must be treated as potentially infectious as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult the BD MAX™ System User's Manual for additional warnings, precautions and procedures.

## 8. Procedure

### 8.1. SAMPLE COLLECTION, STORAGE AND TRANSPORT

The VIASURE SARS-CoV-2, *Flu (A+B)* & RSV Real Time PCR Detection Kit for BD MAX™ System has been validated on nasopharyngeal/ oropharyngeal swab collected in viral transport media (VTM) (Vircell S.L., Spain).

Other types of samples from nasopharyngeal/oropharyngeal swabs in VTM must be validated by the user.

Collection, storage and transport specimens should be maintained per the conditions validated by the user. Overall, respiratory samples should be collected and labelled appropriately in clean containers with or without transport media (depending on sample type) and processed as soon as possible to guarantee the quality of the test. The specimens should be transported at 2 to 8°C for up to 48 hours, following the local and national regulations for the transport of pathogen material. For long term transport (more than 48 hours), we recommend shipping at ≤ -20°C. It is recommended to use fresh specimens for the test. The samples can be stored at 2 to 8°C for up to 48 hours or frozen at -20°C or ideally at -70°C for conservation. Repeated freeze-thaw cycles should be avoided in order to prevent degradation of the sample and nucleic acids.

### 8.2. SAMPLE PREPARATION AND RNA EXTRACTION

Perform the sample preparation according to the recommendations in the instructions for use of extraction kit used, BD MAX™ ExK™ TNA-3. Note that some other samples may require pre-processing. Application-specific extraction preparation procedures should be developed and validated by the user.

1. Pipette 400 µL of nasopharyngeal/ oropharyngeal swab collected in viral transport media (VTM) into a BD MAX™ TNA-3 Sample Buffer Tube and close the tube with a septum cap. Ensure complete mixing by vortexing the sample at high speed for 1 minute. Proceed to BD MAX™ System Operation.

Note: The *Flu A*, *Flu B* & RSV reaction tube has been validated with a sample volume of 200-400 µL and the SARS-CoV-2 (*N1 + N2*) reaction tube with a sample volume of 400-750 µL.

### 8.3. PCR PROTOCOL

Note: Please, refer to the BD MAX™ System User's Manual for detailed instructions.

### 8.3.1. Creating PCR test program for VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System

Note: If you have already created the VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection test, you can skip step 8.3.1 and go directly to 8.3.2.

- 1) On the "Run" screen of the BD MAX™ System, select the "Test Editor" tab.
- 2) Click the "Create" button.
- 3) In the Basic Information tab, within the "Test Name" window, name your test: i.e. VIASURE SARS-CoV-2, Flu (A+B) & RSV (VSARSCoV2,FluA+B,RSV).
- 4) In the "Extraction Type" drop down menu, select "ExK TNA-3".
- 5) In the "Master Mix Format" drop down menu, choose "Dual Master Mix Concentrated Lyophilized MM with Rehydration Buffer (Type 5)".
- 6) In the "Sample extraction parameters" select "User defined" and adjust sample volume to 950 µL.
- 7) In the "Ct Calculation" select "Call Ct at Threshold Crossing".
- 8) If running software version 5.00 or higher and have barcoded foil snap-in tubes, in the "Custom Barcodes" select the following configuration:
  - a. Snap-In 2 Barcode: 1A (concerning Flu A, Flu B & RSV reaction tube)
  - b. Snap-In 3 Barcode: 11 (concerning Rehydration Buffer tube)
  - c. Snap-In 4 Barcode: 1G (concerning SARS-CoV-2 (N1 + N2) reaction tube)
- 9) "PCR Settings" and "Test Steps" must be completed for Snap-In 2 (green) and Snap-In 4 (blue) positions.
- 10) Snap-In 2 (green). In "PCR settings" tab enter the following parameters: "Channel Settings", "Gains" and "Threshold" (Table 2).

| Channel         | Alias | Gain | Threshold | Ct Min | Ct Max |
|-----------------|-------|------|-----------|--------|--------|
| 475/520 (FAM)   | Flu A | 60   | 100       | 0      | 40     |
| 530/565 (HEX)   | IC    | 80   | 300       | 0      | 40     |
| 585/630 (ROX)   | Flu B | 60   | 200       | 0      | 40     |
| 630/665 (Cy5)   | RSV   | 60   | 150       | 0      | 40     |
| 680/715 (Cy5.5) | -     | 0    | 0         | 0      | 0      |

Table 2. PCR settings.

Note: It is recommended to set the minimum threshold values listed above for each channel as a starting point, but the final settings must be determined by the end-user during the result interpretation in order to ensure that thresholds fall within the exponential phase of the fluorescence curves and above any background signal. The threshold value for different instruments may vary due to different signal intensities.

- 11) Snap-In 2 (green). In "PCR settings" tab enter the following parameters "Spectral Cross Talk" (Table 3), as well.

|                    | Channel | False Receiving Channel |         |         |         |         |
|--------------------|---------|-------------------------|---------|---------|---------|---------|
|                    |         | 475/520                 | 530/565 | 585/630 | 630/665 | 680/715 |
| Excitation Channel | 475/520 | -                       | 0.0     | 0.0     | 0.0     | 0.0     |
|                    | 530/565 | 0.0                     | -       | 2.0     | 0.0     | 0.0     |
|                    | 585/630 | 0.0                     | 0.0     | -       | 0.0     | 0.0     |
|                    | 630/665 | 0.0                     | 0.0     | 4.0     | -       | 0.0     |
|                    | 680/715 | 0.0                     | 0.0     | 0.0     | 0.0     | -       |

Table 3. Spectral cross-talk parameters.

12) Snap-In 2 (green). In "Test Steps" tab, enter the PCR protocol (Table 4).

| Step Name                                              | Profile Type  | Cycles | Time (s) | Temperature | Detect |
|--------------------------------------------------------|---------------|--------|----------|-------------|--------|
| Reverse transcription                                  | Hold          | 1      | 900      | 45°C        | -      |
| Initial denaturation                                   | Hold          | 1      | 120      | 98°C        | -      |
| Denaturation and Annealing/Extension (Data collection) | 2-Temperature | 45     | 10       | 95°C        | -      |
|                                                        |               |        | 61.1     | 63°C        | ✓      |

Table 4. PCR protocol.

13) Snap-In 4 (blue). In "PCR settings" tab enter the following parameters: "Channel Settings", "Gains" and "Threshold" (Table 5).

| Channel         | Alias                | Gain | Threshold | Ct Min | Ct Max |
|-----------------|----------------------|------|-----------|--------|--------|
| 475/520 (FAM)   | SARS-CoV-2 N2 target | 80   | 150       | 0      | 40     |
| 530/565 (HEX)   | Endogenous IC        | 80   | 150       | 0      | 35     |
| 585/630 (ROX)   | -                    | 0    | 0         | 0      | 0      |
| 630/665 (Cy5)   | SARS-CoV-2 N1 target | 80   | 150       | 0      | 40     |
| 680/715 (Cy5.5) | -                    | 0    | 0         | 0      | 0      |

Table 5. PCR settings.

Note: It is recommended to set the minimum threshold values listed above for each channel as a starting point, but the final settings must be determined by the end-user during the result interpretation in order to ensure that thresholds fall within the exponential phase of the fluorescence curves and above any background signal. The threshold value for different instruments may vary due to different signal intensities.

14) Snap-In 4 (blue). In "PCR settings" tab enter the following parameters "Spectral Cross Talk" (Table 6), as well.

|                    | Channel | False Receiving Channel |         |         |         |         |
|--------------------|---------|-------------------------|---------|---------|---------|---------|
|                    |         | 475/520                 | 530/565 | 585/630 | 630/665 | 680/715 |
| Excitation Channel | 475/520 | -                       | 3.0     | 0.0     | 0.0     | 0.0     |
|                    | 530/565 | 1.0                     | -       | 0.0     | 0.0     | 0.0     |
|                    | 585/630 | 0.0                     | 0.0     | -       | 0.0     | 0.0     |
|                    | 630/665 | 0.0                     | 0.0     | 0.0     | -       | 0.0     |
|                    | 680/715 | 0.0                     | 0.0     | 0.0     | 0.0     | -       |

Table 6. Spectral cross-talk parameters.

15) Snap-In 4 (blue). In "Test Steps" tab, enter the PCR protocol (Table 7).

| Step Name                                              | Profile Type  | Cycles | Time (s) | Temperature | Detect |
|--------------------------------------------------------|---------------|--------|----------|-------------|--------|
| Reverse transcription                                  | Hold          | 1      | 900      | 45°C        | -      |
| Initial denaturation                                   | Hold          | 1      | 120      | 98°C        | -      |
| Denaturation and Annealing/Extension (Data collection) | 2-Temperature | 45     | 10       | 95°C        | -      |
|                                                        |               |        | 61.1     | 63°C        | ✓      |

Table 7. PCR protocol.

16) Click the "Save Test" button.

### 8.3.2. BD MAX™ Rack set up


- For each sample to be tested, remove one Unitized Reagent Strips from the BD MAX™ ExK TNA-3 kit. Gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes and load on the BD MAX™ System sample racks.
- Remove the required number of BD MAX™ ExK™ TNA Extraction Tubes (B4) (white foil) from their protective pouch. Snap the Extraction Tube(s) (white foil) into its corresponding positions in the TNA strip (Snap position 1, white color coding on the rack. See Figure 1). Remove excess air, and close pouch with the zip seal.
- Determine and separate the appropriate number of *Flu A*, *Flu B* & *RSV* reaction tubes (red or 1A foil) and snap into their corresponding positions in the strip (Snap position 2, green color coding on the rack. See Figure 1).
  - Remove excess air, and close aluminum pouches with the zip seal.
  - In order to carry out a correct rehydration, please make sure that the lyophilized product is in the bottom of the tube and is not adhered to the top area of the tube or to the foil seal. Gently tap each tube on a hard surface to make sure all the product is at the bottom of the tube.
- Remove the required number of Rehydration Buffer tubes (orange or 11 foil) and snap into their corresponding positions in the strip (Snap position 3, non-color coding on the rack. See Figure 1). Remove excess air, and close the pouch with the zip seal.
  - In order to ensure a correct transfer, please make sure that the liquid is in the bottom of the tube and is not adhered to the top area of the tube or to the foil seal. Gently tap each tube on a hard surface to make sure all the product is at the bottom of the tube.

- 5) Determine and separate the appropriate number of SARS-CoV-2 (N1 + N2) reaction tubes (green or 1G foil) and snap into their corresponding positions in the strip (Snap position 4, blue color coding on the rack. See Figure 1).

  - a. Remove excess air, and close aluminum pouches with the zip seal.
  - b. In order to carry out a correct rehydration, please make sure that the lyophilized product is in the bottom of the tube and is not adhered to the top area of the tube or to the foil seal. Gently tap each tube on a hard surface to make sure all the product is at the bottom of the tube.

Figure 1. BD MAX™ TNA Reagent Strip (TNA) from the BD MAX™ ExK TNA-3 kit.

### 8.3.3. BD MAX™ Instrument set up

- 1) Select the "Work List" tab on the "Run" screen of the BD MAX™ System software v4.50A or higher.
- 2) In the "Test" drop down menu, select VSARSCoV2,FluA+B,RSV (if not already created see Section 8.3.1).
- 3) Select the appropriate kit lot number (found on the outer box of extraction kit used) from the pull down menu (optional).
- 4) Enter the Sample Buffer Tube identification number into the Sample tube window of the Worklist, either by scanning the barcode with the scanner or by manual entry.
- 5) Fill the Specimen/Patient ID and/or Accession window of the Worklist and click the "Save" button. Continue until all Sample Buffer Tubes are entered. Ensure that the specimen/patient ID and the Sample Buffer Tubes are accurately matched.
- 6) Place the prepared Sample Buffer Tube into the BD MAX™ Rack(s).
- 7) Load the rack(s) into the BD MAX™ System (Rack A is positioned on the left side of the BD MAX™ System and Rack B on the right side).
- 8) Place the required number of BD MAX™ PCR Cartridge(s) into the BD MAX™ System.
- 9) Close the BD MAX™ System door.
- 10) Click "Start Run" to begin the procedure.

### 8.3.4 BD MAX™ report

- 1) In main menu, click the "Results" button.
- 2) Either double click on your run in the list or press the "view button".
- 3) Click on "Print", select: "Run Details, Test Details and Plot..."
- 4) Click on "Print or Export button" on the "Run Reports" screen.

## 9. Result interpretation

For a detailed description on how to analyze data, refer to the BD MAX™ System User's manual.

The analysis of the data is done by the BD MAX™ software according to the manufacturer's instructions. The BD MAX™ software reports Ct values and amplification curves for each detector channel of each sample tested in the following way:

- Ct value of 0 indicates that there was no Ct value calculated by the software with the specified Threshold (see Table 2). Amplification curve of the sample showing a "0" Ct value must be checked manually.
- Ct value of -1 indicates that no amplification process has occurred.
- Any other Ct value should be interpreted in correlation with the amplification curve and according to the sample interpretation guidelines outlined in Tables 8 and 9.

Check Internal Control signal to verify the correct functioning of the amplification mix. In addition, check that there is no report of BD MAX™ System failure.

Results should be read and analyzed using the following tables:

a. Flu A, Flu B & RSV reaction tube: Snap-In 2

| Flu A<br>(475/520) | Flu B<br>(585/630) | RSV<br>(630/665) | Internal<br>control<br>(530/565) | Interpretation                                                                                                                                                                                                          |
|--------------------|--------------------|------------------|----------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| +                  | +                  | +                | +/-1                             | Flu A, Flu B and RSV RNA Detected <sup>1</sup>                                                                                                                                                                          |
| +                  | -                  | -                | +/-1                             | Flu A RNA Detected, Flu B and RSV RNA Not Detected <sup>1</sup>                                                                                                                                                         |
| +                  | +                  | -                | +/-1                             | Flu A and Flu B RNA Detected, and RSV RNA Not Detected <sup>1</sup>                                                                                                                                                     |
| +                  | -                  | +                | +/-1                             | Flu A and RSV RNA Detected, and Flu B RNA Not Detected <sup>1</sup>                                                                                                                                                     |
| -                  | +                  | -                | +/-1                             | Flu B RNA Detected, Flu A and RSV RNA Not Detected <sup>1</sup>                                                                                                                                                         |
| -                  | +                  | +                | +/-1                             | Flu B and RSV RNA Detected, Flu A RNA Not Detected <sup>1</sup>                                                                                                                                                         |
| -                  | -                  | +                | +/-1                             | RSV RNA Detected, Flu A and Flu B RNA Not Detected <sup>1</sup>                                                                                                                                                         |
| -                  | -                  | -                | + <sup>2</sup>                   | Flu A, Flu B and RSV RNA Not Detected <sup>2</sup>                                                                                                                                                                      |
| -                  | -                  | -                | - <sup>2</sup>                   | Unresolved (UNR) Result obtained in the presence of inhibitors in the PCR reaction or when a general problem (not reported by an error code) with the sample processing and/or amplification steps occurs. <sup>2</sup> |
| IND                | IND                | IND              | IND                              | Indeterminate assay result (IND). Due to BD MAX™ System failure. Assay result displayed in case of an instrument failure linked to an error code.                                                                       |
| INC                | INC                | INC              | INC                              | Incomplete assay result (INC). Due to BD MAX™ System failure. Assay result displayed in case of failure to complete run.                                                                                                |

Table 8. Sample interpretation Flu A, Flu B & RSV reaction tube

+: Amplification occurred  
-: No amplification occurred

**1** A sample is considered positive if the Ct value obtained is less than 40. The internal control may or may not show an amplification signal, because a high copy number of target can cause preferential amplification of target-specific nucleic acids instead of the internal control. In these cases, the detection of the IC is not necessary.

**2** A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive (Ct less than 40). An inhibition of the PCR reaction can be excluded by the amplification of internal control. In case of unresolved results (UNR), absence of internal control signal in negative sample it is recommended to repeat the assay.

b. SARS-CoV-2 (N1 + N2) reaction tube: Snap-In 4

| SARS-CoV-2 (N2 target) (475/520) | Endogenous Internal Control (530/565) | SARS-CoV-2 (N1 target) (630/665) | Interpretation                                                                                                                                                                                                          |
|----------------------------------|---------------------------------------|----------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| +                                | +/- <sup>3</sup>                      | +                                | SARS-CoV-2 N gene RNA Detected <sup>3</sup>                                                                                                                                                                             |
| + <sup>4</sup>                   | +/- <sup>3</sup>                      | -                                | SARS-CoV-2 N gene RNA Detected <sup>3,4</sup>                                                                                                                                                                           |
| -                                | +/- <sup>3</sup>                      | + <sup>4</sup>                   | SARS-CoV-2 N gene RNA Detected <sup>3,4</sup>                                                                                                                                                                           |
| -                                | + <sup>5</sup>                        | -                                | SARS-CoV-2 N gene RNA Not Detected <sup>5</sup>                                                                                                                                                                         |
| -                                | - <sup>5</sup>                        | -                                | Unresolved (UNR) Result obtained in the presence of inhibitors in the PCR reaction or when a general problem (not reported by an error code) with the sample processing and/or amplification steps occurs. <sup>5</sup> |
| IND                              | IND                                   | IND                              | Indeterminate assay result (IND). Due to BD MAX™ System failure. Assay result displayed in case of an instrument failure linked to an error code.                                                                       |
| INC                              | INC                                   | INC                              | Incomplete assay result (INC). Due to BD MAX™ System failure. Assay result displayed in case of failure to complete run.                                                                                                |

Table 9. Sample interpretation SARS-CoV-2 (N1 + N2) reaction tube

+: Amplification occurred

-: No amplification occurred

**3** A sample is considered positive if the Ct value obtained is less than 40. The endogenous Internal Control (IC) may or may not show an amplification signal. Sometimes, the IC detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids.

**4** If only one target site of the N gene amplifies, verify the sigmoid shape of the curve and the intensity of fluorescence. In case of a doubtful interpretation, depending on the available material, it is also recommended to:

- a) re-extract and re-test another aliquot of the same specimen (if possible, increase sample volume to 750 µl) or,
- b) obtain a new specimen and re-test.

**5** In the case of SARS-CoV-2 target sites negative, IC must show an amplification signal with Ct less than 35. The Ct value could be very variable due to the Endogenous Internal Control is a human housekeeping gene that should be present

in all human nucleated cells in the original sample. If there is an absence of signal or Ct value  $\geq 35$  of the endogenous Internal Control, the result is considered as 'Unresolved', and retesting is required.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user; to verify the correct performance of each RT-qPCR steps and review the parameters; and to check the sigmoid shape of the curve and the intensity of fluorescence.

The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

## 10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated with nasopharyngeal/oropharyngeal swab collected in VTM.
- For good test performance, the lyophilized product should be at the bottom of the tube and not adhered to the top area of the tube or the foil seal. Gently tap each tube on a hard surface to make sure all the product is at the bottom of the tube.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or color different from whitish) does not alter the functionality of the test.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from respiratory samples must be extracted.
- This test is a qualitative test and does not provide quantitative values or indicate the number of organisms present.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by SARS-CoV-2, Flu A, Flu B and/or RSV either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.
- The specific primer and probe combinations for detection of conserved regions of *N gene* (SARS-CoV-2) used in VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System have been designed based on the US CDC assay for specific detection of SARS-CoV-2 by amplifying two unique regions of the *N gene*. They do not show significant combined homologies with the human genome, human microflora, SARS-CoV or other coronaviruses, which might result in predictable false positive.
- False Negative results may arise from several factors and their combinations, including:
  - Improper specimens' collection, transport, storage, and/or handling methods.
  - Improper processing procedures (including RNA extraction).
  - Degradation of the viral RNA during sample shipping/storage and/or processing.
  - Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown SARS-CoV-2, Flu and/or RSV variants.
  - A viral load in the specimen below the limit of detection for the assay.

- The presence of RT-qPCR inhibitors or other types of interfering substances.
  - Failure to follow instructions for use and the assay procedure.
- In SARS-CoV-2 (N1 + N2) reaction tube, a single-target site amplification or even random positive results is suggestive of slightly different amplification yield of the target site of the N gene. Samples with low viral load might result in N single target amplification. In case of a doubt, it is recommended referring to a reference laboratory for further testing.
- Some samples (in SARS-CoV-2 (N1 + N2) reaction tube) may fail to exhibit RNase P amplification curves due to low human cell numbers in the original clinical sample. A negative IC signal does not preclude the presence of SARS-CoV-2, Flu and/or RSV RNA in a clinical specimen.
- A positive test result does not necessarily indicate the presence of viable viruses and does not imply that these viruses are infectious or are the causative agents for clinical symptoms. However, a positive result is indicative of the presence of targets viral sequences.
- Negative results do not preclude SARS-CoV-2, Flu and/or RSV infection and should not be used as the sole basis for treatment or other patient management decisions. Optimum specimen types and timing for peak viral levels during infections caused by SARS-CoV-2 and novel Influenza A strain have not been determined. The collection of multiple specimens (types and time points) from the same patient may be necessary to detect the virus.
- If diagnostic tests for other respiratory illnesses are negative and the patient's clinical presentation and epidemiological information suggest that SARS-CoV-2, Flu and/or RSV infection is possible, then a false negative result should be considered, and a re-testing of the patient should be discussed.
- In the case of obtaining Unresolved, Indeterminate or Incomplete results using VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System retesting will be required. Unresolved results may be due to the presence of inhibitors in the sample or an incorrect rehydration of lyophilized reaction mix tube. If there is an instrument failure, Indeterminate or Incomplete results will be obtained.

## 11. Quality control

VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System contains an internal control in each Flu A, Flu B & RSV reaction tube and an endogenous internal control in each SARS-CoV-2 (N1 + N2) reaction tube which confirms the correct performance of the technique.

## 12. Performance characteristics

### 12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System was tested individually in each reaction tube.

The clinical performance of Flu A, Flu B & RSV reaction tube was tested using 344 respiratory specimens (oropharyngeal swabs) from symptomatic patients. These results were compared with those obtained with a molecular detection method (cobas® Influenza A/B & RSV (Roche)).

The results were as follows:

| Flu A, B & RSV reaction tube | cobas® Influenza A/B & RSV (Roche) |     |     |       |
|------------------------------|------------------------------------|-----|-----|-------|
|                              |                                    | +   | -   | Total |
|                              | +                                  | 157 | 2*  | 159   |
|                              | -                                  | 7*  | 178 | 185   |
|                              | Total                              | 164 | 180 | 344   |

Table 10. Comparative results for Flu A.

**Positive percent agreement is >96% and negative percent agreement is >99%.**

\*The low amount of template RNA in this respiratory sample is below the detection limit of the method used.

| Flu A, Flu B & RSV reaction tube | cobas® Influenza A/B & RSV (Roche) |     |     |       |
|----------------------------------|------------------------------------|-----|-----|-------|
|                                  |                                    | +   | -   | Total |
|                                  | +                                  | 99  | 4*  | 103   |
|                                  | -                                  | 1*  | 240 | 241   |
|                                  | Total                              | 100 | 244 | 344   |

Table 11. Comparative results for Flu B.

**Positive percent agreement is >99% and negative percent agreement is >98%.**

\*The low amount of template RNA in this respiratory sample is below the detection limit of the method used.

| Flu A, Flu B & RSV reaction tube | cobas® Influenza A/B & RSV (Roche) |    |     |       |
|----------------------------------|------------------------------------|----|-----|-------|
|                                  |                                    | +  | -   | Total |
|                                  | +                                  | 22 | 4*  | 26    |
|                                  | -                                  | 3* | 315 | 318   |
|                                  | Total                              | 25 | 319 | 344   |

Table 12. Comparative results for RSV.

**Positive percent agreement is >88% and negative percent agreement is >99%.**

\*The low amount of template RNA in this respiratory sample is below the detection limit of the method used.

The clinical performance of SARS-CoV-2 (N1 + N2) reaction tube was tested using 254 respiratory samples (nasopharyngeal swabs in Viracell Transport medium) from patients with clinical suspicion of COVID-19 disease or other similar respiratory diseases. The results were compared with those obtained with the clinical diagnosis performed with Simplexa™ COVID-19 Direct assay with discrepant analysis performed with the Charité protocol.

| SARS-CoV-2 (N1 + N2) reaction tube | Alternative RT-PCR assays |    |     |       |
|------------------------------------|---------------------------|----|-----|-------|
|                                    |                           | +  | -   | Total |
|                                    | +                         | 63 | 2*  | 65    |
|                                    | -                         | 0  | 189 | 189   |
|                                    | Total                     | 63 | 191 | 254   |

Table 13. Comparative results for SARS-CoV-2.

\*Initial diagnose of one of the two samples was invalid and reported to the patient as positive for prevention and quarantine period.

SARS-CoV-2 (N1 + N2) reaction tube detected two positive samples that were not detected using Simplexa™ COVID-19 Direct assay and the Charité protocol.

The Positive Percent Agreement (PPA) and the Negative Percent Agreement (NPA) for SARS-CoV-2 (N1 + N2) reaction tube are >99% and 98%, respectively.

Results show high agreement to detect SARS-CoV-2, Flu A, Flu B and/or RSV viruses using VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System.

## 12.2. Analytical sensitivity

VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System has a detection limit of  $\geq 10$  genome copies per reaction for Flu A,  $\geq 20$  genome copies per reaction for Flu B,  $\geq 2$  genome copies per reaction for RSV and  $\geq 5$  genome copies per reaction for SARS-CoV-2 with a positive rate of  $\geq 95\%$  (Figures 2, 3, 4, 5 and 6).


Figure 2. Dilution series of Flu A ( $2 \times 10^6$ - $2 \times 10^1$  copies per reaction) template run on the BD MAX™ System (475/520 (FAM) channel).

Figure 3. Dilution series of Flu B ( $2 \times 10^6$ - $2 \times 10^1$  copies per reaction) template run on the BD MAX™ System (585/630 (ROX) channel).


Figure 4. Dilution series of RSV ( $2 \times 10^6$ - $2 \times 10^1$  copies per reaction) template run on the BD MAX™ System (630/665 (Cy5) channel).


Figure 5. Dilution series of SARS-CoV-2 ( $N1 + N2$ ) ( $9.9 \times 10^4$ - $9.9 \times 10^0$  and  $5.0 \times 10^0$  genome copies per reaction) template run on the BD MAX™ System (475/520 (FAM) channel).



Figure 6. Dilution series of SARS-CoV-2 (N1 + N2) ( $9.9 \times 10^4$ - $9.9 \times 10^0$  and  $5.0 \times 10^0$  genome copies per reaction) template run on the BD MAX™ System (630/665 (Cy5) channel).

### 12.3. Analytical specificity

The specificity of the SARS-CoV-2, Flu (A+B) & RSV assay was confirmed by testing a panel consisting of different microorganisms representing the most common respiratory pathogens. No cross-reactivity was detected between any of the following microorganisms tested, except the targeted pathogens of each assay:

| Cross-reactivity testing                         |   |                                                              |         |                                                                  |     |
|--------------------------------------------------|---|--------------------------------------------------------------|---------|------------------------------------------------------------------|-----|
| Human Adenovirus types 1-5, 8, 15, 31, 40 and 41 | - | Influenza A/Netherlands/398/2014 (H3N2) virus (clade 3C.3a)  | -<br>/+ | Influenza A/chicken/Hong Kong/G9/1997 x PR8-IBCDC-2 (H9N2) virus | -/+ |
| Bocavirus                                        | - | Influenza A/Netherlands/2393/2015 (H3N2) virus (clade 3C.2a) | -<br>/+ | Influenza A/Chicken/Myanmar/433/2016 (H9N2) virus                | -/+ |
| <i>Bordetella bronchiseptica</i>                 | - | Influenza A/Newcastle/607/2019 (H3N2) virus                  | -<br>/+ | Influenza A/Hong Kong/1073/99 (H9N2) virus                       | -/+ |
| <i>Bordetella holmesii</i>                       | - | Influenza A/New York/39/2012 (H3N2) virus                    | -<br>/+ | Influenza A/Hong Kong/33982/2009 (H9N2) x PR8-IDCDC-RG26 virus   | -/+ |
| <i>Bordetella parapertussis</i>                  | - | Influenza A/Ohio/2/2012 (H3N2) virus                         | -<br>/+ | Influenza B/Brisbane/60/2008 virus                               | -/+ |
| <i>Bordetella pertussis</i>                      | - | Influenza A/Perth/1001/2018 (H3N2) virus                     | -<br>/+ | Influenza B/Colorado/6/2017 virus                                | -/+ |
| <i>Chlamydia caviae</i>                          | - | Influenza A/Singapore/INFIMH-16-0019/2016 (H3N2) virus       | -<br>/+ | Influenza B/Malaysia/2506/2004 virus                             | -/+ |
| <i>Chlamydia psittaci</i> genotype A and C       | - | Influenza A/South Australia/55/2014 (H3N2) virus             | -<br>/+ | Influenza B/Maryland/15/2016 virus                               | -/+ |
| <i>Chlamydophila pneumoniae</i> CM-1             | - | Influenza A/South Australia/55/2014, IVR-175 (H3N2) virus    | -<br>/+ | Influenza B/Netherlands/207/06 virus                             | -/+ |
| Human coronavirus 229E, OC43, NL63 and HKU1      | - | Influenza A/Switzerland/9715293/2013 (H3N2) virus            | -<br>/+ | Influenza B/Netherlands/2518/2016 (clade 1A) virus               | -/+ |

| Cross-reactivity testing                                            |     |                                                                                |         |                                                               |     |
|---------------------------------------------------------------------|-----|--------------------------------------------------------------------------------|---------|---------------------------------------------------------------|-----|
| MERS Coronavirus                                                    | -   | Influenza A/Texas/50/2012 (H3N2) virus                                         | -<br>/+ | Influenza B/Nevada/3/2011 virus                               | -/+ |
| SARS Coronavirus<br>Strain Frankfurt 1                              | -   | Influenza A/Thüringen/5/2017 (H3N2)<br>virus (Clade 3C2a.1)                    | -<br>/+ | Influenza B/New Jersey/1/2012 virus                           | -/+ |
| SARS-CoV-2 strain<br>BetaCoV/Germany/BavPat1/2020<br>p.1            | -/+ | Influenza A/Uruguay/716/2007<br>(H3N2)(NYMC X-175C) virus                      | -<br>/+ | Influenza B/Texas/02/2013 virus                               | -/+ |
| SARS-CoV-2 strain 2019-nCoV/Italy-<br>INMI1                         | -/+ | Influenza A/Victoria/210/2009(H3N2)<br>virus                                   | -<br>/+ | Influenza B/Townsville/8/2016 virus                           | -/+ |
| SARS-CoV-2 isolate<br>Australia/VIC01/2020                          | -/+ | Influenza A/Victoria/361/2011 (H3N2)<br>virus                                  | -<br>/+ | Influenza B/Canberra/11/2016 virus                            | -/+ |
| SARS-CoV-2 isolate Wuhan-Hu-1                                       | -/+ | Influenza A/Victoria/361/2011 IVR-165<br>(H3N2) virus                          | -<br>/+ | Influenza B/Florida/4/2006 virus                              | -/+ |
| SARS-CoV-2 strain<br>2019nCoV/USAWA1/2020                           | -/+ | Influenza A/Anhui/01/2005 (H5N1) virus                                         | -<br>/+ | Influenza B/Florida/07/2004 virus                             | -/+ |
| Enterovirus 68 and 71                                               | -   | Influenza A/Anhui/01/2005 x PR8-IBCDC-<br>RG6 (H5N1) virus                     | -<br>/+ | Influenza B/Guangdong/120/2000<br>virus                       | -/+ |
| Enterovirus Echovirus 11 and 30                                     | -   | Influenza A/chicken/Vietnam/NCVD-<br>016/2008 (H5N1) virus                     | -<br>/+ | Influenza B/Hubei<br>Wujiagang/158/2009 (NYMC BX-39)<br>virus | -/+ |
| Enterovirus Coxsackievirus A24, A9<br>and B3                        | -   | Influenza A/chicken/Vietnam/NCVD-<br>016/2008 x PR8-IDCDC-RG12 (H5N1)<br>virus | -<br>/+ | Influenza B/ Jiangsu/10/2003 virus                            | -/+ |
| <i>Haemophilus influenzae</i> MinnA                                 | -   | Influenza A/chicken/Vietnam/NCVD-<br>03/08 (H5N1) - PR8-IDCDC-RG25a virus      | -<br>/+ | Influenza B/Massachusetts/2/2012<br>virus                     | -/+ |
| Influenza A/Brisbane/02/2018, IVR-<br>190 (H1N1)pdm09 virus         | -/+ | Influenza A/chicken/Yunnan/1251/2003<br>(H5N1) virus                           | -<br>/+ | Influenza B/Netherlands/365/2016<br>(clade 3) virus           | -/+ |
| Influenza<br>A/California/7/2009(H1N1)pdm09<br>virus                | -/+ | Influenza A/common magpie/Hong<br>Kong/645/2006 (H5N1) virus                   | -<br>/+ | Influenza B/Phuket/3073/2013 virus                            | -/+ |
| Influenza A/Dominican<br>Republic/7293/2013 (H1N1)pdm09<br>virus    | -/+ | Influenza A/duck/Hunan/795/2002<br>(H5N1) virus                                | -<br>/+ | Influenza B/Texas/06/2011 virus                               | -/+ |
| Influenza A/Massachusetts/15/2013<br>(H1N1)pdm09 virus              | -/+ | Influenza A/Egypt/321/2007 (H5N1) virus                                        | -<br>/+ | Influenza B/Wisconsin/1/2010 virus                            | -/+ |
| Influenza A/Michigan/45/2015<br>(H1N1)pdm09 virus                   | -/+ | Influenza A/Egypt/321/2007 x PR8-<br>IDCDC-RG11 (H5N1) virus                   | -<br>/+ | Influenza B/Wisconsin/1/2010 BX-41A<br>virus                  | -/+ |
| Influenza A/Netherlands/1250/2016<br>(H1N1)pdm09 virus (clade 6B.1) | -/+ | Influenza A/Egypt/3300-NAMRU3/2008 x<br>PR8-IDCDC-RG13 (H5N1) virus            | -<br>/+ | <i>Legionella bozemanii</i>                                   | -   |
| Influenza A/New<br>Caledonia/20/99(H1N1) virus                      | -/+ | Influenza A/Egypt/N03072/2010 (H5N1) x<br>PR8-IDCDC-RG29 virus                 | -<br>/+ | <i>Legionella dumoffii</i>                                    | -   |
| Influenza A/New York/18/2009<br>(H1N1)pdm09 virus                   | -/+ | Influenza A/Hong Kong/213/2003<br>(H5N1) virus                                 | -<br>/+ | <i>Legionella longbeachae</i>                                 | -   |

| Cross-reactivity testing                                     |     |                                                                    |         |                                                                |     |
|--------------------------------------------------------------|-----|--------------------------------------------------------------------|---------|----------------------------------------------------------------|-----|
| Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus | -/+ | Influenza A/Hubei/1/2010 (H5N1) x PR8-IBCDRCG30 virus              | -<br>/+ | <i>Legionella micdadei</i>                                     | -   |
| Influenza A/Sydney/134/2018 (H1N1)pdm09 virus                | -/+ | Influenza A/India/NIV/2006 xPR8-IBCDC-RG7 (H5N1) virus             | -<br>/+ | <i>Legionella pneumophila</i>                                  | -   |
| Influenza A/Victoria/2040/2018 (H1N1)pdm09 virus             | -/+ | Influenza A/Japanese white eye/Hong Kong/1038/2006 (H5N1) virus    | -<br>/+ | Human metapneumovirus A and B                                  | -   |
| Influenza A/PR/8/34 (H1N1) virus                             | -/+ | Influenza A/Vietnam/1194/2004 (H5N1) virus                         | -<br>/+ | <i>Moraxella catarrhalis</i>                                   | -   |
| Influenza A/Brisbane/117/2018 (H3N2) virus                   | -/+ | Influenza A/Vietnam/1194/2004 (NIBRG-14) (H5N1) virus              | -<br>/+ | <i>Mycoplasma pneumoniae</i>                                   | -   |
| Influenza A/Brisbane/1028/2017 (H3N2) virus                  | -/+ | Influenza A/Vietnam/1203/2004 x PR8-IBCDC-RG (H5N1) virus          | -<br>/+ | <i>Mycobacterium tuberculosis</i> not rifampin resistant       | -   |
| Influenza A/Fujian/411/2002 (H3N2) virus                     | -/+ | Influenza A/Whooper Swan/R65/2006 (H5N1) virus                     | -<br>/+ | Human parainfluenza 1, 2, 3 and 4 viruses                      | -   |
| Influenza A/Hiroshima//52/2005 (IVR-142) (H3N2) virus        | -/+ | Influenza A/pheasant/New Jersey/1355/1998 (H5N2)-PR8-IBCDC-4 virus | -<br>/+ | <i>Pneumocystis jirovecii</i> Type A1 and g885652              | -   |
| Influenza A/Hong Kong/4801/2014 (H3N2) virus                 | -/+ | Influenza A/Duck/Singapore-Q/F119-3/97 (H5N3) virus                | -<br>/+ | Human rhinovirus type C                                        | -   |
| Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus    | -/+ | Influenza A/Duck/Lao/XBY004/2014 (H5N6) (Clade 2.3.4.4) virus      | -<br>/+ | <i>Staphylococcus aureus</i> subsp. <i>aureus</i>              | -   |
| Influenza A/Indiana/8/2011 (H3N2)v virus                     | -/+ | Influenza A/DE-SH/Reiherente/AR8444/2013 (H5N8) virus              | -<br>/+ | <i>Staphylococcus epidermidis</i>                              | -   |
| Influenza A/Indiana/10/2011 (H3N2)v virus                    | -/+ | Influenza A/Turkey/Germany/R2485-86/2014 (H5N8) virus              | -<br>/+ | <i>Streptococcus pneumoniae</i> Z022                           | -   |
| Influenza A/Kansas/14/2017 (H3N2) virus                      | -/+ | Influenza A/turkey/Virginia/2002 x PR8-IBCDC-5 (H7N2) virus        | -<br>/+ | <i>Streptococcus pyogenes</i>                                  | -   |
| Influenza A/Kansas/14/2017, NYMC X-327 (H3N2) virus          | -/+ | Influenza A/Mallard/Netherlands/2/2009 (H7N7) virus                | -<br>/+ | <i>Streptococcus salivarius</i>                                | -   |
| Influenza A/Kumamoto/102/2002 (H3N2) virus                   | -/+ | Influenza A/Mallard/Netherlands/12/2000 (H7N7) - IBCDC-1 virus     | -<br>/+ | Respiratory syncytial virus (RSV) A and B (strain CH93(18)-18) | -/+ |
| Influenza A/Minnesota/11/2010 (H3N2)v virus                  | -/+ | Influenza A/Anhui/1/2013 (H7N9) virus                              | -<br>/+ | Human Respiratory Syncytial Virus strain Long                  | -/+ |
| Influenza A/Minnesota/11/2010 X203 (H3N2)v virus             | -/+ | Influenza A/Guangdong/17SF003/2016 (H7N9) virus                    | -<br>/+ |                                                                |     |

Table 14. Reference pathogenic microorganisms used in this study.

## 12.4. Analytical reactivity

The reactivity of the VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System for **SARS-CoV-2** was evaluated against RNA from Human 2019-nCoV strain BetaCoV/Germany/BavPat1/2020 p.1, Human 2019-nCoV strain 2019-nCoV/Italy-INMI1, SARS-CoV-2 strain 2019nCoV/USA-WA1/2020, synthetic RNA

controls for two variants of the SARS-CoV-2 virus: MT007544.1 (SARS-CoV2 isolate Australia/VIC01/2020) and MN908947.3 (SARS-CoV-2 isolate Wuhan-Hu-1), showing positive result.

The reactivity of the VIASURE SARS-CoV-2, Flu (A+B) & RSV Real Time PCR Detection Kit for BD MAX™ System for **Influenza A** was evaluated against RNA extracted from the following strains: Influenza A/Brisbane/02/2018, IVR-190 (H1N1)pdm09 virus, Influenza A/California/7/2009(H1N1)pdm09 virus, Influenza A/Dominican Republic/7293/2013 (H1N1)pdm09 virus, Influenza A/Massachusetts/15/2013 (H1N1)pdm09 virus, Influenza A/Michigan/45/2015 (H1N1)pdm09 virus, Influenza A/Netherlands/1250/2016 (H1N1)pdm09 virus (clade 6B.1), Influenza A/New Caledonia/20/99(H1N1) virus, Influenza A/New York/18/2009 (H1N1)pdm09 virus, Influenza A/Singapore/GP1908/2015 virus, IVR-180 (H1N1)pdm09 virus, Influenza A/Sydney/134/2018 (H1N1)pdm09 virus, Influenza A/Victoria/2040/2018 (H1N1)pdm09 virus, Influenza A/PR/8/34 (H1N1) virus, Influenza A/Brisbane/117/2018 (H3N2) virus, Influenza A/Brisbane/1028/2017 (H3N2) virus, Influenza A/Fujian/411/2002 (H3N2) virus, Influenza A/Hiroshima//52/2005 (IVR-142) (H3N2) virus, Influenza A/Hong Kong/4801/2014 (H3N2) virus, Influenza A/Hong Kong/4801/2014 NYMC X-263B (H3N2) virus, Influenza A/Indiana/8/2011 (H3N2)v virus, Influenza A/Indiana/10/2011 (H3N2)v virus, Influenza A/Kansas/14/2017 (H3N2) virus, Influenza A/Kansas/14/2017, NYMC X-327 (H3N2) virus, Influenza A/Kumamoto/102/2002 (H3N2) virus, Influenza A/Minnesota/11/2010 (H3N2)v virus, Influenza A/Minnesota/11/2010 X203 (H3N2)v virus, Influenza A/Netherlands/398/2014 (H3N2) virus (clade 3C.3a), Influenza A/Netherlands/2393/2015 (H3N2) virus (clade 3C.2a), Influenza A/Newcastle/607/2019 (H3N2) virus, Influenza A/New York/39/2012 (H3N2) virus, Influenza A/Ohio/2/2012 (H3N2) virus, Influenza A/Perth/1001/2018 (H3N2) virus, Influenza A/Singapore/INFIMH-16-0019/2016 (H3N2) virus, Influenza A/South Australia/55/2014 (H3N2) virus, Influenza A/South Australia/55/2014, IVR-175 (H3N2) virus, Influenza A/Switzerland/9715293/2013 (H3N2) virus, Influenza A/Texas/50/2012 (H3N2) virus, Influenza A/Thüringen/5/2017 (H3N2) virus (Clade 3C2a.1), Influenza A/Uruguay/716/2007 (H3N2)(NYMC X-175C) virus, Influenza A/Victoria/210/2009(H3N2) virus, Influenza A/Victoria/361/2011 (H3N2) virus, Influenza A/Victoria/361/2011 IVR-165 (H3N2) virus, Influenza A/Anhui/01/2005 (H5N1) virus, Influenza A/Anhui/01/2005 x PR8-IBCDC-RG6 (H5N1) virus, Influenza A/chicken/Vietnam/NCVD-016/2008 (H5N1) virus, Influenza A/chicken/Vietnam/NCVD-016/2008 x PR8-IDCDC-RG12 (H5N1) virus, Influenza A/chicken/Vietnam/NCVD-03/08 (H5N1) - PR8-IDCDC-RG25a virus, Influenza A/chicken/Yunnan/1251/2003 (H5N1) virus, Influenza A/common magpie/Hong Kong/645/2006 (H5N1) virus, Influenza A/duck/Hunan/795/2002 (H5N1) virus, Influenza A/Egypt/321/2007 (H5N1) virus, Influenza A/Egypt/321/2007 x PR8-IDCDC-RG11 (H5N1) virus, Influenza A/Egypt/3300-NAMRU3/2008 x PR8-IDCDC-RG13 (H5N1) virus, Influenza A/Egypt/N03072/2010 (H5N1) x PR8-IDCDC-RG29 virus, Influenza A/Hong Kong/213/2003 (H5N1) virus, Influenza A/Hubei/1/2010 (H5N1) x PR8-IDCDCRG30 virus, Influenza A/India/NIV/2006 xPR8-IBCDC-RG7 (H5N1) virus, Influenza A/Japanese white eye/Hong Kong/1038/2006 (H5N1) virus, Influenza A/Vietnam/1194/2004 (H5N1) virus, Influenza A/Vietnam/1194/2004 (NIBRG-14) (H5N1) virus, Influenza A/Vietnam/1203/2004 x PR8-IBCDC-RG (H5N1) virus, Influenza A/Whooper Swan/R65/2006 (H5N1) virus, Influenza A/pheasant/New Jersey/1355/1998 (H5N2)-PR8-IBCDC-4 virus, Influenza A/Duck/Singapore-Q/F119-3/97 (H5N3) virus, Influenza A/Duck/Lao/XBY004/2014 (H5N6) virus (Clade 2.3.4.4), Influenza A/DE-SH/Reiherente/AR8444/2016 (H5N8) virus, Influenza A/Turkey/Germany/R2485-86/2014 (H5N8) virus, Influenza A/turkey/Virginia/2002 x PR8-IBCDC-5 (H7N2) virus, Influenza A/Mallard/Netherlands/2/2009 (H7N7) virus, Influenza A/Mallard/Netherlands/12/2000 (H7N7) - IBCDC-1 virus, Influenza A/Anhui/1/2013 (H7N9) virus, Influenza A/Guangdong/17SF003/2016 (H7N9) virus, Influenza A/Chicken/Hong Kong/G9/1997 x PR8-IBCDC-2 (H9N2) virus, Influenza A/Chicken/Myanmar/433/2016 (H9N2) virus, Influenza A/Hong Kong/1073/99 (H9N2) virus, Influenza A/Hong Kong/33982/2009 (H9N2) x PR8-IDCDC-RG26 virus, showing positive result.

The reactivity of the VIASURE SARS-CoV-2, *Flu (A+B)* & RSV Real Time PCR Detection Kit for BD MAX™ System for **Influenza B** was evaluated against RNA extracted from the following strains: Influenza B/Brisbane/60/2008 virus, Influenza B/Colorado/6/2017 virus, Influenza B/Malaysia/2506/2004 virus, Influenza B/Maryland/15/2016 virus, Influenza B/Netherlands/207/06 virus, Influenza B/Netherlands/2518/2016 (clade 1A) virus, Influenza B/Nevada/3/2011 virus, Influenza B/New Jersey/1/2012 virus, Influenza B/Texas/02/2013 virus, Influenza B/Townsville/8/2016 virus (**B/Victoria lineage**); Influenza B/Canberra/11/2016 virus, Influenza B/Florida/4/2006 virus, Influenza B/Florida/07/2004 virus, Influenza B/Guangdong/120/2000 virus, Influenza B/Hubei Wujiagang/158/2009 (NYMC BX-39) virus, Influenza B/Jiangsu/10/2003 virus, Influenza B/Massachusetts/2/2012 virus, Influenza B/Netherlands/365/2016 (clade 3) virus, Influenza B/Phuket/3073/2013 virus, Influenza B/Texas/06/2011 virus, Influenza B/Wisconsin/1/2010 virus, Influenza B/Wisconsin/1/2010 BX-41A virus (**B/Yamagata lineage**), showing positive result.

The reactivity of the VIASURE SARS-CoV-2, *Flu (A+B)* & RSV Real Time PCR Detection Kit for BD MAX™ System for **RSV** was confirmed against RNA extracted from RSV A and B (strain CH93(18)-18) and Human Respiratory Syncytial Virus strain Long, showing positive result.

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### 13. Bibliography/Bibliografia

1. Huang, C. *et al.* Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. *The Lancet*, 2020. DOI : 10.1016/S0140-6736(20)30183-5.
2. Zhu N. *et al.* A novel coronavirus from patients with pneumonia in China, 2019. *New England Journal of Medicine.*, 2020. DOI : 10.1056/NEJMoA2001017.
3. World Health Organization. MERS situation update. January 2020. Available from <https://applications.emro.who.int/docs/EMCSR254E.pdf?ua=1> Accessed September 2020.
4. Chen N. *et al.* Epidemiological and Clinical Characteristics of 99 Cases of 2019-Novel Coronavirus (2019-nCoV) Pneumonia in Wuhan, China. *The Lancet*, 2020. DOI: 10.1016/S0140-6736(20)30211-7.
5. Lv D.F. *et al.* Dynamic change process of target genes by RT-PCR testing of SARS-Cov-2 during the course of a Coronavirus Disease 2019 patient. *Clinica Chimica Acta* 2020; 506: 172-175.
6. World Health Organization. Clinical management of COVID-19 disease" Interim guidance 27 May 2020. Available from <https://www.who.int/publications/i/item/clinical-management-of-covid-19> Accessed September 2020.
7. Lu R. *et al.* Genomic characterisation and epidemiology of 2019 novel coronavirus: implications for virus origins and receptor binding. *The Lancet*, 2020. DOI : 10.1016/S0140-6736(20)30251-8.
8. Rothe C. *et al.* Transmission of 2019-nCoV Infection from an Asymptomatic Contact in Germany. *New England Journal of Medicine*, 2020. DOI : 10.1056/NEJMc2001468.
9. Centers of Disease Control and Prevention (CDC). Coronavirus Disease 2019 (COVID-19), Symptoms of Coronavirus. Available from <https://www.cdc.gov/coronavirus/2019-ncov/about/symptoms.html> Accessed September 2020.
10. Centers of Disease Control and Prevention (CDC). Coronavirus Disease 2019 (COVID-19), Older Adults. Available from <https://www.cdc.gov/coronavirus/2019-ncov/need-extra-precautions/older-adults.html> Accessed September 2020.
11. World Health Organization. Laboratory testing for 2019 novel coronavirus (2019-nCoV) in suspected human cases. Interim guidance. 19 March 2020. Available from <https://www.who.int/publications-detail/laboratory-testing-for-2019-novel-coronavirus-in-suspected-human-cases-20200117> Accessed September 2020.
12. Yan Y *et al.* Laboratory testing of SARS-CoV, MERS-CoV, and SARS-CoV-2 (2019-nCoV): Current status, challenges, and countermeasures. *Reviews in Medical Virology* 2020; 30(3):e2106.
13. Centers of Disease Control and Prevention (CDC). 2019-Novel Coronavirus (2019-nCoV) Real-time rRT-PCR Panel Primers and Probes. <https://www.cdc.gov/coronavirus/2019-ncov/downloads/rt-pcr-panel-primer-probes.pdf> Accessed September 2020.

14. Chu D.K.W. et al. Molecular Diagnosis of a Novel Coronavirus (2019-nCoV) Causing an Outbreak of Pneumonia. *Clinical Chemistry* 2020;66(4): 549-555.
15. Corman V.M. et al. Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. *European communicable disease bulletin* 2020;25(3).
16. European Centre for Disease Prevention and Control. Novel coronavirus disease 2019 (COVID-19) pandemic: increased transmission in the EU/EEA and the UK – sixth update – 12 March 2020. Stockholm: ECDC; 2020. Available from <https://www.ecdc.europa.eu/sites/default/files/documents/RRA-sixth-update-Outbreak-of-novel-coronavirus-disease-2019-COVID-19.pdf> Accessed September 2020.
17. Lim, Y. X., Ng, Y. L., Tam, J. P., & Liu, D. X. (2016). Human coronaviruses: a review of virus–host interactions. *Diseases*, 4(3), 26.
18. McBride R. et al. The coronavirus nucleocapsid is a multifunctional protein. *Viruses* 2014; 6(8):2991-3018.
19. Sheikh A. et al. Analysis of preferred codon usage in the coronavirus N genes and their implications for genome evolution and vaccine design. *Journal of Virological Methods* 2020; 277:113806.
20. World Health Organization. Laboratory testing strategy recommendations for COVID-19: interim guidance Interim guidance. 21 March 2020. Available from <https://www.who.int/publications/i/item/laboratory-testing-strategy-recommendations-for-covid-19-interim-guidance> Accessed September 2020.
21. G. Neumann et al. Transmission of Influenza A viruses. *Virology* 2015; 234-246.
22. Y. Yang et al. Simultaneous typing and HA/NA subtyping of influenza A and B viruses including the pandemic influenza A/H1N1 2009 by multiplex real-time RT-PCR. *Journal of Virological Methods* 2010; 167(1): 37-44.
23. R.L. Kuo et al. Influenza A/B virus detection and influenza A virus subtyping with emphasis on the novel H7N9 virus by using multiplex real-time RT-PCR. *Journal of Virological Methods* 2014; 208:41-46.
24. World Health Organization. WHO information for molecular diagnosis of influenza virus—update. Available from: [https://www.who.int/influenza/gisrs\\_laboratory/molecular\\_diagnosis/en/](https://www.who.int/influenza/gisrs_laboratory/molecular_diagnosis/en/) . Accessed September 2020.
25. S. Subhash Bawage et al. Recent Advances in Diagnosis, Prevention, and Treatment of Human Respiratory Syncytial Virus. *Advances in Virology* 2013.
26. French, et al. Risk of nosocomial respiratory syncytial virus infection and effectiveness of control measures to prevent transmission events: a systematic review. *Influenza and Other Respiratory Viruses* 2016.
27. X. Yu et al. Human respiratory syncytial virus in children with lower respiratory tract infections or influenza-like illness and its co-infection characteristics with viruses and atypical bacteria in Hangzhou, China. *Journal of Clinical Virology* 2015; 69:1-6.
28. N. Mazur et al. Lower respiratory tract infection caused by respiratory syncytial virus: current management and new therapeutics. *The Lancet Respiratory Medicine* 2015; 3: 888-900.
29. F. de-Paris et al. Optimization of one-step duplex real-time RT-PCR for detection of influenza and respiratory syncytial virus in nasopharyngeal aspirates. *Journal of Virological Methods* 2012; 186(1-2): 189-192.
30. A. Hu et al. Simultaneous detection, subgrouping, and quantitation of respiratory syncytial virus A and B by real-time PCR. *Journal of Clinical Microbiology* 2003; 41(1): 149-154.
31. M. Hindiyeh et al. Evaluation of Simplexa Flu A/B & RSV for direct detection of influenza viruses (A and B) and respiratory syncytial virus in patient respiratory samples. *Journal of Clinical Microbiology* 2013; 51 (7): 2421-2424.

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BD MAX™ is a registered trademark of Becton, Dickinson and Company.

**CerTest Biotec, S.L.**

Pol. Industrial Río Gállego II · Calle J, Nº1  
50840, San Mateo de Gállego, Zaragoza (Spain)  
[www.certest.es](http://www.certest.es)

VIASURE online

F-362 rev01

**VIASURE**

CE

IVD

Real Time PCR Detection Kits

**CerTest**  
BIOTEC