VIASURE

Real Time PCR Detection Kits

by **CerTest**
BIOTEC

SARS-CoV-2 (N1 + N2)

Handbook for the following references/

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ENGLISH

1. Intended use

VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit is a real-time RT-PCR test designed for the qualitative detection of RNA from the SARS-CoV-2 in respiratory samples (nasopharyngeal swab and oropharyngeal swab) and saliva samples from individuals suspected of COVID-19 by their healthcare provider. This test is intended for use as an aid in the diagnosis of SARS-CoV-2 in combination with clinical and epidemiological risk factors. RNA is extracted from respiratory specimens and saliva, amplified using RT-PCR and detected using fluorescent reporter dye probes specific for SARS-CoV-2.

2. Summary and Explanation

Coronavirus are enveloped non-segmented positive-sense RNA viruses and belong to *Coronaviridae* family. There are six coronavirus species known to cause human diseases. Four viruses (229E, OC43, NL63 and HKU1) cause common cold symptoms and the other two (severe acute respiratory syndrome coronavirus (SARS-CoV) and Middle East respiratory syndrome coronavirus (MERS-CoV)) are zoonotic and producing more severe complications. SARS-CoV and MERS-CoV have caused more than 10,000 cumulative cases in the past two decades, with mortality rates of 34% MERS-CoV and 10% SARS-CoV.

In December 2019, some people that worked at or lived around the Huanan seafood market in Wuhan, Hubei Province, China, have presented pneumonia of unknown cause. Deep sequencing analysis of the respiratory samples indicated a novel coronavirus, which was named firstly 2019 novel coronavirus (2019-nCoV) and lately SARS-CoV-2.

Human-to-human transmission of the SARS-CoV-2 has been confirmed, even in the incubation period without symptoms, and the virus could cause severe respiratory illness like those SARS-CoV produced. Although the pneumonia is the principal illness associated, a few patients have developed severe pneumonia, pulmonary edema, acute respiratory distress syndrome, or multiple organ failure and death. Centers of Disease Control and Prevention (CDC) believes that symptoms of SARS-CoV-2 may appear in as few as 2 days or as long as 14 days after exposure, being the most common fever or chills, cough, fatigue, anorexia, myalgia and dyspnea. Less common symptoms are sore throat, nasal congestion, headache, diarrhea, nausea and vomiting. Loss of smell (anosmia) or loss of taste (ageusia) preceding the onset of respiratory symptoms has also been reported. Older adults and people who have severe underlying medical conditions like heart or lung disease or diabetes seem to be at higher risk for developing more serious complications from COVID-19 illness.

Diagnosis of SARS-CoV-2 is performed detecting conventional causes of pneumonia early and detected by next-generation sequencing or real-time RT-PCR methods. Several assays that detect the SARS-CoV-2 have been currently available, such as China CDC (gene targets, *ORF1ab* and *N*), Charité – Germany (gene targets, *RdRP* and *E*) or US CDC (two targets in *N* gene).

CDC recommends upper respiratory tract specimens (nasopharyngeal (NP) swab, oropharyngeal (OP) swabs, nasal mid-turbinate swab, nasal swab, nasopharyngeal wash/aspirate or nasal wash/aspirate (NW) and a saliva specimens collected mainly by a healthcare provider) and/or lower respiratory specimens (sputum, endotracheal

aspirate, or bronchoalveolar lavage in patients with more severe respiratory disease) for the identification of SARS-CoV-2. In addition, other clinical specimens as blood, urine and stool may be collected to monitor the presence of the virus.

3. Principle of the procedure

VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit is designed for the qualitative detection of RNA from the SARS-CoV-2 in respiratory and saliva samples. The detection is done in one step real time RT format where the reverse transcription and the subsequent amplification of specific target sequence occur in the same reaction well. The isolated RNA target is transcribed generating complementary DNA by reverse transcriptase which is followed by the amplification of two conserved regions of N gene (N1 and N2) for SARS-CoV-2 using specific primers and a fluorescent-labeled probe.

VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit is based on the 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bounded to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence can be measured on Real Time PCR platforms.

VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase and retrotranscriptase) in an stabilized format, as well as an **endogenous internal control** to monitor the extraction process and/or discard the inhibition of the polymerase activity. The assay uses a human housekeeping gene as an **endogenous Internal Control (IC)** (human *RNase P* gene). Human housekeeping genes are involved in basic cell maintenance and, therefore, are expected to be present in all nucleated human cells and maintain relatively constant expression levels. N2 target is amplified and detected in FAM channel, N1 target is amplified and detected in Cy5 channel and the endogenous Internal Control which is amplified and detected in HEX channel, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

4. Reagents provided

VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables 1, 2 and 3. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table 1 includes materials and reagents to be used with 8-well strips compatible devices (See Annex 1). Table 2 includes materials and reagents to be used with 96-well plate compatible devices (See Annex 1). Table 3 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips.

Reagent/Material	Description	Colour	Amount
SARS-CoV-2 (N1 + N2) 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and endogenous Internal control in stabilized format	White	6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
SARS-CoV-2 (N1 + N2) Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 x 8-cap strip

Table 1. Reagents and materials provided in VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit with Ref. VS-NCO306L, VS-NCO306H, VS-NCO312L and VS-NCO312H.

Reagent/Material	Description	Color	Amount
SARS-CoV-2 (N1 + N2) 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and endogenous Internal control in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
SARS-CoV-2 (N1 + N2) Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 x 8-cap strip

Table 2. Reagents and materials provided in VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit with Ref VS-NCO313L and VS-NCO313H.

Reagent/Material	Description	Colour	Amount
SARS-CoV-2 (N1 + N2) 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and endogenous Internal control in stabilized format	Transparent	9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
SARS-CoV-2 (N1 + N2) Positive Control	Non-infectious synthetic lyophilized cDNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	9/18 x 4-cap strip

Table 3. Reagents and materials provided in VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit with Ref. VS-NCO336 and VS-NCO372. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- RNA extraction kit.
- Collection and transport system: BD™ Universal Viral Transport System (UTM), Viral Transport Media (VTM) Vircell S.L., Spain) or sterile Virus Transport and preservation Medium (Biocomma®) and equivalent.
- Laboratory freezers: - 30°C to - 10°C and/or ≤ -70°C.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments).

VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTLite Real-Time PCR System, Rotor-Gene® Q (Qiagen) and Roche Molecular Diagnostics Cobas z480 Analyzer. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- The product is intended for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.

- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-NCO313L, VS-NCO313H, VS-NCO336 and VS-NCO372). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For VS-NCO336 and VS-NCO372 (compatible for use with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed. Use separate areas for the preparation of patient samples and controls to prevent false positive results.
- Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink, smoke or apply cosmetic products in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious and/or biohazardous, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, transport, storage, handling, and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose of waste in compliance with local, state, and federal regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

8.1. Specimen collection, transport and storage

The VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection kit has been validated on nasopharyngeal and oropharyngeal specimens collected with synthetic fiber swabs with plastic and placed immediately into a sterile transport tube containing Universal transport medium (UTM) or Viral transport medium (VTM) and saliva samples collected in a sterile clean container (No preservative is required).

Patient samples must be collected, transport and storage according to appropriate laboratory guidelines. For details, refer to the CDC guidelines (Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens for COVID-19) (website <https://www.cdc.gov/coronavirus/2019-ncov/lab/guidelines-clinical-specimens.html>).

8.2. RNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For RNA extraction from respiratory samples you can use your manual or automatic routine optimized system. Also, you can use any commercially available RNA extraction kit and follow the manufacturer's instructions. We have validated the following extraction kits:

- Maxwell® RSC 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® RSC 16 instrument (Promega).
- Total Nucleic Acid Isolation (TNAI) Kit, using COBAS® AmpliPrep (ROCHE).
- MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.).
- QIAamp Viral RNA Mini Kit (QIAGEN).

For RNA extraction from saliva samples you can use your manual or automatic routine optimized system. Also, you can use any commercially available RNA extraction kit and follow the manufacturer's instructions. We have validated the following extraction kits:

- QIAamp Viral RNA Mini Kit (Qiagen).
- Maxwell® RSC 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® RSC 16 instrument (Promega).

8.3. Lyophilized positive control

SARS-CoV-2 (N1 + N2) Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized SARS-CoV-2 (N1 + N2) Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

8.4. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of RNA sample, reconstituted SARS-CoV-2 (N1 + N2) Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps. It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Reverse transcription	15 min	45°C
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 4. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (N2 target), Cy5 (N1 target) and HEX, JOE or VIC (Endogenous Internal Control (IC)). Depending on the equipment used select the proper detection channel (see Annex 2). In Applied Biosystems 7500 Fast Real-Time PCR System and Stratagene Mx3005PTM Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

9. Result interpretation

All the result of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms, and other diagnostic tests. Check Endogenous Internal Control (IC) signal to verify the extraction procedure and/or correct functioning of the amplification mix. The analysis of the controls and samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions. Using the following tables 5 and 6 read and analyze the results.

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal in the positive control well. For a valid diagnostic test run, the following control conditions must be met:

Controls	N2 target (FAM) ¹	N1 target (Cy5) ¹	Endogenous Internal Control (HEX)	Interpretation of Controls
Positive Control (PC)	≤40	≤40	≤40 ²	Valid
Negative Control (NC)	≥40 or no signal	≥40 or no signal	≥40 or no signal	Valid

Table 5. Expected Performance of Controls

1 In cases where either or both of the control assays have failed (an amplification signal is observed in the negative control and/or signals absence in the positive control well for any target channel), all results are reported as 'Invalid' and retesting is required.

2 The positive template control includes human housekeeping *RNase P* gene target; therefore amplification signals are observed in all target channels, including the Endogenous Internal Control.

Assessment of clinical samples test results should be performed after the positive and negative controls have been examined and determined to be valid and acceptable. If one or more controls are not valid, the patient results cannot be interpreted. For interpretation of patient sample results, use the following table:

N2 target (FAM)	N1 target (Cy5)	Endogenous Internal Control (HEX)	Interpretation for patients' samples	
≤40	≤40	≤40 or no signal ¹	Valid	SARS-CoV-2 RNA Detected
≤40 ²	≥40 or no signal	≤40 or no signal ¹	Valid	SARS-CoV-2 RNA Detected ²
≥40 or no signal	≤40 ²	≤40 or no signal ¹	Valid	SARS-CoV-2 RNA Detected ²
≥40 or no signal	≥40 or no signal	≤ 35 ³	Valid	SARS-CoV-2 RNA not Detected ³
≥40 or no signal	≥40 or no signal	≥ 35 or no signal ³	Invalid	Test Failure – Repeat Testing ³

Table 6. Interpretation of patient sample results. Ct values. no signal = no amplification curve.

1 The endogenous Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal). Sometimes, its detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids.

2 If only one target site of the *N* gene amplifies, verify the sigmoid shape of the curve and the intensity of fluorescence. In case of a doubtful interpretation, depending on the available material, it is also recommended to:


- a) re-extract and retest another aliquot of the same specimen or,
- b) repeat RT-qPCR with the same isolated RNA sample, or
- c) obtain a new specimen and retest.

3 In the case of SARS-CoV-2 target sites negative, IC must show an amplification signal with Ct less than 35. The Ct value could be very variable due to the Endogenous Internal Control is a human housekeeping gene that should be present in all human nucleated cells in the original sample. If there is an absence of signal or Ct value ≥ 35 of the endogenous Internal Control, the result is considered as 'Invalid', and retesting is required. It is recommended to repeat the RT-qPCR diluting the RNA sample 1:10 and/or 1:100, or re-extract and retest to check for possible failure in the extraction procedure and/or inhibition issues.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user; to verify the correct performance of each RT-qPCR steps and review the parameters; and to check the sigmoid shape of the curve and the intensity of fluorescence.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated only with RNA extracted respiratory samples (nasopharyngeal swab and oropharyngeal swab) and saliva samples.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by SARS-CoV-2, either the high number of cDNA template copies which contains each SARS-CoV-2 Positive Control vial, samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.
- False Negative results may arise from several factors and their combinations, including:
 - Improper specimens' collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including RNA extraction).
 - Degradation of the viral RNA during sample shipping/storage and/or processing.
 - Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown SARS-CoV-2 variants.
 - A viral load in the specimen below the limit of detection for the assay.
 - The presence of RT-qPCR inhibitors or other types of interfering substances. The impacts of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent COVID-19 or used during the treatment of the infection have not been evaluated.
 - Failure to follow instructions for use and the assay procedure.
- A single-target site amplification or even random positive results is suggestive of slightly different amplification yield of the target site of the *N* gene. Samples with low viral load might result in *N* single target amplification. In case of a doubt, it is recommended referring to a reference laboratory for further testing.

- Some samples may fail to exhibit *RNase P* amplification curves due to low human cell numbers in the original clinical sample. A negative IC signal does not preclude the presence of SARS-CoV-2 RNA in a clinical specimen.
- A positive test result does not necessarily indicate the presence of viable viruses and does not imply that these viruses are infectious or are the causative agents for clinical symptoms. However, a positive result is indicative of the presence of targets viral sequences (*N* genes).
- Negative results do not preclude SARS-CoV-2 infection and should not be used as the sole basis for treatment or other patient management decisions. Optimum specimen types and timing for peak viral levels during infections caused by SARS-CoV-2 have not been determined. The collection of multiple specimens (types and time points) from the same patient may be necessary to detect the virus.
- If diagnostic tests for other respiratory illnesses are negative and the patient's clinical presentation and epidemiological information suggest that SARS-CoV-2 infection is possible, then a false negative result should be considered, and a re-testing of the patient should be discussed.

11. Quality control

VIASURE SARS-CoV-2 (*N1 + N2*) Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the endogenous internal control (IC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE SARS-CoV-2 (*N1 + N2*) Real Time PCR Detection Kit was tested using 173 nasopharyngeal specimens collected in sterile Vircell® transport medium or sterile Virus transport and preservation medium (Biocomma®) from patients with clinical suspicion of COVID-19 disease. These results were compared with those obtained with a molecular detection method (Simplexa™ COVID-19 Direct assay (DiaSorin Molecular LLC, Cypress, CA)). The results were as follows. The results were as follows:

VIASURE SARS-CoV-2 (<i>N1 + N2</i>) Real Time PCR Detection Kit	Simplexa™ COVID-19 Direct assay		
		+	-
+	72	0	72
-	0	101	101
Total	72	101	173

Table 7. Comparative results for SARS-CoV-2.

The Positive Percent Agreement (PPA), Negative Percent Agreement (NPA) and Overall Percent Agreement (OPA) for VIASURE SARS-CoV-2 (*N1 + N2*) Real Time PCR detection kit were calculated in relation to the results from Simplexa™ COVID-19 Direct assay (DiaSorin Molecular LLC, Cypress, CA) as shown in Table 8. PPA, NPA,

OPA, and 95% confidence intervals were calculated using GraphPad QuickCalcs. The 95% confidence intervals were calculated using the modified Wald method.

Microorganism	PPA (%)	NPA (%)	OPA (%)
SARS-CoV-2	100.0 (93.9 to 100.0)	100.0 (95.6 to 100.0)	100.0 (97.4 to 100.0)

Table 8 PPA, NPA and OPA values for VIASURE SARS-CoV-2 Real Time PCR detection kit (95% confidence interval).

Besides, the clinical performance of VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit was tested using 152 nasopharyngeal specimens. Seven positive samples were in Vircell® transport medium and the rest were collected in Virus transport and preservation medium (Biocomma®). These last samples were inactivated before initial clinical diagnosis. These results were compared with those obtained with a molecular detection method (Allplex™ 2019-nCoV Assay (Seegene, Seoul, South Korea)). The results were as follows. The results were as follows:

VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit	Allplex™ 2019-nCoV Assay			
		+	-	Total
	+	74	0	74
	-	3*	75	78
	Total	77	75	152

Table 9. Comparative results for SARS-CoV-2.

* The low amount of template RNA in this respiratory sample is below the detection limit of the method used.

The Positive Percent Agreement (PPA), Negative Percent Agreement (NPA) and Overall Percent Agreement (OPA) for VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR detection kit were calculated in relation to the results from Allplex™ 2019-nCoV Assay (Seegene, Seoul, South Korea) as shown in Table 10. PPA, NPA, OPA, and 95% confidence intervals were calculated using GraphPad QuickCalcs. The 95% confidence intervals were calculated using the modified Wald method.

Microorganism	PPA (%)	NPA (%)	OPA (%)
SARS-CoV-2	96.1 (88.7 to 99.1)	100.0 (94.1 to 100.0)	98.1 (94.1 to 99.6)

Table 10. PPA, NPA and OPA values for VIASURE SARS-CoV-2 Real Time PCR detection kit (95% confidence interval).

Besides, the clinical performance of VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit was conducted testing in parallel 24 anonymized saliva specimens and oropharyngeal swabs from patients with signs and symptoms of COVID-19 infection. The results were as follows:

VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit (saliva specimens)	VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit (oropharyngeal swabs)			
		+	-	Total
	+	11	0	11
	-	3	10	13
Total	14	10	24	

Table 11. Comparative results for different types of clinical specimens.

The values obtained with the oropharyngeal samples were taken as the reference values and it was demonstrated that the results obtained when analyzing the saliva samples in parallel, were as expected.

Results show high agreement to detect SARS-CoV-2 using VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit has a detection limit of 10 genome copies per reaction for SARS-CoV-2 (N2 target) and 20 genome copies per reaction SARS-CoV-2 (N1 target).

Figure 2. Dilution series of N2 target (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (FAM channel).



Figure 3. Dilution series of N1 target (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Cy5 channel).

12.3. Analytical specificity

The specificity of the SARS-CoV-2 assay was confirmed by testing a panel consisting of different microorganisms representing the most common respiratory pathogens. No cross-reactivity was detected between any of the following microorganisms tested.

Cross-reactivity testing					
Human Adenovirus types 1-5, 8, 15, 31, 40 and 41	-	Influenza A/New Caledonia/20/99(H1N1) virus	-	<i>Legionella longbeachae</i>	-
Human Bocavirus	-	Influenza A/California/7/2009(H1N1)pdm09	-	<i>Legionella micdadei</i>	-
<i>Bordetella bronchiseptica</i>	-	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	-	<i>Legionella pneumophila</i>	-
<i>Bordetella holmesii</i>	-	Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus	-	Human metapneumovirus A and B	-
<i>Bordetella parapertussis</i>	-	Influenza A/Victoria/210/2009 (H3N2)	-	<i>Moraxella catarrhalis</i>	-
<i>Bordetella pertussis</i>	-	Influenza A/Thüringen/5/17 (H3N2) virus	-	<i>Mycoplasma pneumoniae</i>	-
<i>Chlamydia caviae</i>	-	Influenza A/Switzerland/9715293/2013 (H3N2) virus	-	<i>Mycobacterium tuberculosis</i> not rifampin resistant	-
<i>Chlamydia psittaci</i> genotype A and C	-	Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus	-	Human parainfluenza 1, 2, 3 and 4 viruses	-
<i>Chlamydophila pneumoniae</i> CM-1	-	Influenza A/South Australia/55/2014, IVR-175 (H3N2) virus	-	<i>Pneumocystis jirovecii</i>	-
Human coronavirus 229E, OC43, NL63 and HKU1	-	Influenza A/DE-SH/Reiherente/AR8444/ 2016 (H5N8) virus	-	Human rhinovirus type C	-
MERS Coronavirus	-	Influenza A/Anhui/1/2013 (H7N9) virus	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-
SARS Coronavirus Strain Frankfurt 1	-	Influenza B/Brisbane/60/2008	-	<i>Staphylococcus epidermidis</i>	-
Enterovirus 68 and 71	-	Influenza B/Florida/04/06 virus	-	<i>Streptococcus pneumoniae</i> Z022	-
Enterovirus Echovirus 11 and 30	-	Influenza B/Phuket/3073/2013 virus	-	<i>Streptococcus pyogenes</i>	-
Enterovirus Coxsackievirus A24, A9 and B3	-	<i>Legionella bozemanii</i>	-	<i>Streptococcus salivarius</i>	-
<i>Haemophilus influenzae</i> MinnA	-	<i>Legionella dumoffii</i>	-	Respiratory syncytial virus (RSV) A and B	-

Table 12. Reference pathogenic microorganisms used in this study

12.4. Analytical reactivity

The reactivity of VIASURE SARS-CoV-2 (N1 + N2) Real Time PCR Detection Kit was evaluated against RNA from Human 2019-nCoV strain BetaCoV/Germany/BavPat1/2020 p.1, Human 2019-nCoV strain 2019-nCoV/Italy-INMI1, SARS-CoV-2 strain 2019nCoV/USA-WA1/2020, synthetic RNA controls for two variants of the SARS-CoV-2 virus: MT007544.1 (SARS-CoV2 isolate Australia/VIC01/2020) and MN908947.3 (SARS-CoV-2 isolate Wuhan-Hu-1), showing positive results.

ANNEX 1

COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

VIASURE Real Time PCR Detection Kits are available in a ready-to-use lyophilized format placed inside wells with different dimensions, low or high profile. Depending on the thermal block of the equipment to be used, one measure or another will fit. Please, consult the table and check the specifications of your equipment. If the equipment does not appear in the list below, please contact with your supplier. This table is for guidance, it is recommended to check the equipment before running the (RT)-qPCR.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast / 7500 Fast Dx Real-Time PCR System ^{(1) (5)}
	QuantStudio™ 12K Flex 96-well Fast
	QuantStudio™ 6 Flex 96-well Fast
	QuantStudio™ 7 Flex 96-well Fast
	QuantStudio™ 3 Fast Real-Time PCR System ⁽²⁾
	QuantStudio™ 5 Fast/ QuantStudio™ 5 Real-Time PCR System
	StepOne Plus™ Real-Time PCR System ⁽²⁾
	StepOne™ Real-Time PCR System ⁽²⁾
	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96 Fast
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
	Mini Opticon™ Real-Time PCR Detection System ⁽³⁾
Bio Molecular Systems	Mic Real Time PCR Cyclers ⁽⁴⁾
Cepheid	SmartCycler® ⁽⁴⁾
Precision System Science Co., Ltd. (PSS)	geneLEAD VIII System ⁽⁴⁾
Qiagen	Rotor-Gene® Q ⁽⁴⁾
Roche	LightCycler® 480 Real-Time PCR System ⁽⁵⁾
	LightCycler® 96 Real-Time PCR System ⁽⁵⁾
	Cobas z480 Analyzer ⁽⁵⁾

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbott m2000 RealTime System ⁽⁵⁾
	7300 Real-Time PCR System ^{(2) (5)}
Applied Biosystems	7500 Real-Time PCR System ⁽⁵⁾
	7900 HT Real-Time PCR System ⁽²⁾
	ABI PRISM 7000 ⁽³⁾
	ABI PRISM 7700 ⁽²⁾
	QuantStudio™ 12K Flex 96-well
	QuantStudio™ 6 Flex 96-well
	QuantStudio™ 7 Flex 96-well
	QuantStudio™ 3 Real-Time PCR System ⁽²⁾
	QuantStudio™ 5 Fast/ QuantStudio™ 5 Real-Time PCR System
	ViiA™ 7 Real-Time PCR System
	Analytik Jena Biometra
qTOWER 2.0	
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
	iCycler iQ™ Real-Time PCR Detection System
	iCycler iQ™5 Real-Time PCR Detection System
	MyiQ™ Real-Time PCR Detection System ⁽³⁾
Bio Molecular Systems	MyiQ™2 Real-Time PCR Detection System ⁽³⁾
	Mic Real Time PCR Cyclers ⁽⁴⁾
Cepheid	SmartCycler® ⁽⁴⁾
DNA-Technology	DTprime Real-time Detection Thermal Cyclers
	DTlite Real-Time PCR System
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽⁴⁾
Precision System Science Co., Ltd. (PSS)	geneLEAD VIII System ⁽⁴⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

(1) Select Ramp Speed "Standard".
 (2) No detection in Cy5 channel.
 (3) Detection in FAM and HEX channels only
 (4) The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Mic, SmartCycler®, Rotor-Gene® Q or geneLEAD VIII System tubes.
 (5) Shell Frame grid plate which fits in these qPCR System is necessary.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

ANNEX 2

DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be "none". Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Lightcycler®480II Roche	FAM	465/510	Colour Compensation is required for Roche Thermocyclers
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Cobas z 480 Roche	FAM	465/510	Colour Compensation is required for Roche Thermocyclers
	HEX	540/580	
	ROX	540/610	
	Cy5	610/670	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene/Agilent Technologies	FAM	FAM	Passive reference option for ROX must be "none"
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range must be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

ANNEX 3

OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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13. Bibliography/Bibliografía


- Huang, C. *et al.* Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. *The Lancet*, 2020. DOI : 10.1016/S0140-6736(20)30183-5.
- Zhu N. *et al.* A novel coronavirus from patients with pneumonia in China, 2019. *New England Journal of Medicine.*, 2020. DOI : 10.1056/NEJMoa2001017.
- World Health Organization. MERS situation update. January 2020. Available from <https://applications.emro.who.int/docs/EMCSR254E.pdf?ua=1> Accessed September 2020.
- Chen N. *et al.* Epidemiological and Clinical Characteristics of 99 Cases of 2019-Novel Coronavirus (2019-nCoV) Pneumonia in Wuhan, China. *The Lancet*, 2020. DOI: 10.1016/S0140-6736(20)30211-7.

5. Lv D.F. *et al.* Dynamic change process of target genes by RT-PCR testing of SARS-Cov-2 during the course of a Coronavirus Disease 2019 patient. *Clinica Chimica Acta* 2020; 506: 172-175.
6. World Health Organization. Clinical management of COVID-19 disease" Interim guidance 27 May 2020. Available from <https://www.who.int/publications/i/item/clinical-management-of-covid-19> Accessed September 2020.
7. Lu R. *et al.* Genomic characterisation and epidemiology of 2019 novel coronavirus: implications for virus origins and receptor binding. *The Lancet*, 2020. DOI : 10.1016/S0140-6736(20)30251-8.
8. Rothe C. *et al.* Transmission of 2019-nCoV Infection from an Asymptomatic Contact in Germany. *New England Journal of Medicine*, 2020. DOI : 10.1056/NEJMc2001468.
9. Centers of Disease Control and Prevention (CDC). Coronavirus Disease 2019 (COVID-19), Symptoms of Coronavirus. Available from <https://www.cdc.gov/coronavirus/2019-ncov/about/symptoms.html> Accessed September 2020.
10. Centers of Disease Control and Prevention (CDC). Coronavirus Disease 2019 (COVID-19), Older Adults. Available from <https://www.cdc.gov/coronavirus/2019-ncov/need-extra-precautions/older-adults.html> Accessed September 2020.
11. World Health Organization. Laboratory testing for 2019 novel coronavirus (2019-nCoV) in suspected human cases. Interim guidance. 19 March 2020. Available from <https://www.who.int/publications-detail/laboratory-testing-for-2019-novel-coronavirus-in-suspected-human-cases-20200117> Accessed September 2020.
12. Yan Y *et al.* Laboratory testing of SARS-CoV, MERS-CoV, and SARS-CoV-2 (2019-nCoV): Current status, challenges, and countermeasures. *Reviews in Medical Virology* 2020; 30(3):e2106.
13. Centers of Disease Control and Prevention (CDC). 2019-Novel Coronavirus (2019-nCoV) Real-time rRT-PCR Panel Primers and Probes. <https://www.cdc.gov/coronavirus/2019-ncov/downloads/rt-pcr-panel-primer-probes.pdf> Accessed September 2020.
14. Chu D.K.W. *et al.* Molecular Diagnosis of a Novel Coronavirus (2019-nCoV) Causing an Outbreak of Pneumonia. *Clinical Chemistry* 2020;66(4): 549-555.
15. Corman V.M. *et al.* Detection of 2019 novel coronavirus (2019-nCoV) by real-time RT-PCR. *European communicable disease bulletin* 2020;25(3).
16. European Centre for Disease Prevention and Control. Novel coronavirus disease 2019 (COVID-19) pandemic: increased transmission in the EU/EEA and the UK – sixth update – 12 March 2020. Stockholm: ECDC; 2020. Available from <https://www.ecdc.europa.eu/sites/default/files/documents/RRA-sixth-update-Outbreak-of-novel-coronavirus-disease-2019-COVID-19.pdf> Accessed September 2020.
17. Lim, Y. X., Ng, Y. L., Tam, J. P., & Liu, D. X. (2016). Human coronaviruses: a review of virus–host interactions. *Diseases*, 4(3), 26.
18. McBride R. *et al.* The coronavirus nucleocapsid is a multifunctional protein. *Viruses* 2014; 6(8):2991-3018.
19. Sheikh A. *et al.* Analysis of preferred codon usage in the coronavirus N genes and their implications for genome evolution and vaccine design. *Journal of Virological Methods* 2020; 277:113806.
20. World Health Organization. Laboratory testing strategy recommendations for COVID-19: interim guidance Interim guidance. 21 March 2020. Available from <https://www.who.int/publications/i/item/laboratory-testing-strategy-recommendations-for-covid-19-interim-guidance> Accessed September 2020.

14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Tabla A.1 TERMOCICLADORES CON BLOQUE DE BAJO PERFIL		Tabla A.2 TERMOCICLADORES CON BLOQUE DE PERFIL ALTO	
Fabricante	Modelo	Fabricante	Modelo
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System	Abbott	Abbott m2000 RealTime System ⁽⁵⁾
Applied Biosystems	7500 Fast / 7500 Fast Dx Real-Time PCR System ^{(1) (5)}	Applied Biosystems	7300 Real-Time PCR System ^{(2) (5)}
	QuantStudio™ 12K Flex 96-well Fast		7500 Real-Time PCR System ⁽⁵⁾
	QuantStudio™ 6 Flex 96-well Fast		7900 HT Real-Time PCR System ⁽²⁾
	QuantStudio™ 7 Flex 96-well Fast		ABI PRISM 7000 ⁽³⁾
	QuantStudio™ 3 Fast Real-Time PCR System ⁽²⁾		ABI PRISM 7700 ⁽²⁾
	QuantStudio™ 5 Fast/ QuantStudio™ 5 Real-Time PCR System		QuantStudio™ 12K Flex 96-well
	StepOne Plus™ Real-Time PCR System ⁽²⁾		QuantStudio™ 6 Flex 96-well
	StepOne™ Real-Time PCR System ⁽²⁾		QuantStudio™ 7 Flex 96-well
	ViiA™ 7 Fast Real-Time PCR System		QuantStudio™ 3 Real-Time PCR System ⁽²⁾
BIONEER	Exicycler™ 96 Fast		QuantStudio™ 5 Fast/ QuantStudio™ 5 Real-Time PCR System
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System		ViiA™ 7 Real-Time PCR System
	Mini Opticon™ Real-Time PCR Detection System ⁽³⁾		TOptical
Bio Molecular Systems	Mic Real Time PCR Cyclers ⁽⁴⁾	Analytik Jena Biometra	qTOWER 2.0
Cepheid	SmartCycler® ⁽⁴⁾	BIONEER	Exicycler™ 96
Precision System Science Co., Ltd. (PSS)	geneLEAD VIII System ⁽⁴⁾	Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Qiagen	Rotor-Gene® Q ⁽⁴⁾		iCycler iQ™ Real-Time PCR Detection System
Roche	LightCycler® 480 Real-Time PCR System ⁽⁵⁾		iCycler iQ™5 Real-Time PCR Detection System
	LightCycler® 96 Real-Time PCR System ⁽⁵⁾		MyiQ™ Real-Time PCR Detection System ⁽³⁾
	Cobas z480 Analyzer ⁽⁵⁾	MyiQ™2 Real-Time PCR Detection System ⁽³⁾	
		Bio Molecular Systems	Mic Real Time PCR Cyclers ⁽⁴⁾
		Cepheid	SmartCycler® ⁽⁴⁾
		DNA-Technology	DTprime Real-time Detection Thermal Cycler
			DTlite Real-Time PCR System
		Eppendorf	Mastercycler™ep realplex
		Qiagen	Rotor-Gene® Q ⁽⁴⁾
		Precision System Science Co., Ltd. (PSS)	geneLEAD VIII System ⁽⁴⁾
		Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
			Mx3005P™ Real Time PCR System

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: December 2020

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
www.certest.es

VIASURE online

F-362 rev01

VIASURE

Real Time PCR Detection Kits

CerTest
BIOTEC