VIASURE

Real Time PCR Detection Kits

by CerTest
BIOTEC

Salmonella, Campylobacter & Shigella/EIEC

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VIASURE <i>Salmonella, Campylobacter & Shigella/EIEC</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-SCS106L
VIASURE <i>Salmonella, Campylobacter & Shigella/EIEC</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-SCS106H
VIASURE <i>Salmonella, Campylobacter & Shigella/EIEC</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-SCS112L
VIASURE <i>Salmonella, Campylobacter & Shigella/EIEC</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-SCS112H
VIASURE <i>Salmonella, Campylobacter & Shigella/EIEC</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-SCS113L
VIASURE <i>Salmonella, Campylobacter & Shigella/EIEC</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-SCS113H
VIASURE <i>Salmonella, Campylobacter & Shigella/EIEC</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-SCS136
VIASURE <i>Salmonella, Campylobacter & Shigella/EIEC</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-SCS172

ENGLISH

1. Intended use

VIASURE *Salmonella*, *Campylobacter* & *Shigella*/EIEC Real Time PCR Detection Kit is designed for specific detection and differentiation of *Salmonella*, *Campylobacter* and *Shigella*/enteroinvasive *Escherichia coli* (EIEC) in human stool samples from patients with signs and symptoms of gastrointestinal infection. This test is intended for use as an aid in the diagnosis of the *Salmonella*, *Campylobacter* and/or *Shigella*/EIEC in combination with clinical and epidemiological risk factors. DNA is extracted from stool specimens, multiplied using Real Time amplification and detected using specific primers and a fluorescent reporter dye probe for *Salmonella*, *Campylobacter* and/or *Shigella*/EIEC.

2. Summary and Explanation

The *Salmonella* genus is divided taxonomically into six *Salmonella enterica* subspecies and *Salmonella bongori*. *Salmonella* causes two types of diseases: enteric (typhoid) fever and acute gastroenteritis commonly referred to as salmonellosis. Transmission of *Salmonella* occurs through contaminated food (poultry, poultry products, beef, pork, eggs, milk, and seafood), water or contact with infected animals. Patients infected with *Salmonella* frequently suffer nausea, vomiting, abdominal cramps/pains, achiness, fever, intense diarrhea and/or headache.

The *Campylobacter* genus consists of 26 species, 2 provisional species, and 9 subspecies. The most common species of *Campylobacter* associated with human illness are *C. jejuni* and *C. coli*, but other species can also cause human infections. Poultry is a major reservoir and source of transmission of *Campylobacter* to humans. Other risk factors include consumption of animal products and water, contact with animals, and even person-to-person transmission (fecal-oral or via fomites). Infection with *Campylobacter* causes gastroenteritis characterised by fever, vomiting, headaches, and abdominal pain with watery or bloody diarrhea, for a median duration of 6 days.

Shigellosis is a dysentery characterized by severely bloody and mucus-containing diarrhea. The disease is caused by any of the four species of *Shigella* (*S. dysenteriae*, *S. flexneri*, *S. boydii* and *S. sonnei*) and by enteroinvasive *Escherichia coli* (EIEC). The differentiation of *Shigella* and enteroinvasive *E. coli* is complicated due to the ability to cause dysentery using of the same method of invasion. *Shigella* is a significant agent of foodborne illness, especially with foods that require hand processing and/or are prepared from raw or previously cooked products without re-heating.

Target genes most frequently employed for rapid detection and identification of *Salmonella*, *Campylobacter* and *Shigella*/EIEC using Real Time PCR assays are the virulence *invA* gene (*Salmonella* invasion protein gene), 16S rRNA gene and *ipaH* gene, respectively.

3. Principle of the procedure

VIASURE *Salmonella*, *Campylobacter* & *Shigella*/EIEC Real Time PCR Detection Kit is designed for the diagnosis of the *Salmonella*, *Campylobacter* and/or *Shigella*/EIEC in human stool samples. After DNA isolation, the identification of *Salmonella*, *Campylobacter* and/or *Shigella*/EIEC is performed by the amplification of a

conserved region of the *invA* gene for *Salmonella*, 16S rRNA gene for *Campylobacter* and *ipaH* gene for *Shigella/EIEC*, using specific primers and a fluorescent-labelled probe.

VIASURE *Salmonella*, *Campylobacter* & *Shigella/EIEC* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE *Salmonella*, *Campylobacter* & *Shigella/EIEC* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase, retrotranscriptase) in a stabilized format, as well as an internal control to discard the inhibition of the polymerase activity. *Salmonella* DNA targets are amplified and detected in FAM channel, *Campylobacter* DNA targets are amplified and detected in Cy5 channel, *Shigella/EIEC* DNA targets are amplified and detected in ROX channel and the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

4. Reagents provided

VIASURE *Salmonella*, *Campylobacter* & *Shigella/EIEC* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables 1, 2 and 3. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table 1 includes materials and reagents to be used with 8-well strips compatible devices (See Annex 1). Table 2 includes materials and reagents to be used with 96-well plate compatible devices (See Annex 1). Table 3 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips.

Reagent/Material	Description	Colour	Amount
<i>Salmonella</i> , <i>Campylobacter</i> & <i>Shigella/EIEC</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Salmonella</i> , <i>Campylobacter</i> & <i>Shigella/EIEC</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE *Salmonella*, *Campylobacter* & *Shigella/EIEC* Real Time PCR Detection Kit with Ref. VS-SCS106L, VS-SCS106H, VS-SCS112L and VS-SCS112H.

Reagent/Material	Description	Color	Amount
<i>Salmonella</i> , <i>Campylobacter</i> & <i>Shigella</i> /EIEC 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Salmonella</i> , <i>Campylobacter</i> & <i>Shigella</i> /EIEC Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 X 8-cap strip

Table 2. Reagents and materials provided in VIASURE *Salmonella*, *Campylobacter* & *Shigella*/EIEC Real Time PCR Detection Kit with Ref VS-SCS113L and VS-SCS113H.

Reagent/Material	Description	Colour	Amount
<i>Salmonella</i> , <i>Campylobacter</i> & <i>Shigella</i> /EIEC 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	Transparent	9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>Salmonella</i> , <i>Campylobacter</i> & <i>Shigella</i> /EIEC Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	9/18 X 4-cap strip

Table 3. Reagents and materials provided in VIASURE *Salmonella*, *Campylobacter* & *Shigella*/EIEC Real Time PCR Detection Kit with Ref. VS-SCS136 and VS-SCS172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *Salmonella*, *Campylobacter* & *Shigella*/EIEC Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortex.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments).

VIASURE *Salmonella*, *Campylobacter* & *Shigella*/EIEC Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR

Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTLite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- The product is indented for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-SCS113L, VS-SCS113H, VS-SCS136 and VS-SCS172). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-SCS136 and VS-SCS172 (compatible with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.

- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

8.1. Sample preparation

Stool samples should be collected in clean containers and processed as soon as possible to guarantee the quality of the test. We recommend to use fresh samples.

For longer storage, the samples must be frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation. Freezing and thawing cycles are not recommended.

Perform the sample preparation according to the recommendations appearing in the instructions for use of extraction kit used.

8.1.1. DNA extraction

For DNA extraction from human stool samples you can use your manual or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- QIAamp DNA Mini kit (QIAGEN).
- QIAamp DNA Stool Mini Kit (QIAGEN).
- Invisorb® Spin Universal Kit (Strattec).
- RIDA® Xtract (r-Biopharm).
- Maxwell®RSC Blood DNA Kit, using the Maxwell® 16 instrument (Promega).

8.2. Lyophilized positive control

Salmonella, *Campylobacter* & *Shigella*/EIEC Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *Salmonella*, *Campylobacter* & *Shigella*/EIEC Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

8.3. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *Salmonella*, *Campylobacter* & *Shigella/EIEC* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® kits).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 5. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (*Salmonella*), Cy5 (*Campylobacter*), ROX (*Shigella/EIEC*) and HEX, JOE or VIC channels (Internal Control (IC)). Depending on the equipment used select the proper detection channel (see Annex 2). In Applied Biosystems 7500 Fast Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *Salmonella*, *Campylobacter* & *Shigella/EIEC* in the positive control well. Check Internal Control signal to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

Using the following table read and analyze the results:

Salmonella (FAM)	Campylobacter (Cy5)	Shigella/EIEC (ROX)	Internal control (HEX)	Negative control	Positive control	Interpretation
+	+	+	+/-	-	+	Salmonella, Campylobacter and Shigella/EIEC Positive
-	-	-	+	-	+	Salmonella, Campylobacter and Shigella/EIEC Negative
+	-	-	+/-	-	+	Salmonella Positive, Campylobacter and Shigella/EIEC Negative
+	+	-	+/-	-	+	Salmonella and Campylobacter Positive and Shigella/EIEC Negative
+	-	+	+/-	-	+	Salmonella and Shigella/EIEC Positive and Campylobacter Negative
-	+	-	+/-	-	+	Campylobacter Positive, Salmonella and Shigella/EIEC Negative
-	+	+	+/-	-	+	Campylobacter and Shigella/EIEC Positive, Salmonella Negative
-	-	+	+/-	-	+	Shigella/EIEC Positive, Salmonella and Campylobacter Negative
+	+	+	+	+	+	Experiment fail
-	-	-	-	-	-	Experiment fail

Table 6. Sample interpretation

+: Amplification curve


-: No amplification curve

A sample is considered positive if the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated only with DNA extracted from human faecal samples and urine samples.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *Salmonella*, *Campylobacter* and/or *Shigella/EIEC*, either samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

11. Quality control

VIASURE *Salmonella*, *Campylobacter* & *Shigella*/EIEC Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

12. Performance characteristics

12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE *Salmonella*, *Campylobacter* & *Shigella*/EIEC Real Time PCR Detection Kit was tested using 400 faecal samples from symptomatic patients. The results were compared with those obtained by commercial Real Time RT-PCR kit ("RIDA®GENE Bacterial Stool Panel" (r-Biopharm) for *Salmonella* and *Campylobacter*).

The results were as follows:

VIASURE <i>Salmonella</i> , <i>Campylobacter</i> & <i>Shigella</i> /EIEC Real Time PCR Detection Kit	RIDA®GENE Bacterial Stool Panel (r-Biopharm)			
		+	-	Total
	+	21	1*#	23
	-	0	378	379
	Total	21	379	400

Table 7. Comparative results for *Salmonella*

* The low amount of template DNA detected in these faecal samples is below the detection limit of the method used and so it may produce random positive results.

Besides, this sample was culture positive for *Salmonella* G.

VIASURE <i>Salmonella</i> , <i>Campylobacter</i> & <i>Shigella</i> /EIEC Real Time PCR Detection Kit	RIDA®GENE Bacterial Stool Panel (r-Biopharm)			
		+	-	Total
	+	58	13*#	71
	-	0	329	329
	Total	58	342	400

Table 8. Comparative results for *Campylobacter*

* The low amount of template DNA detected in these faecal samples is below the detection limit of the method used and so it may produce random positive results.

Besides, these 13 positive samples have been evaluated by an additional commercial Real Time PCR Kit (Mericon-*Campylobacter* spp, QIAGEN), confirming our results.

Overall, a total amount of 400 faecal samples from symptomatic patients were also tested for *Shigella*/EIEC, which was detected in 1 sample. This result was confirmed with that obtained by a commercial Real Time PCR Kit

("RIDA®GENE EHEC/EPEC" (r-Biopharm), which detects and differentiates EHEC, STEC (EHEC), EPEC, *Shigella dysenteriae* type 1 and *Shigella/EIEC*).

The results show a high sensitivity and specificity to detect *Salmonella*, *Campylobacter* and *Shigella/EIEC* using VIASURE *Salmonella*, *Campylobacter* & *Shigella/EIEC* Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE *Salmonella*, *Campylobacter* & *Shigella/EIEC* Real Time PCR Detection Kit has a detection limit of ≥ 10 DNA copies per reaction for *Salmonella*, *Campylobacter* and *Shigella/EIEC* (Figure 2, 3 and 4).

Figure 2. Dilution series of *Salmonella* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (FAM channel).


Figure 3. Dilution series of *Campylobacter* (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (Cy5 channel).



Figure 4. Dilution series of *Shigella*/EIEC (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (ROX channel)

12.3. Analytical specificity

The specificity of the *Salmonella*, *Campylobacter* and *Shigella*/EIEC assays were confirmed by testing a panel consisting of different microorganisms representing the most common enteric pathogens or flora present in the intestine. No cross-reactivity was detected between almost any of the following microorganisms tested, except the targeted pathogens of each assay.

Cross-reactivity testing					
<i>Helicobacter pylori</i>	-	<i>Campylobacter lari</i>	-/+	<i>Listeria monocytogenes</i>	-
<i>Helicobacter hepaticus</i>	-	<i>Campylobacter fetus</i>	-/+	<i>Candida albicans</i>	-
<i>Helicobacter cinaedi</i>	-	<i>Campylobacter coli</i>	-/+	<i>Arcobacter butzleri</i>	-
<i>Helicobacter heilmannii</i>	-	<i>Campylobacter jejuni</i> subsp. <i>jejuni</i>	-/+	<i>Pseudomonas aeruginosa</i>	-
<i>Salmonella typhi</i>	-/+	<i>Campylobacter upsaliensis</i>	-/+	<i>Enterococcus faecalis</i>	-
<i>Salmonella paratyphi</i> A	-/+	<i>Proteus vulgaris</i>	-	<i>Yersinia enterocolitica</i> O:3	-
<i>Salmonella paratyphi</i> B	-/+	<i>Aeromonas hydrophila</i> subsp. <i>hydrophila</i>	-	<i>Yersinia enterocolitica</i> O:9	-
<i>Salmonella typhimurium</i>	-/+	<i>Citrobacter freundii</i>	-	<i>Bacteroides fragilis</i>	-
<i>Salmonella bongori</i>	-/+	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-	<i>Cryptosporidium parvum</i>	-
<i>Salmonella enteritidis</i>	-/+	<i>Serratia liquefaciens</i>	-	<i>Giardia intestinalis</i>	-
<i>Salmonella enterica</i> subsp. <i>enterica</i>	-/+	<i>Vibrio parahaemolyticus</i>	-	<i>Entamoeba histolytica</i>	-
<i>Salmonella pullorum</i>	-/+	<i>Clostridium difficile</i>	-	Adenovirus serotypes 40	-
<i>Salmonella gallinarum</i>	-/+	<i>Clostridium perfringens</i>	-	Adenovirus serotypes 41	-
<i>Shigella flexneri</i>	-/+	Enterotoxigenic <i>Escherichia coli</i>	-	Rotavirus A	-
<i>Shigella dysenteriae</i>	-/+	Enteropathogenic <i>Escherichia coli</i>	-	Norovirus Genotypes I and II	-
<i>Shigella sonnei</i>	-/+	<i>Klebsiella oxytoca</i>	-	Astrovirus Genotype I-VIII	-
<i>Shigella boydii</i>	-/+				

Table 9. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE *Salmonella*, *Campylobacter* & *Shigella*/EIEC Real Time PCR Detection Kit for *Salmonella* was evaluated against *Salmonella enterica* subsp. *enterica*, *Salmonella paratyphi* A, *Salmonella paratyphi* B, *Salmonella typhimurium*, *Salmonella enteritidis*, *Salmonella typhi*, *Salmonella bongori*, *Salmonella pullorum*,

Salmonella gallinarum, *Salmonella mbandaka*, *Salmonella braenderup* and *Salmonella tennessee* showing positive results.

The reactivity of VIASURE *Salmonella*, *Campylobacter* & *Shigella*/EIEC Real Time PCR Detection Kit for *Campylobacter* was evaluated against *Campylobacter jejuni* subsp. *jejuni*, *Campylobacter coli*, *Campylobacter lari*, *Campylobacter upsaliensis*, *Campylobacter fetus*, *Campylobacter concisus*, *Campylobacter hyointestinalis*, *Campylobacter gracilis*, *Campylobacter helveticus*, *Campylobacter curvus* and *Campylobacter rectus* showing positive results.

The reactivity of VIASURE *Salmonella*, *Campylobacter* & *Shigella*/EIEC Real Time PCR Detection Kit for *Shigella*/EIEC was evaluated against *Shigella dysenteriae*, *Shigella flexneri*, *Shigella boydii* and *Shigella sonnei* showing positive results.

ANNEX 1

COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast Real-Time PCR System ⁽¹⁾
Applied Biosystems	7500 Fast Dx Real-Time PCR System ⁽¹⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 3 Fast Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Fast Real-Time PCR System
Applied Biosystems	StepOne Plus™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	StepOne™ Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	Mini Opticon™ Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
Qiagen	Rotor-Gene® Q ⁽³⁾
Roche	LightCycler® 480 Real-Time PCR System ⁽⁴⁾
Roche	LightCycler® 96 Real-Time PCR System ⁽⁴⁾
Roche	Cobas z480 Analyzer ⁽⁴⁾

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

(1) Select Ramp Speed "Standard".

(2) See Annex 3 to check optical measurement exposure setting.

(3) The product should be reconstituted following the appropriate procedure (see Test Procedure, section 8.3) and transferred into the specific tubes designed to perform on Rotor-Gene® Q or SmartCycler® instruments.

(4) Shell Frame grid plate which fits in these Roche qPCR System is necessary.

(5) No detection in Cy5 channel.

(6) Detection in FAM and HEX channels only.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

ANNEX 2

DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation is required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

ANNEX 3

OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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13. Bibliography/Bibliografía

1. J. Hoorfar *et al.* Automated 59 Nuclease PCR Assay for Identification of *Salmonella enterica*. *Journal of Clinical Microbiology* 2000; 38: 3429–3435.
2. D. Rodríguez-Lázaro *et al.* A rapid and direct real time PCR-based method for identification of *Salmonella* spp. *Journal of Microbiological Methods* 2003; 54: 381– 390.
3. F. Barletta *et al.* Multiplex real-time PCR for detection of *Campylobacter*, *Salmonella*, and *Shigella*. *Journal of Clinical Microbiology* 2013; 51: 2822-2829.
4. N.O. Kaakoush *et al.* Global Epidemiology of *Campylobacter* Infection. *Clinical Microbiology Reviews* 2015; 28(3): 687-720.
5. S.M. Man. The clinical importance of emerging *Campylobacter* species. *Nature Reviews Gastroenterology & Hepatology* 2011; 8(12): 669-685.
6. R.F. de Boer *et al.* Detection of *Campylobacter* species and *Arcobacter butzleri* in stool samples by use of real-time multiplex PCR. *Journal of Clinical Microbiology* 2013; 51(1): 253-259.
7. N. Botteldoorn *et al.* Quantification of *Campylobacter* spp. in chicken carcass rinse by real-time PCR. *Journal of Applied Microbiology* 2008; 105(6):1909-1918.
8. A.Ub-Din *et al.* Relationship among *Shigella* spp. and enteroinvasive *Escherichia coli* (EIEC) and their differentiation. *Brazilian Journal of Microbiology* 2014; 45(4): 1131-1138.
9. V. D. Thiem *et al.* Detection of *Shigella* by a PCR assay targeting the *ipaH* gene suggests increased prevalence of shigellosis in Nha Trang, Vietnam. *Journal of Clinical Microbiology* 2004; 42(5): 2031-2035.
10. W. S. Lin *et al.* A quantitative PCR assay for rapid detection of *Shigella* species in fresh produce. *Journal of Food Protection* 2010; 73(2): 221-233.
11. W. Mokhtari *et al.* Real-Time PCR using SYBR Green for the detection of *Shigella* spp. in food and stool samples. *Molecular and Cellular Probes* 2013; 27: 53-59.
12. S. Ghosh *et al.* Genetic characterization of *Shigella* spp. isolated from diarrhoeal and asymptomatic children. *Journal of Medical Microbiology* 2014; 63: 903-910.

14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Tabla A.2 TERMOCICLADORES CON BLOQUE DE PERFIL ALTO	
Fabricante	Modelo
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System ⁽⁵⁾
Applied Biosystems	ABI PRISM 7000 ⁽⁶⁾
Applied Biosystems	ABI PRISM 7700 ⁽⁵⁾
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System ⁽⁶⁾
Bio-Rad	MyiQ™2 Real-Time PCR Detection System ⁽⁶⁾
Cepheid	SmartCycler® ⁽³⁾
DNA-Technology	DTprime Real-time Detection Thermal Cycler ⁽²⁾
DNA-Technology	DTlite Real-Time PCR System ⁽²⁾
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q ⁽³⁾
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System ⁽²⁾
VIASURE	VIASURE 96 Real Time PCR System ⁽²⁾

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: March 2019

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
www.certest.es

VIASURE online

F-362 rev01

VIASURE

CE

IVD

Real Time PCR Detection Kits

CerTest
BIOTEC