

Compatibility guide for the most common real-time PCR equipment

VIASURE Real Time PCR kits are available in a ready-to-use lyophilized format within low-profile or high-profile wells, 2 ml vials or special tubes for Rotorgene.

Depending on the equipment used, it will adjust to one size or another. Please refer to the table and verify the specifications of your equipment.

If the equipment does not appear in the list, contact your supplier. This table is indicative, it is recommended to verify the equipment before executing the (RT) -qPCR.

| Low Profile Cyclers (0,1ml) | |
|-----------------------------|--|
| Manufacturer | Model |
| Agilent Technologies | AriaMx/AriaDx Real-Time PCR System |
| Applied Biosystems | 7500 Fast / 7500 Fast Dx Real-Time PCR System ^{(1) (6)} |
| | QuantStudio™ 12K Flex 96-well Fast |
| | QuantStudio™ 6 Flex 96-well Fast |
| | QuantStudio™ 7 Flex 96-well Fast |
| | QuantStudio™ 3 Fast Real-Time PCR System (3) |
| | QuantStudio™ 5 Fast/ QuantStudio™ 5 Real-Time PCR System |
| | StepOne Plus™ Real-Time PCR System (3) |
| | StepOne™ ⁽²⁾ |
| Azure Biosystems | ViiA™ 7 Fast |
| | Azure Cielo 3 ⁽⁴⁾ |
| | Azure Cielo 6 |
| BIONEER | Exicycler™ 96 Fast |
| Bio-Rad | CFX96TM / CFX96TM IVD Real-Time PCR Detection System |
| | Mini Opticon™ Real-Time PCR Detection System ⁽⁴⁾ |
| Roche | LightCycler®480 Real-Time PCR System ^{(6) (7)} |
| | LightCycler®96 Real-Time PCR System ⁽⁶⁾ |
| | Cobas z480 Analyzer ⁽⁶⁾⁽⁷⁾ |

| Special Formats ⁽⁴⁾ | |
|--|--------------------------|
| Manufacturer | Model |
| Bio Molecular Systems | Mic Real Time PCR Cycler |
| Cepheid | SmartCycler® |
| Qiagen | Rotor-Gene® Q |
| Precision System Science Cp., Ltd. (PSS) | geneLEAD VIII System |

| High Profile Cyclers (0,2ml) | |
|------------------------------|--|
| Manufacturer | Model |
| Abbott | Abbott m2000 ⁽⁶⁾ |
| Agilent | Mx3000P™/ Mx 3005P™ |
| Analytik Jena | qTower |
| | 7300 ^{(3) (6)} |
| | 7500 ⁽⁶⁾ |
| | 7900 ⁽³⁾ |
| | ABI PRISM 7000 ⁽³⁾ |
| | ABI PRISM 7700 ⁽³⁾ |
| | QuantStudio™ 12K Flex 96-well |
| | QuantStudio™ 6 Flex 96-well |
| | QuantStudio™ 7 Flex 96-well |
| | QuantStudio™ 3 Real-Time PCR System ⁽³⁾ |
| Applied Biosystems | QuantStudio™ 5 Fast/ QuantStudio™ 5 Real-Time PCR System |
| | ViiA™ 7 Real-Time PCR System |
| | BIOER |
| BIONEER | Exicycler™ 96 |
| Bio-Rad | CFX96TM Deep Well / CFX96TM Deep Well IVD |
| | iCycler iQTM Real-Time PCR Detection System |
| | iCycler iQTM5 Real-Time PCR Detection System |
| | My iQTM Real-Time PCR Detection System ⁽⁴⁾ |
| | My iQTM2 Real-Time PCR Detection System ⁽⁴⁾ |
| DNA-Technology | DTprime |
| | DTlite |
| Eppendorf | Mastercycler™ ep realplex |
| Qiagen | QIAquant 96 |
| VIASURE | V-Lab96 |

(1) Select Ramp Speed "Standard".

(2) No Cy5 caption.

(3) No ROX caption.

(4) Only FAM and HEX caption.

(5) The product must be reconstituted following the appropriate procedure (see Test procedure) and transferred to the specific tubes for Mic, SmartCycler®, Rotor-Gene® Q or geneLEAD VIII System.

(6) A special grid is needed to fit these real-time PCR kits.

(7) Specific compensation color is required.

Detection channels for the most common real-time PCR systems

| Cyclers | Viasure Channel | Channel Detection | Comments |
|---|-----------------|-------------------|--|
| Bio-Rad CFX96™ | FAM | FAM | Some wells may show abnormal RFU values during the first few cycles of a run with a non-sigmoid ascending line. If so, go to Settings menu, select the Apply Fluorescence Drift Correction for the baseline settings to correct it. |
| | HEX | HEX | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| ABI 7500 Applied Biosystems | FAM | FAM | The passive reference option for ROX should be "none". Some wells may have abnormal RFU values during the first few cycles of a run that show a non-sigmoid ascending line. If you see this effect, modify the baseline by selecting the Cycle Start and Cycle End values so that the baseline ends before significant fluorescence is detected. |
| | HEX | VIC | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| Lightcycler®480II Roche | FAM | 465/510 | For this equipment, a specific color compensation is required with the VS-CCK kit. |
| | HEX | 533/580 | |
| | ROX | 533/610 | |
| | Cy5 | 618/660 | |
| Cobas z 480 Roche | FAM | 465/510 | For this equipment, a specific color compensation is required with the VS-CCK kit. |
| | HEX | 540/580 | |
| | ROX | 540/610 | |
| | Cy5 | 610/670 | |
| Smartcycler® Cepheid | FAM | Channel 1 | |
| | HEX | Channel 2 | |
| | ROX | Channel 3 | |
| | Cy5 | Channel 4 | |
| Abbott m2000rt | FAM | FAM | |
| | HEX | VIC | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| Mx3000PTM Mx 3005PTM Agilent Technologies | FAM | FAM | The passive reference option for ROX should be "none". |
| | HEX | HEX | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| AriaMx Agilent | FAM | FAM | |
| | HEX | HEX | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| Rotor-Gene®Q Qiagen | FAM | Green | In channel settings, click the "Gain Optimization" button and then go to "Optimize Acquiring". The fluorescence target sample range should be between 5 and 10 FI for each channel. Also select the "Perform Optimization Before 1st Acquisition" option. |
| | HEX | Yellow | |
| | ROX | Orange | |
| | Cy5 | Red | |
| Mic Real Time PCR Cycler bms | FAM | Green | In "Run Profile" menu, introduce correct parameters "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and thermal protocol. Go to "Cycling", select "Acquire on" option for all the channels. Use gain default values "Gain" for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10) |
| | HEX | Yellow | |
| | ROX | Orange | |
| | Cy5 | Red | |
| Exicycler™ 96 BIONEER | FAM | FAM | |
| | HEX | JOE | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |