
**CerTest**  
BIOTEC

**CerTest Turbilatex**

**PLUS CRP TURBILATEX® Combo**

**PLUS CRP latex turbidimetric assay**

Cat. TL-022CRP200ED, TL-022CRP800ED

**IU-TL -022CRP200ED-800ED en-es v.0522**

## ENGLISH

### INTENDED USE

---

PLUS CRP Turbilatex® is a latex turbidimetric assay **only for the quantitative detection of C-reactive protein in human serum samples** (not to be used for body fluid as whole blood or plasma). This product may be used to quantify samples of High Sensitivity CRP and C-Reactive Protein (CRP).

This assay is indicated to evaluate the amount of C-reactive protein in serum samples.

This assay is simple and widely applicable. This product is optimized for several automated analyser.

For professional *in vitro* diagnostic use only.

### SYNTHESIS

---

C-reactive protein (CRP) is a suitable marker of inflammation process thanks to its long half-life, independence from age, long stability, and the availability of a wide range of number of diagnostic methods available, is involved in innate immunity as the first defence system against infections.

Cardiovascular accidents are the main cause of death in developed countries. It has been established that inflammation is a fundamental process in the pathogenesis of arteriosclerosis in all its stages, initiation, growth, and complication of arteriosclerotic plaque, considered an inflammatory response to the damage.

HS-CRP is a very promising novel biochemical marker for the prediction of future or recurrent coronary events. American and European prospective studies have been highly consistent regarding the ability of HS-CRP to predict future Coronary Heart Disease (CHD) risk in both men and women.

In acute inflammatory processes caused by bacterial infections, postoperative processes or tissue injury, CRP increases up to 500 mg/L after 6 hours and peaks after 48 hours.

CRP determination is a useful test for the detection of acute infections as well as the monitoring of inflammatory processes, such as acute rheumatic diseases and gastrointestinal diseases. The determination of CRP has more advantages than erythrocyte sedimentation rate (ESR) and leukocyte count. It is more sensitive; the rise occurs earlier, and concentrations decrease more rapidly to baseline values after healing.

### PRINCIPLE

---

PLUS CRP Turbilatex® is based on the principle of a quantitative immunoturbidimetric assay for the quantitative detection of C-Reactive Protein from serum samples.

The intended use of the test is **exclusively to determine the amount of C-Reactive Protein in serum samples.**

PLUS CRP Turbilatex® is based on agglutination reactions. These involve *in vitro* aggregation of microscopic latex particles. This aggregation consists in the specific reaction between antigen and antibodies. Polystyrene latex particles are coated with specific IgG antibodies against CRP. The sample is mixed with a suspension containing anti-CRP antibodies bound to latex particles, if CRP is present in the sample, it will react and form an aggregate. If no antigen is present in the sample the mixture will keep its appearance as a smooth suspension. Such turbidity is measured as an increase in absorbance at determinate wavelength and is proportional to the quantity of CRP contained in the sample.

## PRECAUTIONS

---

- For professional *in vitro* diagnostic use only.
- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in Turbidimetric technique, autoanalyzer use and *in vitro* diagnostic procedures.
- The analyser should be ready before performing any assay.
- Read the instructions for use carefully before use the kit.
- Do not use after expiration date.
- Do not use the reagents if pack is damaged or opened.
- Components provided in the kit are approved for use in PLUS CRP Turbilatex®. Do not use with other commercial kit components (from another supplier).
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink, or smoke in the working area.
- All the specimens should be considered potentially hazardous as well as all the materials that have been exposed to the samples and they must be handled according to the local or national safety regulations. Must be handled in the same manner as an infectious agent. Use proper infection control practices. These practices should include, but are not limited to, personal protective equipment (PPE), such as laboratory coat, surgical or appropriate mask, or face shield, disposable gloves, and eye protection. Take necessary precautions during the collection, transport, storage, handling, and disposal of the samples.
- The solutions should be discarded in a proper container after testing following local regulations.

## STORAGE AND STABILITY

---

Kit components must be stored at refrigerated temperature (2-8°C). Do not freeze.

Reagents are stable up to the expiration date printed on the label, always considering that reagent containers must be properly closed to avoid any contamination, kept away from the sunlight, and conserved at 2-8°C.

## REAGENTS

---

### Materials provided

<u>References</u>	<u>Cat.:</u> <u>TL-022CR200ED</u>	<u>Cat.:</u> <u>TL-022CR800ED</u>
- TL-022CR01 CRP Turbilatex® Reagent 1	2 x 22mL	8 x 22mL
- TL-022CR02/1 CRP Turbilatex® Reagent 2	1 x 12mL	
- TL-022CR02 Turbilatex® Reagent 2		2 x 22mL
- TL-022CR70/ TL-022CR71/ TL-022CR72 / TL-022CR73 CRP Turbilatex® Calibrators		4 x 0,30mL
- TL-022CR08 CRP Turbilatex® Control 1		1 x 0,5mL
- TL-022CR09 CRP Turbilatex® Control 2		1 x 0,5mL

### Materials and instruments that could be required but no provided.

- Specimen collection container tube.
- Disposable gloves and laboratory equipment.
- Automated analyser.
- Extraction needle for sample collection.

## SPECIMEN COLLECTION AND PREPARATION

---

Collect enough quantity of human blood samples. These samples should be collected in clean and dry normal extraction tubes (no preservatives or additives). The samples must be centrifugated to remove blood cells and plasma and get the serum. Serum samples can be directly analysed or stored in the refrigerator (2-8°C) for 7 days prior to testing.

Extraction tube can be directly introduced in the analyser after centrifugation.

## ANALYTICAL PROCEDURE

---

R1 and R2 are ready to use.

### Preparation of the calibration curve


For calibration use only the following materials:

TL-022CR70, TL-022CR71, TL-022CR72, TL-022CR73, C-Reactive Protein Turbilatex® Calibrator vials: liquid calibrator containing CRP antigen at different concentrations.

Concentration is indicated on the label of each vial.

Calibrate the system at least once a month is extremely recommended. Recalibrate the system when reagent lot is changed or when the controls are out of the assigned range given on the control label or certificate of analysis provided with the kit.

C-Reactive Protein Turbilatex® Calibrator vials (calibrator 0 to 3) are ready to use.

Picture 1. Prepared calibrators for performing the calibration curve (2 points and blank)

## Quality control

For quality control only use the following materials:

TL-022CR08/ TL-022CR09 PLUS CRP Turbilatex® Control 1 and 2: liquid controls. Concentration is indicated on the label of the vial. The use of control material is recommended to verify test precision.

Controls should be assayed every day before running patient serum samples extract to validate the calibration curve.

If the obtained results are out of the tolerance range, the equipment, the reagents, or the technique, must be reviewed.

PLUS CRP Turbilatex® Control vials are ready to use.

**Analytical procedure, see annex information.**


## INTERPRETATION OF RESULTS

Positive results: higher or equal than the cut-off fixed by the clinical lab.

Recommended: 3 µg of C-Protein reactive /ml of serum for diagnostic procedures and 10 µg C-Reactive Protein/ml serum for CRP diagnostic procedures.

Positive results determine the abnormal presence of C-Reactive Protein in human serum samples.

### Control vial

Picture 2. Prepared controls 1 and 2 (check control concentration on the label)

HS-CRP Turbilatex® cut-off value:

C-Protein Reactive concentration values lower than 3 µg of C-Protein reactive/ml of serum are considered normal values and that is not indicative of an inflammation and values higher than 3 µg of C-Protein reactive/ml of serum are indicative of a severe inflammation.

CRP Turbilatex® cut-off value:

C-Reactive Protein concentration values less than 10 µg C-Reactive Protein/ml serum are considered normal and are not indicative of acute inflammation or infection and values above this value of 10 µg C-Reactive Protein /ml serum is indicative of the onset of an inflammatory reaction or infection.

## PERFORMANCE CHARACTERISTICS

### Analytical sensitivity

Limit of Quantification (LOQ): Values below 1.0 µg/mL give in poorly reproducible results.

### Prozone

Studies have been made up to a concentration of 640 µg /mL of serum and no false negative results have been observed. Studies using higher concentrations have not been carried out.

Data obtained by the analyser Biolis 24i (Tokyo Boeki).

### Clinical sensitivity and specificity

An evaluation was performed comparing a turbidimetric assay (HS-CRP Turbilatex®, CerTest) and another commercial test (CRP-Ultra Spinreact). The results were as follows:

#### HS-CRP Turbilatex®

	CRP-Ultra Spinreact			Total
		+	-	
HS-CRP Turbilatex®	+	58	3	61
	-	0	29	29
Total	58	32	90	

HS-CRP Turbilatex® vs CRP-Ultra Spinreact		
	Mean Value	95% confidence interval
<b>Sensitivity</b>	100.0%	93.8-100.0%
<b>Specificity</b>	90.6%	75.-98.0%
<b>PPV</b>	95.1%	86.3-99.0%
<b>NPV</b>	100.0%	88.1-100.0%

**CRP Turbilatex®**

		CRP-Ultra Spinreact		
		+	-	Total
CRP Turbilatex®	+	52	1	53
	-	4	43	47
	Total	56	44	100

CRP Turbilatex® vs CRP-Ultra Spinreact		
	Mean Value	95% confidence interval
<b>Sensitivity</b>	92.9%	82.7-98.0
<b>Specificity</b>	97.7%	88.0-99.9%
<b>PPV</b>	98.1%	89.9-100.0%
<b>NPV</b>	91.5%	79.6-97.6%

The results showed a high sensitivity and specificity to detect C-Reactive Protein using HS-CRP Turbilatex® and C-Reactive Protein (CRP) Turbilatex®.

## Interferences

An evaluation was performed to determine the interferences of PLUS CRP Turbilatex®. No interferences were founded against other substances tested: BSA, Ascorbic acid, Sucrose, Octanoic acid and Butyric acid.

## Cross reactivity

An evaluation was performed to determine the cross reactivity; no cross reactivity was founded against the antigens tested, such as: paracetamol, human transferrin, human haemoglobin, urea and bilirubin.

## LIMITATIONS

1. PLUS CRP Turbilatex® should be used only with serum samples.
2. Positive results determine the presence of C-reactive protein in serum samples. The clinical diagnosis should not be based on the results of a single test as both clinical and laboratory data should be evaluated and reviewed by the specialist.

## EXPECTED VALUES

### High Sensitivity (HS-CRP)

Concentration values below 3 µg/ml are considered low cardiovascular risk.

Concentrations of 3 to 5 µg/ml are considered medium risk.

Concentrations above 5 µg/ml are considered high cardiovascular risk.

### **C-Reactive Protein**

Concentrations below 10.0 µg/ml indicate mild inflammation, mild infections such as gingivitis, flu or cold.

Concentrations of 10.0 to 40.0 µg/ml generally indicate more severe infections and moderate infections, such as chicken pox, COVID-19 or other respiratory infection.

Concentrations above 40 µg/ml generally indicate bacterial infection.

Concentrations above 200.0 µg/ml may indicate septicaemia, a life-threatening situation.

**\* Each laboratory should establish its own reference range.**

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1. Libby P, Ridker PM. Novel inflammatory markers of coronary risk. Theory versus practice. *Circulation* 2003; 107: 499-511
2. Ridker PM. High-sensitivity C-reactive protein: Potential for global risk assessment in the primary prevention of cardiovascular disease. *Circulation* 2001; 103: 1813-18
3. Kindmark C-O. The concentration of C-Reactive Protein in sera from healthy individuals. *Scand J Clin Lab Invest* 1972; 29: 407-411.
4. B. Barceló, I.Llompарт, E. Fuentespina. Comparación de tres métodos automatizaos para la medida con gran detectabilidad de la proteína C reactiva. *Química Clínica* 2005; 24 (6) 468-473.
5. Luis Alonso González Naranjo, José Fernando Molina Restrepo. Evaluación de la inflamación en el laboratorio. *Revista Colombiana de Reumatología*. Vol.17 N° 1, Marzo 2010, pp. 35-47.
6. Nader Rifai and Paul M. Ridker. High-Sensitivity C-Reactive Protein: A Novel and Promising Marker of Coronary Heart Disease. *Clinical Chemistry* 47:3 403-411 (2001)
7. Oliver Chenillot, Joseph Henny, Josiane Steinmetz, Bernard Herbeth, Carola Wagner and Gérard Siest. High Sensitivity C-Reactive Protein: Biological Variations and References Limits. *Clin Chem Lab Med* 38(10): 1003 (2000).

## SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO *IN VITRO*

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The application notes for use our turbidimetric products with the different equipments could be consulted in our website:  
**[www.certest.es](http://www.certest.es)**

Note/Nota: Turbilatex® registered brand name/ Turbilatex® marca registrada.  
Nº: CTM-010179299; community (OAMI)/ N.º: CTM-010179299; comunitaria (OAMI)  
Nº: 4,650,387; EEUU (UPSTO)

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