CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test)
Reference kits: SC8010-SCE : 1 Test/Kit, SC8030-SCE : 3 Test/Kit
SC8200-SCE: 20Test/kit, SC8500-SCE: 50 Test/kit
Visual Reading of the results

INTENDED USE

CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test) is a test for self-testing for lay users. CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test) is a coloured chromatographic immunoassay (one step assay) for the qualitative detection of nucleoprotein (N) antigen of SARS-CoV-2 from human nasal swab samples from patients suspected of COVID-19 infection.

CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test) offers a simple and a highly sensitive screening assay as an aid in the diagnosis of SARS-CoV-2. The result of this test should not be the only data for the diagnosis and should be followed by others clinical information and even by additional confirmatory techniques.

The test is for *in vitro* diagnostic use only. This test can be used for non-prescription home use (self-testing) with nasal swab samples.

Read carefully the instructions for use provided into the kit.

INTRODUCTION

An unseen pneumonia spread swiftly in Wuhan, China, in December 2019. Scientists sequenced and recognized a novel β -coronavirus. The disease, caused by the Severe Acute Respiratory Syndrome Corona Virus 2 (SARS-CoV-2), is called as Coronavirus disease 2019 (COVID-19).

Clinical forms of this disease comprise from mild to very severe symptoms: pneumonia, fever and respiratory symptoms are the most frequent. Other symptoms of viral infection include a sore throat, aches, pain, and shortness of breath. In some cases, runny nose, nausea, and diarrhoea also occur.

The outbreak of newly discovered coronavirus begins mainly via discharge from the nose or droplets of saliva, once an infected individual sneezes or coughs, and the estimated time for incubation is within 2 weeks.

TEST PRINCIPLE

CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test) is based on the principle of a qualitative immunochromatographic assay for the qualitative detection of nucleoprotein (N) antigen of SARS-CoV-2 from human nasal swab samples from patients suspected of COVID-19 infection.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against SARS-CoV-2 and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-SARS-CoV-2) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-SARS-CoV-2 monoclonal antibodies-red polystyrene microspheres), which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-SARS-CoV-2 antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no SARS-CoV-2 antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-SARS-CoV-2 antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appears. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For *in vitro* diagnostic use only.
- The test has been designed for lay users. Children or very elder people should be assisted by an adult.
- Read carefully and completely the instructions for use before performing the test.
- Do not use after expiration date. Do not use if pouch (aluminium bag) is damage or the seal of the pouch is broken. Do not use the test if the silica gel bag is not included into the pouch.
- The test device, swab and sample extraction vial are single use. Do not reuse.
- Keep test kit and materials out of the reach of children and pets before and after use.
- All the specimens should be considered potentially hazardous (especially samples from people suspected of SARS-CoV-2 infection) and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.

- Do not eat or smoke while handling specimens.
- Use of personal protection materials such as gloves are recommended.

- The tests (swabs, devices, vials) should be discarded in urban waste container after testing in accordance with the applicable local regulations.

- Clean up spills thoroughly using an appropriate disinfectant.

- Before and after testing, wash or disinfect your hands properly.

- Reagent (sample extraction vial) contain preservatives (<0.1% sodium azide). Do not ingest and avoid any contact with skin or mucous membrane. In case of accidental contact rinse the affected area with clear water. If irritation persists go to a doctor. Safety data sheets are available on request.

- Components provide in the kit are approved for use in the CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test). Do not use any other commercial kit component or components of the test from other lots. Use the material provided into the kit.

- Sterile swabs provided in the kit should be only used for the nasal sample collection. They cannot be reuse.

- Do not touch the head of the sterile swabs provided when opening the primary packaging to avoid contamination.

- Do not dip the swab into buffer (sample extraction vial) or other liquid before inserting the swab into the nose. Open the primary container of the swab and take the sample.

- Only use the liquid from the sample extraction vial provided in the kit. The use of other liquids will conduct to inaccurate results.

- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.

- All positive results should be processed following local or national laws and regulations.

- The way to interpret the results is through the visual reading of some coloured lines (green / red) that will appear after performing the test (depending on the obtained results), the interpretation must be made by people who do not have problems with the correct identification of these colours.

MATERIAL PROVIDED. PRODUCT REFERENCE: SC8010-SCE

- CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test) (1 test)
- Sample extraction vial with sample diluent buffer (1 vial)
- Sterile swabs (1 swab)*
- Instructions for use (IFU) + Quick Guide (1 IFU + 1 Quick Guide)

MATERIAL PROVIDED. PRODUCT REFERENCE: SC8030-SCE

- CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test) (3 test)
- Sample extraction vial with sample diluent buffer (3 vials)
- Sterile swabs (3 swabs)*
- Instructions for use (IFU) + Quick Guide (1 IFU + 1 Quick Guide)

MATERIAL PROVIDED. REF. SC8200-SCE

- CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test) (20 test)
- Sample extraction vial with sample diluent (20 vials)
- Sterile swabs (20 swabs)*
- Plastic bags (20 units)
- Instructions for use (IFU)+ Quick Guide (1IFU +20 Quick Guide)

MATERIAL PROVIDED. REF. SC8500-SCE

- CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test) (50 test)
- Sample extraction vial with sample diluent (50 vials)
- Sterile swabs (50 swab)*
- Plastic bags (50 units)
- Instructions for use (IFU)+ Quick Guide (1IFU +50 Quick Guide)

		Jiangsu Changfeng Medical Industry Co., Ltd Touqiao Town, Guangling District, Yangzhou, Jiangsu 225109 China		Ulms Service & Consulting GmbH Obere Seengasse 34/2, 69124 Heidelberg, Germany. Email: info@ulms-service.com
--	--	---	--	--

MATERIAL REQUIRED BUT NO PROVIDED

- Timer/watch

PREPARATION

Prepare a clean, clear area with adequate lighting to perform the test.

The kit components (tests, vials with sample diluent, sterile swabs) must be conditioned to room temperature before use.

Before using the product, wash your hands carefully and/or use disinfectant gel.

Review and identify the material supplied in the kit, ref: SC8(010/030/200/500)-SCE (see section MATERIALS PROVIDED) and prepare the non-supplied but necessary (timer/ watch).

Read carefully and completely the instructions for use (IFU) and quick guide.


TEST PROCEDURE

Step 1. Take the nasal sample with the sterile swab following the instructions indicated in the section SPECIMEN COLLECTION.

SPECIMEN COLLECTION

Nasal swab (see illustration):

- Before using the product, wash your hands carefully and/or use disinfectant gel. Make sure the hands are clean and dry before starting the procedure.
- Remove the sterile swab from its packing. To open the packing, follow the instructions indicated on it, open in the area enabled for it and never touch the cottony part of the swab (this part is only for collecting the specimens).
- Insert the swab into the nostril (approx. 2cm), rotating against the nasal wall several times (**5 times**) to ensure swab contains cells as well as mucus. Do not apply any pressure. Remove the swab from the nostril carefully and repeat the procedure with the other nostril.
- Proceed immediately with the next section. (Step 2). Samples should be taken and tested immediately.


Specimen collection may feel uncomfortable. Do not insert the swab any deeper if you feel strong resistance or pain.

If you are swabbing others please wear a face mask.

For children (for very young children) probably do not need or cannot insert the swab far into the nostril. Do not continue the test if child feel any pain.

Failure when swabbing properly may cause false negative result.

Step 2. Take the sample extraction vial, open the blue cap of the sample collection vial (prefilled with sample diluent) and immediately put the swab with the sample into the sample extraction vial. Mix and homogenize the solution by rotating the swab forcefully against the side of the extraction vial to allow the solution to mix for **15 seconds**. Best results are obtained when it is mixed to extract as much liquid as possible from the swab avoiding splashes and aerosols. Rotate the swab against the side of the sample extraction vial as the swab is withdrawn and squeeze the vial with your fingers to extract as much liquid as possible. Discard the swab. Close the extraction vial with its blue cap, please be sure that both caps (blue and white-translucent caps) are properly closed.

If the sample diluent is in contact with the skin, wash the affected area with abundant water. In case of irritation contact your doctor. Product safety data sheets are available upon request.

Step 3 . Remove the CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test) from its sealed bag just before using it. For opening the aluminum bag, take advantage of any of the two slits of the bag. Once opened place the test on a horizontal and well-lit surface.


Image of the aluminum bag containing the test (card) (in the right side of the bag can see the two slits for aperture of the bag).

Image of the test (card).

A silica gel bag is included inside the primary container (pouch) next to the card (test). It basically allows to guarantee the correct maintenance of the product until the opening of the pouch. The silica gel should not be used for anything else, it can be thrown away after the pouch squeeze. Do not swallow.

Step 4. Take vial with the sample extracted and open the white-translucent cap and add **3 full drops** (approx. 100-120µL) in the circular window (S). After adding the drops into the sample window, start the timer and close the white-translucent cap.

You can continue with the test even if you accidentally add 4 drops onto the circular window (S)

Step 5. Wait 10 minutes (use the timer/watch) after adding the drops into the sample window (S), leave the card it in a horizontal position without moving during this period. **Read the results at 10 minutes (see section INTERPRETATION OF RESULTS).** Do not read the test result later than 10 minutes. Any result read after more than 10 minutes will be considered invalid.

Step 6. Discard the material used that has been in contact with the sample in urban waste container . Other materials, as IFU, kit box can be recycled in the proper container. Follow applicable local regulations.

The Sample diluent buffer provided into the vials has SARS-CoV-2 inactivation capability (99.98% after 10 minutes). Anyway, please, handle the sample carefully to avoid possible infection of other people.

INTERPRETATION OF RESULTS (please refer to the illustration below)

SARS-CoV-2	INTERPRETATION OF THE RESULTS
1. POSITIVE GREEN-RED 	Read the results at 10 minutes: POSITIVE: In addition to the GREEN control line (C), a RED line also appears in the site marked with the letter T (test line) in the results window. There is SARS-CoV-2 antigen presence detected in the sample. The result would be considered as positive independent on the intensity of the RED test line. Any result read after more than 10 minutes will be considered invalid. What to do now? Isolate yourself immediately, inform your doctor about the situation and result of the test and he will indicate you the steps to follow. To confirm the results, a qPCR test should be performed. There is a small possibility of the test to give a positive result that was a wrong result (false positive result).
2. NEGATIVE GREEN 	Read the results at 10 minutes: NEGATIVE: Only one GREEN line appears across the results window in the site marked with the letter C (control line). This line must be ALWAYS appears and indicates the test is working correctly (sample volume and followed procedure is adequate) There is no SARS-CoV-2 antigen presence detected in the sample. The intensity of the GREEN control line could vary, but band must be clear. Any result read after more than 10 minutes will be considered invalid. What to do now? Very low possibility that it is currently infectious, continue with the hygiene and prevention measures recommended by the health authorities. Consider that the result of the test is valid for the time you perform the test, but it does not guarantee the result or immunity after it. In case of doubt, after other possible exposure to the virus (SARS-CoV-2) after performing the test it would be necessary to repeat the test. You should contact your doctor in case you develop symptoms or become more severe. Could be possible for this test to give a negative result that was an incorrect result (false negative) in some cases of people with COVID-19 infection. If this is your case, the healthcare professional will consider the results with other aspects (symptoms, duration of the illness, exposures to COVID-19 and others considered aspects) and inform you of the following steps. Until you have an answer for your healthcare professional (doctor), you should isolate yourself. You could repeat the test after 24 or 48 hours, COVID-19 infection could not be detected with precision in the first phases of the infection.
3. INVALID Any other result 	Read the results at 10 minutes: INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor. Invalid result, we recommend repeating the assay using the same sample with another test. Any result read after more than 10 minutes will be considered invalid. What to do now? If the result is not meaningful or invalid, the test must be repeated with a new vial extraction and a new sample. Please, make sure of correctly reading the instructions for use and follow all the steps. If your test result is still invalid, please contact with your doctor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigen presents in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

LIMITATIONS

- The test must be carried immediately after opening the sealed bag.
- CerTest SARS-CoV-2 (Ag Nasal sample Test-Self Testing) should be used only with human nasal swab samples. Follow the instructions indicated in the section SPECIMEN COLLECTION. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper nasal specimens must be obtained. False negative results could be obtained if sample was collected incorrectly.
- Samples collected after 6 days from the onset of symptoms may conduct to false negative results.
- The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test. A red coloured band appearing in the test line (T) in the results window with a green line appearing in the control line (C) in the results window, the results must be considered as a positive.
- Positive results determine the presence of SARS-CoV-2 antigens. A diagnosis of SARS-CoV-2 infection should only be made by a physician after all clinical and labora-

tory findings have been evaluated and must be based on the correlation of the results with further clinical observations.

- Positive results do not dismiss co-infections with other pathogens.
- A negative result is not meaningful because it is possible the antigen concentration in the nasal samples was lower than the detection limit value. If the symptoms or situation persist, it is recommended that all negative results undergo confirmatory testing using other method and/or virus identification by cell culture, PCR or other technique. Please contact with your doctor to inform of the results and follow the doctor's indications.
- CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test) has been used with nasal samples from asymptomatic patients, detecting those with high viral loads.
- The result of the test is valid for the time you perform the test, and it must be repeated in case of doubts after exposure to situation of possible infection.
- The immune response cannot be evaluated using this test. Other test methods are required for this purpose.
- Due to very low cross reactivity with high concentrations of SARS-CoV and MERS a false positive result may occur in the case of infection with SARS-CoV or MERS (very low possibility).
- To obtain accurate results, do not use bloody nasal samples.
- The nose should be wiped before the specimen collection.
- Very viscous nasal samples may cause non-specific reaction in the test.
- Information about cross reactivity and interference studies are available under request.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity
CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test) has been shown in clinical evaluations, performed by health care professionals (990 samples were evaluated).

Results:
Sensitivity: 82.7% of the positive nasal samples correctly identify with a confidence interval (75.8-88.3%).
Specificity: 99.2% of the negatives nasal samples correctly identify with a confidence interval (98.3%-99.7%).
 At higher virus concentrations, which corresponds with low Ct (values of Ct \leq 30), CerTest SARS-CoV-2 showed a sensitivity of 92.2% (confidence interval: 85.8% to 96.4%) and a specificity of 99.2 (confidence interval: 98.3%-99.7%) detecting SARS-CoV-2
 A feasibility study (with lay users) has also been performed, being the conclusions of this study:
 98.8% (173/175) of lay users carried out the test without requiring assistance.
 99.4% (174/175) of the different types of the results (positive, negative or invalid) were interpreted correctly.

Cross reactivity and microbial interference

52 different microorganism, possible cross-reactants, and microbial interference, as well as a human nasal pooled wash (representative of normal respiratory microbial flora) had been evaluated with no impact on the performance of CerTest SARS-CoV-2. These were the following:

Adenovirus, Astrovirus, Bocavirus, Coronavirus Strain 229E, Coronavirus Strain HKU1, Coronavirus Strain NL63, Coronavirus Strain OC43, Enterovirus, Influenza A, Influenza B, Metapneumovirus human (hMPV), Norovirus GI, Norovirus GII, Parainfluenza virus, Respiratory Syncytial Virus, Rhinovirus, Rotavirus, Bordetella pertussis, Campylobacter jejuni, Chlamydia pneumoniae, Clostridium difficile antigen GDH, Clostridium difficile Toxin A, Clostridium difficile Toxin B, Escherichia coli O:157, Haemophilus influenzae, Helicobacter pylori, Legionella pneumophila, Listeria Monocytogenes, Mycobacterium tuberculosis, Mycoplasma pneumoniae, Pneumocystis jirovecii, Salmonella enteritidis, Salmonella typhi, Salmonella typhimurium, Shigella boydii, Shigella dysenteriae, Shigella flexneri, Shigella sonnei, Staphylococcus aureus, Staphylococcus epidermidis, Streptococcus A, Streptococcus pneumoniae, Streptococcus pneumoniae, Streptococcus pyogenes, Yersinia enterocolitica O:3, Yersinia enterocolitica O:9, Cryptosporidium parvum, Giardia lamblia, Entamoeba histolytica, pig haemoglobin, and bovine haemoglobin.

Interferences

47 potentially exogenous interfering substances had been evaluated with no impact on the performance of CerTest SARS-CoV-2. These were the following:
 Metronidazole, Ampicillin, Oseltamivir, Amantadine, Ribavirin, Codeine (Toseina), Benzocaine (Angileptol), Cloperastine (Flutax), Carbocisteine (Niston mucolitico), Loratadine, Dexchloropheniramine (Polarumine), Ebastine (Ebastel), Acetyl Salicylic (Adiro), Ibuprofen (Espidifen), Paracetamol (Dolocetyl), Metamizole (Nolotil), Prednisone, Omeprazole, Naso GEL, CVS Nasal Spray (Cromolyn), Afrin (Oxymetazoline), CVS Nasal Drops (Phenylephrine), ZICAM, Homeopathic, Loperamide hydrochloride (Fortasec), Heparin (Hibor), Almagato (Almag), Fosfamyacin (Monurof), Acetylcysteine (Fluimucil), Dextketo-profen trometamol (Enantyum), Levofloxacin, Ciprofloxacin, Rifampicin (Rifaldin), Phenoxymethylpenicillin potassium, Ambroxol hydrochloride (Mucosan), Macrogol 3350 (Movicol), Lysine Carbocysteinate (Pectox), Hydroxyzine dihydrochloride, Lorazepam, Amoxicillin, Mercaptopurine , Biotine, Sore Throat Phenol spray, Tobramycin, Mupirocin, Fluticasone Propionate and Chloraseptic (Menthol/Benzocaine).

6 potentially endogenous substances have been also evaluated with neither no impact: Human haemoglobin, human transferrin, human calprotectin, human lactoferrin, mucine and human blood.

In addition, HAMAs (Human anti-mouse antibodies) were evaluated with no impact.

FOR TECHNICAL SUPPORT

Please contact: self.covid19@certest.es

SYMBOLS FOR IVD COMPONENTS AND REAGENTS

	In vitro diagnostic device		Use by		Do not re-sterilize
	Consult instructions for use		Contains sufficient for <->- test		Unique Device Identifier (UDI)
	Do not re-use		Manufacturer		Authorized representative in the European Community
	Keep dry		DIL Sample diluent		Sterilized using ethylene
	Temperature limitation		LOT Batch number		Do not use if package is damaged
	This symbol is used for warning and cautions		REF Catalogue number		CE mark

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MATERIAL SUMINISTRADO. REF. SC8010-SCE

- CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test) (1 test)

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MATERIAL SUMINISTRADO. REF. SC8030-SCE

- CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test) (3 test)

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MATERIAL SUMINISTRADO. REF. SC8200-SCE

-CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test) (20 test)

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MATERIAL SUMINISTRADO. REF. SC8500-SCE

-CerTest SARS-CoV-2 (Ag-Nasal Sample-Self Test) (50 test)

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CE 0197		Jiangsu Changfeng Medical Industry Co., Ltd. Touqian Town, Guangling District, Yangzhou, Jiangsu 225109 China	CE-IMP	Ulms Service & Consulting GmbH Seegasse 34/2, 69124 Heidelberg, Germany. Email: info@ulms-service.com
---------	--	--	--------	---

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