**VIASURE**

**Viral Positive Control Kit**

**SARS-CoV-2 Delta (B.1.617.2)**

**CE IVD**

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Table A 1. References / *Referencias.*

## Content

1. Intended use .....	3
2. Summary and Explanation .....	3
3. Principle of the procedure .....	4
4. Reagents provided .....	4
5. Reagents and equipment to be supplied by the user .....	4
6. Transport and storage conditions .....	5
7. Warnings and Precautions .....	5
8. Test procedure .....	6
8.1 Reconstitution of the Positive Control .....	6
8.2 Nucleic Acid Extraction .....	6
8.3 Real-Time PCR .....	7
9. Result interpretation .....	8
10. Limitations of the test .....	8
11. Quality control .....	9
12. Performance characteristics .....	9

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---

### 1. Intended use

VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit is intended for monitoring the whole reverse transcription quantitative real-time PCR (RT-qPCR) process, from nucleic acid extraction to reverse transcription and amplification steps.

VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit contains non-replicative and non-infectious SARS-CoV-2 Delta variant (B.1.617.2) viral particles that are expected to be processed along with respiratory samples from individuals with clinical suspicion of SARS-CoV-2 infection.

This Kit is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in molecular techniques and *in vitro* diagnostic procedures.

### 2. Summary and Explanation

The SARS-CoV-2 virus was detected for the first time in Wuhan (China) on December 31<sup>st</sup>, 2019. During 2020, the virus spread progressively throughout the world and on March 11<sup>th</sup>, the WHO decreed the state of pandemic.<sup>1</sup>

The appearance of mutations is a natural and expected event within the evolution process of the virus. In fact, some specific mutations define the viral genetic groups that are currently circulating globally. The mutations identified to date remain within the expected patterns for a coronavirus. Thanks to the genetic sequencing of the pathogen worldwide, it has been possible to establish patterns of dispersal and evolution of the virus.<sup>2,3</sup>

Delta variant (B.1.617.2 lineage), firstly reported in India in December 2020, was closely associated with a huge COVID-19 increase in this country during Spring 2021. It has spread over 60 countries due to its ability to invade the host's immune system. Delta variant has 23 mutations compared to the first identified SARS-CoV-2 strain. Twelve of those mutations are in the Spike protein and facilitate the viral fusion.<sup>4</sup> On May 2021, the European Centre of Disease Control and Prevention (ECDC) upgraded B.1.617.2 to a variant of concern (VOC), after new information became available on increased transmissibility and some impact on vaccine effectiveness.<sup>5</sup>

The appearance of variants that increase the transmissibility of the virus, its virulence or that escape the action of the neutralizing antibodies generated after natural infection or the vaccine, constitute a first-order public health problem that can have an important impact on control of the pandemic. Diagnosis of SARS-CoV-2 variants can be problematic. Among the diagnostic methods, the real-time PCR technique has stood out for its sensitivity, specificity, and its ability to process large numbers of samples in a reduced turnaround time.<sup>6,7</sup>

VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit improves the diagnosis of SARS-CoV-2 Beta variant infection by identifying false negative or false positive results that may occur because of a failure during the nucleic acid extraction process or the RT-qPCR reaction. Many factors can lead to invalid results, such as low RNA yields after nucleic acid extraction, low reverse transcription efficiency or cross-contamination.

### 3. Principle of the procedure

VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit contains non-replicative and non-infectious recombinant viral particles in a lyophilized format, which have been created using the lentiviral system<sup>8</sup>. This system is based on the ability of lentivirus to package foreign RNA inside cells. The cells are co-transfected with two different lentiviral vectors: the packaging vector, which encodes for lentiviral structural proteins, lentiviral enzymes and lentiviral gene regulatory proteins, and the transfer vector, which contains the viral genome sequence that is completely customizable. This viral genome is bounded to lentiviral enzymes and enclosed by a lentiviral capsid that is surrounded by a matrix. This, in turn, is surrounded by a lipid bilayer (envelope) taken from the membrane of a cell when a newly formed viral particle buds from the cell. This viral particle is non-infectious as it lacks the envelope glycoprotein that allows the virus to enter the cell and, in addition, it is heat-treated<sup>9</sup>. This viral particle is also non-replicative as it is a defective virus without the ability to produce new viral particles.

VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit contains viral particles of SARS-CoV-2 Delta (B.1.617.2) variant, which contain the whole genome of this variant of the virus.

VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit is designed for monitoring the performance of all RT-qPCR assays whose target regions belong to SARS-CoV-2 Delta variant genome. The RNA sequences are homologous to the SARS-CoV-2 strain SARS-CoV-2/human/JPN/SARS-CoV-2, B.1.617.2 lineage, Delta variant/2021 (GenBank: OK091006.1).

### 4. Reagents provided

VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit includes the following materials and reagents detailed in Table 1.

Reagent/Material	Description	Colour	Amount
Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control	Non-infectious lyophilized viral particle	Red	4/1 vials
Viral Rehydration Buffer	Solution to reconstitute the lyophilized product	Transparent vial with green cap	1 vial x 1 mL
Viral Negative Control	Non template control	Transparent vial with violet cap	4/1 vials x 1 mL

Table 1. Reagents and materials provided in VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit.

### 5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit:

- Real-Time PCR instrument (thermocycler).
- Nucleic acid extraction kit.
- Real-Time PCR assay.
- Ribonuclease (RNase)/ Deoxyribonuclease (DNase) free tubes.
- Vortex.
- Micropipettes.
- Filter tips.
- Powder-free disposable gloves.

## 6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- The **Viral Rehydration Buffer** is intended for **four uses**. It is recommended to separate it into four aliquots (see section 8 "Test procedure") and store it at 2-40°C until the expiration date which is stated on the label.
- Each vial of **Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control** is intended for **one use**. Once this Positive Control has been reconstituted (see section 8 "Test procedure"), it must be subjected to nucleic acid extraction and RT-qPCR. If the reconstituted vial cannot be processed immediately after rehydration, store at 4°C or 25°C for up to 24 hours or freeze at -20°C or -80°C for a longer period (do not refreeze).
- Each vial of **Viral Negative Control** is intended for **one use**. Store at 2-40°C until use.
- It is recommended to separate the RNA extracted from the Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control into aliquots to minimize freeze-thaw cycles. The RNA has been validated to remain stable after 10 freeze-thaw cycles at either -20°C or -80°C.

## 7. Warnings and Precautions

- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and *in vitro* diagnostic procedures (including training on the real-time PCR instrument (thermocycler) and nucleic acid extraction system).
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
- Protect the Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control against from humidity. Prolonged exposure to humidity may affect product performance.
- An appearance of the Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or colour different from whitish) does not alter the functionality of the Positive Control.
- Design a unidirectional workflow. It should begin with the aliquots of the Viral Rehydration Buffer in a clean area, and then move to the extraction area, where the Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control is rehydrated and extracted along with the Viral Negative Control and the samples with clinical suspicion of SARS-CoV-2 infection. Finally, move on to the nucleic acid amplification and detection area. Do not return samples, equipment, and reagents to the area in which the previous step was performed.
- Always avoid microbial and RNase/DNase contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles, and mask. Do not eat, drink, or smoke or apply cosmetic products in the working area. Once you finish the test wash your hands.
- The Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control and the Viral Negative Control should be treated as clinical specimens and they should undergo the same process. Specimens must be treated as potentially infectious and/or biohazardous, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, transport, storage, handling, and disposal of samples.

- The viral particles that constitute the Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control are non-replicative since they are not capable of producing new viral particles. In addition, they are non-infectious, as they lack the envelope glycoprotein, and they are heat-treated. The viral particles confer resistance to puromycin, so uninfected cells are naturally sensitive to this antibiotic, while infected cells are resistant. The non-infectivity of the virus has been proven by the absence of colonies resistant to puromycin after infection with the Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose of waste in compliance with local, state, and federal regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- In accordance with Regulation (EC) No 1907/2006 (REACH), VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit does not require Material Safety Data Sheets on account of their classification as non-hazardous to health and the environment because it does not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP) or which are in concentrations higher than the value established in the mentioned regulation for their declaration.
- Consult each nucleic acid extraction kit, automated extraction system, RT-qPCR kit or real-time PCR instrument's reference manual for additional warnings, precautions, and procedures.

## 8. Test procedure

### 8.1 Reconstitution of the Positive Control

The Viral Rehydration Buffer is intended for four uses, so it is recommended to separate it into **four aliquots (240 µl/each** in RNase/DNase-free tubes) to avoid contamination between uses. This step should be performed in a clean area separate from the Extraction Area and before handling the Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control. The aliquots can be stored at 2-40°C until the expiration date which is stated on the label.

Reconstitute the lyophilized Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control (red vial) by adding **220 µl** of the supplied **Viral Rehydration Buffer** (transparent vial with green cap previously aliquoted) and vortex thoroughly. It is recommended to perform this step in the Extraction Area and immediately processing it along with the Viral Negative Control (transparent vial with violet cap) and the clinical specimens. If the reconstituted vial cannot be processed immediately after rehydration, see section 6 "Transport and storage conditions".

In case the Viral Rehydration Buffer (supplied in the kit) is replaced by another rehydration buffer, this must be validated by the final user.

### 8.2 Nucleic Acid Extraction

The Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control and the Viral Negative Control must be subjected to the same nucleic acid extraction process than clinical samples. Any manual or automatic optimized system from human clinical specimens can be used for nucleic acid extraction, following the manufacturer's instructions or the in-house protocol.

The following commercially available nucleic acid extraction kits included in table 2 have been validated with VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit.

Nucleic acid extraction kit	Equipment
MagDEA® Dx SV (Precision System Science Co.)	MagLEAD® 12gC Instrument
QIAamp® Viral RNA Mini Kit (Qiagen)	N/A
Invisorb® Spin Universal Kit (Invitek)	N/A
Maxwell® RSC Viral Total Nucleic Acid Purification Kit (Promega)	Maxwell® RSC Instrument
MagMAX™ Viral/Pathogen Nucleic Acid Isolation Kit (Thermo Fisher Scientific)	KingFisher Flex System
MagCore® Viral Nucleic Acid Extraction Kit (RBC Bioscience)	MagCore®HF16 automated Nucleic Acid Extractor System
DNA/RNA Extraction Kit (Nanjing Vazyme Medical Technology Co.)	Vazyme VNP-32YL

Table 2. Nucleic acid extraction kits compatible with VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit.

### 8.3 Real-Time PCR

The RNA obtained in the nucleic acid extraction process can be amplified using any in-house or commercially available RT-qPCR kit whose targets belong to SARS-CoV-2 Delta variant genome. Manufacturer's instructions for the Real-Time PCR kits or the in-house protocol must be followed.

These controls must not be substituted for the positive and negative controls provided with the RT-qPCR kits.

VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit has also been validated with the following commercial RT-qPCR kits, following the manufacturer's instructions.

- VIASURE SARS-CoV-2 Variant I Real Time PCR Detection Kit.
- VIASURE *Flu A, Flu B* & SARS-CoV-2 Real Time PCR Detection Kit.
- VIASURE SARS-CoV-2, *Flu* & *RSV* Real Time PCR Detection Kit.
- VIASURE SARS-CoV-2 *S gene* Real Time PCR Detection Kit.
- VIASURE SARS-CoV-2 Real Time PCR Detection Kit.
- VIASURE SARS-CoV-2 (*N1 + N2*) Real Time PCR Detection Kit.
- VIASURE SARS-CoV-2 *Triplex (ORF1ab, E & N genes)* Real Time PCR Detection Kit.
- VIASURE SARS-CoV-2 *del 69/70, ORF1ab & N genes* Real Time PCR Detection Kit.
- VIASURE SARS-CoV-2 (*HV 69/70 del, L452R, ORF1ab & N gene*) Real Time PCR Detection Kit.

VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit must follow the same diagnostic procedure as samples with clinical suspicion of SARS-CoV-2 infection. If the nucleic acid extraction or RT-qPCR conditions differ from those validated here, the procedure must be validated by the final user.

## 9. Result interpretation

The Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control and the Viral Negative Control must be included in each run of nucleic acid extraction and RT-qPCR. The Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control is used to monitor the performance of the whole process, whereas the Viral Negative Control is used to monitor the cross-contamination that may occur during the whole process. Both controls are required to validate the reaction and accurately interpret the results of the tested samples.

A valid Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control result is considered when an amplification curve is observed in the RT-qPCR assay. The Cq value obtained is dependent on the performance of the nucleic acid extraction method and the sensitivity of the RT-qPCR assay used.

A valid Viral Negative result is considered when no amplification curve is observed in the RT-qPCR assay or when an amplification curve is observed at a Cq value greater than the cut-off Cq indicated in the manufacturer's instructions of the RT-qPCR assay used.

If one or more controls are not valid (see Table 3), the sample results cannot be interpreted.

Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Result	Viral Negative Control Result	Interpretation of Controls
Positive	Negative	<b>Valid</b>
Positive	Positive	<b>Test failure *</b>
Negative	Negative	<b>Test failure **</b>
Negative	Positive	<b>Test failure ***</b>

Table 3. Interpretation of the results obtained with VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit.

\*Possible cross-contamination problem.

\*\*Possible failure of nucleic acid extraction or RT-qPCR reaction.

\*\*\*Possible cross-contamination problem, nucleic acid extraction or RT-qPCR failure or exchange of controls.

All samples that have been processed along with a Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control or/and Viral Negative Control with a test failure result, must be retested from the nucleic acid extraction.

## 10. Limitations of the test

- The Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control and the Viral Negative Control must not be substituted for the positive and negative controls provided with the RT-qPCR kits.
- Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control may contain traces of plasmid DNA.
- Changes in the nucleotide sequences of primer or probes to detect new or unknown SARS-CoV-2 variants can lead to a failure in the identification of the Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control.
- Since the entire genome of SARS-CoV-2 Delta variant has been included in the design, it can rarely happen that primer or probes used do not hybridize.

- Performance of VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit depends on the efficiency of the nucleic acid extraction method and the RT-qPCR assay used by each laboratory. Each test procedure must be validated by the final user.
- Reconstitution of the Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control with a reagent different from the Viral Rehydration Buffer (supplied with the kit) must be validated by the final user.
- Test failure results may arise from several factors and their combinations, including:
  - Improper transport, storage, and/or handling of the kit.
  - Improper processing procedures (reconstitution of the positive control, nucleic acid extraction, RT-qPCR assay).
  - Degradation of the viral RNA during storage and/or processing.
  - Failure to follow instructions for use and the assay procedure.

Consult section 9 "Result interpretation" in the present document for more information.

- VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit should not be used as calibrator. The concentration of Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control is obtained by digital PCR following a specific workflow. Deviations from this workflow can lead to concentration variations.

## 11. Quality control

The expected values of VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit depends on the diagnostic procedure (nucleic acid extraction method and RT-qPCR assay), so each test procedure must be validated by the final user.

The concentration of Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control is stated in the Certificate of analysis of each batch, available upon request.

## 12. Performance characteristics

The performance characteristics of VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit have been validated using the nucleic acid extraction systems and RT-qPCR kits mentioned in section 8 of the present document. The results obtained by testing the VIASURE Viral SARS-CoV-2 Delta (B.1.617.2) Positive Control Kit using different workflows, operators, and reagents lots show similar amplification curves, demonstrating good performance, reproducibility, and repeatability.

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

## Bibliography/Bibliografía

1. N, Z. *et al.* A Novel Coronavirus from Patients with Pneumonia in China, 2019. *N. Engl. J. Med.* **382**, 727–733 (2020).
2. Centers of Disease Control and Prevention (CDC). Emerging SARS-CoV-2 Variants | CDC. <https://www.cdc.gov/coronavirus/2019-ncov/science/science-briefs/scientific-brief-emerging-variants.html>.
3. Centers of Disease Control and Prevention (CDC). SARS-CoV-2 Variant Classifications and Definitions. <https://www.cdc.gov/coronavirus/2019-ncov/cases-updates/variant-surveillance/variant-info.html>.
4. Shiehzadegan, S., Alaghemand, N., Fox, M. & Venketaraman, V. Analysis of the Delta Variant B.1.617.2 COVID-19. *Clin. Pract.* 2021, Vol. 11, Pages 778-784 **11**, 778–784 (2021).
5. Threat Assessment Brief: Emergence of SARS-CoV-2 B.1.617 variants in India and situation in the EU/EEA | ECDC |. <https://www.ecdc.europa.eu/en/publications-data/threat-assessment-emergence-sars-cov-2-b1617-variants>.
6. K, W. *et al.* mRNA-1273 vaccine induces neutralizing antibodies against spike mutants from global SARS-CoV-2 variants. *bioRxiv Prepr. Serv. Biol.* (2021) doi:10.1101/2021.01.25.427948.
7. Y, W. *et al.* Escape from neutralizing antibodies by SARS-CoV-2 spike protein variants. *Elife* **9**, 1 (2020).
8. T, S., MA, B. & Y, I. Lentiviral vectors: basic to translational. *Biochem. J.* **443**, 603–618 (2012).
9. F, K., L, S., S, C., L, K. & HM, B. Efficient inactivation of pseudotyped HIV-based lentiviral vectors and infectious HIV. *J. Virol. Methods* **276**, (2020).

## Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico in vitro

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## Trademarks

Modification rights reserved. All rights reserved. © CerTest Biotec, S.L.

All other trademarks that may appear in this package insert are the property of their respective owners.

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