VIASURE

Real Time PCR Detection Kit

Bordetella
for BD MAX™ System

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Note: The user should notify the manufacturer and the competent authority of the Member State in which he is established as a user and/or patient of any serious incident related to the product.

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Content

1.	Intended purpose	5
2.	Summary and Explanation	5
3.	Principle of the procedure	5
4.	Reagents provided	6
5.	Reagents and equipment to be supplied by the user	6
6.	Transport and storage conditions.....	7
7.	Precautions for users	7
8.	Test procedure	8
8.1.	Sample collection, transport, and storage.....	8
8.2.	Sample preparation and DNA extraction	9
8.3.	PCR protocol	9
8.3.1.	Creating PCR test program for VIASURE Bordetella Real Time PCR Detection Kit for BD MAX™ System	9
8.3.2.	BD MAX™ Rack set up	11
8.3.3.	BD MAX™ Instrument set up	12
8.3.4.	BD MAX™ report	12
9.	Result interpretation	13
10.	Limitations of the test	14
11.	Quality control.....	15
12.	Analytical Performance characteristics	16
12.1.	Analytical Linearity	16
12.2.	Analytical sensitivity. Limit of Detection (LoD)	17
12.3.	Measuring range	17
12.4.	Accuracy	17
12.4.1.	Trueness (Veracity).....	17
12.4.2.	Precision.....	18
12.5.	Analytical specificity and reactivity	19
12.5.1.	Analytical Specificity	20
12.5.1.1.	Cross-reactivity and exclusivity assay.....	20
12.5.2.	Analytical reactivity	22
12.6.	Metrological traceability	22
13.	Clinical performance characteristics.....	22

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ENGLISH

1. Intended purpose

VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System is an automated qPCR test designed for the qualitative detection and differentiation of DNA from *Bordetella pertussis*, *Bordetella parapertussis* and/or *Bordetella holmesii* in respiratory samples (nasopharyngeal swabs and nasopharyngeal aspirates) from patients suspected of respiratory infection by their healthcare professional (HCP). This test is intended to be used as an aid in the identification of *Bordetella pertussis*, *Bordetella parapertussis* and/or *Bordetella holmesii*, in combination with patient's clinical signs and symptoms and epidemiological risk factors. The assay uses the BD MAX™ System for automated extraction of DNA and subsequent qPCR, employing the reagents provided combined with universal reagents and disposables for the BD MAX™ System. DNA is extracted from respiratory samples, amplified using real-time PCR, and detected with fluorescent reporter dye probes specific for *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii*.

2. Summary and Explanation

The genus *Bordetella* is comprised of 8 species, 4 of which are known to infect humans: *B. pertussis*, *B. parapertussis*, *B. holmesii*, and *B. bronchiseptica*. The most important cause for whooping cough (pertussis) is *B. pertussis* infection, followed by *B. parapertussis*. *B. holmesii* has been isolated from patients with a serious underlying disease, whereas *B. bronchiseptica* is usually restricted to animals but occasionally has also been isolated from immunocompromised patients.

Whooping cough is a very contagious disease which spreads from person to person usually by coughing or sneezing, or by close contact with an infected person in a common breathing space. The clinical course of the illness is divided into three stages, which include the following clinical features: catarrhal (coryza, low-grade fever, mild and occasional cough), paroxysmal (paroxysms of numerous and rapid coughs, cyanosis, vomiting and exhaustion) and convalescent (gradual recovery and less persistent paroxysmal coughs).

Despite vaccination, pertussis remains endemic in most areas of the world. Reliable diagnosis is required to start appropriate treatment and prophylaxis of contacts if needed, particularly unvaccinated infants in whom pertussis might present as a life-threatening disease. Nucleic acid amplification tests, including PCR and more recently real-time PCR, overcome some of the limitations of culture and serological methods for the diagnosis of *Bordetella* infections. Most of the PCR tests are based on detection of insertion sequences (IS) present in multiple copies per genome, increasing the sensitivity of PCR tests.

3. Principle of the procedure

VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System is designed for the simultaneous qualitative detection of DNA from *Bordetella pertussis*, *Bordetella parapertussis* and/or *Bordetella holmesii* in respiratory samples. After DNA isolation, the identification of these strains of *Bordetella* is performed by the amplification of a conserved region of the insertion sequences IS481 for *Bordetella pertussis/holmesii*, pIS1001 for *Bordetella parapertussis* and hIS1001 for *Bordetella holmesii*, using specific primers and a fluorescent-labelled probe.

VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence is measured on the BD MAX™ System.

VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System contains in each tube all the components necessary for real-time PCR assay (specific primers/probes, dNTPs, buffer, polymerase) in a stabilized format, as well as an **internal control** to verify if the amplification reaction works properly. If the PCR reaction is inhibited due to the presence of an artifact/inhibitor or because the reaction mix is incorrectly rehydrated, this control may become negative.

Target	Channel	Gene
<i>B. pertussis</i> / <i>B. holmesii</i>	475/520	IS481 sequence
<i>B. holmesii</i>	585/630	hIS1001 sequence
<i>B. paraptussis</i>	630/665	pIS1001 sequence
Internal control (IC)	530/565	-

Table 1. Target, channel and genes.

4. Reagents provided

VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System includes the following materials and reagents detailed in Table 2:

Reagent/Material	Description	Barcode	Amount
<i>Bordetella</i> reaction tube	A mix of enzymes, primers-probes, buffer, dNTPs, stabilizers and internal control in stabilized format	1C foil	2 pouches of 12 transparent tubes
Rehydration Buffer tube	Solution to reconstitute the stabilized product	11 foil	1 pouch of 24 transparent tubes

Table 2. Reagents and materials provided in VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System with Cat. N°. VS-BDT124 (444204).

5. Reagents and equipment to be supplied by the user

The following list includes the materials and equipment that are required for use but not included in the VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System.

- Real-time PCR instrument: BD MAX™ System (Ref: 441916).
- BD MAX™ ExK™ TNA-3 (Ref:442827 or 442828).
- BD MAX™ PCR Cartridges (Ref: 437519).
- Vortex.
- Micropipettes (accurate between 2 and 1000 µL).
- Nuclease-free water.
- Filter tips.
- Powder-free disposable gloves.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- After opening the aluminum pouches which contain the reaction tubes, the product can be used up to 28 days. Keep the vial away from light.

The following table summarises the conditions for transport, storage and use of the overall kit and each component:

Component	Transport Conditions	Storage Conditions	In-use conditions
Entire VIASURE <i>Bordetella</i> Real Time PCR Detection Kit for BD MAX™ System	2-40°C during the shelf life stated in the label.	Before use: 2-40°C during the shelf life stated in the label.	* See <i>in-use conditions</i> of each component.
<i>Bordetella</i> reaction tube (1C foil)		Before use: 2-40°C during the shelf life stated in the label. Once pouch is opened with the silica gel: 25°C±5°C for up to 28 days.	Immediate use at Room temperature.

Table 3. Summary of the conditions for transport, storage and use of the VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System and each component.

7. Precautions for users

- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and *in vitro* diagnostic procedures.
- For *in vitro* diagnostic use.
- Instructions for use of the VIASURE product and the BD MAX™ System User's Manual must be read carefully before using the VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System. Do not perform the assay until the information about procedures, safety precautions and limitations described therein have been understood.
- Do not use expired reagents and/or materials.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective box is open or broken upon arrival.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches.
- Close protective pouches of reagents promptly with the zip seal after each use to protect the master mix from sunlight. Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and/or kits and/or lots.
- Protect reagents from humidity. Prolonged exposure to humidity may affect product performance.
- An appearance of the reaction mixture in stabilized format, normally found at the bottom of the tube, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or color different from whitish) does not alter the functionality of the test.
- Make sure reaction tube and rehydration buffer tube are snapped into place securely during the BD MAX™ rack set up.

- In cases where other PCR tests are conducted in the same general area of the laboratory, care must be taken to ensure that the VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System, BD MAX™ ExK™ TNA-3 extraction kit, or any additional reagents required for testing and the BD MAX™ System are not contaminated. Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents and cartridges (BD MAX™ PCR Cartridge).
- To avoid contamination of the environment by amplicons, do not break apart the BD MAX™ PCR Cartridge after use. The seals of the BD MAX™ PCR Cartridge are designed to prevent contamination.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment, and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles, and mask. Do not eat, drink, smoke or apply cosmetic products in the working area. Wash your hands after finishing the test.
- Samples must be treated as potentially infectious and/or biohazardous, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, transport, storage, handling, and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose of waste in compliance with local and state regulations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- In accordance with Regulation (EC) No 1907/2006 (REACH), VIASURE Real Time PCR Detection Kits do not require Material Safety Data Sheets on account of their classification as non-hazardous to health and the environment, because they do not contain substances and/or mixtures which meet the hazard classification criteria available in Regulation (EC) No 1272/2008 (CLP), or which are in concentrations higher than the value established in the mentioned regulation for their declaration.
- Consult the BD MAX™ System User's Manual for additional warnings, precautions and procedures.
- The certificate of analysis is not included with the device; however, it could be downloaded from Certest Biotec S.L. website (www.certest.es) in case of need.

8. Test procedure

8.1. Sample collection, transport, and storage

The VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System has been tested on respiratory samples (nasopharyngeal swabs and aspirates). Other types of samples must be validated by the user.

Collection, storage, and transport of specimens should be maintained per the conditions validated by the user. Overall, respiratory samples should be collected and labelled appropriately in clean containers with or without transport media (depending on sample type) and processed as soon as possible to guarantee the quality of the test. The nasopharyngeal aspirates should be transported at 4°C for up to 72 hours, following the local and national

regulations for the transport of pathogen material. For long term transport (more than 72 hours), we recommend shipping at -20°C or lower. However, nasopharyngeal swabs in Universal Transport media (UTM) should be transported frozen at -20°C or lower. It is recommended to use fresh specimens for the test. The samples can be stored at 4°C for up to 72 hours (nasopharyngeal aspirates) or frozen at -20°C or ideally at -70°C (in UTM) for conservation. Repeated freeze-thaw cycles should be avoided in order to prevent degradation of the sample and nucleic acids.

The respiratory samples must be collected, transported and stored according to appropriate laboratory guidelines. For details, refer to the CDC's Infectious diseases laboratory test directory (2023). Website <https://www.cdc.gov/laboratory/specimen-submission/cdc-lab-tests.pdf>.

Please note: Specimen collection, transport, and storage conditions indicated above are suggested based on the recommendations for respiratory samples collected in transport medium and intended to be used for *Bordetella pertussis*, *Bordetella parapertussis* and/or *Bordetella holmesii* DNA detection, which appear in the referenced CDC's Infectious diseases laboratory test directory. Nevertheless, we recommend following laboratory guidelines, and/or microbiology laboratory policy manual for proper transport and preservation of samples.

8.2. Sample preparation and DNA extraction

Perform the sample preparation according to the recommendations in the instructions for use of extraction kit used, BD MAX™ ExK™ TNA-3.

1. Pipette 200 µL of sample into a BD MAX™ ExK™ TNA-3 Sample Buffer Tube and close the tube with a septum cap. Ensure complete mixing by vortexing the sample at high speed for 1 minute. Proceed to BD MAX™ System Operation.

Note that application-specific extraction preparation procedures should be developed and validated by the user and that some other samples may require pre-processing.

8.3. PCR protocol

Note: Please, refer to the BD MAX™ System User's Manual for detailed instructions.

8.3.1. Creating PCR test program for VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System

Note: If you have already created the test for the VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System, you can skip step 8.3.1 and go directly to 8.3.2.

- 1) On the "Run" screen of the BD MAX™ System, select the "Test Editor" tab.
- 2) Click the "Create" button.
- 3) In the Basic Information tab, within the "Test Name" window, name your test: i.e. VIASURE *Bordetella*.
- 4) In the "Extraction Type" drop down menu, select "ExK TNA-3".
- 5) In the "Master Mix Format" drop down menu, choose "Type 5".

- a. Note: Product may be used in combination with an additional VIASURE for BD MAX™ test, then select “Dual Master Mix Concentrated Lyophilized MM with Rehydration Buffer (Type 5)”.
- 6) In the “Sample extraction parameters” select “User defined” and adjust sample volume to 350 µL.
- 7) In the “Ct Calculation” select “Call Ct at Threshold Crossing”.
- 8) If running software version 5.00 or higher and have barcoded foil snap-in tubes, in the “Custom Barcodes” select the following configuration:
- Snap-In 2 Barcode: 1C (concerning *Bordetella* reaction tube).
 - Snap-In 3 Barcode: 11 (concerning Rehydration Buffer tube).
 - Snap-In 4 Barcode: another VIASURE reaction tube (different foil) if you choose the format “Dual Master Mix Concentrated Lyophilized MM with Rehydration Buffer (Type 5)” (Section 8.3.1).
- 9) In “PCR settings” tab enter the following parameters: “Channel Settings”, “Gains” and “Threshold” (Table 4).
- a. Note: Product may be used in combination with an additional VIASURE for BD MAX™ test, PCR Settings and Test Steps should be completed for snap 2 (green) and snap 4 (blue) positions.

Channel	Alias	Gain	Threshold	Ct Min	Ct Max
475/520 (FAM)	<i>B. pertussis/B. holmesii</i>	50	200	0	40
530/565 (HEX)	IC	80	200	0	35
585/630 (ROX)	<i>B. holmesii</i>	50	200	0	40
630/665 (Cy5)	<i>B. parapertussis</i>	50	200	0	40
680/715 (Cy5.5)	-	-	-	-	-

Table 4. PCR settings.

Note: It is recommended to set the minimum threshold values listed above for each channel as a starting point, but the final settings must be determined by the end-user during the result interpretation, in order to ensure that thresholds fall within the exponential phase of the fluorescence curves and above any background signal. The threshold value for different instruments may vary due to different signal intensities.

- 10) In “PCR settings” tab enter the following parameters “Spectral Cross Talk” (Table 5), as well.

		False Receiving Channel				
	Channel	475/520	530/565	585/630	630/665	680/715
Excitation Channel	475/520	-	0.0	0.0	0.0	0.0
	530/565	1.0	-	0.0	0.0	0.0
	585/630	0.0	0.0	-	2.0	0.0
	630/665	0.0	0.0	3.0	-	0.0
	680/715	0.0	0.0	0.0	0.0	-

Table 5. Spectral cross-talk parameters.

- 11) In “Test Steps” tab, enter the PCR protocol (Table 6).

Step Name	Profile Type	Cycles	Time (s)	Temperature	Detect
Initial denaturation	Hold	1	120	98°C	-
Denaturation and Annealing/Extension (Data collection)	2-Temperature	45	10	95°C	-
			41	63°C	✓


Table 6. PCR protocol.

- 12) Click the "Save Test" button.

8.3.2. BD MAX™ Rack set up

- For each sample to be tested, remove one Unitized Reagent Strips from the BD MAX™ ExK™ TNA-3 kit. Gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes and load on the BD MAX™ System sample racks.
- Remove the required number of BD MAX™ ExK™ TNA Extraction Tubes (B4) (white foil) from their protective pouch. Snap the Extraction Tube(s) (white foil) into its corresponding positions in the TNA strip (Snap position 1, white color coding on the rack. See Figure 1). Remove excess air, and close pouch with the zip seal.
- Determine and separate the appropriate number of *Bordetella* reaction tubes (1C foil) and snap into their corresponding positions in the strip (Snap position 2, green color coding on the rack. See Figure 1).
 - Remove excess air, and close aluminum pouches with the zip seal.
 - In order to carry out a correct rehydration, please make sure that the lyophilized product is in the bottom of the tube and is not adhered to the top area of the tube or to the foil seal. Gently tap each tube on a hard surface to make sure all the product is at the bottom of the tube.
 - Note: If you choose the format "Dual Master Mix Concentrated Lyophilized MM with Rehydration Buffer (Type 5)" (Section 8.3.1), determine, and separate the appropriate number of additional VIASURE reaction tubes (different foil) and snap into their corresponding positions in the strip (Snap position 4, blue color coding on the rack. See Figure 1). Remove excess air, and close aluminum pouches with the zip seal.
- Remove the required number of Rehydration Buffer tubes (11 foil) and snap into their corresponding positions in the strip (Snap position 3, non-color coding on the rack. See Figure 1).
 - Remove excess air, and close the pouch with the zip seal.
 - In order to ensure a correct transfer, please make sure that the liquid is in the bottom of the tube and is not adhered to the top area of the tube or to the foil seal. Gently tap each tube on a hard surface to make sure all the buffer is at the bottom of the tube.

Figure 1. BD MAX™ TNA Reagent Strip (TNA) from the BD MAX™ ExK™ TNA-3 kit.

8.3.3. BD MAX™ Instrument set up

- 1) Select the "Work List" tab on the "Run" screen of the BD MAX™ System software v4.50A or higher.
- 2) In the "Test" drop down menu, select VIASURE *Bordetella* (if not already created see Section 8.3.1).
- 3) Select the appropriate kit lot number (found on the outer box of extraction kit used) from the pull-down menu (optional).
- 4) Enter the Sample Buffer Tube identification number into the Sample tube window of the Work List, either by scanning the barcode with the scanner or by manual entry.
- 5) Fill the Specimen/Patient ID and/or Accession window of the Work List and click the "Save" button. Continue until all Sample Buffer Tubes are entered. Ensure that the specimen/patient ID and the Sample Buffer Tubes are accurately matched.
- 6) Place the prepared Sample Buffer Tube into the BD MAX™ Rack(s).
- 7) Load the rack(s) into the BD MAX™ System (Rack A is positioned on the left side of the BD MAX™ System and Rack B on the right side).
- 8) Place the required number of BD MAX™ PCR Cartridge(s) into the BD MAX™ System.
- 9) Close the BD MAX™ System door.
- 10) Click "Start Run" to begin the procedure.

8.3.4. BD MAX™ report

- 1) In main menu, click the "Results" button.
- 2) Either double click on your run in the list or press the "view button".
- 3) Click on "Print", select: "Run Details, Test Details and Plot...".
- 4) Click on "Print or Export button" on the "Run Reports" screen.

9. Result interpretation

For a detailed description on how to analyze data, refer to the BD MAX™ System User's manual.

The analysis of the data is done by the BD MAX™ software according to the manufacturer's instructions. The BD MAX™ software reports Ct values and amplification curves for each detector channel of each sample tested in the following way:

- Ct value of 0 indicates that there was no Ct value calculated by the software with the specified Threshold (see Table 4). Amplification curve of the sample showing a "0" Ct value must be checked manually.
- Ct value of -1 indicates that no amplification process has occurred.
- Any other Ct value should be interpreted in correlation with the amplification curve and according to the sample interpretation guidelines outlined in Table 7.

Check Internal Control signal to verify the correct functioning of the amplification mix. In addition, check that there is no report of BD MAX™ System failure.

Results should be read and analyzed using the following table:

<i>B. pertussis</i> / <i>B. holmesii</i> (475/520)	<i>B. holmesii</i> (585/630)	<i>B. parapertussis</i> (630/665)	Internal Control (530/565)	Interpretation for patient's individual samples
+	+	+	+/- ¹	<i>B. pertussis</i> and/or <i>B. holmesii</i> and <i>B. parapertussis</i> DNA Detected ¹
-	-	-	+ ²	<i>B. pertussis</i> , <i>B. holmesii</i> and <i>B. parapertussis</i> DNA Not Detected ²
+	-	-	+/- ¹	<i>B. pertussis</i> DNA Detected, <i>B. holmesii</i> and <i>B. parapertussis</i> DNA Not Detected ¹
+	+	-	+/- ¹	<i>B. pertussis</i> and/or <i>B. holmesii</i> DNA Detected, <i>B. parapertussis</i> DNA Not Detected ¹
+	-	+	+/- ¹	<i>B. pertussis</i> and <i>B. parapertussis</i> DNA Detected, <i>B. holmesii</i> DNA Not Detected ¹
-	+ ³	-	+/- ¹	<i>B. holmesii</i> DNA Detected, <i>B. pertussis</i> and <i>B. parapertussis</i> DNA Not Detected ¹
-	+ ³	+	+/- ¹	<i>B. holmesii</i> and <i>B. parapertussis</i> DNA Detected, <i>B. pertussis</i> DNA Not Detected ¹
-	-	+	+/- ¹	<i>B. parapertussis</i> DNA Detected, <i>B. pertussis</i> and <i>B. holmesii</i> DNA Not Detected ¹
-	-	-	- ²	Unresolved (UNR) Result obtained in the presence of inhibitors in the PCR reaction or when a general problem (not reported by an error code) with the sample processing and/or amplification steps occurs. ²
IND	IND	IND	IND	Indeterminate assay result (IND). Due to BD MAX™ System failure. Assay result displayed in case of an instrument failure linked to an error code.
INC	INC	INC	INC	Incomplete assay result (INC). Due to BD MAX™ System failure. Assay result displayed in case of failure to complete run.

Table 7. Sample interpretation.

+: Amplification occurred.

-: No amplification occurred.

1 A sample is considered positive if the Ct value obtained is ≤ 40 . The Internal Control (IC) may or may not show an amplification signal. Sometimes, the IC detection is not necessary because a high copy number of the target can cause preferential amplification of target-specific nucleic acids.

2 A sample is considered negative if the sample shows no amplification signal in the detection system but the internal control is positive ($Ct \leq 35$). An inhibition of the PCR reaction can be excluded by the amplification of internal control. In case of unresolved results (UNR), absence of internal control signal in negative sample it is recommended to repeat the assay following the indications below.

3 *B. holmesei* DNA is detected in two different channels, 475/520 (FAM) and 585/630 (ROX), but due to higher sensitivity of 585/630 (ROX) in positive samples around the LoD it is possible to obtain signal only in that channel.

In case of a continued ambiguous result, it is recommended to review the instructions for use, the extraction process used by the user; to verify the correct performance of each PCR steps and review the parameters; and to check the sigmoid shape of the curve and the intensity of fluorescence.

NOTE: Sufficient volume is available for one repeat test from the Sample Buffer Tube. For prepared BD MAX™ Sample Buffer Tubes stored at 2–8 °C or 25°C, retesting must be performed within 24 hours.

The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms, and other diagnostic tests.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms, and other diagnostic tests.
- Although this assay can be used with other types of samples, it has been validated with nasopharyngeal swabs and nasopharyngeal aspirates.
- For good test performance, the lyophilized product should be at the bottom of the tube and not adhered to the top area of the tube or the foil seal. Gently tap each tube on a hard surface to make sure all the product is at the bottom of the tube.
- Possible crosstalk might be observed in empty channels of the BD MAX™ System if there is no target to be detected, so it is required to select only the channels where these are amplified when interpretation of results is performed. Please contact to viasuresupport@certest.es for any queries.
- The quality of the test depends on the quality of the sample; nucleic acid must be properly extracted from clinical samples.
- This test is a qualitative test and does not provide quantitative values or indicate the number of organisms present. It is not possible to correlate the Ct values obtained by PCR with the concentration of the sample as they depend on the thermal cycler used and the run itself.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- Please note the intended measurement range of the assay, as samples with concentrations above or below this range may give erroneous results.
- There is a possibility of false positive results due to cross-contamination by *Bordetella* suspicious samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

- The specific primer and probe combinations for detection of the IS481, hIS1001 and pIS1001 sequences, used in VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System, do not show significant combined homologies with the human genome, human microflora, or other respiratory microorganisms, which might result in predictable false positive.
- False Negative results may arise from several factors and their combinations, including:
 - Improper specimens' collection, transport, storage, and/or handling methods.
 - Improper processing procedures (including DNA extraction).
 - Degradation of the DNA during sample shipping/storage and/or processing.
 - Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown *Bordetella* strains.
 - Organism levels in the specimen below the limit of detection for the assay.
 - The presence of qPCR inhibitors or other types of interfering substances. The impact of antibiotics used to prevent the infection or used during the treatment of the infection have not been evaluated.
 - The effect of interfering substances has only been evaluated for those indicated in section 12.5.1.1 of this instructions for use. Please, see this section to check the most common endogenous and exogenous substances that induce a total or partial interference of qPCR reaction. Other substances not indicated in this part could lead to erroneous results.
 - Failure to follow instructions for use and the assay procedure.
- A positive test result does not necessarily indicate the presence of viable bacteria and does not imply that these bacteria are infectious or are the causative agents for clinical symptoms. However, a positive result is indicative of the presence of *Bordetella* targets sequences.
- A negative result does not preclude the presence of *B. pertussis*, *B. parapertussis* and/or *B. holmesii* DNA in a clinical specimen and should not be used as the sole basis for treatment or other patient management decisions. Optimum specimen types and timing for peak bacterial levels during infections caused by *B. pertussis*, *B. parapertussis* and/or *B. holmesii* have not been determined. The collection of multiple specimens (types and time points) from the same patient may be necessary to detect the pathogen.
- If diagnostic tests for other respiratory illnesses are negative and clinical observations, patient history and epidemiological information *B. pertussis*, *B. parapertussis* and/or *B. holmesii* infection, then a false negative result should be considered, and a re-testing of the patient should be discussed.
- In the case of obtaining Unresolved, Indeterminate or Incomplete results using VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System, retesting will be required. Unresolved results may be due to the presence of inhibitors in the sample or an incorrect rehydration of lyophilized reaction mix tube. If there is an instrument failure, Indeterminate or Incomplete results will be obtained.

11. Quality control

VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System contains an Internal Control (IC) in each reaction tube which confirms the correct performance of the technique.

12. Analytical Performance characteristics

12.1. Analytical Linearity

The linearity of the assay was determined by testing nasopharyngeal swab samples (ESwab™ -Copan-) spiked with serial dilutions containing a known concentration of quantified suspension of the targeted *Bordetella* strains (*B. pertussis*, *B. parapertussis* and/or *B. holmesii*). The concentrations of serial dilutions of each *Bordetella* strain were 1.2×10^7 CFU/mL to 1.2×10^{-1} CFU/mL for *Bordetella pertussis*; 1.2×10^7 CFU/mL to 12 CFU/mL for *Bordetella parapertussis*; and 4.1×10^5 CFU/mL to 4.1×10^{-3} CFU/mL for *Bordetella holmesii*.

Example of the amplification plot resulting from an assay is included below:

Figure 1. Dilution series of *Bordetella pertussis* (7.07×10^4 - 7.07×10^{-3} CFU per reaction) template run on the BD MAX™ System (475/520 (FAM) channel).


Figure 2. Dilution series of *Bordetella holmesii* (2.42×10^3 - 2.42×10^{-4} CFU per reaction) template run on the BD MAX™ System (475/520 (FAM) channel).


Figure 3. Dilution series of *Bordetella holmesii* (2.42×10^3 - 2.42×10^{-4} CFU per reaction) template run on the BD MAX™ System (585/630 (ROX) channel).



Figure 4. Dilution series of *Bordetella parapertussis* (7.07×10^4 - 7.07×10^{-2} CFU per reaction) template run on the BD MAX™ System (630/665 (Cy5) channel).

12.2. Analytical sensitivity. Limit of Detection (LoD)

Analytical sensitivity or limit of detection (LoD) of VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System was analysed using nasopharyngeal swabs samples in ESwab™ (Copan), and in VTM (Viracell®). The results obtained showed that VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System has a detection limit of 1.2 CFU/mL of sample for *B. pertussis*, 4.1×10^{-2} CFU/mL of sample for *B. holmesii*, and 12 CFU/mL of sample for *B. parapertussis*, with a positive rate of $\geq 95\%$, on nasopharyngeal swab samples (same LoD with both types of transport medium).

12.3. Measuring range

The measuring range of the assay was determined by testing a series of ten-fold dilutions containing a known concentration specific DNA belonging to the targeted *Bordetella* strains (*B. pertussis*, *B. parapertussis*, and/or *B. holmesii*). Results allowed to confirm the correct detection of the three targets from 1.2×10^7 CFU/mL for *B. pertussis* and *B. parapertussis*, and 4.1×10^5 CFU/mL para *B. holmesii* to the specific LoD of each target.

12.4. Accuracy

12.4.1. Trueness (Veracity)

The veracity of VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System was evaluated by testing reference material listed below:

Veracity of the test for *Bordetella pertussis* was evaluated against the following strains with a positive result:

- *Bordetella pertussis* (NCTC 10739).
- *Bordetella pertussis* (Bergey et al. 1923) Moreno-Lopez 1952, strain BORPER 79936 (CECT 7974).

Veracity of the test for *Bordetella parapertussis* was evaluated against the following strains with a positive result:

- *Bordetella parapertussis* (NCTC 5952).

Veracity of the test for *Bordetella holmesii* was evaluated against the following strains with a positive result:

- *Bordetella holmesii* (NCTC 12912).
- *Bordetella holmesii* Weyant et al. 1995 type strain (DSM 13416).

12.4.2. Precision

To determine the precision of VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System, intra-assay (repeatability), and inter-assay, inter-batch and inter-thermocycler assays (reproducibility) were performed.

Intra-assay, Repeatability

The repeatability was tested by analysing replicates of all samples in the same run using VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System. A summary of results is shown in the table below.

Sample	Target	Channel	$\bar{x}$ (Ct)	σ	CV%
10xLoD Positive samples	<i>B. pertussis</i>	475/520 (FAM)	30.58	0.16	0.52
	<i>B. holmesii</i>	475/520 (FAM)	34.62	0.92	2.65
	<i>B. holmesii</i>	585/630 (ROX)	32.10	0.54	1.69
	<i>B. parapertussis</i>	630/665 (Cy5)	33.10	0.71	2.15
2-3xLoD Positive samples	<i>B. pertussis</i>	475/520 (FAM)	32.50	0.84	2.57
	<i>B. holmesii</i>	475/520 (FAM)	37.13	0.88	2.36
	<i>B. holmesii</i>	585/630 (ROX)	33.38	0.48	1.42
	<i>B. parapertussis</i>	630/665 (Cy5)	34.43	0.84	2.44
Negative sample	<i>B. pertussis</i> <i>B. holmesii</i> <i>B. parapertussis</i>	475/520 (FAM)/ 585/630 (ROX)/ 630/665 (Cy5)	Neg	n.a.	n.a.
	Internal Control	530/565 (HEX)	25.08	0.31	1.22

Table 8. Intra-assay repeatability results of VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System. (Ct) = threshold cycle. ($\bar{x}$) = arithmetic mean Ct value, (σ) = standard deviation, (CV %) = coefficient of variation, Neg = negative, n.a.= not applicable.

Inter-assay, Reproducibility

The inter-assay reproducibility was determined by testing replicates of the different samples on three different days by three different operators, using the VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System. A summary of results is shown in the table below.

Sample	Target	Channel	$\bar{x}$ (Ct)	σ	CV%
10xLoD Positive samples	<i>B. pertussis</i>	475/520 (FAM)	30.83	0.40	1.30
	<i>B. holmesii</i>	475/520 (FAM)	35.30	1.29	3.64
	<i>B. holmesii</i>	585/630 (ROX)	32.35	0.78	2.41
	<i>B. parapertussis</i>	630/665 (Cy5)	33.02	0.85	2.57
2-3xLoD Positive samples	<i>B. pertussis</i>	475/520 (FAM)	32.38	0.42	1.29
	<i>B. holmesii</i>	475/520 (FAM)	36.61	1.38	3.77
	<i>B. holmesii</i>	585/630 (ROX)	33.42	1.04	3.12
	<i>B. parapertussis</i>	630/665 (Cy5)	34.93	0.80	2.28
Negative sample	<i>B. pertussis</i> <i>B. holmesii</i> <i>B. parapertussis</i>	475/520 (FAM)/ 585/630 (ROX)/ 630/665 (Cy5)	Neg	n.a.	n.a.
	Internal Control	530/565 (HEX)	25.18	0.68	2.71

Table 9. Inter-assay reproducibility results of VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System. (Ct) = threshold cycle. ($\bar{x}$) = arithmetic mean Ct value, (σ) = standard deviation, (CV %) = coefficient of variation, Neg = negative, n.a.= not applicable.

Inter-batch, Intermediate precision

The inter-batch intermediate precision values were determined with replicates of the different samples by using three batches of VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System. A summary of results is shown in the table below.

Sample	Target	Channel	$\bar{x}$ (Ct)	σ	CV%
10xLoD Positive samples	<i>B. pertussis</i>	475/520 (FAM)	31.08	0.38	1.21
	<i>B. holmesii</i>	475/520 (FAM)	34.49	1.17	3.39
	<i>B. holmesii</i>	585/630 (ROX)	31.86	0.63	1.98
	<i>B. parapertussis</i>	630/665 (Cy5)	32.48	0.39	1.20
2-3xLoD Positive samples	<i>B. pertussis</i>	475/520 (FAM)	32.89	0.37	1.13
	<i>B. holmesii</i>	475/520 (FAM)	37.09	1.36	3.67
	<i>B. holmesii</i>	585/630 (ROX)	33.93	0.57	1.69
	<i>B. parapertussis</i>	630/665 (Cy5)	34.54	0.65	1.88
Negative sample	<i>B. pertussis</i> <i>B. holmesii</i> <i>B. parapertussis</i>	475/520 (FAM)/ 585/630 (ROX)/ 630/665 (Cy5)	Neg	n.a.	n.a.
	Internal Control	530/565 (HEX)	24.01	0.48	1.99

Table 10. Inter-batch reproducibility results of VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System. (Ct) = threshold cycle. ($\bar{x}$) = arithmetic mean Ct value, (σ) = standard deviation, (CV %) = coefficient of variation, Neg = negative, n.a.= not applicable.

Inter-thermocycler (site-to-site)

The inter-thermocycler reproducibility values were determined with replicates of the same samples used for intra-assay, inter-assay and inter-batch, using the VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System. These assays were run at two laboratory sites with two different BD MAX™ System (Model: 441916– Serial Number: CT0871 and Model: BD MAX CLINICAL – Serial Number: MX0040). A summary of results is shown in the table below.

Sample	Target	Channel	$\bar{x}$ (Ct)	σ	CV%
10xLoD Positive samples	<i>B. pertussis</i>	475/520 (FAM)	30.82	0.32	1.04
	<i>B. holmesii</i>	475/520 (FAM)	34.37	0.84	2.44
	<i>B. holmesii</i>	585/630 (ROX)	32.16	0.48	1.49
	<i>B. parapertussis</i>	630/665 (Cy5)	32.95	0.52	1.57
2-3xLoD Positive samples	<i>B. pertussis</i>	475/520 (FAM)	32.78	0.77	2.35
	<i>B. holmesii</i>	475/520 (FAM)	36.94	1.07	2.90
	<i>B. holmesii</i>	585/630 (ROX)	33.50	0.62	1.85
	<i>B. parapertussis</i>	630/665 (Cy5)	34.53	0.69	2.00
Negative sample	<i>B. pertussis</i> <i>B. holmesii</i> <i>B. parapertussis</i>	475/520 (FAM)/ 585/630 (ROX)/ 630/665 (Cy5)	Neg	n.a.	n.a.
	Internal Control	530/565 (HEX)	25.26	0.42	1.65

Table 11. Overall site-to-site reproducibility results of VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System. (Ct) = threshold cycle. ($\bar{x}$) = arithmetic mean Ct value, (σ) = standard deviation, (CV %) = coefficient of variation, Neg = negative, n.a.= not applicable.

12.5. Analytical specificity and reactivity

The analytical specificity and reactivity were evaluated for the VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System *in silico* and using different starting material such as certified reference strains, certified reference RNA/DNAs and material from the EQAs programmes.

12.5.1. Analytical Specificity

12.5.1.1. Cross-reactivity and exclusivity assay

Cross-reactivity *in silico* analysis

The *in silico* Analytical Specificity was assessed by using reference sequences of different pathogens from NCBI Genbank (<https://www.ncbi.nlm.nih.gov/genbank/>), as well as bioinformatic tools such as BLAST (<http://blast.ncbi.nlm.nih.gov/Blast.cgi>) and an in-house bioinformatic analysis software. Results described below only indicate sequences with a homology percentage higher than 80%. Aligned sequences with a percentage of alignment less than 80% of homology were considered unlikely to be detected. Results obtained demonstrated that **all analyzed sequences were below 80% of homology with *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii* primers and probe set**. Therefore, the VIASURE *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii* target designs of VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System should not cause false positives in detecting *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii* when other organisms are present.

Cross-reactivity wet testing

The specificity of the *Bordetella* assay was confirmed by testing a panel consisting of different microorganisms associated to respiratory diseases. No cross-reactivity was detected between any of the following microorganisms tested:

Cross-reactivity testing					
Human Adenovirus types 1-5, 8, 15, 31, 40 and 41	-	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	-	<i>Legionella micdadei</i>	-
Bocavirus	-	Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus	-	<i>Legionella pneumophila</i> serogroup 1	-
<i>Bordetella bronchiseptica</i>	-	Influenza A/ South Australia/55/2014, IVR-175 virus	-	MERS Coronavirus	-
<i>Chlamydia caviae</i>	-	Influenza A/Switzerland/9715293/2013 (H3N2) virus	-	Human metapneumovirus A and B	-
<i>Chlamydia psittaci</i> genotype A and C	-	Influenza A/Turkey/Germany R2485+86/2014 (H5N8) virus	-	<i>Moraxella catarrhalis</i>	-
<i>Chlamydophila pneumoniae</i> strain CM-1	-	Influenza A/Thüringen/5/2017 (H3N2) virus (Clade 3C2a.1)	-	<i>Mycoplasma pneumoniae</i>	-
Human coronavirus 229E, OC43 and NL63	-	Influenza A/Victoria/210/2009 (H3N2) virus	-	Human parainfluenza 1, 2, 3 and 4 viruses	-
<i>Haemophilus influenzae</i> MinnA	-	Influenza B/Brisbane/60/2008 virus	-	<i>Pneumocystis jirovecii</i> Type A1 and g885652	-
Influenza A/Anhui/1/2013 (H7N9) virus	-	Influenza B/Florida/04/06 virus	-	Respiratory syncytial virus (RSV) types A and B	-
Influenza A/New Caledonia/20/99(H1N1) virus	-	Influenza B/Phuket/3073/2013 virus	-	Human rhinovirus	-
Influenza A/California/7/2009(H1N1)pdm09	-	<i>Legionella bozemanii</i> serovar 1	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-
Influenza A/DE-SH/Reiherente/AR8444/ 2016 (H5N8) virus	-	<i>Legionella dumoffii</i>	-	<i>Streptococcus pneumoniae</i> 2022	-
Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus	-	<i>Legionella longbeachae</i>	-		

Table 12. Reference pathogenic microorganisms used in this study.

Interference and inhibitors of PCR

A study of interfering substances was performed to test the possible interfering effect of endogenous and exogenous substances on VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System. A total of 20 potentially interfering substances were evaluated in triplicate with VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System. All prepared samples were extracted using the BD MAX™ ExK TNA-3 on the BD MAX™ System in DUAL format, so that 6 replicates of each condition were analyzed. The potentially interfering substances were added to the negative nasopharyngeal pool matrix enriched with *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii* at a concentration of 3x LoD.

Positive Matrix Control and Negative Matrix Control (PMC and NMC, respectively) were used for the interpretation of results. PMC corresponds to the negative nasopharyngeal matrix pool spiked with *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii* at a concentration of 3x LoD without any interfering substance, whereas NMC corresponds to the negative nasopharyngeal matrix pool without spiking with any reference microorganisms nor interfering substance. The following results were obtained:

Substance Name	Concentration tested	Result
PMC	-	-
NMC	-	-
Azithromycin	1.1x10 ⁻² mg/mL	N.I
Benzocaine	3 mg/mL	N.I
Carbocysteine	5 mg/mL	N.I
Phenylephrine	0.33 mg/mL	N.I
Fluticasone	1.26x10 ⁻⁶ mg/mL	N.I
Galphimia glauca, luffa operculata, Sabbadilla	12.5 mg/mL	N.I
Genomic DNA	3.5x10 ⁻³ mg/mL	N.I
Human mucus (Sputum)	1 % (v/v)	N.I
Albumin	10 mg/mL	N.I
Mucin	2.5 mg/mL	N.I
Mupirocin	1.5x10 ⁻³ mg/mL	N.I
N-acetylcysteine (NAC)	0.15 mg/mL	N.I
Oseltamivir	3.99x10 ⁻⁴ mg/mL	N.I
Oxymetazoline hydrochloride	0.1 mg/mL	N.I
Nicotine	0.03 mg/mL	I
	7.5x10 ⁻³ mg/mL	N.I
	3.75x10 ⁻³ mg/mL	N.I
Sodium chloride	0.09 mg/mL	N.I
Tobramycin	3.3x10 ⁻² mg/mL	N.I
Triglycerides	15 mg/mL	N.I
Whole Blood	1 % (v/v)	N.I
Zanamivir	3.3 mg/mL	N.I

Table 13. Potential interference substances. N.I: No reportable interfere / I: Interfere.

Different potentially interfering substances, both endogenous and exogenous, were tested on VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System. The results obtained lead to conclude that, at the concentrations tested, no interference of any the evaluated substances is observed at the concentrations indicated for for each condition.

12.5.2. Analytical reactivity

Analytical Reactivity in silico assay

The *in silico* Analytical Reactivity (Inclusivity) of VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System was assessed by using a search and/or alignment tools as BLAST (<http://blast.ncbi.nlm.nih.gov/Blast.cgi>) and an in-house bioinformatic analysis software. Only complete sequences have been used for the inclusivity analysis. Results are shown in the table below:

Microorganism	Gene	Number of analyzed sequences	Sequences detected of the <i>Bordetella</i> strain analyzed without mismatches	Fw mismatches	Rv mismatches	Probe mismatches
<i>Bordetella pertussis</i>	<i>IS481</i>	125 (up to 01/12/2022)	100% (125/125)	-	-	-
<i>Bordetella parapertussis</i>	<i>pIS1001</i>	186 (up to 01/12/2022)	100% (186/186)	-	-	-
<i>Bordetella holmesii</i>	<i>hIS1001</i>	140 (up to 01/12/2022)	100% (140/140)	-	-	-

Table 14. Number of mismatches and percentage of *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii* sequences with or without mismatches as result of the alignment for the Inclusivity analysis.

To sum up, the inclusivity analysis shown a correct detection of *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii* with the VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System.

Analytical Reactivity wet testing

The reactivity of VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System was evaluated against DNA extracted from *B. pertussis*, *B. parapertussis*, and *B. holmesii* as templates, showing positive results.

12.6. Metrological traceability

This assay is not designed for measuring purposes.

13. Clinical performance characteristics

The clinical performance of VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System was tested using clinical nasopharyngeal swabs, and nasopharyngeal aspirates, already characterized as positive or negative for *B. pertussis*, *B. parapertussis* and/or *B. holmesii*. The results were as follows:

	Site	Sample type	Workflow	Target
1	CerTest Biotec (Zaragoza, Spain)	Nasopharyngeal swabs (Cerba Xpert)	BD MAX™ ExK™ TNA-3 kit + BD MAX™ System	<i>B. pertussis</i>
2		Nasopharyngeal aspirates (Cerba Xpert)		<i>B. pertussis</i> <i>B. parapertussis</i> <i>B. holmesii</i>
3		Nasopharyngeal swabs and aspirates (Cerba Xpert)		<i>B. pertussis</i> <i>B. parapertussis</i> <i>B. holmesii</i>
4		Nasopharyngeal swabs (Eurofins)		<i>B. pertussis</i>
5		Nasopharyngeal swabs (Cerba Xpert + Eurofins + simulated samples)		<i>B. pertussis</i> <i>B. parapertussis</i> <i>B. holmesii</i>

Table 15. Site, sample type, workflow and target.

True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, NPV and likelihood ratio (LR) for VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System were calculated in relation to each comparator assay as shown in the following table:

Site	Comparator assay	Target	TP	TN	FP	FN	Sensitivity	Specificity	PPV	NPV	LR+	LR-
1	SmartCycler Bordetella® pertussis/parapertussis assay (Cepheid)	<i>B. pertussis</i>	33	100	0	1	97% (84% – 99%)	100% (96% – 100%)	1 (0.89 – 1)	0.99 (0.95 – 0.99)	193.3 (12.16 – 3072.7)	0.04 (0.01 – 0.13)
2	SmartCycler Bordetella® pertussis/parapertussis assay (Cepheid)	<i>B. pertussis</i>	43	34	1	0	100% (91% – 100%)	97% (85% – 99%)	0.97 (0.88 – 0.99)	1 (0.90 – 1)	23.73 (4.95 – 113.70)	0.01 (0.001 – 0.187)
		<i>B. parapertussis</i>	28	48	0	2	93% (77% – 99%)	100% (92% – 100%)	1 (0.88 – 1)	0.96 (0.86 – 0.99)	90.10 (5.71 – 1423)	0.081 (0.03 – 0.268)
	SmartCycler Bordetella® pertussis/parapertussis assay (Cepheid) + RIDA®GENE Bordetella (R- biopharm)	<i>B. holmesii</i>	1	76	0	1	50% (10% – 98%)	100% (95% – 100%)	1 (0.03 – 1)	0.99 (0.93 – 0.99)	77 (3.89 – 1524)	0.5 (0.162 – 1.561)
3	SmartCycler Bordetella® pertussis/parapertussis assay (Cepheid)	<i>B. pertussis</i>	76	134	1	1	98% (93% – 100%)	99% (95% – 100%)	0.99 (0.93 – 0.99)	0.99 (0.96 – 0.99)	133.25 (18.9 – 939.24)	0.013 (0.002 – 0.092)
	SmartCycler Bordetella® pertussis/parapertussis assay (Cepheid)	<i>B. parapertussis</i>	28	182	0	2	93% (77% – 99%)	100% (98% – 100%)	1 (0.87 – 1)	0.99 (0.96 – 0.99)	336.48 (21.08 – 5369.8)	0.081 (0.025 – 0.265)
	SmartCycler Bordetella® pertussis/parapertussis assay (Cepheid) + RIDA®GENE Bordetella (R- biopharm)	<i>B. holmesii</i>	1	210	0	1	50% (10% – 98.7%)	100% (98% – 100%)	1 (0.17 – 1)	0.99 (0.97 – 0.99)	211 (10.6 – 4199.5)	0.5 (0.162 – 1.554)
4	Validated in-house mix pertussis/parapertussis assay	<i>B. pertussis</i>	43	54	0	3	93% (82.1% – 98.6%)	100% (96.6% – 100%)	1 (0.966 – 1)	1 (0.92 – 1)	0.95 (0.85 – 0.99)	101.81 (6.44 – 1609.2)
5	Initial characterization*	<i>B. pertussis</i>	98	172	0	4	96% (90.3% – 98.9%)	100% (97.9% – 100%)	1 (0.85 – 1)	1 (0.81 – 1)	37.17 (2.41 – 573.29)	0.022 (0.001 – 0.347)
		<i>B. parapertussis</i>	12	262	0	0	100% (73.5% – 100%)	100% (98.7% – 100%)	1 (0.73 – 1)	1 (0.87 – 1)	55.77 (3.56 – 872.35)	0.039 (0.003 – 0.593)
		<i>B. holmesii</i>	20	254	0	0	100% (83.2% – 100%)	100% (98.6% – 100%)	1 (0.83 – 1)	1 (0.83 – 1)	41 (2.65 – 634.61)	0.024 (0.002 – 0.378)

Table 16. True positive and negative values, false positive and negative values, sensitivity, specificity, predictive positive values (PPV), Predicted negative values (NPV) and the likelihood ratios (LR) for VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System.

* Initial characterization includes SmartCycler Bordetella® pertussis/parapertussis assay (Cepheid) and validated in-house mix pertussis/parapertussis assay for *B. pertussis* and *B. parapertussis* detection, and SmartCycler Bordetella® pertussis/parapertussis assay (Cepheid) + RIDA®GENE Bordetella (R-biopharm) for *B. holmesii* detection, in Cerba Xpert samples.

Results show high agreement to detect *B. pertussis*, *B. parapertussis* and/or *B. holmesii* using VIASURE *Bordetella* Real Time PCR Detection Kit for BD MAX™ System.

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Centro	Ensayo comparador	Diana	TP	TN	FP	FN	Sensibilidad	Especificidad	PPV	NPV	LR+	LR-
1	SmartCycler Bordetella® pertussis/parapertussis assay (Cepheid)	<i>B. pertussis</i>	33	100	0	1	97% (84% – 99%)	100% (96% – 100%)	1 (0.89 – 1)	0.99 (0.95–0.99)	193.3 (12.16–3072.7)	0.04 (0.01–0.13)
2	SmartCycler Bordetella® pertussis/parapertussis assay (Cepheid)	<i>B. pertussis</i>	43	34	1	0	100% (91% – 100%)	97% (85% – 99%)	0.97 (0.88 – 0.99)	1 (0.90 – 1)	23.73 (4.95 – 113.70)	0.01 (0.001 – 0.187)
		<i>B. parapertussis</i>	28	48	0	2	93% (77%–99%)	100% (92%–100%)	1 (0.88–1)	0.96 (0.86 – 0.99)	90.10 (5.71 – 1423)	0.081 (0.03 – 0.268)
	SmartCycler Bordetella® pertussis/parapertussis assay (Cepheid) + RIDA®GENE Bordetella (R-biopharm)	<i>B. holmesii</i>	1	76	0	1	50% (10%–98%)	100% (95%–100%)	1 (0.03 – 1)	0.99 (0.93 – 0.99)	77 (3.89 – 1524)	0.5 (0.162 – 1.561)
3	SmartCycler Bordetella® pertussis/parapertussis assay (Cepheid)	<i>B. pertussis</i>	76	134	1	1	98% (93%–100%)	99% (95%–100%)	0.99 (0.93–0.99)	0.99 (0.96–0.99)	133.25 (18.9–939.24)	0.013 (0.002–0.092)
	SmartCycler Bordetella® pertussis/parapertussis assay (Cepheid)	<i>B. parapertussis</i>	28	182	0	2	93% (77%–99%)	100% (98%–100%)	1 (0.87–1)	0.99 (0.96–0.99)	336.48 (21.08–5369.8)	0.081 (0.025–0.265)
	SmartCycler Bordetella® pertussis/parapertussis assay (Cepheid) + RIDA®GENE Bordetella (R-biopharm)	<i>B. holmesii</i>	1	210	0	1	50% (10%–98.7%)	100% (98%–100%)	1 (0.17–1)	0.99 (0.97–0.99)	211 (10.6–4199.5)	0.5 (0.162–1.554)
4	Validated in-house mix pertussis/parapertussis assay	<i>B. pertussis</i>	43	54	0	3	93% (82.1%–98.6%)	100% (96.6%–100%)	1 (0.966–1)	1 (0.92–1)	0.95 (0.85–0.99)	101.81 (6.44–1609.2)
5	Caracterización inicial*	<i>B. pertussis</i>	98	172	0	4	96% (90.3%–98.9%)	100% (97.9%–100%)	1 (0.85–1)	1 (0.81–1)	37.17 (2.41–573.29)	0.022 (0.001–0.347)
		<i>B. parapertussis</i>	12	262	0	0	100% (73.5%–100%)	100% (98.7%–100%)	1 (0.73–1)	1 (0.87–1)	55.77 (3.56–872.35)	0.039 (0.003–0.593)
		<i>B. holmesii</i>	20	254	0	0	100% (83.2%–100%)	100% (98.6%–100%)	1 (0.83–1)	1 (0.83–1)	41 (2.65–634.61)	0.024 (0.002–0.378)

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Bibliography/Bibliografía

1. K. Kösters *et al.* Real-time LightCycler PCR for detection and discrimination of *Bordetella pertussis* and *Bordetella parapertussis*. *Journal of Clinical Microbiology* 2002; 40(5): 1719-1722.
2. K.M. Tatti *et al.* Novel multitarget real-time PCR assay for rapid detection of *Bordetella* species in clinical specimens. *Journal of Clinical Microbiology* 2011; 49(12): 4059-4066.
3. L. Roorda *et al.* A real-time PCR assay with improved specificity for detection and discrimination of all clinically relevant *Bordetella* species by the presence and distribution of three Insertion Sequence elements. *BMC Research Notes* 2011; 4:11.
4. A. Tizolova *et al.* Development of real-time PCR assay for differential detection of *Bordetella bronchiseptica* and *Bordetella parapertussis*. *Diagnostic Microbiology and Infectious Disease* 2014; 78(4):347-351.
5. V. Kolodkina *et al.* Multiplex real-time PCR assay for detection and differentiation of *Bordetella pertussis* and *Bordetella parapertussis*. *Iranian Journal of Microbiology* 2014; 6(3): 140-148.
6. Centers for Disease Control and Prevention. Pertussis (Whooping Cough) (<https://www.cdc.gov/pertussis/>).
7. World Health Organization. Pertussis (<http://www.who.int/topics/pertussis/en/>).

Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Trademarks

BD MAX™ is a registered trademark of Becton, Dickinson and Company.

Modification rights reserved. All rights reserved. © CerTest Biotec, S.L.

All other trademarks that may appear in this package insert are the property of their respective owners.

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Australian sponsor information /

Becton Dickinson Pty Ltd., 66 Waterloo Road.

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Macquarie Park NSW 2113, Australia

New Zealand sponsor information /

Becton Dickinson Limited, 14B George Bourke Drive.

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