# VIASURE

## Real Time PCR Detection Kits

by CerTest  
BIOTEC

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## ENGLISH

### 1. Intended use

VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit is designed for the specific identification and differentiation of *Haemophilus influenzae*, *Neisseria meningitidis* and/or *Streptococcus pneumoniae* in cerebrospinal fluid or blood samples from patients with signs and symptoms of meningitis. This test is intended to be used as an aid in the diagnosis of *H. influenzae*, *N. meningitidis* and/or *S. pneumoniae* in combination with clinical and epidemiological risk factors. DNA is extracted from specimens, multiplied using Real Time amplification and detected using fluorescent reporter dye probes specific for *H. influenzae*, *N. meningitidis* and *S. pneumoniae*.

### 2. Summary and Explanation

The most common causes of bacterial meningitis in adults are *Neisseria meningitidis*, *Streptococcus pneumoniae*, and *Haemophilus influenzae*. These organisms are spread from person to person by close contact with respiratory secretions. Once acquired, each species can colonize the mucosa of the nasopharynx and oropharynx, which is known as pharyngeal carriage. From there, they may cross the mucosa and enter the blood. Once in the blood, they can reach the meninges, causing meningitis, or other body sites causing other syndromes. The term "meningitis" describes inflammation of the membranes (meninges) that surrounds and protects the brain and spinal cord. Meningitis symptoms include sudden onset of fever, headache, and stiff neck. There are often other symptoms, such as nausea, vomiting, photophobia and altered mental status. Bacterial meningitis is a lifethreatening condition that requires prompt recognition and treatment.

*Haemophilus influenzae*, a pleomorphic gram-negative coccobacillus. Is a common commensal microorganism of the upper respiratory tract. It is a human-only pathogen that can cause severe invasive disease, including meningitis, pneumonia, and septicemia. *H. influenzae* strains are divided based on the presence or absence of a polysaccharide capsule; there are 6 encapsulated serotypes (Hia–Hif) and nonencapsulated, nontypeable *H. influenzae* (NTHi) strains. Among them, Hib strains are considered the most prevalent.

*Neisseria meningitidis* may either be encapsulated or unencapsulated. However, nearly all invasive *N. meningitidis* organisms are encapsulated, or surrounded by a polysaccharide capsule. This capsular polysaccharide is used to classify *N. meningitidis* into 12 serogroups. Six of these serogroups cause the great majority of infections in people: A, B, C, W135, X, and Y. Meningitis epidemics are generally caused by serogroup A, although outbreaks have also been caused by serogroups C, W135, and X.

*Streptococcus pneumoniae* cause severe infections like meningitis, community-acquired pneumonia (CAP), bacteremia, bronchitis, sinusitis, and otitis media. More than 90 different serotypes of *S. pneumoniae* have been identified to date on the basis of the biochemical structure of the capsular polysaccharide that is a major virulence factor. The distribution of serotypes can vary with age, geography and time.

Cerebrospinal fluid culture is considered the diagnostic reference standard for bacterial meningitis, and bacterial isolation is important for antimicrobial susceptibility testing and molecular epidemiology. However, CSF culture requires at least a day or more, and has limited sensitivity. Real-Time Polymerase Chain Reaction (qPCR) of CSF

and blood has been suggested as a rapid diagnostic test for bacterial meningitis, and amplification of DNA from non-viable bacteria could potentially facilitate diagnosis in culture negative cases. Using multiple analytic approaches, we found that qPCR assays on CSF and blood specimens were highly accurate for diagnosis of *S. pneumoniae*, *N. meningitidis*, and *H. influenzae* meningitis.

### 3. Principle of the procedure

VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit is designed for the diagnosis of *H. influenzae*, *N. meningitidis* and/or *S. pneumoniae* in clinical samples. After DNA isolation, the identification of *H. influenzae*, *N. meningitidis* and *S. pneumoniae* is performed by the amplification of a conserved region of the *hdp* gene for *Haemophilus influenzae*, *lytA* and *piaA* genes for *Streptococcus pneumoniae* and *ctrA* gene for *Neisseria meningitidis*, using specific primers and a fluorescent-labelled probe.

VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit is based on the 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bounded to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of target template. This fluorescence can be measured on Real Time PCR platforms.

VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase) in an stabilized format, as well as an internal control to monitor PCR inhibition. *Haemophilus influenzae* DNA targets are amplified and detected in the FAM channel, *Streptococcus pneumoniae* DNA targets are amplified and detected in the ROX channel, *Neisseria meningitidis* DNA targets are amplified and detected in the Cy5 channel and the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

### 4. Reagents provided

VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables 1, 2 and 3. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table 1 includes materials and reagents to be used with 8-well strips compatible devices (See Annex 1). Table 2 includes materials and reagents to be used with 96-well plate compatible devices (See Annex 1). Table 3 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips.

Reagent/Material	Description	Colour	Amount
<i>H. influenzae</i> , <i>N. meningitidis</i> & <i>S. pneumoniae</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>H. influenzae</i> , <i>N. meningitidis</i> & <i>S. pneumoniae</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNAse free	RNAse/DNAse free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit with Ref. VS-HNS106L, VS-HNS106H, VS-HNS112L and VS-HNS112H.

Reagent/Material	Description	Color	Amount
<i>H. influenzae</i> , <i>N. meningitidis</i> & <i>S. pneumoniae</i> 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>H. influenzae</i> , <i>N. meningitidis</i> & <i>S. pneumoniae</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNAse free	RNAse/DNAse free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 X 8-cap strip

Table 2. Reagents and materials provided in VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit with Ref. VS-HNS113L and VS-HNS113H.

Reagent/Material	Description	Colour	Amount
<i>H. influenzae</i> , <i>N. meningitidis</i> & <i>S. pneumoniae</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	Transparent	9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>H. influenzae</i> , <i>N. meningitidis</i> & <i>S. pneumoniae</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNAse free	RNAse/DNAse free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	9/18 X 4-cap strip

Table 3. Reagents and materials provided in VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit with Ref. VS-HNS136 and VS-HNS172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

## 5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortexer.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments).

VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTlite Real-time Detection Thermal Cycler, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).


To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

## 6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

## 7. Precautions for users

- The product is intended for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-HNS113L, VS-HNS113H, VS-HNS136 and VS-HNS172). Remove any excess air in the pouches prior to sealing.

- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-HNS136 and VS-HNS172 (compatible with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

## 8. Test procedure

### 8.1. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of the extraction kit used.

For DNA extraction from cerebrospinal fluid or blood samples you can use your manual or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- Maxwell® 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® 16 instrument (Promega).
- ZP02006 MagPurix Bacterial DNA Extraction Kit, using the MagPurix 12A instrument (Zinexts Life Science Corp.).
- Invisorb® Spin Universal Kit (Stratec).

### 8.2. Lyophilized positive control

*H. influenzae*, *N. meningitidis* & *S. pneumoniae* Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNAse free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

### 8.3. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 4. PCR protocol

Fluorogenic data should be collected during the extension step (\*) through the FAM (*H. influenzae*), ROX (*S. pneumoniae*), Cy5 (*N. meningitidis*) and HEX, JOE or VIC channels (CI). Depending on the equipment used select the proper detection channel (see Annex 2). In the Applied Biosystems 7500 Fast Real-Time PCR System and Stratagene Mx3005P™ Real Time PCR System check that the passive reference option for ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

## 9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *H. influenzae*, *N. meningitidis* & *S. pneumoniae* in the positive control well. Check Internal Control signal to verify the correct functioning of the amplification mix. The

analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

Using the following table read and analyze the results:

<i>H. influenzae</i> (FAM)	<i>S. pneumoniae</i> (ROX)	<i>N. meningitidis</i> (Cy5)	Internal control (HEX)	Negative Control	Positive Control	Interpretation
+	+	+	+/-	-	+	<i>H. influenzae, S. pneumoniae and N. meningitidis Positive</i>
-	-	-	+	-	+	<i>H. influenzae, S. pneumoniae and N. meningitidis Negative</i>
+	-	-	+/-	-	+	<i>H. influenzae Positive, S. pneumoniae and N. meningitidis Negative</i>
+	+	-	+/-	-	+	<i>H. influenzae and S. pneumoniae Positive, and N. meningitidis Negative</i>
+	-	+	+/-	-	+	<i>H. influenzae and N. meningitidis Positive, and S. pneumoniae Negative</i>
-	+	-	+/-	-	+	<i>S. pneumoniae Positive, H. influenzae and N. meningitidis Negative</i>
-	+	+	+/-	-	+	<i>S. pneumoniae and N. meningitidis Positive, H. influenzae Negative</i>
-	-	+	+/-	-	+	<i>N. meningitidis Positive, H. influenzae and S. pneumoniae Negative</i>
-	-	-	-	-	+	<b>Experiment fail</b>
+	+	+	+	+	-	<b>Experiment fail</b>

Table 5. Sample interpretation

+: Amplification curve


-: No amplification curve

A sample is considered positive if the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

## 10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples, it has been validated with DNA extracted from cerebrospinal fluid and blood samples.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *H. influenzae*, *N. meningitidis* and *S. pneumoniae* either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

## 11. Quality control

VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

## 12. Performance characteristics

### 12.1. Clinical sensitivity and specificity

VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit was evaluated with different EQA panels (QCMD programmes). These panels consist of 18 clinical specimens dissolved in transport medium and saline samples. The results were compared with the EQA programme final reports. All *H. influenza* (3/18), *N. meningitidis* (2/18) and *S. pneumoniae* (3/18) samples could be detected. In addition, *Streptococcus agalactiae* sample could be confirmed as a negative.

In conclusion, the results show a high sensitivity and specificity to detect *H. influenzae*, *N. meningitidis* and *S. pneumoniae* using VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit.

### 12.2. Analytical sensitivity

VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit has a detection limit of  $\geq 10$  DNA copies per reaction for *H. influenzae*, *N. meningitidis* and *S. pneumoniae* (Figure 2, 3 and 4).

Figure 2. Dilution series of *H. influenzae* ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel FAM).


Figure 3. Dilution series of *S. pneumoniae* ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel ROX).



Figure 4. Dilution series of *N. meningitidis* ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel Cy5).

### 12.3. Analytical specificity

The specificity of the *H. influenzae*, *N. meningitidis* and *S. pneumoniae* assay was confirmed by testing a panel consisting of different microorganisms representing the most common meningitis and encephalitis pathogens. No cross-reactivity was detected between almost any of the following microorganisms tested, except the targeted pathogens of each assay.

Cross-reactivity testing					
<i>E. coli</i> 0:1285;O18:H7:K1	-	<i>Cryptococcus gattii</i>	-	HSV-1	-
<i>Enterococcus durans</i>	-	Dengue 1 virus strain Hawaii A	-	HSV-2	-
<i>Enterococcus faecalis</i>	-	Dengue 2 virus strain New Guinea C	-	HHV 6 Type A and B	-
<i>Enterococcus faecium</i> VZA1	-	Dengue 3 virus strain H87	-	Varicella-Zoster Virus	-
<i>Enterococcus faecium</i> Serotype 11	-	Dengue 4 virus strain H241	-	Coxsackievirus A9	-
<i>Haemophilus influenzae</i>	-	Chikungunya virus strain S27 Petersfield	-	Coxsackievirus A24	-
<i>Listeria innocua</i>	-	Yellow Fever Virus	-	Coxsackievirus B3	-
<i>Listeria ivanovii</i>	-	St Louis Encephalitis Virus virus strain 17D	-	Echovirus 11	-
<i>Listeria monocytogenes</i> Serotype 1/2b	-	West Nile Virus Heja	-	Echovirus 30	-
<i>Listeria monocytogenes</i> Serovar 4b/Strain CIP 59.53	-	West Nile Virus Ug37	-	Enterovirus 68	-
<i>Neisseria meningitidis</i> Serogroup A	-	West Nile Virus strain H160/99	-	Enterovirus 71	-
<i>Streptococcus agalactiae</i>	-	Tick-Borne Encephalitis virus	-	Parechovirus Type 3	-
<i>Streptococcus pneumoniae</i>	-	Epstein-Barr virus	-	Citomegalovirus	-
<i>Staphylococcus aureus</i>	-				

Table 6. Reference pathogenic microorganisms used in this study.

## 12.4. Analytical reactivity

The reactivity of VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit for *H. influenzae* was evaluated against *Haemophilus influenzae* strain MinnA showing positive result.

The reactivity of VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit for *S. pneumoniae* was evaluated against *Streptococcus pneumoniae* strain Z022 showing positive result.

The reactivity of VIASURE *H. influenzae*, *N. meningitidis* & *S. pneumoniae* Real Time PCR Detection Kit for *N. meningitidis* was evaluated against *N. meningitidis* Serogroup A showing positive result.

## ANNEX 1

**COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT**

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast Real-Time PCR System <sup>(1)</sup>
Applied Biosystems	7500 Fast Dx Real-Time PCR System <sup>(1)</sup>
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 3 Fast Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	QuantStudio™ 5 Fast Real-Time PCR System
Applied Biosystems	StepOne Plus™ Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	StepOne™ Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	Mini Opticon™ Real-Time PCR Detection System <sup>(6)</sup>
Cepheid	SmartCycler® <sup>(3)</sup>
Qiagen	Rotor-Gene® Q <sup>(3)</sup>
Roche	LightCycler ®480 Real-Time PCR System <sup>(4)</sup>
Roche	LightCycler ®96 Real-Time PCR System <sup>(4)</sup>
Roche	Cobas z480 Analyzer <sup>(4)</sup>

(1)Select Ramp Speed "**Standard**".

(2)See Annex 3 to check optical measurement exposure setting.

(3)The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Rotor-Gene® Q or SmartCycler® tubes.


(4)Shell Frame grid plate which fits in these Roche qPCR System is necessary.

(5)No detection in Cy5 channel.

(6)Detection in FAM and HEX channels only.

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbott m2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	ABI PRISM 7000 <sup>(6)</sup>
Applied Biosystems	ABI PRISM 7700 <sup>(5)</sup>
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	TOptical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System <sup>(6)</sup>
Bio-Rad	MyiQ™2 Real-Time PCR Detection System <sup>(6)</sup>
Cepheid	SmartCycler® <sup>(3)</sup>
DNA-Technology	DTprime Real-time Detection Thermal Cycler <sup>(2)</sup>
DNA-Technology	DTlite Real-Time PCR System <sup>(2)</sup>
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q <sup>(3)</sup>
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005PT™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System <sup>(2)</sup>
VIASURE	VIASURE 96 Real Time PCR System <sup>(2)</sup>

Table A1/A2. Compatible low and high profile Real Time PCR systems.

## ANNEX 2

**DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT**

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480 II	FAM	465/510	Colour Compensation is required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 Fl for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 µl) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

## ANNEX 3

**OPTICAL MEASUREMENT EXPOSURE SETTING**

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -500\*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

\*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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## 13. Bibliography/Bibliografía

1. V. Thors et al. Population density profiles of nasopharyngeal carriage of 5 bacterial species in pre-school children measured using quantitative PCR offer potential insights into the dynamics of transmission. *Human Vaccines & Immunotherapeutics* 2016;12:2, 375—382.
2. M. G. S. Carvalho et al. Evaluation and Improvement of Real-Time PCR Assays Targeting *lytA*, *ply*, and *psaA* Genes for Detection of Pneumococcal DNA. *Journal of clinical microbiology* 2007; 2460-2466.
3. C.E. Corless et al. Simultaneous Detection of *Neisseria meningitidis*, *Haemophilus influenzae*, and *Streptococcus pneumoniae* in Suspected Cases of Meningitis and Septocemia Using Real-Time PCR. *Journal of clinical microbiology* 2001; 1553-1558.
4. K. Trzcinski et al. Superiority of Trans-Oral over Trans-Nasal Sampling in Detecting *Streptococcus pneumoniae* Colonization in Adults. *Plos One* 2013; 8(3): e60520.

## 14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico in vitro

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- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: June 2019

**CerTest Biotec, S.L.**

Pol. Industrial Río Gállego II · Calle J, N°1  
50840, San Mateo de Gállego, Zaragoza (Spain)  
[www.certest.es](http://www.certest.es)

VIASURE online

F-362 rev01

**VIASURE**

Real Time PCR Detection Kits

**CerTest**  
BIOTEC