CERTEST

E. coli

O157+VT1+VT2

ONE STEP

E. coli O157+VT1+VT2 COMBO

CARD TEST

CERTEST BIOTEC S.L.

For information purposes only

ENGLISH

INTENDED USE

CerTest *E. coli* O157+VT1+VT2 one step combo card test is a coloured chromatographic immunoassay for the simultaneous qualitative detection of *Escherichia coli* O157:H7 and Verotoxins 1 and 2 (VT1 and VT2) produced by *E. coli* in stool samples and *E. coli* suspected colonies in stool culture.

CerTest *E. coli* O157+VT1+VT2 combo card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of enterohaemorrhagic *E. coli* (EHEC) infection.

INTRODUCTION

Infection with *Escherichia coli* O157:H7 (Enterohemorrhagic *Escherichia coli*, EHEC) presents with a wide spectrum of clinical manifestations, including asymptomatic carriage, nonbloody diarrhoea, haemorrhagic colitis, the haemolytic-uremic syndrome (HUS), and thrombotic thrombocytopenic purpura (TTP). Not only is *Escherichia coli* O157:H7 an important agent for haemorrhagic colitis, it is also one of the leading causes of bacterial diarrhoea.

Transmission of *Escherichia coli* (EHEC) is primarily food-borne. Undercooked meat is the most common culprit, dairy products and secondary person-to-person spread is also important. The organism produces at least two Verotoxins (Shiga-like toxins). Verotoxins can be classified into two main categories: Verotoxin 1 (VT1) and Verotoxin 2 (VT2). EHEC strains may produce VT1 or VT2 only or both VT1 and VT2 simultaneously. These toxins are thought to have direct pathogenic significance in *Escherichia coli* O157:H7 infection. This infection is usually diagnosed from a positive stool culture, from the presence of Shiga-like toxins, or both. Timely collection (within 7 days of illness onset) of a stool sample for culture is imperative for a high recovery rate.

EHEC are also known as Shiga toxin-producing *E. coli* (STEC) or Verotoxin-producing *E. coli*. Many different STEC serotypes have been isolated from human cases of HC and HUS. O157:H7 is epidemiologically the most frequent serotype; however the improvement of microbiological detection methods has reported an increasing frequency of other serotypes such as O145:H2, O111:H2, O26:H11, O103:H2, and O121:H19. Regardless of serotype, VT1 and VT2 are uniformly expressed by virtually all human STEC. This universal expression of these toxins, in addition to their important role in the life-threatening complications of STEC infections, makes these antigens attractive immunoprophylactic targets.

TEST PRINCIPLE

CerTest *E. coli* O157+VT1+VT2 is based on the principle of a qualitative immunochromatographic assay for the determination of *Escherichia coli* O157:H7, VT1 and VT2 produced by *E. coli* in stool samples and *E. coli* suspected colonies in stool culture.

Strip A consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Escherichia coli* O157:H7 and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Escherichia coli* O157:H7) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

Strip B consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against VT1 and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-VT1) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

Strip C consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against VT2 and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-VT2) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming two coloured conjugate complexes.

If the sample is *Escherichia coli* O157:H7 positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*Escherichia coli* O157:H7 monoclonal antibodies-red polystyrene microspheres) in the strip A, if the sample is VT1 positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-VT1 monoclonal antibodies-red polystyrene microspheres) in the strip B and if the sample is VT2 positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-VT2 monoclonal antibodies-red polystyrene microspheres) in the strip C, which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*Escherichia coli* O157:H7 antibodies present on the membrane of strip A (test line), the anti-VT1 antibodies present on the membrane of strip B (test line) and the anti-VT2 antibodies present on the membrane of strip C (test line) capture the coloured conjugate and the red lines will be visible in the strips. These bands are used to interpret the result.

If the sample is negative, there is no *Escherichia coli* O157:H7, VT1 and VT2 presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-*Escherichia coli* O157:H7, anti-VT1 and anti-VT2 antibodies present on the membranes (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in the three strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on three membranes will capture control green-conjugate complex and three control lines will always appear. The presence of these green lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest *E. coli* O157+VT1+VT2 combo card test. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.


SPECIMEN COLLECTION AND PREPARATION

Stool samples

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 24 hours prior to testing. For longer storage, maximum 2 weeks, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Ensure only the amount needed is thawed because of freezing and defrosting cycles are not recommended. Homogenise stool samples as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

1. Take out the cap of the collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2), to collect faecal sample (approx. 125mg) and add it to the collection tube. For liquid samples, add approx. 125µL in the collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

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

Suspected EHEC colonies in stool culture

The Sorbitol MacConkey Agar (SMACA) without tellurite or cefixime is the method of choice for the isolation of EHEC (aerobic atmosphere, 24 hours/37°C). After 24 hours of incubation in Sorbitol MacConkey the EHEC colonies will be appear.

1. Take out the cap of the collection tube and add 50µg of Polymixin B (1).
2. Examine Sorbitol MacConkey plates after 1 day incubation. Select EHEC typical colonies. Use the stick or an inoculating needle to pick up 3 or 4 suspected EHEC colonies (2) and add them to the collection tube. Incubate for 30min/35°C (3).
3. Close the tube with the diluent and suspected colonies. Shake the tube in order to assure good sample dispersion (4).

E. coli O157:H7 colonies in SMACA

MATERIALS

MATERIALS PROVIDED

- CerTest *E. coli* O157+VT1+VT2 combo card tests
- Instructions for use
- Collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer
- SMAC agar
- Polymixin B

TEST PROCEDURE

Allow tests, samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove CerTest *E. coli* O157+VT1+VT2 combo card test from its sealed bag just before using it.
3. Take the stool collection tube, cut the end of the cap (4) and dispense 4 drops in the circular window marked with the letter A (5), 4 drops, using the same tube, in the circular window marked with the letter B (6), and 4 drops, using the same tube, in the circular window marked with the letter C (7). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of extraction diluent until seeing the liquid running through the reaction zone.

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



VT2 Strip C-procedure

Add 4 drops in C-circular window

INTERPRETATION OF THE RESULTS (please refer to the illustration below)


- A: Green→Negative *E. coli* O157:H7
B: Green→Negative VT1
C: Green→Negative VT2
- A: Green/Red→Positive *E. coli* O157:H7
B: Green/Red→Positive VT1
C: Green/Red→Positive VT2
- A: Green/Red→Positive *E. coli* O157:H7
B: Green/Red→Positive VT1
C: Green→Negative VT2
- A: Green/Red→Positive *E. coli* O157:H7
B: Green→Negative VT1
C: Green/Red→Positive VT2
- A: Green/Red→Positive *E. coli* O157:H7
B: Green→Negative VT1
C: Green→Negative VT2
- A: Green→Negative *E. coli* O157:H7
B: Green/Red→Positive VT1
C: Green/Red→Positive VT2
- A: Green→Negative *E. coli* O157:H7
B: Green/Red→Positive VT1
C: Green→Negative VT2
- A: Green→Negative *E. coli* O157:H7
B: Green→Negative VT1
C: Green/Red→Positive VT2

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	A (<i>E. coli</i> O157)	B (VT1)	C (VT2)	Interpretation of the results
1.	- GREEN	- GREEN	- GREEN	There are no <i>Escherichia coli</i> O157:H7, VT1 or VT2 presence. No infection caused by EHEC.
2.	+ GREEN-RED	+ GREEN-RED	+ GREEN-RED	There are <i>Escherichia coli</i> O157:H7, VT1 and VT2 presence. Infection caused by EHEC.
3.	+ GREEN-RED	+ GREEN-RED	- GREEN	There is <i>Escherichia coli</i> O157:H7 and VT1 presence. Infection caused by EHEC.
4.	+ GREEN-RED	- GREEN	+ GREEN-RED	There is <i>Escherichia coli</i> O157:H7 and VT2 presence. Infection caused by EHEC.
5.	+ GREEN-RED	- GREEN	- GREEN	There is <i>Escherichia coli</i> O157:H7 presence. Infection caused by EHEC.
6.	- GREEN	+ GREEN-RED	+ GREEN-RED	There is VT1 and VT2 presence. Infection caused by EHEC.
7.	- GREEN	+ GREEN-RED	- GREEN	There is VT1 presence. Infection caused by EHEC.
8.	- GREEN	- GREEN	+ GREEN-RED	There is VT2 presence. Infection caused by EHEC.
9.	Any other result			Invalid result either A, B or C, we recommend repeating the assay using the same sample with another test.

INVALID: Total absence of any control coloured line (GREEN) regardless the appearance or not of the test lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

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QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. Only fresh or fresh-frozen unpreserved and unfixed stool samples can be tested.
3. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
4. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
5. After one week of infection, the number of bacteria in faeces is decreasing, making the sample less reactive. Stool samples should be collected within one week of the onset of symptoms
6. CerTest *E. coli* O157+VT1+VT2 should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
7. Positive results determine the presence of *Escherichia coli* O157:H7, VT1 and VT2 of EHEC (including non-O157 STEC) in faecal samples, nevertheless, a positive result should be followed up with additional laboratory techniques (biochemical methods or by PCR) to confirm the results and the strain. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
8. A negative result is not meaningful because of it is possible the antigens concentration in the stool samples is lower than the detection limit value. If the symptoms or situation still persist, EHEC determination should be carried out on a sample from an enrichment culture.
9. Verotoxin 1 (Shiga-like toxin 1) produced by EHEC is very similar to one produced by *Shigella dysenteriae* type 1 (Shiga toxin). Verotoxin 1 differs from true Shiga toxin by one to seven aminoacids differences. CerTest *E. coli* O157+VT1+VT2 may give a positive result with Shiga toxin produced by *Shigella dysenteriae* type 1. EHEC and *Shigella dysenteriae* type 1 determination should be carried out on selective culture and with biochemical analysis.

EXPECTED VALUES

Escherichia coli O157:H7 causes 73,000 illnesses in the United States annually. That means 8,598 cases, 17% require hospitalizations, 4% hemolytic uremic syndrome cases and 0.5% deaths. The main transmissions routes are: 52% foodborne, 21% unknown, 14% person to person, 31% waterborne, 3% animal contact and 0.3% laboratory related.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed using an immunochromatographic test (CerTest *E. coli* O157+VT1+VT2, CerTest) and these results were confirmed using culture (positive results). The results were as follows:

		Culture		
		+	-	Total
IC test: CerTest <i>E. coli</i> O157+VT1+VT2 (<i>E. coli</i> O157:H7)	+	5	0	5
	-	0	60	60
	Total	5	60	65

	Sensitivity	Specificity	PPV	NPV
CerTest <i>E. coli</i> O157+VT1+VT2 (<i>E. coli</i> O157:H7) vs culture	>99%	>99%	>99%	>99%

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An evaluation, with faecal samples, was performed using an immunochromatographic test (CerTest *E. coli* O157+VT1+VT2, CerTest) and these results were compared with a commercial available test (Shiga Toxin Quik Chek, TechLab®) for Verotoxins 1 and 2. The results were as follows:

		Shiga Toxin Quik Chek, TechLab®		
CerTest <i>E. coli</i> O157+VT1+VT2 (VT1 and VT2)		+	-	Total
	+	4	0	4
	-	0	30	30
	Total	4	30	34

	Sensitivity	Specificity	PPV	NPV
CerTest <i>E. coli</i> O157+VT1+VT2 vs Shiga Toxin Quik Chek	>99%	>99%	>99%	>99%

The results showed a high sensitivity and specificity to detect *Escherichia coli* O157:H7 and Verotoxins 1 and 2 of *Escherichia coli* using CerTest *E. coli* O157+VT1+VT2.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *E. coli* O157+VT1+VT2; no cross reactivity against gastrointestinal pathogens occasionally present in faeces:

<i>Campylobacter coli</i>	<i>Listeria monocytogenes</i>	<i>Shigella flexneri</i>
<i>Campylobacter jejuni</i>	<i>Salmonella enteritidis</i>	<i>Shigella sonnei</i>
<i>Clostridium difficile</i>	<i>Salmonella paratyphi</i>	<i>Staphylococcus aureus</i>
<i>E. coli</i> O157:H7 (Strip B and C)	<i>Salmonella typhi</i>	Verotoxin 1 EHEC (Strip A and C)
<i>E. coli</i> O26	<i>Salmonella typhimurium</i>	Verotoxin 2 EHEC (Strip A and B)
<i>E. coli</i> O111	<i>Shigella boydii</i>	<i>Yersinia enterocolitica</i>
<i>Helicobacter pylori</i>	<i>Shigella dysenteriae</i>	

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