CERTEST

Salmonella paratyphi

ONE STEP *Salmonella paratyphi* CARD
TEST

CERTEST BIOTEC S.L.

CerTest
BIOTEC

The logo consists of a dark blue circle containing the text 'CerTest' in a white, sans-serif font, with 'BIOTEC' in a smaller, white, sans-serif font below it.

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ENGLISH

INTENDED USE

CerTest *Salmonella paratyphi* one step card test is a coloured chromatographic immunoassay for the qualitative detection of *Salmonella paratyphi* in stool samples and *Salmonella paratyphi* suspected colonies in stool culture.

CerTest *Salmonella paratyphi* card test offers a simple and a highly sensitive screening assay to make a presumptive diagnosis of *Salmonella paratyphi* infection and it could be used to identify of suspected isolates of *Salmonella paratyphi* from selective media (stool culture).

INTRODUCTION

Enteric fever (typhoid and paratyphoid fever) is a major human bacterial infection. Although the disease is not common in industrialised countries, it remains an important and persistent health problem in developing nations. That enteric fever is a major public health problem in those countries, with *Salmonella enterica serovar typhi* (*Salmonella typhi*) the most common aetiologic agent but with an apparently increasing number of cases due to *Salmonella paratyphi*. Because risk factors such as poor sanitation, lack of a safe drinking water supply and low socio economic conditions in resource-poor countries are amplified by the evolution of multidrug resistant *salmonellae* with reduced susceptibility to fluoroquinolone, which is associated with increased mortality and morbidity.

In Europe, *Salmonella typhi* and *Salmonella paratyphi* infections occurs among travellers returning from disease endemic areas.

Enteric fever caused by *Salmonella paratyphi* is indistinguishable from that caused by *Salmonella typhi*. This fever usually develops one to three weeks after exposure and varies in severity. Symptoms include high fever, weakness, lethargy, muscle pain, headache, loss of appetite and diarrhoea or constipation. Pink spots appear on the chest; examinations will usually reveal enlargement of the liver and spleen. In severe cases, symptoms of altered mental status and meningitis (fever, stiff neck, seizures) have been reported.

TEST PRINCIPLE

CerTest *Salmonella paratyphi* is based on the principle of a qualitative immunochromatographic assay for the determination of *Salmonella paratyphi* in stool samples and *Salmonella paratyphi* suspected colonies in stool culture.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Salmonella paratyphi* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Salmonella paratyphi*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*Salmonella paratyphi* monoclonal antibodies-red polystyrene microspheres), which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*Salmonella paratyphi* antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no *Salmonella paratyphi* antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-*Salmonella paratyphi* antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appears. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.

- Components provide in the kit are approved for use in the CerTest *Salmonella paratyphi* test. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.


SPECIMEN COLLECTION AND PREPARATION

Stool samples

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Ensure only the amount needed is thawed because of freezing and defrosting cycles are not recommended. Homogenise stool samples as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

1. Take out the cap of the collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2), to collect faecal sample (approx. 125mg) and add it to the collection tube. For liquid samples, add approx. 125µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample. Shake the tube in order to assure good sample dispersion (3).

Suspected *Salmonella paratyphi* colonies in stool culture

The Hektoen Enteric Agar (HEA) is the method of choice for the isolation of *Salmonella paratyphi* (aerobic atmosphere, 24 hours/37°C). After 24 hours of incubation in Hektoen the typical *Salmonella paratyphi* colonies will be green.


1. Examine Hektoen plates after 1 day incubation. Select *Salmonella paratyphi* typical colonies. Take out the cap of the collection tube (1). Use the stick or an inoculating needle to pick up 3 or 4 suspected *Salmonella paratyphi* colonies (2) and add them to the collection tube.

Salmonella paratyphi colonies in HEA

2. Close the tube with the diluent and suspected colonies. Shake the tube in order to assure good sample dispersion (3).

Note: Because it is possible for the *Salmonella paratyphi* content in the stool to be too small, a negative result cannot be considered to be the same as a negative antigen finding. It is recommended to carry out enrichment in Rappaport-Vassiliadis medium (aerobic atmosphere, 24 hours/37°C). After 24 hours of incubation, spread in Hektoen Enteric Agar and incubate (aerobic atmosphere, 24h/37°C).


Salmonella Shigella Agar (SS) is another method for the isolation of *Salmonella paratyphi* (aerobic atmosphere, 24-48 hours/35-37°C). After the incubation period in SS the typical *Salmonella paratyphi* colonies will be transparent without any black center. The SS agar is a very inhibitor that sometimes prevents the growth of certain species of *Salmonella*. We recommend using this medium unless accompanied by another inhibitor HEA.

Salmonella paratyphi colonies in SS

1. Examine SS plates after incubation period. Select *Salmonella paratyphi* typical colonies. Take out the cap of the collection tube (1). Use the stick or an inoculating needle to pick up 3 or 4 suspected *Salmonella* colonies (2) and add them to the collection tube.
2. Close the tube with the diluent and suspected colonies. Shake the tube in order to assure good sample dispersion (3).

MATERIALS

MATERIALS PROVIDED

- CerTest *Salmonella paratyphi* card tests
- Instructions for use
- Collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer
- Hektoen Enteric Agar (HEA)/*Salmonella Shigella Agar (SS)*

TEST PROCEDURE


Allow tests, samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the collection tube in order to assure good sample dispersion.
2. Remove CerTest *Salmonella paratyphi* card test from its sealed bag just before using it.
3. Take the collection tube, cut the end of the cap (4) and dispense 4 drops in the circular window marked with the letter S (5). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF RESULTS (please refer to the illustration below)

	<i>Salmonella paratyphi</i>	Interpretation of the results
1.	- GREEN	There is no <i>Salmonella paratyphi</i> presence. No infection caused by <i>Salmonella paratyphi</i> .
2.	+ GREEN-RED	There is <i>Salmonella paratyphi</i> presence. <i>Salmonella paratyphi</i> infection, which might mean high fever, weakness, lethargy, muscle pain, headache, loss of appetite and diarrhoea or constipation.
3.	Any other result	Invalid result, we recommend repeating the assay using the same sample with another test.

NEGATIVE: Only the GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new card test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

- The test must be carried out within 2 hours after opening the sealed bag.
- An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the buffer and repeat the test.
- The intensity of test line may vary from very strong at high antigens concentrations to faint when the antigens concentration is close to the detection limit value of the test.
- CerTest *Salmonella paratyphi* should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
- Positive results determine the presence of *Salmonella paratyphi* in faecal samples, nevertheless, a positive result should be followed up with additional laboratory techniques (biochemical and serological methods or by PCR) to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
- A negative result is not meaningful because of it is possible the *antigens* concentration in the stool sample is lower than the detection limit value. If the symptoms or situation still persist, a *Salmonella paratyphi* determination should be carried out on a sample from an enrichment culture.

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EXPECTED VALUES

An estimated 22 million cases of typhoid fever and 200.000 related deaths occur worldwide each year; and additional 6 million cases of paratyphoid fever are estimated to occur annually.

Approximately 400 cases of typhoid fever and 150 cases of paratyphoid fever are reported to CDC each year among persons with onset of illness in the United States, most of whom are recent travellers.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit)

Detection limit value of CerTest *Salmonella paratyphi* is 3.125×10^7 CFU/mL.

Clinical sensitivity and specificity

An evaluation, with *Salmonella paratyphi* culture, was performed using an immunochromatographic test (CerTest *Salmonella paratyphi*, CerTest) and these results were confirmed with a commercial available immunochromatographic test (Singlepath@*Salmonella*, Merck). The results were as follows:

		IC test: Singlepath@ <i>Salmonella</i>		
		+	-	Total
IC test: CerTest <i>Salmonella paratyphi</i>	+	10	0	10
	-	0	30	30
Total		10	30	40

	Sensitivity	Specificity	PPV	NPV
CerTest <i>Salmonella paratyphi</i> vs Singlepath@ <i>Salmonella</i>	>99%	>99%	>99%	>99%

The results showed a high sensitivity and specificity to detect *Salmonella paratyphi* using CerTest *Salmonella paratyphi*.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *Salmonella paratyphi*; no cross reactivity against gastrointestinal pathogens occasionally present in faeces:

<i>Campylobacter coli</i>	<i>Helicobacter pylori</i>	<i>Shigella flexneri</i>
<i>Campylobacter jejuni</i>	<i>Listeria monocytogenes</i>	<i>Shigella sonnei</i>
<i>Clostridium difficile</i>	<i>Shigella boydii</i>	<i>Staphylococcus aureus</i>
<i>Escherichia coli</i> O157:H7	<i>Shigella dysenteriae</i>	<i>Yersinia enterocolitica</i>

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1. PFEIFER Y, MATTEN J, RABSCH W. *Salmonella enterica* serovar *typhi* with CTX-M β -lactamase, Germany. Emerg Infect Dis. 2009 Sep. <http://www.cdc.gov/EID/content/15/9/1533.htm>.
2. ABUBAKAR, I. "A systematic review of the clinical, public health and cost-effectiveness of rapid diagnostic tests for the detection and identification of bacterial intestinal pathogens in faeces and food". Health Technology Assessment 2007; Vol. 11: No. 36.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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The logo for CerTest BIOTEC features a light blue circle. Inside the circle, the word "CerTest" is written in a white, sans-serif font, with the "C" and "T" being larger and bolder than the other letters. Below "CerTest", the word "BIOTEC" is written in a smaller, white, all-caps, sans-serif font.

CerTest
BIOTEC

For Informational Purposes
Only

For Information Only Purposes

For Information Only Purposes

CerTest
BIOTEC S.L.

Pol. Industrial Río Gállego II, Calle J, Nº 1,
50840, San Mateo de Gállego, Zaragoza (SPAIN)
www.certest.es