CERTEST

Strep A

ONE STEP Group A
Streptococcal CARD TEST

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ENGLISH

INTENDED USE

CerTest *Strep A* one step card test is a coloured chromatographic immunoassay for the qualitative detection of Group A *Streptococcal* from throat swabs and suspected Group A *Streptococcal* colonies recovered from culture.

CerTest *Strep A* card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of Group A *Streptococcal* respiratory infection.

INTRODUCTION

Group A *Streptococcal* (Strep) infections are caused by Group A *Streptococcus*, a bacterium responsible for a variety of health problems. These infections can range from mild skin infection or sore throat to severe, life-threatening conditions such as toxic shock syndrome (multi-organ failures) and necrotizing fasciitis (soft tissue disease), commonly known as flesh eating disease. Most people are familiar with strep throat, which along with minor skin infection, is the most common form of the disease. Health experts estimate that more than 10 million mild infections (throat and skin) like these occur every year.

Conventional identification procedures for Group A *Streptococcus* from throat swabs involve the isolation and subsequent identification of viable pathogens by techniques that require 24 to 48 hours or longer. CerTest *Strep A* detects either viable or nonviable organisms directly from a throat swab, providing results within 5 minutes.

TEST PRINCIPLE

CerTest *Strep A* is based on the principle of a qualitative immunochromatographic assay for the determination of Group A *Streptococcal* from throat swabs or and suspected Group A *Streptococcal* colonies recovered from culture.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against Group A *Streptococcus* and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-Group A *Streptococcus*) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to blue polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*Strep A* monoclonal antibodies-red polystyrene microspheres) which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*Strep A* antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no Group A *Streptococcus* antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-*Strep A* antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control blue-conjugate complex and the control line will always appears. The presence of this blue line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest *Strep A* test. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.


SPECIMEN COLLECTION AND PREPARATION

Samples should be process as soon as possible after collection. If this is not possible, the samples can be stored in the refrigerator (2-8 °C) for 8 hours prior to testing.

Specimen preparation (see illustration):

Throat swab method:

1. Remove the swab from its packing.
2. Collect specimen with a sterile swab from the tonsils and/or the back of the throat. Take care to avoid the teeth, gums, tongue or cheek surfaces.
3. Insert the swab into the tonsils and/or the back of the throat.
4. Process the swab as soon as possible after collecting the specimen.

Specimen collection

Suspected Group A Streptococcus colonies:

The Blood Agar is the method of choice for the isolation of Group A *Streptococcus* (anaerobic atmosphere, 72 hours/37°C). After 72 hours of incubation in blood the typical Group A *Streptococcus* colonies will be translucent white showing beta haemolysis.

1. Examine Blood plates after incubation period. Select *Streptococcus* typical colonies. Use a sterile swab to pick up 1 or 3 colonies suspected (showing characteristic beta haemolysis).
2. Follow the instructions in the test procedure section to test the swab.

-hemolytic activity of *Streptococcus pyogenes* on blood agar

MATERIALS

MATERIALS PROVIDED

- CerTest *Strep A* card tests
- Reagent A (2M Sodium Nitrite)
- Reagent B (0.15M Acetic Acid)
- Swabs
- Disposable pipettes
- Testing tubes
- Instructions for use
- CerTest *Strep A* Positive Control swab + Instructions for use

MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer
- Blood Agar


TEST PROCEDURE

Allow tests, samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.


1. Add 5 drops (1) Reagent A (light pink) and 5 drops (2) Reagent B in a testing tube (the solution should turn light yellow (colourless)) immediately put the swab into the tube.
2. Mix the solution by rotating the swab forcefully against the side of the tube at least 1 minute. Best results are obtained when the specimen is vigorously extracted in the solution (3). Extract as much liquid as possible from the swab, squeezing the sides of the tube or rotating the swab against the side of the tube as the swab is withdrawn. Discard the swab.
3. Remove CerTest *Strep A* card test from its sealed bag just before using it.
4. Use a separate pipette and test for each sample or control. Dispense exactly 4 drops from the testing tube, into the circular window marked with the letter S (4).
5. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

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INTERPRETATION OF RESULTS (please refer to the illustration below)

	Strep A	Interpretation of the results
1.	- BLUE	There is no Group A <i>Streptococcus</i> (<i>Strep A</i>) presence. No infection caused by Group A <i>Streptococcal</i> .
2.	+ BLUE-RED	There is Group A <i>Streptococcus</i> (<i>Strep A</i>) presence. Infection caused by Group A <i>Streptococcal</i> .
3.	Any other result	Invalid result, we recommend repeating the assay using the same sample with another test.

NEGATIVE: Only the BLUE line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the BLUE control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

INVALID: Total absence of the control coloured line (BLUE) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are the mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the result window will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A blue line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

The colour of the liquid changes from light pink to light yellow (colourless) as you adds Reagent B to Reagent A. This is an internal extraction reagent control. The colour change means that you mixed the extraction reagent control, that you mixed the extraction reagent properly and that the reagents are functioning properly.

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LIMITATIONS

- 1.The test must be carried out within 2 hours after opening the sealed bag.
- 2.This test does not differentiate between carriers and acute infection. Pharyngitis may be caused by organisms other than Group A *Streptococcus*.
- 3.CerTest *Strep A* should be used only with throat swabs or colonies taken directly from a plate. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper throat swab specimens must be obtained.
- 4.The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
- 5.Positive results determine the presence of Group A *Streptococcus* infection. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
- 6.A negative result is not meaningful because of it is possible the antigens concentration in the throat samples is lower than the detection limit value. If the symptoms or situation still persist, a Group A *Streptococcus* determination should be carried out, on a sample from an enrichment culture.

EXPECTED VALUES

Although “strep throat” can occur at any age and at any time of the year, it mainly affects school-age children during the winter and spring.

PERFORMANCE CHARACTERISTICS

Clinical sensitivity and specificity

An evaluation, with throat specimens from patients with pharyngitis symptoms, was performed comparing the results obtained by an immunochromatographic test (CerTest *Strep A*, CerTest) and another rapid test commercial available (OSOM® *Strep A* Test, Genzyme Diagnostics). The results were as follows:

IC test: CerTest <i>Strep A</i>	IC test: OSOM® <i>Strep A</i> Test		
	+	-	Total
+	4	0	4
-	0	20	20
Total	4	20	24

	Sensitivity	Specificity	PPV	NPV
CerTest <i>Strep A</i> vs OSOM® <i>Strep A</i> Test	>99%	>99%	>99%	>99%

The results showed a high sensitivity and specificity to detect Group A *Streptococcus* using CerTest *Strep A*.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *Strep A*; no cross reactivity against organisms that cause other respiratory infections:

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Adenovirus

Streptococcus Grupo D: Enterococcus

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1. Dos Santos, Ana Gabriele P., Berezin, Fitan N.; Comparative analysis of clinical and laboratory methods for diagnosing Streptococcal sore throat. J. Pediatr. (Rio J.) 2005; 81 (1): 23-8.
2. Borbeau, Paul P. and Heither, Barbara J.; Use of swab without Transport media for the Gen-Probe Group A Strep Direct Test. Journal of Clinical Microbiology, July 2004, p. 3207-3211.
3. Petts, D.N.; Evaluation of a modified nitrous acid extraction latex agglutination kit for grouping beta-hemolytic Streptococci and Enterococci. Journal of Clinical Microbiology, Apr. 1995, p.1016-1018.
4. Sheeler, Robert P. et al., Accuracy of Rapid Strep testing in patients who have had recent Streptococcal Pharyngitis. JABFP, July-August 2002, Vol 15 No 4, p. 261-265.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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