CERTEST
Yersinia
enterocolitica
0:3

ONE STEP
Yersinia enterocolitica 0:3
CARD TEST

CERTEST BIOTEC S.L.

For Information Purposes Only

CERTEST *Yersinia enterocolitica* O:3

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ENGLISH

INTENDED USE

CerTest *Yersinia enterocolitica* O:3 one step card test is a coloured chromatographic immunoassay for the qualitative detection of *Yersinia enterocolitica* serotype O:3 in stool samples and *Yersinia* suspected colonies in stool culture.

CerTest *Yersinia enterocolitica* O:3 card test offers a simple and a highly sensitive screening assay to make a presumptive diagnosis of *Yersinia enterocolitica* O:3 infection and it could be used to serotype of suspected isolates of *Yersinia enterocolitica* O:3 from selective media (stool culture).

INTRODUCTION

Yersinia enterocolitica (*Y. enterocolitica*), a Gram-negative, oxidase-negative, and facultatively anaerobic species, is highly heterogeneous and can be divided into several bioserotypes, only a few of which are known to associate with human disease. Six different biotypes (biotype 1A, 1B, 2-5) and numerous serotypes of *Y. enterocolitica* have been described. Eleven of those serotypes have frequently been associated with infections in humans. Most *Y. enterocolitica* strains associates with human yersiniosis belong to bioserotypes 1B/O:8, 2/O:5,27, 2/O:9, 3/O:3 and 4/O:3.

Yersinia enterocolitica is a common pathogen for both humans and animals. Infection with *Y. enterocolitica* can cause a variety of symptoms depending on the age of the person infected. The predominant clinical feature is abdominal pain and diarrhoea with or without fever. Arthritis and erythema nodosum are common complications. Infection with *Y. enterocolitica* occurs most often in young children. Common symptoms in children are fever, abdominal pain, and diarrhea, which is often bloody. Symptoms typically develop 4 to 7 days after exposure and may last 1 to 3 weeks or longer. There is indirect evidence that foods particularly pork products, are important sources of human infections.

The epidemiology of *Y. enterocolitica* infections is complex and poorly understood. Most cases of yersiniosis occur sporadically without an apparent source. *Y. enterocolitica* is thought to be a significant food-borne pathogen, although pathogenic isolates have seldom been isolated from foods, except from edible pig offal. In case control studies, a correlation has been demonstrated between the consumption of raw or undercooked pork and yersiniosis. Genotypes of *Y. enterocolitica* strains found in humans, further supporting the association between yersiniosis and consumption of pork.

TEST PRINCIPLE

CerTest *Yersinia enterocolitica* O:3 is based on the principle of a qualitative immunochromatographic assay for the determination of *Yersinia enterocolitica* serotype O:3 in stool samples and *Yersinia* suspected colonies in stool culture.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against *Y. enterocolitica* serotype O:3 and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-*Y. enterocolitica* serotype O:3) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-*Y. enterocolitica* O:3 monoclonal antibodies-red polystyrene microspheres) which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-*Y. enterocolitica* O:3 antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no *Y. enterocolitica* O:3 presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-*Y. enterocolitica* antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appears. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provided in the kit are approved for use in the CerTest *Yersinia enterocolitica* O:3 test. Do not use any other commercial kit component.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.

SPECIMEN COLLECTION AND PREPARATION

Stool samples

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 1-2 days prior to testing. For longer storage, maximum 1 year, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

Specimen preparation (see illustration):

1. Take out the cap of the collection tube (1) and use the stick to pick up sufficient sample quantity. Then, introduce the stick once into 4 different parts of the stool sample (2), to collect faecal sample (approx. 125mg) and add it to the collection tube. For liquid samples, add approx. 125µL in the collection tube using a micropipette.
2. Close the collection tube with the diluent and stool sample. Shake the collection tube in order to assure good sample dispersion (3).


Suspected *Yersinia* colonies in stool culture

The CIN media (faecal culture) is the method of choice for the isolation of *Yersinia spp.* (aerobes, 48 hours/30°C). After incubation in CIN Agar (Cefsulodin-Irgasan-Novobiocin Agar) the typical *Yersinia enterocolitica* colonies will have deep-red centres surrounded by a transparent, pale border giving the appearance of a "bull's-eye".

1. Examine CIN plates after 1 day incubation. Select small (1-2mm diameter) colonies having deep red centre with sharp border surrounded by clear colourless zone with entire edge. Take out the cap of the collection tube (1). Use the stick or an inoculating needle to pick up 3 or 4 suspected *Yersinia* colonies (2) and add them to the collection tube.
2. Close the collection tube with the diluent and suspected colonies. Shake the collection tube in order to assure good sample dispersion (3).

Yersinia colonies in CIN Agar

Note: After 42 to 48 hours of incubation, the colonies are often completely pink.

MATERIALS

MATERIALS PROVIDED

- CerTest *Yersinia enterocolitica* O:3 card tests
- Instructions for use
- Collection tubes with diluent

MATERIALS REQUIRED BUT NO PROVIDED


- Specimen collection container
- Disposable gloves
- Timer
- CIN media (Agar)

TEST PROCEDURE


Allow tests, samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the collection tube in order to assure good sample dispersion.
2. Remove CerTest *Yersinia enterocolitica* O:3 card test from its sealed bag just before using it.
3. Take the collection tube, cut the end of the cap (4) and dispense 4 drops in the circular window marked with the letter S (5). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of diluent until seeing the liquid running through the reaction zone.

INTERPRETATION OF THE RESULTS (please refer to the illustration below)

	<i>Yersinia enterocolitica</i> O:3	Interpretation of the results
1.	- GREEN	There is no <i>Yersinia enterocolitica</i> O:3 presence. No infection caused by <i>Yersinia enterocolitica</i> O:3.
2.	+ GREEN-RED	There is <i>Yersinia enterocolitica</i> O:3 presence. Infection caused by <i>Yersinia enterocolitica</i> O:3.
3.	Any other result	Invalid result, we recommend repeating the assay using the same sample with another test.

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NEGATIVE: Only the GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation persist, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigens present in the specimen. However, neither the quantitative value nor the rate of increase in antigens can be determined by this qualitative test.

QUALITY CONTROL

Internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. An excess of sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
3. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
4. CerTest *Yersinia enterocolitica* O:3 should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
5. Positive results determine the presence of *Yersinia enterocolitica* O:3 in faecal samples; nevertheless, a positive result should be followed up with additional laboratory techniques (biochemical and serological methods or by PCR) to confirm the results. A confirmed infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
6. A negative result is not meaningful because of it is possible the antigens concentration in the stool sample is lower than the detection limit value. If the symptoms or situation persist, a *Yersinia enterocolitica* O:3 determination should be carried out on a sample from an enrichment culture.

EXPECTED VALUES

In Europe and the United States, most of the human pathogenic *Y. enterocolitica* strains are classified as biotype 4, serotype O:3. After campylobacteriosis and salmonellosis, yersiniosis ranks third among the notifiable bacterial zoonoses in the EU.

Y. enterocolitica is a relatively infrequent cause of diarrhea and abdominal pain. Based on data from the Foodborne Diseases Active Surveillance Network (FoodNet), which measures the burden and sources of specific diseases over time, approximately one culture-confirmed *Y. enterocolitica* infection per 100,000 persons occurs each year. Children are infected more often than adults, and the infection is more common in winter.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit)

Detection limit value of CerTest *Yersinia enterocolitica* O:3 is 3×10^5 CFU/mL.

Clinical sensitivity and specificity

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest *Yersinia enterocolitica* O:3, CerTest) and another commercial available agglutination test (*Yersinia enterocolitica* Agglutination Kit, Progen). The results were as follows:

		AG test: <i>Yersinia enterocolitica</i> Agglutination Kit		
		+	-	Total
IC test: CerTest <i>Yersinia enterocolitica</i> O:3	+	7	0	7
	-	0	56	56
Total		7	56	63

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Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest *Yersinia enterocolitica* O:3; no cross reactivity against gastrointestinal pathogens occasionally present in faeces:

<i>Campylobacter jejuni</i>	<i>Salmonella enteritidis</i>	<i>Shigella flexneri</i>
<i>Campylobacter coli</i>	<i>Salmonella paratyphi</i>	<i>Shigella sonnei</i>
<i>Clostridium difficile</i>	<i>Salmonella typhi</i>	<i>Staphylococcus aureus</i>
<i>Escherichia coli</i> O157:H7	<i>Salmonella typhimurium</i>	<i>Yersinia enterocolitica</i> O:9
<i>Helicobacter pylori</i>	<i>Shigella boydii</i>	
<i>Listeria monocytogenes</i>	<i>Shigella dysenteriae</i>	

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1. WEYNANTS V., et al. "Detection of *Yersinia enterocolitica* Serogroup O:3 by a PCR Method". *Journal of clinical microbiology*, May 1996. p. 1224-1227, Vol. 34, No.5.
2. FREDRIKSSON-AHOMAA M, et al. "Low occurrence of Pathogenic *Yersinia enterocolitica* in Clinical, Food, and Environmental Samples: a Methodological Problem". *Clinical microbiology reviews*, Apr. 2003. p. 22-229, Vol. 16, No. 2.
3. ROSNER B., et al. "Epidemiology of reported *Yersinia enterocolitica* infections in Germany, 2001-2008". *BMC Public Health*, 2010,10:337.

SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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