# VIASURE

## Real Time PCR Detection Kits

by CerTest  
BIOTEC

### *M. tuberculosis complex*

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VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 6 x 8-well strips, low profile	VS-MTC106L
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 6 x 8-well strips, high profile	VS-MTC106H
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 12 x 8-well strips, low profile	VS-MTC112L
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 12 x 8-well strips, high profile	VS-MTC112H
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 96-well plate, low profile	VS-MTC113L
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 96-well plate, high profile	VS-MTC113H
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 9 x 4-well strips, Rotor-Gene®	VS-MTC136
VIASURE <i>M. tuberculosis complex</i> Real Time PCR Detection Kit 18 x 4-well strips, Rotor-Gene®	VS-MTC172

For information purposes only

## ENGLISH

---

### 1. Intended use

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit is designed for specific detection of *M. tuberculosis* complex species in respiratory samples from patients with signs and symptoms of Tuberculosis infection. This test is intended for use as an aid in the diagnosis of the Tuberculosis in combination with clinical and epidemiological risk factors. DNA is extracted from specimens, multiplied using real time amplification and detected using fluorescent reporter dye probes specific for *M. tuberculosis* complex strains.

### 2. Summary and Explanation

Tuberculosis is a contagious, chronic and granulomatous disease caused by *M. tuberculosis*. This disease was declared in 1993 as a "global health emergency" due to its magnitude as a public health problem. Infection occurs through the inhalation of aerosols that contain the pathogen and are transmitted by people with active pulmonary tuberculosis. After inhalation, the bacteria deposit in the alveoli and spread through the lymphatic circulation. Further dissemination to other parts of the lung and occasionally to other organs is achieved through hematogenous circulation. The most common form of the disease is pulmonary tuberculosis, although tuberculous meningitis, miliary tuberculosis (disseminated), intestinal tuberculosis, lymphadenitis, osteomyelitis and Pott's disease (affected bones) also occur. Primary infection leads to active disease in approximately 10% of infected people and in 80% of cases in the period of two years. In the remaining 90%, the immune system controls the infection and the individual is not infectious or asymptomatic. It is estimated that one third of the world population has latent tuberculosis; that is, these people are infected with the bacillus, but (still) have not become ill or can transmit the infection. In this clinical state, TB bacilli can remain inactive for years (latent TB). However, when the immune system weakens, the latent infection can be reactivated. In a person infected with HIV, the risk of reactivation of latent TB is more than 10% per year, compared to a lifetime risk of 10-20% for HIV negative people.

When the active form of the disease occurs, the symptoms (cough, fever, night sweats, weight loss, etc.) can be mild for many months. As a result, patients sometimes take time to seek medical attention and transmit the bacteria to others. Over a year, a tubercular patient can infect 10 to 15 people through close contact. If they do not receive the appropriate treatment, up to two thirds of the tubercular patients die.

The application of molecular techniques for diagnosis and typing, more sensitive, specific and quick than traditional tests, allow to improve the knowledge of the epidemiology of the infection and facilitate decisions for its control.

### 3. Principle of the procedure

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit is designed for the diagnosis of Tuberculosis caused by *M. tuberculosis* complex strains in clinical samples. After DNA isolation, the detection of *M. tuberculosis* complex strains is performed by the amplification of regions of the insertion sequence IS6110 and IS1081 using specific primers and fluorescent-labeled probes.

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase) in a stabilized format, as well as an internal control to monitor PCR inhibition. The Insertion sequences IS6110 and IS1081 are amplified and detected in FAM channel and the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

#### 4. Reagents provided

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables 1, 2 and 3. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table 1 includes materials and reagents to be used with 8-well strips compatible devices (See Annex 1). Table 2 includes materials and reagents to be used with 96-well plate compatible devices (See Annex 1). Table 3 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips.

Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis complex</i> 8-well strips	A mix of enzymes, primers probes, buffer, dNTPs and stabilizers in stabilized format	White	6/12 x 8-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis complex</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	6/12 X 8-cap strip

Table 1. Reagents and materials provided in VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit with Ref. VS-MTC106L, VS-MTC106H, VS-MTC112L and VS-MTC112H.

Reagent/Material	Description	Color	Amount
<i>M. tuberculosis complex</i> 96-well plate	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	White	1 plate
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis complex</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
Tear-off 8-cap strips	Optical caps for sealing plate during thermal cycling	Transparent	12 X 8-cap strip

Table 2. Reagents and materials provided in VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit with Ref. VS-MTC113L and VS-MTC113H.

Reagent/Material	Description	Colour	Amount
<i>M. tuberculosis complex</i> 4-well strips	A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format	Transparent	9/18 x 4-well strip
Rehydration Buffer	Solution to reconstitute the stabilized product	Blue	1 vial x 1.8 mL
<i>M. tuberculosis complex</i> Positive Control	Non-infectious synthetic lyophilized DNA	Red	1 vial
Negative control	Non template control	Violet	1 vial x 1 mL
Water RNase/DNase free	RNase/DNase free water	White	1 vial x 1 mL
4-cap strips	Optical caps for sealing wells during thermal cycling	Transparent	9/18 X 4-cap strip

Table 3. Reagents and materials provided in VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit with Ref. VS-MTC136 and VS-MTC172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

## 5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).
- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortexer.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with and Qiagen/Corbett Rotor-Gene® instruments).

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-

Technology DTLite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

## 6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

## 7. Precautions for users

- The product is indented for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-MTC113L, VS-MTC113H, VS-MTC136 and VS-MTC172). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-MTC136 and VS-MTC172 (compatible with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.
- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

## 8. Test procedure

### 8.1. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of extraction kit used.

For DNA extraction from clinical samples (sputums, BAL and others) you can use your manually or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions for use. We have validated the following extraction kits:

- Viasure RNA-DNA Extraction kit (VIASURE), recommended.
- MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.).

### 8.2. Lyophilized positive control

*M. tuberculosis complex* Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *M. tuberculosis complex* Positive Control (red vial) by adding 100 µL of the supplied Water RNase/DNase free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

### 8.3. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *M. tuberculosis complex* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene®).

Load the plate or the strips in the thermocycler.

3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

Cycles	Step	Time	Temperature
1	Initial denaturation	2 min	95°C
45	Denaturation	10 seg	95°C
	Annealing/Extension (Data collection*)	50 seg	60°C

Table 5. PCR protocol

Fluorogenic data should be collected during the extension step (\*) through the FAM (Insertion sequences IS6110 and IS1081) and HEX, JOE or VIC channels (Internal Control (IC)). Depending on the equipment used select the proper detection channel (see Annex 2). In Applied Biosystems 7500 Fast Real-Time PCR System, Applied Biosystems StepOne Plus™ Real-Time PCR System and Stratagene Mx3005P™ Real-Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

## 9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *M. tuberculosis complex* in the positive control well. Check Internal Control signal to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

Using the following table read and analyze the results:

Insertion sequences IS6110 and IS1081 (FAM)	Internal control (HEX)	Negative control	Positive control	Interpretation
+	+/-	-	+	<i>M. tuberculosis complex</i> strains Positive
	+	-	+	<i>M. tuberculosis complex</i> strains Negative
+	+	+	+	Experiment fail
-	-	-	-	Experiment fail

Table 6. Sample interpretation

+: Amplification curve


-: No amplification curve

A sample is considered positive if the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.

The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

## 10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated with DNA extracted from sputum and BAL samples.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by *M. tuberculosis complex* Positive Control either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.

## 11. Quality control

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

## 12. Performance characteristics

### 12.1. Clinical sensitivity and specificity

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit was tested using 70 positive DNAs (MTBC and NTM strains) from clinical samples (sputums and BAL) from symptomatic patients. The obtained results were compared with those obtained using traditional culture and characterized by conventional PCR follow by sequencing.

Selected samples were positive to MTBC and NTM species and the strain distribution of sublineages was as follow: L2 sublineage Beijing, L4 Euro-American sublineages T, Haarlem, LAM, X and *Mycobacterium africanum* L5 and L6. In addition, 5 DNAs belong to NTM species (*M. avium*, *M. smegmatis*, *M. fortuitum* and *M. abscessus*).

The results were as follows:

VIASURE <i>Mycobacterium tuberculosis complex</i> Real Time PCR Detection Kit	PCR conventional and sequencing			
		+	-	Total
+		65	0	65
-		0	5	5
Total		65	5	70

Table 7. Comparative results.


The results show a high sensitivity and specificity to detect *M. tuberculosis complex* using VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit.

### 12.2. Analytical sensitivity

VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit has a detection limit of  $\geq 10$  DNA copies per reaction (Figure 2).

Figure 2. Dilution series of Insertion sequences *IS6110* and *IS1081* ( $10^7$ - $10^1$  copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel FAM).

### 12.3. Analytical specificity

The specificity of the *M. tuberculosis complex* assay was confirmed by testing a panel consisting of different microorganisms representing the most common respiratory pathogens. No cross-reactivity was detected between almost any of the following microorganisms tested, except the targeted pathogens of each assay.

Cross-reactivity testing					
<i>Bordetella pertussis</i>	-	<i>Moraxella catarrhalis</i>	-	Influenza A/Anhui/1/2013 (H7N9) virus	-
<i>Bordetella parapertussis</i>	-	<i>Haemophilus influenzae</i> MinnA	-	Influenza B/Brisbane/60/2008-like virus	-
<i>Bordetella holmesii</i>	-	<i>Chlamydia pneumoniae</i>	-	Influenza B/Florida/04/06 virus	-
<i>Bordetella bronchiseptica</i>	-	<i>Chlamydia caviae</i>	-	Influenza B/Phuket/3073/2013 virus	-
<i>Legionella bozemanii</i>	-	<i>Chlamydia psittaci</i> genotypes A and C	-	Human parainfluenza 1, 2, 3 and 4 viruses	-
<i>Legionella micdadei</i>	-	Influenza A/New Caledonia/20/99(H1N1) virus	-	Human metapneumovirus A and B	-
<i>Legionella dumoffii</i>	-	Influenza A/California/7/2009(H1N1) virus	-	Human coronavirus 229E, OC43 and NL63	-
<i>Legionella pneumophila</i>	-	Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	-	MERS Coronavirus	-
<i>Legionella longbeache</i>	-	Influenza A/Perth/16/2009(H3N2) virus	-	Human rhinovirus	-
<i>Mycoplasma pneumoniae</i>	-	Influenza A/Thüringen/5/17 (H3N2) virus	-	Human Adenovirus	-
<i>Streptococcus pneumoniae</i>	-	Influenza A/Switzerland/9715293/2013 (H3N2) virus	-	Respiratory Syncytial virus (RSV A and RSV B)	-
<i>Staphylococcus aureus subsp. aureus</i>	-	Influenza A/Turkey/Germany R2485+86/2014 virus (H5N8) virus	-	Human Bocavirus	-
<i>Methicillin-resistant Staphylococcus aureus</i>	-	Influenza A/DE-SH/Reihenente/AR8444/2016 (H5N8) virus	-	<i>Mycobacterium smegmatis</i>	-
<i>Mycobacterium peregrinum</i>	-	<i>Mycobacterium marinum</i>	-	<i>Mycobacterium kansasii</i>	-
<i>Mycobacterium abscessus</i>	-	<i>Mycobacterium avium subsp. avium</i>	-	<i>Mycobacterium gordonae</i>	-

Table 8. Reference pathogenic microorganisms used in this study.

#### 12.4. Analytical reactivity

The reactivity of VIASURE *M. tuberculosis complex* Real Time PCR Detection Kit was evaluated against *Mycobacterium bovis* AF 2122/97, *Mycobacterium tuberculosis* H37Rv, *Mycobacterium tuberculosis* TMC 331, *Mycobacterium tuberculosis* X004439, *Mycobacterium africanum* and *Mycobacterium microti* strains showing positive results.

For information purposes only

ANNEX 1

**COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT**

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

Table A.1 LOW PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Agilent Technologies	AriaMx/AriaDx Real-Time PCR System
Applied Biosystems	7500 Fast Real-Time PCR System <sup>(1)</sup>
Applied Biosystems	7500 Fast Dx Real-Time PCR System <sup>(1)</sup>
Applied Biosystems	QuantStudio™ 12K Flex 96-well Fast
Applied Biosystems	QuantStudio™ 6 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 7 Flex 96-well Fast
Applied Biosystems	QuantStudio™ 3 Fast Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	QuantStudio™ 5 Fast Real-Time PCR System
Applied Biosystems	StepOne Plus™ Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	StepOne™ Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	ViiA™ 7 Fast Real-Time PCR System
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ / CFX96™ IVD Real-Time PCR Detection System
Bio-Rad	Mini Opticon™ Real-Time PCR Detection System <sup>(4)</sup>
Cepheid	SmartCycler® <sup>(3)</sup>
Qiagen	Rotor-Gene® Q <sup>(3)</sup>
Roche	LightCycler® 480 Real-Time PCR System <sup>(4)</sup>
Roche	LightCycler® 96 Real-Time PCR System <sup>(4)</sup>
Roche	Cobas z480 Analyzer <sup>(4)</sup>

Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS	
Manufacturer	Model
Abbott	Abbottm2000 RealTime System
Applied Biosystems	7300 Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	7500 Real-Time PCR System
Applied Biosystems	7900 HT Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	ABI PRISM 7000 <sup>(6)</sup>
Applied Biosystems	ABI PRISM 7700 <sup>(5)</sup>
Applied Biosystems	QuantStudio™ 12K Flex 96-well
Applied Biosystems	QuantStudio™ 6 Flex 96-well
Applied Biosystems	QuantStudio™ 7 Flex 96-well
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System <sup>(5)</sup>
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System
Applied Biosystems	ViiA™ 7 Real-Time PCR System
Analytik Jena Biometra	Optical
Analytik Jena Biometra	qTOWER 2.0
BIONEER	Exicycler™ 96
Bio-Rad	CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System
Bio-Rad	iCycler iQ™ Real-Time PCR Detection System
Bio-Rad	iCycler iQ™5 Real-Time PCR Detection System
Bio-Rad	MyiQ™ Real-Time PCR Detection System <sup>(6)</sup>
Bio-Rad	MyiQ™2 Real-Time PCR Detection System <sup>(6)</sup>
Cepheid	SmartCycler® <sup>(3)</sup>
DNA-Technology	DTprime Real-time Detection Thermal Cycler <sup>(2)</sup>
DNA-Technology	DTlite Real-Time PCR System <sup>(2)</sup>
Eppendorf	Mastercycler™ep realplex
Qiagen	Rotor-Gene® Q <sup>(3)</sup>
Stratagene / Agilent Technologies	Mx3000P™ Real Time PCR System
Stratagene / Agilent Technologies	Mx3005P™ Real Time PCR System
VIASURE	VIASURE 48 Real Time PCR System <sup>(2)</sup>
VIASURE	VIASURE 96 Real Time PCR System <sup>(2)</sup>

(1) Select Ramp Speed "Standard".  
 (2) See Annex 3 to check optical measurement exposure setting.  
 (3) The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Rotor-Gene® Q or SmartCycler® tubes.  
 (4) Shell Frame grid plate which fits in these Roche qPCR System is necessary.  
 (5) No detection in Cy5 channel.  
 (6) Detection in FAM and HEX channels only.

Table A1/A2. Compatible low and high profile Real Time PCR systems.

ANNEX 2

**DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT**

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

REAL-TIME PCR THERMOCYCLER	VIASURE CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation is required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option for ROX must be none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene®Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Mic Real Time PCR Cycler bms	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10)
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Table A3: Detection fluorescence channels of different Real Time PCR systems.

## ANNEX 3

**OPTICAL MEASUREMENT EXPOSURE SETTING**

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -500\*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel -1000.
- DTLite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

\*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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
## 13. Bibliography/Bibliografía

1. Selim *et al.* Duplex real-time PCR assay targeting insertion elements IS1081 and IS6110 for detection of *Mycobacterium bovis* in lymph nodes of cattle. *Biotechnology in Animal Husbandry* 30 (1), p 45-59, 2014.
2. Jesús Gonzalo-Asensio *et al.* New insights into the transposition mechanisms of IS6110 and its Dynamic distribution between *Mycobacterium tuberculosis* Complex lineages. *Plos Genetics*, 14(4) April 2018.
3. Ramadass Balamurugan *et al.* PCR Amplification of the IS6110 Insertion Element of *Mycobacterium tuberculosis* in Fecal Samples from Patients with Intestinal Tuberculosis. *Journal of Clinical Microbiology*. 2006 May; 44(5): 1884–1886.
4. Paul H.M. Savelkoul *et al.*, Detection of *Mycobacterium tuberculosis* complex with Real Time PCR: Comparison of different prime-probe sets base on the IS6110 element. *Journal of Microbiological Methods*. Volume 66, Issue 1, July 2006, Pages 177-180
5. L.A.S. Lira *et al.* Evaluation of a IS6110-Taqman real-time PCR assay to detect *Mycobacterium tuberculosis* in sputum samples of patients with pulmonary TB. *Journal of Applied Microbiology* 2012.
6. Francesco Broccolo *et al.* Rapid Diagnosis of Mycobacterial Infections and Quantitation of *Mycobacterium tuberculosis* Load by two real-time calibrated PCR assays. *Journal of clinical microbiology*. Oct. 2003, p.4565-4572.
7. Santos FCF *et al.* Performance of the IS6110-TaqMan assay in the diagnosis of extrapulmonary tuberculosis from different biological samples. *Revista da Sociedade Brasileira de Medicina Tropical* 51 (3):331-337, may-june, 2018.
8. Dilia Fontalvo Rivera y Doris Gómez Camargo. Genes del *Mycobacterium tuberculosis* involucrados en la patogenicidad y resistencia a antibióticos durante la tuberculosis pulmonar y extrapulmonar. *Revista Médicas UIS*. 2015;28 (1):39-51.
9. I. Dorronsoro y L. Torroba. Microbiología de la tuberculosis. An. Sist. Navar. 2007 Vol 30, Suplemento 2.
10. Fernando Alcaide Fernández de Vega. Nuevos métodos de identificación de micobacterias. Control Calidad SEIMC.
11. Bergmann JS, Keating WE, Woods GL. Clinical evaluation of the BDProbe Tc ET system for rapid detection of *Mycobacterium tuberculosis*. *Journal of Clinical Microbiology* 2000; 38:863-865.
12. Gianna M. *et al.*, Evaluation of the BDProbeTc ET System for Direct Detection of *Mycobacterium tuberculosis* in Pulmonary and Extrapulmonary Samples: a Multicenter Study. *Journal of Clinical Microbiology*. Apr. 2003, p. 1779-1782.
13. Gülnur Tarhan. Molecular Diagnostic Test Used in the Diagnosis of Tuberculosis. *Cohesive Journal of Microbiology and Infectious Disease*.
14. Emilia Cercenado y Rafael Cantón. Procedimientos en Microbiología Clínica. *Recomendaciones de la Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica*. 9ª. Micobacterias.

15. <https://www.tbfacts.org/tb-tests/>
16. Onno W. Akkerman *et al.* Comparison of 14 molecular assays for detection of *Mycobacterium tuberculosis* complex in bronchoalveolar lavage fluid. *Journal of Clinical Microbiology*. 2013 Nov; 51(11): 3505–3511.

#### 14. Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico in vitro

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For information purposes only

- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: March 2019

For information purposes only

For information purposes only

**CerTest Biotec, S.L.**  
Pol. Industrial Río Gállego II · Calle J, Nº1  
50840, San Mateo de Gállego, Zaragoza (Spain)  
[www.cerTest.es](http://www.cerTest.es)

VIASURE online

F-362 rev01

**VIASURE**

Real Time PCR Detection Kits

**CerTest**  
BIOTEC