VIASURE

Real Time PCR Detection Kits

by CerTest
BIOTEC

M. tuberculosis complex + non-tuberculous mycobacteria

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ENGLISH

1. Intended use

VIASURE M. tuberculosis complex + non-tuberculous mycobacteria Real Time PCR Detection Kit is designed for specific detection of DNA of the genus mycobacterium, the M. tuberculosis complex (MTBC) and the specie M. tuberculosis in mycobacterial culture, clinical isolates of mycobacteria and sputum (with positive or negative acid-fast bacilli (AFB) smear) from patients with signs and symptoms of Tuberculosis infection. This test is intended for use as an aid in the diagnosis of M. tuberculosis complex + non-tuberculous mycobacteria in combination with clinical and epidemiological risk factors. DNA is extracted from culture, strains or sputum, multiplied using real time amplification and detected using fluorescent reporter dye probes specific for the selected targets.

2. Summary and Explanation

Tuberculosis is a contagious, chronic and granulomatous disease caused mainly by the specie M. tuberculosis. This disease was declared in 1993 as a "global health emergency" due to its magnitude as a public health problem. Tuberculosis usually affects the lungs, but it can also affect other parts of the body such as the brain, kidneys or spine. Only people who suffer from pulmonary tuberculosis transmit the infection. Primary infection leads to active disease in approximately 10% of infected people and in 80% of cases in the period of two years. In the remaining 90%, the immune system controls the infection and the individual is not infectious or asymptomatic. It is estimated that one third of the world population has latent tuberculosis; that is, these people are infected with the bacillus, but (still) have not become ill or can transmit the infection. In this clinical state, TB bacilli can remain inactive for years (latent TB). However, when the immune system weakens, the latent infection can be reactivated. In a person infected with HIV, the risk of reactivation of latent TB is more than 10% per year, compared to a lifetime risk of 10-20% for HIV negative people.

Mycobacteria are a group of microorganisms that constitute one of the most serious health problems worldwide. Three groups can be defined within the genus Mycobacterium: 1) Tuberculosis complex that produces tuberculosis and is formed by the species M. tuberculosis, M. bovis (including M. bovis BCG), M. africanum, M. microti, M. caprae, M. canetti, M. pinnipedii, M. mungi and M. suricattae; 2) M. leprae that produces leprosy; 3) Other non-tuberculous mycobacteria (NTM) of the Mycobacteria other than tuberculosis that are opportunistic and produce non-tuberculous symptoms with less pathogenic power. The isolation of NTM is increasingly frequent and its differentiation from the M. tuberculosis complex is of great clinical and public health importance, since it defines the isolation of patients in special rooms of health centres and the study of patient contacts.

The symptoms of active tuberculosis will depend on where mycobacteria develop within the body. The general symptoms of tuberculosis include fatigue, weight loss, fever and night sweats. Active pulmonary tuberculosis can cause breathing difficulties, chest pain and bloody expectorations. The symptoms of tuberculosis in other parts of the body depend on the area or organ affected. Over a year, a tubercular patient can infect 10 to 15 people through close contact. If they do not receive the appropriate treatment, up to two thirds of the tubercular patients die.

The early diagnosis of tuberculosis is essential for the control of the disease in order to interrupt the chain of transmission. The application of molecular techniques for diagnosis and typing, more sensitive, specific and quick than traditional tests, allow to improve the knowledge of the epidemiology of the infection and facilitate decisions for its control.

3. Principle of the procedure

VIASURE M. tuberculosis complex + non-tuberculous mycobacteria Real Time PCR Detection Kit is designed for the diagnosis of pulmonary Tuberculosis caused by *M. tuberculosis* or *M. tuberculosis* complex species and NTM infections caused by strains classified as non-tuberculous mycobacteria in clinical samples. After DNA isolation, the detection of the genus mycobacterium is performed by the amplification of a region of the 16S rRNA, the detection of *M. tuberculosis* complex is performed by the amplification of the insertion sequences IS6110 and IS1081 and the detection and differentiation of *M. tuberculosis* specie is performed by the amplification of a fragment of the TbD1 region using specific primers and fluorescent-labeled probes.

VIASURE M. tuberculosis complex + non-tuberculous mycobacteria Real Time PCR Detection Kit is based on 5' exonuclease activity of DNA polymerase. During DNA amplification, this enzyme cleaves the probe bound to the complementary DNA sequence, separating the quencher dye from the reporter. This reaction generates an increase in the fluorescent signal which is proportional to the quantity of the target template. This fluorescence could be measured on Real Time PCR platforms.

VIASURE M. tuberculosis complex + non-tuberculous mycobacteria Real Time PCR Detection Kit contains in each well all the components necessary for real time PCR assay (specific primers/probes, dNTPS, buffer, polymerase) in a stabilized format, as well as an internal control to monitor PCR inhibition. Two regions of the Insertion sequences IS6110 and IS1081 are amplified and detected in FAM channel, a fragment of the 16S rRNA is detected in ROX channel, a fragment of TbD1 region is amplified in Cy5 channel and the internal control (IC) is detected in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

4. Reagents provided

VIASURE M. tuberculosis complex + non-tuberculous mycobacteria Real Time PCR Detection Kit includes the following materials and reagents detailed in Tables 1, 2 and 3. Based on the commercial presentation and the Real Time PCR platform used, the stabilized PCR reaction mix could be placed inside different wells and could be marketed on multiple formats. Table 1 includes materials and reagents to be used with 8-well strips compatible devices (See Annex 1). Table 2 includes materials and reagents to be used with 96-well plate compatible devices (See Annex 1). Table 3 includes materials and reagents for use with Qiagen/Corbett Rotor-Gene® instruments for 4-well strips.

| Reagent/Material | Description | Colour | Amount |
|--|--|-------------|---------------------|
| <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria 8-well strips | A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format | White | 6/12 x 8-well strip |
| Rehydration Buffer | Solution to reconstitute the stabilized product | Blue | 1 vial x 1.8 mL |
| <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Positive Control | Non-infectious synthetic lyophilized DNA | Red | 1 vial |
| Negative control | Non template control | Violet | 1 vial x 1 mL |
| Water RNase/DNAse free | RNAse/DNAse free water | White | 1 vial x 1 mL |
| Tear-off 8-cap strips | Optical caps for sealing wells during thermal cycling | Transparent | 6/12 X 8-cap strip |

Table 1. Reagents and materials provided in VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit with Ref. VS-MTD106L, VS-MTD106H, VS-MTD112L and VS-MTD112H.

| Reagent/Material | Description | Color | Amount |
|--|--|-------------|------------------|
| <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria 96-well plate | A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format | White | 1 plate |
| Rehydration Buffer | Solution to reconstitute the stabilized product | Blue | 1 vial x 1.8 mL |
| <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Positive Control | Non-infectious synthetic lyophilized DNA | Red | 1 vial |
| Negative control | Non template control | Violet | 1 vial x 1 mL |
| Water RNase/DNAse free | RNAse/DNAse free water | White | 1 vial x 1 mL |
| Tear-off 8-cap strips | Optical caps for sealing plate during thermal cycling | Transparent | 12 X 8-cap strip |

Table 2. Reagents and materials provided in VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit with Ref. VS-MTD113L and VS-MTD113H.

| Reagent/Material | Description | Colour | Amount |
|--|--|-------------|---------------------|
| <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria 4-well strips | A mix of enzymes, primers probes, buffer, dNTPs, stabilizers and Internal control in stabilized format | Transparent | 9/18 x 4-well strip |
| Rehydration Buffer | Solution to reconstitute the stabilized product | Blue | 1 vial x 1.8 mL |
| <i>M. tuberculosis</i> complex + non-tuberculous mycobacteria Positive Control | Non-infectious synthetic lyophilized DNA | Red | 1 vial |
| Negative control | Non template control | Violet | 1 vial x 1 mL |
| Water RNase/DNAse free | RNAse/DNAse free water | White | 1 vial x 1 mL |
| 4-cap strips | Optical caps for sealing wells during thermal cycling | Transparent | 9/18 X 4-cap strip |

Table 3. Reagents and materials provided in VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit with Ref. VS-MTD136 and VS-MTD172. For use with Qiagen/Corbett Rotor-Gene® instruments and compatible accessories with strips of 4 tubes 0.1 ml (72-Well Rotor and Locking Ring 72-Well Rotor).

5. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit.

- Real Time PCR instrument (thermocycler).

- DNA extraction kit.
- Centrifuge for 1.5 mL tubes and PCR-well strips or 96-well plate (if available).
- Vortexer.
- Micropipettes (0.5-20 µL, 20-200 µL).
- Filter tips.
- Powder-free disposable gloves.
- Loading block (for use with Qiagen/Corbett Rotor-Gene® instruments).

VIASURE M. tuberculosis complex + non-tuberculous mycobacteria Real Time PCR Detection Kit has been validated on the following equipments: Applied Biosystems 7500 Fast Real-Time PCR System, Bio-Rad CFX96™ Real-Time PCR Detection System, Agilent Technologies AriaMx Real-Time PCR System, DNA-Technology DTprime Real-time Detection Thermal Cycler, DNA-Technology DTlite Real-Time PCR System, Rotor-Gene® Q (Qiagen), SmartCycler® (Cepheid), Roche Molecular Diagnostics Cobas z480 Analyzer, VIASURE 48 Real Time PCR System and VIASURE 96 Real Time PCR System. When using the Applied Biosystems 7500 Fast with strips it is recommend to place a plate holder to reduce the risk of crushed tube (Ref. PN 4388506).

To check thermocycler compatibility, see Annex 1, to check most common detection channels see Annex 2 and to check optical measurement exposure setting see Annex 3.

6. Transport and storage conditions

- The kits can be shipped and stored at 2-40°C until the expiration date which is stated on the label.
- Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles. Positive control has been validated as still being stable after 6 freeze-thaw cycles.
- Keep components away from sunlight.

7. Precautions for users

- The product is intended for use by professional users only, such as laboratory or health professionals and technicians, trained in molecular biological techniques.
- Do not use past expiration date.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Do not use reagents if desiccant is not present or broken inside reagent pouches.
- Do not remove desiccant from reagent pouches once is open.
- Close protective pouches of reagents promptly with the zip seal after each use (if available, Ref. VS-MTD113L, VS-MTD113H, VS-MTD136, VS-MTD172). Remove any excess air in the pouches prior to sealing.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different envelopes and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity. Prolonged exposure to humidity may affect product performance.
- For references VS-MTD136 and VS-MTD172 (compatible with Qiagen/Corbett Rotor-Gene® instruments) use the loading block to pipette reagents and samples into each tube and to help with fitting caps properly and avoid cross contamination.

- Design a unidirectional workflow. It should begin in the Extraction Area and then move to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles and mask. Do not eat, drink or smoke in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Consult safety data sheets, upon request.
- Consult each Real Time PCR instrument's reference manual for additional warnings, precautions and procedures.

8. Test procedure

8.1. DNA extraction

Perform the sample preparation according to the recommendations appearing in the instructions for use of extraction kit used.

For DNA extraction from mycobacterial culture, clinical isolates and sputum*, you can use your manually or automatic routine optimized system. Also, you can use any commercially available DNA extraction kit and follow the manufacturer's instructions for use. We have validated the following extraction kits:

- Total Nucleic Acid Isolation (TNAI) Kit, using COBAS® AmpliPrep (ROCHE).
- Maxwell® 16 Viral Total Nucleic Acid Purification Kit, using the Maxwell® 16 instrument (Promega).
- GenoLyse® (Hain Lifescience GmbH).
- MagDEA Dx SV kit, using the magLEAD® 12gC instrument (Precision System Science Co.).

**In order to improve the yield and quality of bacterial DNA from sputum specimens, it is recommended a sample pre-treatment with N-acetyl-L-cysteine-sodium hydroxide (NALC-NaOH). After decontamination, the cell pellets should be resuspended in a maximum final volume of 1 mL of phosphate buffer. In addition, usage of small elution volumes (50-100 µL) may raise the DNA concentration.*

8.2. Lyophilized positive control

M. tuberculosis complex + non-tuberculous mycobacteria Positive Control contains high copies of the template, the recommendation is to open and manipulate it in a separate laboratory area away from the other components. Reconstitute the lyophilized *M. tuberculosis complex + non-tuberculous mycobacteria Positive Control* (red vial) by adding 100 µL of the supplied Water RNase/DNAse free (white vial) and vortex thoroughly.

Once the positive control has been re-suspended, store it at -20°C. We recommend to separate it in aliquots to minimize freeze and thaw cycles.

8.3. PCR protocol

Determine and separate the number of required reactions including samples and controls. One positive and negative control must be included in each run for each assay. Peel off protective aluminium seal from plates or strips.

- 1) Reconstitute the number of wells you need.

Add 15 µL of Rehydration Buffer (blue vial) into each well.

- 2) Adding samples and controls.

Add 5 µL of DNA sample, reconstituted *M. tuberculosis complex + non-tuberculous mycobacteria* Positive Control (red vial) or Negative Control (violet vial) in different wells and close them with the provided caps.

It is recommended to briefly centrifuge the 8-well strips or 96-well plate, or gently tap each strip onto a hard surface to ensure that all the liquids are at the bottom of the tubes (for Qiagen/Corbett Rotor-Gene® Real Time PCR kit).

Load the plate or the strips in the thermocycler.

- 3) Set up the thermocycler (to check compatibility see Annex 1).

Program the thermocycler following the conditions listed below and start the run:

| Cycles | Step | Time | Temperature |
|--------|--|--------|-------------|
| 1 | Initial denaturation | 2 min | 95°C |
| 45 | Denaturation | 10 seg | 95°C |
| | Annealing/Extension (Data collection*) | 50 seg | 60°C |

Table 4. PCR protocol

Fluorogenic data should be collected during the extension step (*) through the FAM (Insertion sequences IS6110 and IS1081), ROX (16S rRNA), Cy5 (TbD1 region) and HEX, JOE or VIC channels (Internal Control (IC)). Depending on the equipment used select the proper detection channel (see Annex 2). In Applied Biosystems 7500 Fast Real-Time PCR System and Stratagene Mx3005PTM Real Time PCR System check that passive reference option ROX is none. In the Applied Biosystems 7500 Fast Real-Time PCR System select Ramp Speed Standard in Select New Experiment/Advanced Setup/Experiment Properties.

9. Result interpretation

The use of positive and negative controls in each run, validate the reaction by checking the absence of signal in the negative control well and the presence of signal for *M. tuberculosis complex + non-tuberculous mycobacteria* in the positive control well. Check Internal Control signal to verify the correct functioning of the amplification mix. The analysis of the samples is done by the software of the used real time PCR equipment itself according to manufacturer's instructions.

Using the following table read and analyze the results:

| Insertion sequences IS6110 and IS1081 (FAM) | 16S rRNA (ROX) | TbD1 region (Cy5) | Internal control (HEX) | Negative control | Positive control | Interpretation |
|---|----------------|-------------------|------------------------|------------------|------------------|--|
| + | + * | + | +/- | - | + | <i>M. tuberculosis</i> specie Positive or <i>M. tuberculosis</i> specie and NTM Positives (Ct value ROX channel < Ct value FAM and/or Cy5 channels) |
| + | - # | + | +/- | - | + | <i>M. tuberculosis</i> specie Positive |
| + | - # | - # | +/- | - | + | MTBC strains different from <i>M. tuberculosis</i> Positive |
| + | + * | - | +/- | - | + | MTBC strains different from <i>M. tuberculosis</i> Positive (Ct value FAM channel < Ct value ROX channel) or MTBC strains different from <i>M. tuberculosis</i> and NTM Positives (Ct value ROX channel < Ct value FAM channel) |
| - | + † | - | +/- | - | + | NTM positive |
| - | - | - | + | - | + | MTBC and NTM negative |
| - | - | - | - | - | + | Experiment fail |
| + | + | + | + | + | - | Experiment fail |

Table 5. Sample interpretation

+: Amplification curve

-: No amplification curve

*: A sample is considered NTM positive if Ct value obtained in ROX channel (16S rRNA) is lesser than Ct values obtained in FAM (Insertion sequences IS6110 and IS1081) and/or Cy5 (TbD1 region) channels. Please be aware that amplification curves with later Ct may appear in FAM and/or Cy5 channels due to cross-contamination and/or environmental contamination.

#: Because of the amplified targets, FAM signal is stronger than ROX or ROX and Cy5 signals. It could be observed preferential amplification curve in this channel and absence in ROX or ROX and Cy5 channels. This result also indicates the presence of MTBC strain different from *M. tuberculosis* or *M. tuberculosis* specie.


†: Please be aware that amplification curves with later Ct may appear in ROX channel due to environmental NTM.

A sample is considered positive if the Ct value obtained is less than 40 and the internal control shows or not an amplification signal. Sometimes, the detection of internal control is not necessary because a high copy number of target can cause preferential amplification of target-specific nucleic acids.

A sample is considered negative, if the sample shows no amplification signal in the detection system but the internal control is positive. An inhibition of the PCR reaction can be excluded by the amplification of internal control.

Figure 1. Correct run of negative and positive control run on the Bio-Rad CFX96™ Real-Time PCR Detection System.


The result is considered invalid if there is signal of amplification in negative control or absence of signal in the positive well. We recommend to repeat the assay again.

In case of absence of internal control signal in sample wells we recommend to repeat the assay diluting the sample 1:10 or to repeat the extraction to check for possible problems of inhibition.

In case of a doubtful interpretation result, it is recommended to verify the correct performance of each of the steps and review the parameters and the sigmoid shape of the curve. If the situation is not solved, it is recommended to repeat the assay, preferably in duplicate. The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.

10. Limitations of the test

- The results of the test should be evaluated by a health care professional in the context of medical history, clinical symptoms and other diagnostic tests.
- Although this assay can be used with other types of samples it has been validated only with DNA extracted from mycobacterial culture, clinical isolates of mycobacteria and sputum (with positive or negative acid-fast bacilli (AFB) smear). In order to avoid cross-contamination, use dilutions of a culture suspension and/or clinical strains.
- The quality of the test depends on the quality of the sample; proper extracted nucleic acid from clinical samples must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination *M. tuberculosis complex* + non-tuberculous mycobacteria Positive Control either by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions. We recommend use dilutions of a culture suspension and strains.

11. Quality control

VIASURE M. tuberculosis complex + non-tuberculous mycobacteria Real Time PCR Detection Kit contains a positive and a negative control that must be included in each run to correctly interpret the results. Also, the internal control (IC) in each well confirms the correct performance of the technique.

12. Performance characteristics


12.1. Clinical sensitivity and specificity

The clinical performance of VIASURE M. tuberculosis complex + non-tuberculous mycobacteria Real Time PCR Detection Kit was tested using DNA from 52 clinical strains (12 non-tuberculous mycobacteria (NTM) strains, 13 Mycobacterium tuberculosis complex (MTBC) strains, 25 M. tuberculosis strains from different lineages and sublineages, and 2 samples with mixed strains) and 130 decontaminated sputum specimens (49 specimens from patients with a respiratory bacterial infection, 30 specimens from patients with a respiratory NTM infection, and 51 specimens from patients with a Human tuberculosis (TB) diagnosis). These results were compared with those obtained with the reference methods in diagnosing active TB (culture identification). Besides, in patients with a respiratory NTM infection, NTM species were identified with INNO-LiPA® MYCOBACTERIA v2 (Fujirebio). In patients with TB diagnosis, MTBC identification was performed with the SD BIOLINE TB Ag MPT64 Rapid test (Abbot).

Sensitivity, specificity and agreement values for VIASURE M. tuberculosis complex + non-tuberculous mycobacteria Real Time PCR Detection Kit compared with culture identification are shown in next table:

| | Sensitivity (%) (95% CI) | Specificity (%) (95% CI) | Agreement (%) | Kappa |
|-----------------------------------|-----------------------------|--------------------------------|----------------|-------|
| Clinical strains (n=50) | | | | |
| NTM | 12/12 (100) (69.9-100) | 38/38 (100) (88.6-100) | 49/50 (98.0) | 0.968 |
| MTBC | 13/13 (100) (71.7-100) | 36/37 (97.3) (84.2-99.9) | | |
| M. tuberculosis | 24/25 (96.0) (77.7-99.8) | 25/25 (100) (83.4-100) | | |
| Clinical specimens (n=130) | | | | |
| NTM | 23/30 (76.7) (57.3-89.4) | 91/100 (91.0)* (83.2-95.5) | 106/130 (81.5) | 0.717 |
| MTBC/M. tuberculosis | 50/51 (98.0) (88.2-99.9) | 67/79 (84.8)* (74.6-91.6) | | |
| Overall (n=182) | | | | |
| NTM | 35/42 (83.3) (68.0-92.5) | 129/138 (93.5)* (87.6-96.8) | 155/180 (86.1) | 0.819 |
| MTBC/ M. tuberculosis | 88/89 (98.9) (93.0-99.9) | 79/91 (86.8)* (77.7-92.7) | | |

Table 6. Sensitivity, specificity and agreement values between VIASURE M. tuberculosis complex + non-tuberculous mycobacteria Real Time PCR Detection Kit and culture identification. CI: confidence interval. *Contamination in initial steps of sputum processing cannot be ruled out.

In conclusion, the results show a high sensitivity and specificity to identify NTM DNA, *M. tuberculosis* complex DNA and the differentiation of the specie *M. tuberculosis* DNA using VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit.

12.2. Analytical sensitivity

VIASURE *M. tuberculosis* complex + non-tuberculous mycobacteria Real Time PCR Detection Kit has a detection limit of ≥ 10 DNA copies per reaction for the insertion sequences IS6110 and IS1081, the 16S rRNA and the TbD1 region (Figures 2, 3 and 4).

Figure 2. Dilution series of Insertion sequences IS6110 and IS1081 (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel FAM).


Figure 3. Dilution series of 16S rRNA (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel ROX).



Figure 4. Dilution series of TbD1 region (10^7 - 10^1 copies/rxn) template run on the Bio-Rad CFX96™ Real-Time PCR Detection System (channel Cy5)

12.3. Analytical specificity

The specificity of the *M. tuberculosis* complex + non-tuberculous mycobacteria assay was confirmed by testing a panel consisting of different microorganisms representing the most common respiratory pathogens. No cross-reactivity was detected between almost any of the following microorganisms tested, except the targeted pathogens of each assay.

| Cross-reactivity testing | | | | | |
|---|-----|---|-----|---|---|
| <i>Bordetella pertussis</i> | - | <i>Mycobacterium bovis</i> | -/+ | Influenza A/California/7/2009(H1N1) virus | - |
| <i>Bordetella parapertussis</i> | - | <i>Mycobacterium celatum</i> | -/+ | Influenza A/Michigan/45/2015 (H1N1)pdm09 virus | - |
| <i>Bordetella holmesii</i> | - | <i>Mycobacterium chelonae</i> | -/+ | Influenza A/Perth/16/2009(H3N2) virus | - |
| <i>Bordetella bronchiseptica</i> | - | <i>Mycobacterium chimaera-intracellulare</i> | -/+ | Influenza A/Thüringen/5/17 (H3N2) virus | - |
| <i>Chlamydia caviae</i> | - | <i>Mycobacterium fortuitum</i> | -/+ | Influenza A/Switzerland/9715293/2013 (H3N2) virus | - |
| <i>Chlamydophila pneumoniae</i> | - | <i>Mycobacterium genavense</i> | -/+ | Influenza A/Turkey/Germany R2485+86/2014 virus (H5N8) virus | - |
| <i>Chlamydia psittaci</i> genotypes A and C | - | <i>Mycobacterium gordoneae</i> | -/+ | Influenza A/DE-SH/Reiherente/AR8444/ 2016 (H5N8) virus | - |
| <i>Haemophilus influenzae</i> MinnA | - | <i>Mycobacterium intracellulare</i> | -/+ | Influenza A/Anhui/1/2013 (H7N9) virus | - |
| <i>Legionella bozemani</i> | - | <i>Mycobacterium kansasii</i> | -/+ | Influenza B/Brisbane/60/2008-like virus | - |
| <i>Legionella micdadei</i> | - | <i>Mycobacterium marinum</i> | -/+ | Influenza B/Florida/04/06 virus | - |
| <i>Legionella dumoffii</i> | - | <i>Mycobacterium microti</i> | -/+ | Influenza B/Phuket/3073/2013 virus | - |
| <i>Legionella pneumophila</i> | - | <i>Mycobacterium peregrinum</i> | -/+ | <i>Streptococcus pneumoniae</i> | - |
| <i>Legionella longbeache</i> | - | <i>Mycobacterium smegmatis</i> | -/+ | <i>Staphylococcus aureus</i> subsp. <i>aureus</i> | - |
| <i>Moraxella catarrhalis</i> | - | <i>Mycobacterium scrofulaceum</i> | -/+ | Human Adenovirus 1, 2, 3, 4, 5, 6, 7, 8, 15, 31, 40 and 41 | - |
| <i>Mycobacterium abscessus</i> | -/+ | <i>Mycobacterium tuberculosis</i> H37Rv | -/+ | Human Bocavirus | - |
| <i>Mycobacterium africanum</i> | -/+ | <i>Mycobacterium tuberculosis</i> TMC 331 | -/+ | Human coronavirus 229E, OC43 and NL63 | - |
| <i>Mycobacterium africanum</i> lineage L5 and L6 | -/+ | <i>Mycobacterium tuberculosis</i> X004439 | -/+ | Human metapneumovirus A and B | - |
| <i>Mycobacterium avium</i> subsp. <i>avium</i> | -/+ | <i>Mycobacterium tuberculosis</i> lineage 1, 2, 3, and 4 (T sublineage, LAM sublineage, and X sublineage) | -/+ | Human parainfluenza 1, 2, 3 and 4 viruses | - |
| <i>Mycobacterium avium</i> complex (MAC) | -/+ | <i>Mycobacterium xenopi</i> | -/+ | Human rhinovirus | - |
| <i>Mycobacterium avium</i> subsp. <i>paratuberculosis</i> | -/+ | <i>Mycoplasma pneumoniae</i> | - | MERS Coronavirus | - |
| <i>Mycobacterium avium</i> subsp. <i>silvaticum</i> | -/+ | Influenza A/New Caledonia/20/99(H1N1) virus | - | Respiratory Syncytial virus (RSV) | - |
| <i>Mycobacterium bovis</i> BCG | -/+ | | | | |

Table 7. Reference pathogenic microorganisms used in this study.

12.4. Analytical reactivity

The reactivity of VIASURE M. tuberculosis complex + non-tuberculous mycobacteria Real Time PCR Detection Kit for the specie *M. tuberculosis* was evaluated against *Mycobacterium tuberculosis* H37Rv, *Mycobacterium tuberculosis* TMC 331 and *Mycobacterium tuberculosis* X004439 strains, and strains of *M. tuberculosis* lineage 1, 2, 3, and 4 (T sublineage, LAM sublineage, and X sublineage), showing positive results.

The reactivity of VIASURE M. tuberculosis complex + non-tuberculous mycobacteria Real Time PCR Detection kit for the mycobacterium genus was evaluated against the following strains: *Mycobacterium bovis* AF 2122/97, *Mycobacterium bovis* BCG, *Mycobacterium tuberculosis* H37Rv, *Mycobacterium tuberculosis* TMC 331, *Mycobacterium tuberculosis* X004439, *Mycobacterium tuberculosis* lineage 1, 2, 3, and 4 (T sublineage, LAM sublineage, and X sublineage), *Mycobacterium africanum*, *Mycobacterium africanum* lineage L5 and L6, *Mycobacterium microti*, *Mycobacterium abscessus*, *Mycobacterium avium* complex (MAC), *Mycobacterium avium* subsp. *avium*, *Mycobacterium avium* subsp. *paratuberculosis*, *Mycobacterium avium* subsp. *silvaticum*, *Mycobacterium celatum*, *Mycobacterium chelonae*, *Mycobacterium chelonae* complex (group III- *M. abscessus*), *Mycobacterium chimaera-intracellulare*, *Mycobacterium fortuitum*, *Mycobacterium genavense*, *Mycobacterium gordoneae*, *Mycobacterium intracellulare*, *Mycobacterium kansasii*, *Mycobacterium marinum*, *Mycobacterium peregrinum*, *Mycobacterium smegmatis*, *Mycobacterium scrofulaceum*, and *Mycobacterium xenopi*, showing positive results.

The reactivity of VIASURE M. tuberculosis complex + non-tuberculous mycobacteria Real Time PCR Detection kit for MTBC strains was evaluated against the following strains: *Mycobacterium bovis* AF 2122/97, *Mycobacterium bovis* BCG, *Mycobacterium tuberculosis* H37Rv, *Mycobacterium tuberculosis* TMC 331, *Mycobacterium tuberculosis* X004439, *Mycobacterium tuberculosis* lineage 1, 2, 3, and 4 (T sublineage, LAM sublineage, and X sublineage), *Mycobacterium africanum*, *Mycobacterium africanum* lineage L5 and L6, and *Mycobacterium microti*, showing positive results.

ANNEX 1

COMPATIBILITY WITH THE MOST COMMON REAL TIME PCR EQUIPMENT

Low profile strips can be used in all PCR thermocyclers equipped with a low profile block, like the systems listed in table A.1. High profile strips can be used in all PCR thermocyclers equipped with a high or regular profile block, like the systems listed in table A.2. If you do not find your thermocycler in the list below, please contact with your supplier.

| Table A.1 LOW PROFILE BLOCK THERMOCYCLERS | |
|---|---|
| Manufacturer | Model |
| Agilent Technologies | AriaMx/AriaDx Real-Time PCR System |
| Applied Biosystems | 7500 Fast Real-Time PCR System ⁽¹⁾ |
| Applied Biosystems | 7500 Fast Dx Real-Time PCR System ⁽¹⁾ |
| Applied Biosystems | QuantStudio™ 12K Flex 96-well Fast |
| Applied Biosystems | QuantStudio™ 6 Flex 96-well Fast |
| Applied Biosystems | QuantStudio™ 7 Flex 96-well Fast |
| Applied Biosystems | QuantStudio™ 3 Fast Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | QuantStudio™ 5 Fast Real-Time PCR System |
| Applied Biosystems | StepOne Plus™ Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | StepOne™ Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | ViiA™ 7 Fast Real-Time PCR System |
| BIONEER | Exicycler™ 96 |
| Bio-Rad | CFX96™ / CFX96™ IVD Real-Time PCR Detection System |
| Bio-Rad | Mini Opticon™ Real-Time PCR Detection System ⁽⁶⁾ |
| Cepheid | SmartCycler® ⁽³⁾ |
| Qiagen | Rotor-Gene® Q ⁽³⁾ |
| Roche | LightCycler ®480 Real-Time PCR System ⁽⁴⁾ |
| Roche | LightCycler ®96 Real-Time PCR System ⁽⁴⁾ |
| Roche | Cobas z480 Analyzer ⁽⁴⁾ |

(1)Select Ramp Speed "Standard".

(2)See Annex 3 to check optical measurement exposure setting.

(3)The product should be reconstituted following the appropriate procedure (see Test Procedure) and transferred into the specific Rotor-Gene® Q or SmartCycler® tubes.

(4)Shell Frame grid plate which fits in these Roche qPCR System is necessary.

(5)No detection in Cy5 channel.

(6)Detection in FAM and HEX channels only.

| Table A.2 HIGH PROFILE BLOCK THERMOCYCLERS | |
|--|--|
| Manufacturer | Model |
| Abbott | Abbott m2000 RealTime System |
| Applied Biosystems | 7300 Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | 7500 Real-Time PCR System |
| Applied Biosystems | 7900 HT Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | ABI PRISM 7000 ⁽⁶⁾ |
| Applied Biosystems | ABI PRISM 7700 ⁽⁵⁾ |
| Applied Biosystems | QuantStudio™ 12K Flex 96-well |
| Applied Biosystems | QuantStudio™ 6 Flex 96-well |
| Applied Biosystems | QuantStudio™ 7 Flex 96-well |
| Applied Biosystems | QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | QuantStudio™ 5 Real-Time PCR System |
| Applied Biosystems | ViiA™ 7 Real-Time PCR System |
| Analytik Jena Biometra | TOptical |
| Analytik Jena Biometra | qTOWER 2.0 |
| BIONEER | Exicycler™ 96 |
| Bio-Rad | CFX96™ Deep Well / CFX96™ Deep Well IVD Real-Time PCR Detection System |
| Bio-Rad | iCycler iQ™ Real-Time PCR Detection System |
| Bio-Rad | iCycler iQ™5 Real-Time PCR Detection System |
| Bio-Rad | MyIQ™ Real-Time PCR Detection System ⁽⁶⁾ |
| Bio-Rad | MyIQ™2 Real-Time PCR Detection System ⁽⁶⁾ |
| Cepheid | SmartCycler® ⁽³⁾ |
| DNA-Technology | DTprime Real-time Detection Thermal Cycler ⁽²⁾ |
| DNA-Technology | DTlite Real-Time PCR System ⁽²⁾ |
| Eppendorf | Mastercycler™ep realplex |
| Qiagen | Rotor-Gene® Q ⁽³⁾ |
| Stratagene / Agilent Technologies | Mx3000PT™ Real Time PCR System |
| Stratagene / Agilent Technologies | Mx3005PT™ Real Time PCR System |
| VIASURE | VIASURE 48 Real Time PCR System ⁽²⁾ |
| VIASURE | VIASURE 96 Real Time PCR System ⁽²⁾ |

Table A1/A2. Compatible low and high profile Real Time PCR systems.

ANNEX 2

DETECTION CHANNELS FOR THE MOST COMMON REAL TIME PCR EQUIPMENT

The fluorescence detection channels for some of most common Real Time PCR Thermocyclers are specified in Table A3.

| REAL-TIME PCR THERMOCYCLER | VIASURE CHANNEL | DETECTION CHANNEL | OBSERVATIONS |
|-------------------------------|-----------------|-------------------|---|
| Bio-Rad CFX96™ | FAM | FAM | Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, in the Settings menu, select the option Apply Fluorescence Drift Correction for Baseline Settings to correct it. |
| | HEX | HEX | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| ABI 7500 Applied Biosystems | FAM | FAM | Passive reference option for ROX must be none. Some wells may have abnormally drifting RFU values during the initial few cycles of a run showing a non-sigmoidal ascendant line. If you see this effect, please modify the baseline: Select the Start Cycle and End Cycle values so that the baseline ends before significant fluorescence is detected. |
| | HEX | VIC | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| Roche Lightcycler®480II | FAM | 465/510 | Colour Compensation is required |
| | HEX | 533/580 | |
| | ROX | 533/610 | |
| | Cy5 | 618/660 | |
| Smartcycler® Cepheid | FAM | Channel 1 | |
| | HEX | Channel 2 | |
| | ROX | Channel 3 | |
| | Cy5 | Channel 4 | |
| Abbott m2000rt | FAM | FAM | |
| | HEX | VIC | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| Mx3000P™ Mx 3005P™ Stratagene | FAM | FAM | Passive reference option for ROX must be none |
| | HEX | VIC | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| AriaMX Agilent | FAM | FAM | |
| | HEX | HEX | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |
| Rotor-Gene®Q Qiagen | FAM | Green | In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition". |
| | HEX | Yellow | |
| | ROX | Orange | |
| | Cy5 | Red | |
| Mic Real Time PCR Cycler bms | FAM | Green | In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the appropriate thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10) |
| | HEX | Yellow | |
| | ROX | Orange | |
| | Cy5 | Red | |
| Exicycler™ 96 BIONEER | FAM | FAM | |
| | HEX | JOE | |
| | ROX | ROX | |
| | Cy5 | Cy5 | |

Table A3: Detection fluorescence channels of different Real Time PCR systems.

ANNEX 3

OPTICAL MEASUREMENT EXPOSURE SETTING

Optical measurement parameters of some thermocyclers must be adjusted to be suitable for operation with "VIASURE Real Time PCR Detection Kits". This assay has been validated with the following set exposition values:

- DTprime Real-time Detection Thermal Cycler (DNA-Technology) and VIASURE 96 Real Time PCR System (CerTest Biotec S.L.): FAM channel -500*, HEX channel – 1000, ROX channel – 1000 and Cy5 channel - 1000.
- DTlite Real-Time PCR System (DNA-Technology) and VIASURE 48 Real Time PCR System (CerTest Biotec S.L.): FAM channel - 500, HEX channel - 500, ROX channel – 500 and Cy5 channel - 500.

*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

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13. Bibliography/Bibliografía

1. Selim et al. Duplex real-time PCR assay targeting insertion elements IS1081 and IS6110 for detection of *Mycobacterium bovis* in lymph nodes of cattle. *Biotechnology in Animal Husbandry* 30 (1), p 45-59, 2014.
2. Jesús Gonzalo-Asensio et al. New insights into the transposition mechanisms of IS6110 and its Dynamic distribution between *Mycobacterium tuberculosis* Complex lineages. *Plos Genetics*, 14(4) April 2018.
3. Ramadass Balamurugan et al. PCR Amplification of the IS6110 Insertion Element of *Mycobacterium tuberculosis* in Fecal Samples from Patients with Intestinal Tuberculosis. *Journal of Clinical Microbiology*. 2006 May; 44(5): 1884-1886.
4. Paul H.M. Savelkoul et al., Detection of *Mycobacterium tuberculosis* complex with Real Time PCR: Comparison of different prime-probe sets base on the IS6110 element. *Journal of Microbiological Methods*. Volume 66, Issue 1, July 2006, Pages 177-180
5. L.A.S. Lira et al. Evaluation of a IS6110-Taqman real-time PCR assay to detect *Mycobacterium tuberculosis* in sputum samples of patients with pulmonary TB. *Journal of Applied Microbiology* 2012.
6. Francesco Broccolo et al. Rapid Diagnosis of Mycobacterial Infections and Quantitation of *Mycobacterium tuberculosis* Load by two real-time calibrated PCR assays. *Journal of clinical microbiology*. Oct. 2003, p.4565-4572.
7. Santos FCF et al. Performance of the IS6110-TaqMan assay in the diagnosis of extrapulmonary tuberculosis from different biological samples. *Revista da Sociedade Brasileira de Medicina Tropical* 51 (3):331-337, may-june, 2018.
8. Dilia Fontalvo Rivera y Doris Gómez Camargo. Genes del *Mycobacterium tuberculosis* involucrados en la patogenicidad y resistencia a antibióticos durante la tuberculosis pulmonar y extrapulmonar. *Revista Médicas UIS*. 2015;28 (1):39-51.
9. I. Dorronsoro y L. Torroba. Microbiología de la tuberculosis. An. Sist. Navar. 2007 Vol 30, Suplemento 2.

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(5) No lectura en canal Cy5.

(6) Lectura solo en canales FAM y HEX.

| Tabla A.2 TERMOCICLADORES CON BLOQUE DE PERFIL ALTO | |
|---|---|
| Fabricante | Modelo |
| Abbott | Abbott m2000 RealTime System |
| Applied Biosystems | 7300 Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | 7500 Real-Time PCR System |
| Applied Biosystems | 7900 HT Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | ABI PRISM 7000 ⁽⁶⁾ |
| Applied Biosystems | ABI PRISM 7700 ⁽⁵⁾ |
| Applied Biosystems | QuantStudio™ 12K Flex 96-well |
| Applied Biosystems | QuantStudio™ 6 Flex 96-well |
| Applied Biosystems | QuantStudio™ 7 Flex 96-well |
| Applied Biosystems | QuantStudio™ 3 Real-Time PCR System ⁽⁵⁾ |
| Applied Biosystems | QuantStudio™ 5 Real-Time PCR System |
| Applied Biosystems | ViiA™ 7 Real-Time PCR System |
| Analytik Jena Biometra | TOptical |
| Analytik Jena Biometra | qTOWER 2.0 |
| BIONEER | Exicycler™ 96 |
| Bio-Rad | CFX96™ / CFX96™ IVD Real-Time PCR Detection System |
| Bio-Rad | iCycler iQ™ Real-Time PCR Detection System |
| Bio-Rad | iCycler iQ™5 Real-Time PCR Detection System |
| Bio-Rad | MyIQ™ Real-Time PCR Detection System ⁽⁶⁾ |
| Bio-Rad | MyIQ™2 Real-Time PCR Detection System ⁽⁶⁾ |
| Cepheid | SmartCycler® ⁽³⁾ |
| DNA-Technology | DTprime Real-time Detection Thermal Cycler ⁽²⁾ |
| DNA-Technology | DTlite Real-Time PCR System ⁽²⁾ |
| Eppendorf | Mastercycler™ep realplex |
| Qiagen | Rotor-Gene® Q ⁽³⁾ |
| Stratagene / Agilent Technologies | Mx3000PT™ Real Time PCR System |
| Stratagene / Agilent Technologies | Mx3005PT™ Real Time PCR System |
| VIASURE | VIASURE 48 Real Time PCR System ⁽²⁾ |
| VIASURE | VIASURE 96 Real Time PCR System ⁽²⁾ |

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
- CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.
- ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.
- LightCycler® is a registered trademark of Roche.
- Mx3000P™, Mx3005™ and AriaMx are registered trademarks of Agilent Technologies.
- Mastercycler™ is a registered trademark of Eppendorf.
- Rotor-Gene®Q is a registered trademark of Qiagen.
- SmartCycler® is a registered trademark of Cepheid.

Revision: January 2020

CerTest Biotec, S.L.

Pol. Industrial Río Gállego II · Calle J, Nº1
50840, San Mateo de Gállego, Zaragoza (Spain)
www.certest.es

VIASURE online

F-362 rev01

VIASURE

Real Time PCR Detection Kits

CerTest
BIOTEC