

Rotavirus and Adenovirus are most common and major causes of severe gastroenteritis in infants and young children. Pattern also observed in adults. They are transmitted by faecal-oral contact. The main symptoms of viral gastroenteritis are watery diarrhoea and vomiting.

The affected person may also have headache, fever and abdominal cramps ("stomach ache"). In general,


the symptoms begin 1 to 2 days following infection with a virus that causes gastroenteritis and may last for 1 to 10 days, depending on which virus causes the illness (*Rotavirus* 3 days, and *Adenovirus* 5-8 days).

CerTest Rotavirus + Adenovirus one step combo card test offers a simple and highly sensitive screening assay to make a presumptive diagnosis of Rotavirus and/or Adenovirus infection.



Rotavirus + Adenovirus detection kit

Test procedure



Interpretation of the results



ANY OTHER RESULT: INVALID

Clinical Sensitivity and Specificity

	Sensitivity	Specificity	PPV	NPV		Sensitivity	Specificity	PPV	NPV
CerTest <i>Rotavirus+Adenovirus</i> (<i>Rotavirus</i>) vs. Ridascreen® <i>Rotavirus</i> Elisa Test	>99%	>98%	>94%	>99%	CerTest Rotavirus+Adenovirus (Adenovirus) vs. PCR	>99%	>99%	>99%	>99%



CSQ MED ISO 13485