



TECHNICAL SHEET

PRODUCT NAME: Inactivated E. coli O157 antigen (native extract)

Code: MT-28EC7U

Physical State: Liquid

Source: E. coli O157 culture

Description: Inactivated E. coli O157 antigen from culture.

Storage Conditions: Short term storage (<1 week) at 2-8°C. Long term storage at -20°C. Avoid multiple/freeze thaw cycles by storing multiple aliquots at -20°C.

Health & Safety Information: This product contains sodium azide as preservative. Although the amount is lower than 0.1%, this substance is poisonous and hazardous and it should be handled by trained staff only. Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.

Presentation: PBS, <0.1% Sodium Azide (NaN₃)

THIS PRODUCT IS INTENDED FOR *IN VITRO* USE ONLY.

Date: 20/01/2016