



# TECHNICAL SHEET

**PRODUCT NAME:** C. difficile Toxin A recombinant protein (fragment without toxin activity)

**Code:** MT-24TXA, MT-24TXAs1

**Physical State:** Liquid

**Source:** *Escherichia coli*

**Description:** Soluble extract containing recombinant Toxin A inactivated protein

**Storage Conditions:** Storage 2-8°C.

**Health & Safety Information:** This product contains sodium azide as preservative. Although the amount is 0.09%, this substance is poisonous and hazardous and it should be handled by trained staff only.

Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.

**Presentation:** Phosphate buffer with additives and salts, pH 7.5, 0.09% Sodium Azide (NaN<sub>3</sub>).

THIS PRODUCT IS INTENDED FOR *IN VITRO* USE ONLY.

**Date:** 17/08/2018