



# TECHNICAL SHEET

**PRODUCT NAME:** **Campylobacter jejuni recombinant outer membrane protein**

**Code:** MT-25CEP, MT-25CEPs1

**Physical State:** Liquid

**Source:** *Escherichia coli*

**Description:** *Campylobacter jejuni* recombinant outer membrane protein expressed and purified from *Escherichia coli*.

**Appearance:** Clear

**Purity:** >90 %

**Storage Conditions:** Storage at 2-8°C.

**Health & Safety Information:** This product contains sodium azide as preservative. Although the amount is 0.09%, this substance is poisonous and hazardous and it should be handled by trained staff only. Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.

**Presentation:** TRIS Buffer with additives and salts, pH 8, 0.09% Sodium Azide (NaN<sub>3</sub>)

THIS PRODUCT IS INTENDED FOR *IN VITRO* USE ONLY.

**Date:** 17/08/2018