



# TECHNICAL SHEET

**PRODUCT NAME:** Inactivated *Campylobacter jejuni* antigen (native extract)

**Code:** MT-28CJU, MT-28CJU51

**Physical State:** Liquid

**Source:** *Campylobacter jejuni* culture

**Description:** Inactivated *Campylobacter jejuni* antigen from culture.

**Storage Conditions:** Short term storage at 2-8°C. Long term storage at -20°C. Avoid multiple/freeze thaw cycles by storing multiple aliquots at -20°C.

**Health & Safety Information:** This product contains sodium azide as preservative. Although the amount 0.09%, this substance is poisonous and hazardous and it should be handled by trained staff only. Good Laboratory Practices should be followed when handling this material. The end user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with the respective national regulations.

**Presentation:** Phosphate buffer with salts, pH 7.4, 0.09% Sodium Azide (NaN<sub>3</sub>).

THIS PRODUCT IS INTENDED FOR *IN VITRO* USE ONLY.

**Date:** 17/08/2018